



RYPO

Thinking about your safety

INDEX

1. THREE-LAYER HYGIENIC MASK (3PLY)
2. SURGICAL MASK II-R
3. FFP2 MASKS
4. FFP3 MASKS
5. FFP2 MASK FOR CHILDREN
6. INFRARED THERMOMETER CONTACTLESS
7. DISPOSABLE PROTECTIVE SUIT
8. MEDICAL GOWN
9. DISPOSABLE HATS
10. SHOE COVERS
11. HIGH RISK BOOT COVERS
12. TEST DYNAMIKER
13. BIOCREDIT COVID-19 Rapid Diagnostic Test
14. NEWGENE COVID-19 Rapid Diagnostic Test
15. NITRILE GLOVES
16. VINYL GLOVES
17. OXYGEN THERAPY ACCELERATOR
18. HYDROALCOHOLIC SPRAY
19. HYDROALCOHOLIC GEL VITALCURE
20. HYDROALCOHOLIC GEL PRADY
21. FACE SHIELD
22. SYRINGES
23. NEEDLES

THREE-LAYER HYGIENIC MASK (3PLY) Bag of 10 or 50 Units.

Structure: It consists of a central body with a nose clip and a mask made of a fine folded interwoven fabric.

Particles blocked: Bacteria, dust, sneezing, excretion, etc.

> 95%

Material: Material grade BFE99 non-allergic filter cloth.

Instructions for use: Place the elastic band over the ears and adjust the coverage of the bridge of the nose comfortably by pulling the mask until it covers the chin. Discard after use.

Test Report: Aitex.



SURGICAL MASK II-R Bag of 50 or 10 Units.

Structure: Consists of a central body with a nose clip and a mask made of a fine folded interwoven fabric.

Particles blocked: Bacteria, dust, sneezing, expectorations, etc.

> 98%

Material: Material grade BFE99 non-allergic filter cloth.

Instructions for use: Place the elastic band over the ears and adjust the coverage of the bridge of the nose comfortably by pulling the mask to cover the chin. Discard after use.

Test Report: Aitex.



FFP2 MASK.

Bag of 10 Units and Boxes of 20 Individual Units

Class FFP2 masks offer maximum protection against air pollution. Total leakage can be up to 5% and 99% of all particles down to 0.6 microns are filtered out. This type of mask also filters out poisonous, oncogenic and radioactive particles.



FFP3 MASKS 10 Individual Units Box

Class FFP3 protection masks offer maximum protection against air pollution. Total leakage can be up to 5% and 99% of all particles down to 0.6 microns are filtered out. This type of mask also filters out poisonous, oncogenic and radioactive particles.
 Instructions for use: Place the elastic band over the ears and adjust the coverage of the bridge of the nose comfortably by pulling the mask until it covers the chin.



FFP2 MASK FOR CHILDREN Box of 10 Individual Units

Class FFP2 protective masks offer maximum protection against air pollution. Total leakage can be up to 5% and 99% of all particles down to 0.6 microns are filtered out. This type of mask also filters out poisonous, oncogenic and radioactive particles.



INFRARED THERMOMETER CONTACTLESS 1 Unit Box

High precision and reliability non-contact infrared digital thermometer. It takes only 1 second to display the data of the temperature obtained.

Its dual mode allows the measurement of both body temperature and that of different surfaces, such as bath water, food (baby's milk), etc.



FHT-1

Infrarrojo sin contacto

Termómetro

CATEGORY 3 AND 4 DISPOSABLE PROTECTIVE SUIT

Diver or Tyvek type coverall, with impervious protection to solid particles and splashes of chemical products of limited intensity.

This diver made of SMS material is resistant. Its comfort makes it a pleasant companion in many areas of work.

Adjustable hood, elasticated waist, sleeves and legs and zipper flap.



MEDICAL GOWN

Made of lightweight and breathable polypropylene fabric, to tie at the waist, with trikot cuffs for a secure fit on the arms, suitable for activities in the medical, care and hygiene sector.

20 gr 25 gr 30 gr

Sterile and non-sterile



DISPOSABLE HATS

Made from lightweight and breathable spun polypropylene

Latex free

With corrugated elastic Size M (Ø approx. 52 cm)

White, green, blue, yellow, red and orange color

Sterile and non-sterile



SHOE COVERS

Made of chlorinated polyethylene
With latex elastic around the opening perimeter Hammered
15 x 42

They are used in clean spaces, in cleaning surfaces as temporary protection; in contaminated environments of all kinds (for example in construction sites, agriculture, gardening etc), when entering and stepping on hygienic areas and laboratories.

PPE Category I according to standard 89/686 / CEE and regulation EU 2016/425



HIGH RISK ANTI-SLIP BOOT COVERS TYPE 5B / 6B CAT. III

Knee-high boot with non-slip sole.
Stitched and taped seams for greater strength and waterproofing.
Protection against a range of commonly used chemicals.

EN ISO 13688:2013
EN ISO 14605:2005+A1:2009
EN ISO 14126:2003/AC:2004




TEST DYNAMIKER Antigens / Antibodies / PCR



BIOCREDIT COVID-19 Rapid Diagnostic Test

Biocredit Antigen Detection Kit Allows early detection of infected patients without symptoms or who have not yet developed the disease.



NEWGENE NEWGENE COVID-19 ANTIGEN DETECTION KIT

Newgene Antigen Detection Kit.

It allows the early detection of infected patients without symptoms or who have not yet developed the disease.

1. Use saliva as a sample
2. Rapid detection, results in 15 minutes
3. High precision.
4. Easy to use.





NITRILE GLOVES 100/200 Units per Box

Nitrile gloves are the perfect alternative to conventional latex gloves. They do not contain latex, they reliably protect against external influences, respecting more the skin and achieving fewer reactions on the user's skin. These types of gloves are more durable than latex gloves.



In addition, thanks to the rolled edge that adheres securely, nitrile gloves can be quickly put on and taken off, fitting the hand perfectly.

Nitrile is a synthetic material recommended for those who have allergies or have latex sensitive skin. Nitrile gloves, being more resistant than latex, are recommended for those who need to wear them for hours.

They are the most used in hospitals, pharmacies and laboratories.



NITRILE GLOVES 100/200 Units per Box



EN ISO 374-5: 2016





93/42/EEC 2007/47/EEC

VINYL GLOVES 100/200 Units. Per Box

Pack of 100/200 units of transparent powder-free vinyl gloves. They offer great sensitivity to the touch, with a smooth surface on the fist, palm and fingers.

Ideal for use in the healthcare field as they have EN-455 certification and 374-5 as well as in industries such as food, electronics or cleaning

In addition, thanks to the rolled edge that sticks securely, the vinyl gloves can be quickly put on and taken off, fitting the hand perfectly.



OXYGEN THERAPY ACCELERATOR

DEVICE FOR CONTINUOUS PRESSURE IN THE RESPIRATORY TRACT

Portable and field system for automatic and autonomous airway care, ventilator with coupling mouth to mask or intubated patient, with oxygen inlet. This device is autonomous. DOES NOT NEED ELECTRICITY. Provides respiratory support through an endotracheal tube or mask using a continuous flow source of gas.

It works with at least two 1 + 1 O₂ cylinders and compressed air.

It works and responds equally with an O₂ cylinder (under medical supervision). Saturation of 82% to 97% of O₂ to the patient in a very short time from its application. Composition of components customized to criteria, resources, components, etc. The medical results are the same as any conventional respirator. Helps to optimize the "Helmet" system. It has an auxiliary valve for administration of treatments



HYDROALCOHOLIC SPRAY 500ml / 1L / 5L

Hydroalcoholic gel for direct application to the hands without the need for water or rinsing. Antiseptic product for personnel in all areas, including healthcare and the food industry. Bactericidal, fungicidal, virucidal, mycobactericidal.

How to use: Dose a bit of gel in the palm of your hand and rub until completely evaporated.

Do not rinse.

Meets the standards: UNE-EN 1500, UNE-EN 12791, UNE-pr EN 12054, UNE-EN 1650, UNE-EN 14476, UNE-EN 14348.



HYDROALCOHOLIC GEL VITALCURE 200 mL

Gel designed for the exhaustive and deep cleaning of hands with aloe vera and hyaluronic acid for maximum skin care and comfort. Alcohol content = 70% v / v.

How to use: Dose a bit of gel in the palm of your hand and rub until completely evaporated. Do not rinse.

Precautions: Highly flammable. irritates the eyes. Inhalation of vapors may cause drowsiness and dizziness. Keep out of reach of children. Keep away from any flame or source of sparks. In case of contact with the eyes, rinse immediately and abundantly with water and consult your doctor. Do not eat.



HYDROALCOHOLIC GEL PRADY 100ml. / 500ml / 5L

Gel designed for the exhaustive and deep cleaning of hands with aloe vera and hyaluronic acid for maximum skin care and comfort. Alcohol content = 70% v / v.

How to use: Dose a bit of gel in the palm of your hand and rub until completely evaporated. Do not rinse.

Precautions: Highly flammable. irritates the eyes. Inhalation of vapors may cause drowsiness and dizziness. Keep out of reach of children. Keep away from any flame or source of sparks. In case of contact with the eyes, rinse immediately and abundantly with water and consult your doctor. Do not eat.



FACE SHIELD

Full face shield: it is transparent, lightweight, effective in protecting the eyes, mouth and nose from debris, drops, aerosols, aerosols and splashes.

Breathable and comfortable: each side is equipped with an elastic band and a sponge headband, ideal for long-time wear.



SYRINGES

All the following products have passed the TÜV Rheinland LGA Products GmbH Certificates

1. EN ISO 13485:2016
2. EC Certificate (Directive 93/42/EEC Annex V Production Quality Assurance Medical Devices)

Products:

- Infusion Sets for Single Use
- Scalp Vein Sets for Single Use
- Transfusion Sets for Single Use
- Sterile Hypodermic Syringes for Single Use
- Sterile Hypodermic Needles for Single Use
- Infusion Sets with Precision Filters for Single Use
- Irrigation Syringes
- Disposable Transfer Sets
- Disposable Insulin Syringes



Two part Syringe



Syringe without needle



Syringe with Needle (Luer Slip)



Syringe with Needle (Luer Lock)



Syringe with Precision Filter



Irrigation/Feeding Syringe

SUBCUTANEOUS HYPODERMIC NEEDLE WITH SECURITY SYSTEM LENGTH (25mm Ø: 0.5)

Certificates of TÜV Rheinland LGA Products GmbH:

EN ISO 13485: 2016

The hypodermic needle is used to inject substances or medications into the body.

It is made up of a hollow needle normally used with a two or three body syringe.

