

**LIVE
DEMO!**

Accurate Liquid Distribution

Demonstrating 3 different flow patterns




DP-11

Switching Delivery Pump

A Variety of Flow Patterns!


1-in 5-out or free port selection by programming

Accurate Flow Volume!



Perfect distribution

Dispensing
1 solvent into 5 bottles
by only one pump!




Setup flow rate : 5 mL/min
Setup flow volume : 300 μ L

1

Consistent liquid volume

Look at the exact same
volume in each bottle!




Setup flow rate : 25 mL/min
Setup flow volume : 60 μ L x 5 times

2

A wide flow rate range

You can control even
a subtle amount of liquid!



Setup flow rate : 2.5 mL/min
Setup flow volume : Equal to a droplet

3