



23

Restorative

Restorative

Acid	p.411	Fluoride Releasing Composite	p.419 - p.421
Activator & Primer	p.411	Glass Ionomers	p.421 - p.424
Bonding System	p.412 - p.414	GI Conditioner & Primer	p.424
Universal Composite	p.415 - p.416	Base & Liner	p.424 - p.425
Posterior Composites	p.417	Amalgam	p.425
Bulk Fill Posterior Composite	p.417	Opaquers	p.426
Bulk Fill Flowable Composite	p.417	Temporary Fillings	p.426
Flowable Composite	p.418	Accessories	p.426 - p.427
Direct Composite Veneer	p.419	Others	p.427



Restorative

Acid

Scotchbond™ Etchant Gel / 3M ESPE

✓ *Scotchbond™ Etchant Is A 35% Phosphoric Etchant Which Is Used To Etch Enamel And/ Or Dentin As Indicated In The Instructions Provided With The Various Dental Adhesive Systems*

Low Viscosity



■ Bottle Refill

Volume	Item no.
9 mL/bottle	3M-3007

Scotchbond™ Universal Etchant - Etchant Gel / 3M ESPE



Advantages

- 32% phosphoric acid with a pH of approximately 0.1
- Low-viscosity, gel-like consistency
- Gel has a distinct blue color to assist in complete rinsing away of the etchant

Indications

- Acid etching of enamel and dentin

Content	Item no.
- 2 x Etchant Syringes (3 mL) - 50 x Dispensing Tips	3M-41263

Ultra-Etch® / Ultradent

- Self-limiting
- Penetrates smallest fissures but won't migrate on a vertical surface
- Increases bond strength and reduces sensitivity
- Rinses cleanly — leaves no residue
- Places precisely



■ Ultra-Etch Kit

Content	Item no.
- 4 x Syringes (1.2 mL) - 20 x Blue Micro Tips	UL-163

■ Ultra-Etch Econo Kit

Content	Item no.
- 20 x Syringes (1.2 mL) - 40 x Blue Micro Tips	UL-167

■ Ultra-Etch Refill

Content	Item no.
- 4 x Syringes (1.2 mL)	UL-164

■ Ultra-Etch Econo Refill

Content	Item no.
- 20 x Syringes (1.2 mL)	UL-168



■ Ultra-Etch IndiSpense Syringe

Content	Item no.
- 1 x Syringe (30 mL)	UL-685

■ Ultra-Etch IndiSpense Kit

Content	Item no.
- 1 x Syringe (30 mL) - 20 x Blue Micro Tips - 20 x Empty Syringes (1.2 mL)	UL-383

Activator & Primer

Scotchbond™ Universal DCA Dual Cure Activator / 3M ESPE

- Consists of sodium toluene sulfinate and ethanol. When the activator is mixed with Scotchbond Universal adhesive, it allows for bonding to self- or dual-cure composite, core build-up and cement materials that rely on the self-cure mechanism for polymerizing the material



Volume	Item no.
5 mL/bottle	3M-41264

Adper™ Scotchbond™ Multi-Purpose Primer / 3M ESPE

- Adper™ Scotchbond™ Multi-Purpose Primer is one part of our light cure dental adhesive system that is designed to handle difficult bonding situations, is fast and easy to apply and offers excellent bond strengths to a variety of surfaces



Volume	Item no.
8 mL/bottle	3M-3008

RelyX™ Ceramic Primer / 3M ESPE

- A prehydrolyzed silane-based primer designed to enhance bonds to porcelain, ceramic, pre-cured composite and metal surfaces



Volume	Item no.
5 mL/bottle	3M-2721

Ketac™ N100 Light-Curing Nano-Ionomer Restorative Primer / 3M ESPE



Volume	Item no.
6.5 mL/bottle	3M-3527P



Shofu AZ Primer / Shofu

✓ *For Adhesive Bonding To alumina And Zirconia Porcelains*

- The AZ-Primer has been exclusively developed for modern materials such as alumina and zirconia and guarantees outstanding bonding strengths between the restoration, the luting composite and the tooth structure

Volume	Item no.
5 mL/bottle	SF-3232



Bonding Systems

Total-Etch

Adper™ Single Bond 2 Adhesive / 3M ESPE

✓ A Fast, Easy And Convenient Total-Etch, Single-Component Bonding Agent Offering Exceptional Bond Strength

Advantages

- Excellent shear bond strength provides exceptional bonding performance
- Nanofiller is stable and will not settle out of dispersion
- Unique, easy-dispensing bottle permits visual inspection of exactly how much adhesive remains in the bottle
- Cap snaps securely shut minimising evaporation
- No shaking required - use straight out of the vial
- Ethanol/water-based adhesive maintains its shear bond strength over time to yield an outstanding shelf life

Indications

- All direct, light-cured bonding procedures
- Indirect bonding procedures when combined with 3M ESPE RelyX™ ARC Adhesive Resin Cement
- Root surface desensitisation
- Porcelain and composite repair
- Bonding veneers in combination with RelyX™ Veneer Cement from 3M ESPE



■ Adper Single Bond 2 Vial Refill

Content	Item no.
6 mL/bottle	3M-51202

Self-Etch

Adper™ Prompt™ Self-Etch Adhesive / 3M ESPE

✓ Offers Easy Activation And High Performance
 ✓ The Unique L-Pop™ Delivery System Is Fast And Easy To Use, Allowing You To Tech, Prime And Bond In A Single Step

Advantages

- Bonds to dentin, cut and uncut enamel without a separate etch step
- Clinically proven
- Very low incidence of post-operative sensitivity
- Desensitizes hypersensitive root surfaces
- Simplifies sealant technique
- Cures with any light source – halogen, laser, plasma or LED
- Ensures hygienic delivery of the adhesive material
- Extra small applicator works well with very small preparations and in hard-to-reach areas of the mouth

Indications

- Bonding between dentin/enamel and direct, light-cured composite and compomer filling materials
- Bonding light-cured pit and fissure sealants
- Root surface desensitization



■ Adhesive Unit-Dose Applicators

Packing	Item no.
40's/pack	3M-41925

Self-Etch

Adper™ Easy Bond Self-Etch Adhesive / 3M ESPE

✓ A One-Step, One-Coat Self-Etch Adhesive That Provides A Fast And Easy Application, Along With Strong Bonds

Advantages

- One-Bottle, one-coat convenience
- High bond strength with proven chemistry and filler
- Reduces risk of post-operative sensitivity
- High marginal integrity and low microleak age
- Moisture tolerant - works well on wet or dry enamel and dentin
- Desensitizes hypersensitive root surface
- Uniquely suited to the selective etch technique

Indications

- All classes of fillings with light-curing composite or compomer
- Core build-ups made of light-curing composite
- Repair of composite or compomer fillings
- Root surface desensitization
- Intraoral repair of existing composite, porcelain-fused-to metal, and all-ceramic restorations
- Cementation of indirect restorations made of composite or compomer, ceramic and metal (after light curing and when combined with RelyX™ ARC Adhesive Resin Cement)



■ Adper Easy Bond Self-Etch Adhesive Intro Kit

Content	Item no.
- 1 x Adhesive Vial (5 mL) - 50 x Disposable Mixing Wells - 50 x Disposable Microbrushes	3M-41245

■ Adper Easy Bond Self-Etch Adhesive Vial Refill

Volume	Item no.
5 mL/bottle	3M-41246



Restorative

Total-Etch

Adper™ Scotchbond™ Multi-Purpose Adhesive / 3M ESPE

✓ *A Light Cure Dental Adhesive System That Is Designed To Handle Difficult Bonding Situations, Is Fast And Easy To Apply And Offers Excellent Bond Strengths To A Variety Of Surfaces*



Advantages

- Versatile adhesive system designed to handle difficult bonding situations such as bonding to sclerotic dentin and bonding in moist environments
- Reduces microleakage that can contribute to post-operative sensitivity for increased patient comfort and doctor confidence
- Offers excellent bond strengths to a variety of surfaces for reliable bonding performance
- The system offers complete dental bonding with dual-cure capabilities when used with Adper Scotchbond Multi-Purpose Plus Activator and Adper Scotchbond Multi-Purpose Plus Catalyst
- Easy-to-follow technique guides and color-coded, numbered vials make the system efficient and simple to use

Indications

- Light-cured composite restorations
- Pocerlain and composite repairs
- Bonding of veneers in combination with RelyX™ Veneer Cement from 3M ESPE

Intro Kit

Content	Item no.
- 2 x Phosphoric Etchant Syringes (3 mL) - 1 x Primer Vial (8 mL) - 1 x Adhesive Vial (8 mL) - 3 x Well Mixing Well - Brush Tips/Color-Coded Handles	3M-7540S

Refills

Item	Packing	Item no.
Primer	8 mL/bottle	3M-3008
Adhesive		3M-3009



Scotchbond™ Universal Adhesive / 3M ESPE

- ✓ *This Single-Bottle Adhesive Solution Offers High Bond Strength With One Simple Application Technique And A Quick 35-Second Adhesive Application Time*
- ✓ *Can Be Used In All Application Techniques, Including Total-Etch, Self-Etch And Selective-Etch*
- ✓ *Ideal For Both Direct And Indirect Applications With Virtually No Post-Operative Sensitivity And Can be Used On All Surfaces Without Any Extra Primer*

Advantages

- Combined Total-Etch, Self-Etch and Selective-Etch adhesive
- Consistent bond strength to both moist and dry etched dentin
- Virtually no post-op sensitivity
- Effective on all surfaces
- High bond strength to all indirect surfaces (zirconia, alumina, glass ceramics and metals) without a separate primer
- Excellent marginal integrity for highly esthetic restorations
- One-step, one-coat adhesive application
- Self- and dual-cure compatible with (Scotchbond™ Universal DCA Dual Cure Activator)
- No refrigeration needed
- Store at room temperature for two years
- Flip-top vial and unit-dose delivery
- Improved, easy-to-place etchant

Indications

Direct Placement Indications:

- Bonding light-cured composite or compomer for all classes of direct restorations
- Root surface desensitization
- Sealing of dentin prior to amalgam restorations
- Protective varnish for glass ionomer restorative materials
- Repair of composite or compomer restorations
- Bonding sealants

Indirect Placement Indications:

- Bonding veneers with RelyX™ Veneer Cement
- Adhesive/primer for chemical adhesion when bonding zirconia, alumina, metal or glass ceramic restorations
- Intraoral repair of existing indirect restorations
- Bonding self- or dual-cure core build-up materials and resin cements (with Scotchbond™ Universal DCA Dual Cure Activator)
- Sealing of dentin prior to temporization for indirect restoration placement

Intro Kit

Content	Item no.
- 1 x Adhesive Vial (5 mL) - 50 x Disposable Mixing Wells - 50 x Disposable Microbrushes - 1 x Etchant Syringe (3 mL)	3M-41254

Intro Kit - Unit Dose 40

Content	Item no.
- 40 x Unit Dose - 1 x Etchant Syringe (3 mL)	3M-41255

Refill - Adhesive Vial

Content	Item no.
- 1 x Adhesive Vial (5 mL)	3M-41258



Bonding Systems

Peak Universal Bond / Ultradent

- ✓ Peak Universal Bond's Versatile Formulation Is Ideal For Indirect And Direct Bonding, Including Post And Core Procedures
- ✓ With 7.5% Filler, Its Viscosity Has Been Optimized For Minimal Film Thickness (2µm) And Superior Strength; The Addition Of Chlorhexidine (0.2%) In Peak Universal Bond May Provide The Added Benefit Of Improving Long-Term Bond Strengths



Features

- Universal adhesive can be used in all bonding procedures, direct and indirect
- Added chlorhexidine can ensure long-lasting bonds
- Proven superior bond strengths
- Available in syringe or bottle delivery

Total-Etch Bottle Kit

Content	Item no.
- 1 x Peak Universal Bond bottle (4 mL) - 4 x Ultra-Etch syringes (1.2 mL) - 40 x Blue Micro tips - 50 x Mixing wells - 50 x MicroApplicator brushes	UL-4542

Self-Etch Bottle Kit

Content	Item no.
- 1 x Peak Universal Bond bottle (4 mL) - 4 x Peak SE Primer syringes (1 mL) - 40 x Black Mini Brush tips - 50 x Mixing Wells - 50 x MicroApplicator brushes	UL-4541

4 mL Bottle Refill

Content	Item no.
- 1 x Peak Universal Bond Bottle (4 mL)	UL-4543

Anaesthetics
Articulating
Bur and Abrasives
Cosmetics
Diagnostic
Disposables
Endodontics
Equipment
Handpieces
Implant
Infection Control
Instruments
Isolation
Laboratory
Oral Surgery
Orthodontics
Patient Education and Gift
Pedodontics
Periodontal
Post and Core
Preventive
Prosthodontics
Restorative
Uniform



BeutiBond® / Shofu

- ✓ New Seventh-Generation Bonding Agent Has An Exclusive Chemistry With Unique Dual Adhesive Monomers That Deliver Equal Bond Strength To Enamel And Dentin

Advantages & Features

- Bond strengths comparable to leading multi-component brands
- Single coat for shorter working time
- Efficient 30 second application
- 5µm film thickness
- HEMA-free
- Convenient unit dose, no mixing



BeutiBond Unit Dose

Content	Item no.
- 50 x Unit Doses (0.1 mL) - 50 x Microbrushes - Instructions	SF-1782



BeutiBond Bottle

Volume	Item no.
6 mL/bottle	SF-1781



LuxaBond / DMG

- ✓ Total-Etch Dual-Curing Adhesive System



Advantages & Features

- Self-curing: Hardens even without light
- No premature polymerization of bonding pools
- Compatible with self-curing and dual-curing materials: Suitable for all applications with LuxaCore
- Pre-Bond facilitates quicker curing of the bond mix in the dentine for superior retention
- Easy root canal application with special Endobrushes
- Does not contain acetone



LuxaBond Intro Kit

Content	Item no.
- 1 x Pre-Bond Bottle (5 mL) - 1 x Bond A Bottle (5 mL) - 1 x Bond B Bottle (5 mL) - 1 x Etching Gel-Medium Viscosity Syringe (2 mL) - 1 x Mixing Pad	-25 x Endobrushes (Red) - 25 x Endobrushes (Black) - 3 x Luer-Lock Tips DMG-212384



Restorative

Universal Composites



HENRY SCHEIN® Natural Elegance NH Composite

/ Henry Schein

Made In Switzerland

- Features**
- Contouring natural morphology easily
 - Non-sticky, condensable & sculptable
 - Excellent marginal adaption and contact point
 - Superior blend-in optical properties
 - High wear resistance
 - Less polishing time with high gloss polish
 - Duo shade – save time & money
 - Low polymerization shrinkage



■ Syringe Intro Kit

Content	Item no.
- 1 x Syringe (4 g) <i>(1 each of shade A1/D2, A2, B2, A3/D3, A3.5/B3)</i>	HSPL-900-3856
- Shade Guide	

■ Syringe Refills

Shade	Packing	Item no. (HSPL-)
A1/D2	1 x Syringe (4 g)	900-3538
A2/B2		900-3539
A3/D3		900-3540
A3.5/B3		900-3541
Bleach		900-3543
Transparent		900-3544



■ Capsule Refills

Shade	Item no. (HSPL-)	Shade	Item no. (HSPL-)
A1/D2	900-3545	A3.5/B3	900-3548
A2/B2	900-3546	A4/M5	900-3549
A3/D3	900-3547	C2/C3	900-3811

Filtek™ Z350 XT Universal Restorative

/ 3M ESPE



Advantages

- One simple shading system with endless possibilities
- Wide range of shades and opacities
- Enhanced fluorescence
- Excellent aesthetic results, whether you use a single-dual or multi-shade technique

Indications

- Direct anterior and posterior restorations (Class I, II, III, IV and V)
- Minimally Invasive Dentistry (MID)
- Sandwich technique with glass ionomer resin material
- Indirect anterior and posterior restorations including inlays, onlays and veneers
- Cusp buildup
- Core buildup
- Splinting

■ Filtek Z350 XT Universal Restorative Refills

- Body Shades (Capsules)

Content: 20 x Capsules (0.2 g)

Shade	Item no.	Shade	Item no.
A1B	3M-7019-A1B	B1B	3M-7019-B1B
A2B	3M-7019-A2B	B2B	3M-7019-B2B
A3B	3M-7019-A3B	C3B	3M-7019-C2B
A3.5B	3M-7019-A3.5B	C3B	3M-7019-C3B
A4B	3M-7019-A4B		

■ Filtek Z350 XT Universal Restorative Refills

- Body Shades (Syringes)

Content: 1 x Syringes (4 g)

Shade	Item no.	Shade	Item no.
A1B	3M-7018-A1B	B5B	3M-7018-B5B
A2B	3M-7018-A2B	C1B	3M-7018-C1B
A3B	3M-7018-A3B	C2B	3M-7018-C2B
A3.5B	3M-7018-A3.5B	C3B	3M-7018-C3B
A4B	3M-7018-A4B	D2B	3M-7018-D2B
A6B	3M-7018-A6B	D3B	3M-7018-D3B
B1B	3M-7018-B1B	W	3M-7018-WB
B2B	3M-7018-B2B	XW	3M-7018-XWB
B3B	3M-7018-B3B		

- Dentin Shades

Content: 1 x Syringe (4 g)

Shade	Item no.	Shade	Item no.
A1D	3M-7018-A1D	B3D	3M-7018-B3D
A2D	3M-7018-A2D	C4D	3M-7018-C4D
A3D	3M-7018-A3D	White	3M-7018-WD
A4D	3M-7018-A4D		

- Enamel Shades

Content: 1 x Syringe (4 g)

Shade	Item no.	Shade	Item no.
A1E	3M-7018-A1E	D2E	3M-7018-D2E
A2E	3M-7018-A2E	White	3M-7018-WE
A3E	3M-7018-A3E	Extra Withe	3M-7018-XWE
B1E	3M-7018-B1E		
B2E	3M-7018-B2E		

- Translucent Shades

Content: 1 x Syringe (4 g)

Shade	Item no.	Shade	Item no.
Clear T	3M-7018-CT	Grey T	3M-7018-GT
Blue T	3M-7018-BT	Amber T	3M-7018-AT



Universal Composites

Filtek™ Z250 Universal Restorative

/ 3M ESPE

✓ *Speed And Convenience In A Universal Dental Composite Resin That Provides Exceptional Handling, Long-Lasting Esthetic Results*



Advantages

- Low polymerization shrinkage for less stress on the tooth and less potential for patient sensitivity
- Speed and convenience: easy to use, fast placement technique, and easy to finish to a smooth surface
- Excellent handling with exceptional resistance to slumping
- Strong, long-lasting and natural-looking results
- Outstanding versatility for any area of the mouth

Indications

- Direct anterior and posterior restorations
- Indirect restorations including inlays, onlays and veneers
- Core buildups
- Splinting

Filtek Z250 Syringes

- Introductory Syringe Kit

Content	Item no.
- 3 x Syringes (4 g) (1 each of shade A2, A3, A3.5) - 1 x Scotchbond™ Etchant (4 mL) - 1 x Single Bond 2 (6 g) - Accessories	3M-1370

- Syringe Refills

Content: 1 x Syringe (4 g)

Shade	Item no.	Shade	Item no.
A1	3M-1370-A1	B3	3M-1370-B3
A2	3M-1370-A2	C2	3M-1370-C2
A3	3M-1370-A3	D3	3M-1370-D3
A3.5	3M-1370-A3.5	Universal Dentine	3M-1370-UD
A4	3M-1370-A4	Incisal	3M-1370-I
B1	3M-1370-B1		
B2	3M-1370-B2		

Filtek Z250 Capsules

- Capsule Refills

Content: 20 x Capsules (0.2 g)

Shade	Item no.	Shade	Item no.
A1	3M-6021-A1	B3	3M-6021-B3
A2	3M-6021-A2	C2	3M-6021-C2
A3	3M-6021-A3	D3	3M-6021-D3
A3.5	3M-6021-A3.5	Universal Dentine	3M-6021-UD
A4	3M-6021-A4	Incisal	3M-6021-I
B1	3M-6021-B1		
B2	3M-6021-B2		

Filtek™ Z250 XT Nano Hybrid Universal Restorative

/ 3M ESPE

✓ *Filtek Z250XT Nano Hybrid Restorative Is Formulated With A Broad Range Of Particle Sizes To Allow A Higher Filler Loading. The Result Is High Flexural Strength, Diametral Tensile Strength And Wear Resistance*

Advantages

- Excellent handling
- High strength
- Good wear resistance
- 12 popular shades
- Predictable esthetics
- Indicated for all direct restorations
- Easy to polish

Syringe Refills

Content: 1 x Syringe (3 g)

Shade	Item no.	Shade	Item no.
A1	3M-1470-A1	B2	3M-1470-B2
A2	3M-1470-A2	B3	3M-1470-B3
A3	3M-1470-A3	C2	3M-1470-C2
A3.5	3M-1470-A3.5	D3	3M-1470-D3
A4	3M-1470-A4	OA2	3M-1470-OA2
B1	3M-1470-B1	OA3	3M-1470-OA3



Capsule Refills

Content: 20 x Capsules (0.2 g)

Shade	Item no.	Shade	Item no.
A1	3M-1471-A1	B2	3M-1471-B2
A2	3M-1471-A2	B3	3M-1471-B3
A3	3M-1471-A3	C2	3M-1471-C2
A3.5	3M-1471-A3.5	D3	3M-1471-D3
A4	3M-1471-A4	OA2	3M-1471-OA2
B1	3M-1471-B1	OA3	3M-1471-OA3





Restorative

Posterior Composites

Filtek™ P60 Posterior Restorative Syringe

/ 3M ESPE

✓ *With A Rapid Cure Time, Easy Handling And Fast Placement, You Can Rely On Filtek™ P60 Posterior Restorative For Excellent Results*

- Advantages**
- Superior packability--packs well without sticking to instruments
 - Easy sculpting with exceptional resistance to slumping
 - Excellent marginal adaptation and interproximal contact
 - Fast placement technique takes only 20 seconds to cure an increment of 2.5mm
 - Low polymerization shrinkage means less stress on the tooth and less potential for patient sensitivity
 - Excellent wear resistance and strength means you can be confident of long-lasting posterior restorations

- Indications**
- Direct anterior and posterior restorations
 - Indirect restorations including inlays, onlays and veneers
 - Core buildups
 - Splinting



■ Filtek P60 Posterior Restorative Syringe Refills

Item	Packing	Item no.
A3	1 x Syringe (4 g)	3M-8100-A3
B2		3M-8100-B2
C2		3M-8100-C2

Bulk Fill Posterior Composites

Filtek™ Bulk Fill Posterior Composites

/ 3M ESPE

✓ *A Visible, Light-Activated Restorative Composite Optimized To Create Posterior Restorations Simpler And Faster*

- Advantages**
- The bulk fill material provides excellent strength and low wear for durability
 - The shades are semi-translucent and low-stress curing, enabling up to a 5 mm depth of cure
 - With excellent polish retention, Filtek Bulk Fill Posterior Restorative is also suitable for anterior restorations that call for a semi-translucent shade
 - All shades are radiopaque

- Indications**
- Direct anterior and posterior restorations (including occlusal surfaces)
 - Base/liner under direct restorations
 - Core build-ups
 - Splinting
 - Indirect restorations including inlays, onlays and veneers
 - Restorations of deciduous teeth
 - Extended fissure sealing in molars and premolars
 - Repair of defects in porcelain restorations, enamel and temporaries



■ Syringe Refills

Content: 1 x Syringes (3 g)

Shade	Item no.
A1	3M-4863-A1
A2	3M-4863-A2
A3	3M-4863-A3
B1	3M-4863-B1
C2	3M-4863-C2

■ Capsule Refills

Content: 20 x Capsules (0.2 g)

Shade	Item no.
A1	3M-4864-A1
A2	3M-4864-A2
A3	3M-4864-A3
B1	3M-4864-B1
C2	3M-4864-C2

Bulk Fill Flowable Composites

Filtek™ Bulk Fill Flowable Composites / 3M ESPE

✓ *Offers A Faster, Easier Solution*
 ✓ *With Less Shrinkage Stress Than Many Traditional Flowable And Universal Composites, It Can Be Placed In One 4 mm Increment*

Advantages

- Up to 50% stronger than the leading bulk fill flowable
- Nearly twice the wear resistance of leading bulk fill flowable
- Choice of deliveries: easy-to-use capsules or syringes

Indications

- Base under Class I and II direct restorations
- Liner under direct restorative materials
- Pit and fissure sealant
- Restoration of minimally invasive cavity preparations (including small, non stress-bearing occlusal restorations)
- Class III and V restorations
- Undercut blackout
- Repair of small enamel defects
- Repair of small defects in esthetic indirect restorations
- Repair of resin and acrylic temporary materials
- As a core build-up where at least half the coronal tooth structure is remaining to provide structural support for the crown



■ Syringe Refills

Content:
 - 2 x Syringes (2 g)
 - 20 x Dispensing Tips

Shade	Item no.
A1	3M-4862-A1
A2	3M-4862-A2
A3	3M-4862-A3
Universal	3M-4862-U

■ Capsule Refills

Content:
 - 15 x Capsules (0.2 g)

Shade	Item no.
A1	3M-4861-A1
A2	3M-4861-A2
A3	3M-4861-A3
Universal	3M-4861-U



Flowable Composites

Filtek™ Z350 XT Flowable Restorative

/ 3M ESPE

✓ *Filtek™ Z350 XT Flowable Restorative Offers Flow-On-Demand Handling From Either Capsules Or Syringes*

- Advantages**
- Excellent aesthetics
 - Unsurpassed polish retention
 - Low wear

Indications

- Base/liner under direct restorations
- Undercut blockout
- Class III and V restorations
- Restoration of minimally invasive cavity preparations (incl. small, nonstress-bearing occlusal restorations)
- Repair of small defects in aesthetic indirect restorations
- Pit and fissure sealant
- Repair of resin and acrylic temporary materials



■ Syringe Refills

Content: 2 x Syringes (2 g), 20 x Dispensing Tips

Shade	Item no.	Shade	Item no.
A1	3M-7032-A1	B2	3M-7032-B2
A2	3M-7032-A2	C2	3M-7032-C2
A3	3M-7032-A3	OA3	3M-7032-OA3
A3.5	3M-7032-A3.5	D2	3M-7032-D2
A4	3M-7032-A4	W	3M-7032-W
B1	3M-7032-B1	XW	3M-7032-XW

LuxaFlow® Star

/ DMG

✓ *This Flowable, Light-Curing Composite Is Well-Suited For Esthetic Repairs And Add-Ons Of Luxatemp (And Other Bisacrylate) Temporaries. LuxaFlow Star Was Tailored To Luxatemp Star And Its Shades, Thus Ensuring A Natural Appearance Of Temporaries Even After A Repair.*

Advantages

- Natural fluorescence
- Natural appearance in all light conditions
- Easy to polish
- Fluid during and rigid after the application
- Radiopaque

- Indications**
- Small class III, IV and V restorations
 - Minimal invasive restorations and underfillings



■ Assortment Pack

Content	Item no.
- 4 x Syringes (1.5 g) (1 each of shade A1, A2, A3.5, B1)	DMG-214000
- 20 x Luer-Lock-Tips	

■ Refills

Content: 2 x Syringes (1.5 g)

Shade	Item no.	Shade	Item no.
A1	DMG-214001	A3.5	DMG-214004
A2	DMG-214002	B1	DMG-214005
A3	DMG-214003		

EMBRACE™ WetBond™ Class V / Pulpdent

✓ *Light-Cure, Fluoride-Releasing, Radiopaque Cervical Restorative Resin*

Advantages

- Contains No Bisphenol A, No Bis-GMA, No BPA Derivatives
- Embrace Class V restorative is moisture tolerant, which is particularly advantageous for restoring abfraction lesions and cervical areas, where moisture control is most difficult
- Bonds in a slightly moist field
- Self-etching to dentin
- Self-adhesive properties – bonding agents are optional
- Treats and desensitizes abfractions and exposed roots
- Excellent marginal adaptation, adhesion and retention



No Bonding Needed!



Shade	Content	Item no.
A2	- 1 x Syringe (1.2 g) - Accessories	PULP-EMVA2
A3.5		PULP-EMVA3.5
B2		PULP-EMVB2





Restorative

Direct Composite Veneer

Edelweiss™ Prefabricated Composite Veneers

/ Ultradent

Edelweiss composite veneers correct anatomical deformities, repair damage caused by bruxism, change the shape and length of teeth close a diastema, fix tooth discoloration, or treat caries. At a fraction of the cost for most lab-based systems, Edelweiss unlocks what was once out of reach for many people: a long-lasting, esthetic smile.

Indications

- Laser Sintering heats up and melts the surface to form a highly resistant, very durable finish
- Extremely thin surface (0.2-0.6mm) gives very high translucency and enables less tooth preparation
- Contoured mamelons create a completely natural look
- Economical price allows you to offer esthetic smiles to more patients
- One-appointment procedure—start to finish & No lab time
- Minimally invasive tooth preparation



The First And Only Laser-Sintered Veneer



Edelweiss Kit



Content	Item no.
- 1 x Set of Small Uppers (#6 - #11)	UL-3922
- 1 x Set of Medium Uppers (#6 - #11)	
- 1 x Set of Large Uppers (#6 - #11)	
- 1 x Ultra-Etch Syringe (1.2 mL)	
- 1 x Peak Universal Bond Syringe (1.2 mL)	
- 7 x Jiffy Polishing Cups	
- 5 x Composite Compules (0.3 g) of each of Amelogen Plus Shade (A1, A2, A3, A3.5, N, SL)	
- 20 x Inspirial Bush Tips	
- 10 x Blue Micro Tips	
- 1 x Shade Guide	

Refills

Content: 6's/pack

Items	Item no.
Small Uppers (#6-#11)	UL-39228
Small Uppers (#22-#27)	UL-39231
Medium Loweres (#6-#11)	UL-39229
Medium Loweres (#22-#27)	UL-39232
Large Uppers (#6-#11)	UL-39230

Fluoride Releasing Composites

Universal Composites

Beautiful® Injectable

/ Shofu

Advantages

- The only universal restorative that comes in an injectable syringe, packed with exceptional strength and sculptable beauty at your fingertips
- Formulated with a higher filler load and a uniquely strengthened resin matrix, Beautiful Injectable exhibits a greater degree of polymerization to achieve remarkable mechanical properties



- The filler-resin structure has been optimized to control light diffusion within to mimic the optical characteristics of natural teeth. Inherent natural aesthetics, fluorescence and opalescence enable the clinician to create indiscernible restorations

Indications

- Class I through Class V restorations
- Direct veneers
- Core build-up
- Shape & shade modification
- Minimally invasive anterior and posterior restorations
- Restoration of occlusal surface
- Repair of composite
- Ceramic restorations

Syringe Refills

Shade	Item no.	Shade	Item no.
A0.5	SF-2050	A20	SF-2055
A0.50	SF-2051	A3	SF-2056
A1	SF-2052	A3.5	SF-2057
A10	SF-2053	A4	SF-2058
A2	SF-2054	BW	SF-2059
Universal A2*	SF-2054B	INC	SF-2060

*Universal A2: 1g/syringe

Beautiful® II / Shofu

✓ A Fluoride Releasing, Highly Aesthetic Dental Material For All Classes Of Restorations



Advantages & Features

- Eight-year University of Florida clinical study published in JADA showed no secondary caries, no failures, no post-operative sensitivity and high retention of color match and luster
- S-PRG (Surface Pre-reacted Glass Ionomer) technology – ability to recharge fluoride
- Appropriate for patients with a high caries index
- Excellent fluorescence and light diversifying aesthetics
- Eliminates shade variations associated with polymerization
- High level of radiopacity Enamel-like surface hardness

Syringe Refills

Shade	Item no.	Shade	Item no.
A1	SF-1401	B3	SF-1410
A2	SF-1402	C2	SF-1412
A3	SF-1403	C3	SF-1413
A3.5	SF-1404	Incisal	SF-1415
A4	SF-1405	Bleach	SF-1416
B1	SF-1408	A20	SF-1419
B2	SF-1409	A30	SF-1420



Fluoride Releasing Composites

Flowable Composites

Beautifil® Flow / Shofu

✓ *Designed For A Wide Range Of Applications: Minimal Lesions Such As Class I, III And V Cavities (minimal intervention), Pit And Fissure Sealant, Base And Liner Under Direct Restorations*



- A light-cure fluoride releasing aesthetic dental restorative material developed through the Giomer concept using the S-PRG technology (Surface Pre-Reacted Glass Ionomer) in two distinct viscosities **F02 (Low Flow)** and **F10 (High Flow)**
- Excellent aesthetics and high surface polish
- Unique Anti Plaque Effect
- Fluoride release and recharge
- High elasticity and stable stress relaxation
- Remarkable Radiopacity
- User friendly syringe design; dispense on demand with no oozing of excess paste, 360° Rotary Finger Grip and multi-purpose rubber shield cap

F02 (Low Flow): Effective to use where paste flows moderately and shape should be retained, such as cervical areas

F10 (High Flow): Effective to use where smooth flow is necessary, such as uneven or narrow spaces, and also can be used as a liner

■ Syringe Refills Content: 1 x Syringe (2 g), 5 x Needle Tips

F02 (Low Flow)		F10 (High Flow)	
Shade	Item no.	Shade	Item no.
A1	SF-1431	A1	SF-1461
A2	SF-1432	A2	SF-1462
A3	SF-1433	A3	SF-1463
A3.5	SF-1434	A3.5	SF-1464
A4	SF-1435	A4	SF-1465
A30	SF-1449	A30	SF-1479
Gum	SF-1451	A3T	SF-1482

Beautifil® Flow Plus / Shofu

✓ *Beautifil Flow Plus Is The Next Step In The Evolution Of Restorative Materials, Combining The Convenience Of A Flowable Delivery System, With Strength And Functionality That Surpass Leading Hybrid Composites*

Features

- Physical properties rival leading composites
- Precise handling
- Superior adaptation
- Self-leveling, bubble free consistency
- Radiopacity greater than enamel
- Fluoride release and recharge

Advantages

- Approved for all classes!
- Stacks without slumping and just stays put
- Ensures a tight marginal seal, decreasing voids
- Reduces instrumentation/polishing time
- Allows clear visibility during follow-ups
- Ideal for high caries index patients

■ Syringe Refills Content: 1 x Syringe (2.2 g)

F00 (Zero Flow)



(SF-)			
Shade	Item no.	Shade	Item no.
A0.5	2000	A0.50	2006
A1	2001	A10	2007
A2	2002	A20	2008
A3	2003	Incisal	2010
A3.5	2004	Bleach White	2011
A4	2005		

F03 (Low Flow)



(SF-)			
Shade	Item no.	Shade	Item no.
A0.5	2013	A30	2021
A1	2014	Incisal	2023
A2	2015	Bleach White	2024
A3	2016		
A3.5	2017	Milky White	2026
A4	2018		
A20	2020	Cervical	2027

Bulk Fill Posterior Composites

Beautifil®-Bulk Restorative / Shofu

✓ *Provides Exceptional Condensability And Sculptability Ideal For Posterior Restorations*

✓ *Delivering Superb Shade Stability, Low Shrinkage Stress And Full Polymerization At 4mm*



Advantages

- High fill ratio (87.0wt% 74.5vol%) reduces polymerization shrinkage and shrinkage stress while increasing compressive and flexural strength
- Ideal balance of light diffusion and transmission ensures complete polymerization at 4mm depth of cure
- Outstanding Vickers hardness value
- Optimum translucency creates esthetic shades unaffected by surrounding intra-oral color
- Sustained fluoride release and rechargeability
- High level of radiopacity

■ Syringe Refills

Shade	Packing	Item no.
A	1 x Syringe (4.5 g)	SF-2035
Universal		SF-2034

Bulk Fill Flowable Composites

Beautifil®-Bulk Flowable / Shofu

✓ *Designed To Deliver Aesthetic Results, Low Shrinkage Stress And Complete Polymerization At 4mm Depth Of Cure*



Advantages & Features

- High filler ratio (73wt%) helps to reduce volumetric shrinkage (3.5%) and shrinkage stress while increasing compressive and flexural strength
- Vickers Hardness test results (48 Hv) show above 95% polymerization at 4mm depth of cure
- High flexural strength (119 Mpa) with a low flexural modulus (9 Gpa)
- Natural Opaque shade stability to match surrounding tooth structure
- Superior fluoride release and rechargeability
- Self-leveling allows easy adaptation
- 10 second cure time
- High level of radiopacity, 2.5 Al:mm

■ Syringe Refills

Shade	Packing	Item no.
Universal	1 x Syringe (2.4 g)	SF-2030
Dentin		SF-2031



Anaesthetics
 Articulating
 Bur and Abrasives
 Cosmetics
 Diagnostic
 Disposables
 Endodontics
 Equipment
 Handpieces
 Implant
 Infection Control
 Instruments
 Isolation
 Laboratory
 Oral Surgery
 Orthodontics
 Patient Education and Gift
 Pedodontics
 Periodontal
 Post and Core
 Preventive
 Prosthodontics
 Restorative
 Uniform

Restorative

Fluoride Releasing Composites

MiCD (Minimally Invasive Cosmetic Dentistry)

/ Shofu



✓ *A holistic approach that integrates minimally invasive treatment techniques with aesthetic dentistry to enhance the smile while taking into consideration the psychology, health, function and aesthetics of the patients*

✓ *Using MiCD procedures that preserve maximum natural tooth for better oral health while still achieve dramatic aesthetics result*

The MiCD Concept

- Non Invasive Treatment Options

- Smile Exercise
- Tooth Whitening
- Sectional Orthodontics
- Non-prep Veneers
- Enamel Augmentations
- Adhesive Pontics as long term temporary restoration

- Minimally Invasive Treatment Options

- Selective Contouring of Teeth & Gums
- Direct Restorations with minimal tooth preparation
- Partial Veneers, Inlays, Onlays
- Minimally Invasive implants
- Air Abrasive
- Fissurotomy
- Adhesive bridges

Content	Item no.
- 3 x Beautifil II Individual Syringes (4.5 g) <i>(1 each of Shade A2, A3, A3.5)</i>	SF-MICD KIT
- 3 x Beautifil Flow (2 g) <i>(1 each of F00: A2, A3 and F03: A3)</i>	
- 1 x BeautiBond Set (SF-1781) <i>(1 x BeautiBond - 6 mL Bottle, 1 x Microtip - with 1 handle, 50 x Application Tips, 25 x V-Dish)</i>	
- 1 x Direct Cosmetic Restoration Kit (SF-F272) <i>(8 x Diamond Points FG - #340S, #MI-2P, #121, #102R, #F216, #SF104R, #SF215, #SF219)</i>	
- 1 x Super-Snap Mini Kit, CA (SF-0505) <i>(48 x Super-Snap Disks, 1 x Composite Fine-Point, 1 x Super-Snap Mandrel)</i>	
- 1 x Silicone One Gloss Midi (box) <i>(25 x Silicone Midi Points, 1 x Mandrel CA)</i>	
- 2 x Packs Super-Snap Polystrips, Assorted <i>(Coarse, Medium, Fine & Superfine)</i>	

Glass Ionomers

Light-Cure GI



Ketac™ Nano Light Curing Glass Ionomer Restorative / 3M ESPE

✓ *The First Paste/Paste, Resin-Modified Glass Ionomer Developed With Nanotechnology*

Quick Mix Restorative Trial Kit

Content	Item no.
- 10 x A2 Quick Mix Capsules	3M-3305TK
- 10 x A3 Quick Mix Capsules	
- 1 x Primer (3 mL)	

Advantages & Features

- The first Paste/paste RMGI restorative based on bonded nanofiller technology
- Esthetic glass ionomer restorative
- High fluoride release
- Improved wear resistance
- Light cure on demand
- Radiopaque
- Automatic mixing for fast dispensing from the unit-dose Quick Mix Capsule

Cap슐e Refills Content: 20 x Capsules

Shade	Item no.	Shade	Item no.
A1	3M-3305-A1	A3.5	3M-3305-A3.5
A2	3M-3305-A2	B2	3M-3305-B2
A3	3M-3305-A3		

Ketac™ N100 Light Curing Nano-Ionomer Restorative / 3M ESPE



✓ *Ketac™ N100 Light-Curing Nano-Ionomer Restorative Paste/Paste Form Makes An Ideal Alternative Esthetic Glass Ionomer Solution For Patients Ranging From Pediatric Through Geriatric*

Advantages

- Faster, easier mixing and dispensing - with accurate, reproducible dosing from the award-winning Chicker Dispenser for 3M ESPE
- The first paste/paste RMGI restorative based on bonded nanofiller technology
- Esthetic glass ionomer restorative
- High fluoride release
- Improved wear resistance

Indications

- Primary teeth restorations
- Small Class I restorations
- Class III and V restorations
- Transitional restorations
- Sandwich restorations
- Core build-ups

Introductory Kit

Content	Item no.
- 1 x Clicker™ Dispenser A2 (12 g)	3M-3527I
- 1 x Clicker™ Dispenser A3 (12 g)	
- 1 x Primer (6.5 mL)	
- 50 x Delivery Tips with Pistons	
- 1 x Mixing Pad	
- 60 x Fiber Tips	
- 1 x Applicator Handle	
- 48 x Disposable Mixing Wells	
- Shade Guide, Instructions, Technique Guide	

Trial Kit

Content	Item no.
- 1 x Clicker™ Dispenser A3 (12 g)	3M-3527TK
- 1 x Primer (6.5 mL)	
- 10 x Delivery Tips with Pistons	
- Instructions, Technique Guide	

Refills Content: 1 x Clicker™ Dispenser (12 g)

Shade	Item no.	Shade	Item no.
A1	3M-3527-A1	A3.5	3M-3527-A3.5
A2	3M-3527-A2	B2	3M-3527-B2
A3	3M-3527-A3		



Glass Ionomers

Dual-Cure GI

Photac™ Fil Quick Aplicap™ Glass Ionomer Restorative / 3M ESPE



- ✓ *Photac™ Fil Quick Material Is A Fluoride-Releasing, Resin-Modified Glass Ionomer*
- ✓ *It's A Self-Adhesive Property Eliminates The Need For Bonding And Cavity Conditioning*
- ✓ *The Result Is A Fast And Easy Solution For Treating Anxious Pediatric And Geriatric Patients*

■ Capsule Refills Content: 50 x Capsule

Shade	Item no.	Shade	Item no.
A2	3M-61020	B2	3M-61050
A3	3M-61030	B3	3M-61040
A3.5	3M-61060		

Advantages

- Dual cure provides light-cure control and the assurance of a complete cure
- High fluoride release
- Excellent mechanical properties, plus it's radiopaque and biocompatible
- Convenience, simplicity and accuracy of the Aplicap delivery system
- Conservative cavity preparation, fluoride release and fast, easy technique make it ideal for geriatric and pediatric patients

Indications

- Class III and Class IV restorations
- Small Class I restorations
- Primary tooth restorations
- Temporary restorations

Self-Cure GI

Ketac™ Fil Quick Aplicap™ Glass Ionomer Restorative / 3M ESPE



- ✓ *The Radiopaque Conventional Glass Ionomer Filling Material In The Unique Aplicap Delivery System Offers High Fluoride Release And Helps You Achieve Virtually Invisible Anterior Restorations*

■ Capsule Refills Content: 50 x Capsule

Shade	Item no.	Shade	Item no.
A1	3M-55020	A3.5	3M-55050
A2	3M-55030	A4	3M-55060
A3	3M-55040		

Advantages

- Provides molecular bond to enamel and dentin to guarantee a tight marginal seal
- Offers high fluoride release
- Radiopaque
- Consistent, easy and unit-dose dispensing with high visibility during placement
- Minimal preparation
- Available in eight esthetic shades to better match existing tooth structure

Indications

- Class III and Class IV restorations
- Wedge-shaped defects
- Small Class I fillings
- Core buildups

Ketac™ Molar Quick Aplicap™ / 3M ESPE

- ✓ *Ketac™ Molar Quick Material Is A Bulk-Fill, Packable And Fast-Settling Conventional Glass Ionomer*
- ✓ *High Compressive Strength And Marginal Integrity Make It A Glass Ionomer Of Choice For Posterior Restorations*



Advantages

- High compressive strength
- Minimal abrasion
- Excellent handling characteristics; plastic contouring and condensing properties similar to amalgam
- Aplicap capsule provides consistent quality, simplified handling and a long nozzle for direct application of the material

Indications

- Primary teeth fillings
- Single surface fillings in non-occlusal areas
- Cervical fillings in posterior region
- Base material for single- and multiple-surface composite fillings
- Core buildups prior to crown placement
- Minimally Invasive Dentistry (MID)

■ Capsule Refills

Content: 20 x Capsules

Shade	Item no.
A3	3M-56360
A2	3M-56361
B2	3M-56362

Content: 50 x Capsules

Shade	Item no.
A1	3M-56366
A2	3M-56363
A3	3M-56364

Ketac™ Silver Glass Ionomer Restorative / 3M ESPE



- ✓ *This Silver Glass Ionomer Contains No Amalgam And Offers High, Long-Term Fluoride Release As Well As Adhesion To Enamel And Dentin To Ensure Sealed, Lasting Restorations Without Need For Undercut*

Advantages

- Silver particles are sintered into the glass using a 3M ESPE process, which makes Ketac Silver glass ionomer restorative an excellent thermal and electrical insulator
- Offers high, sustained fluoride release
- Contains no amalgam
- Aplicap™ and Maxicap™ delivery systems provide accurate, consistent mixing and easy dispensing
- Maxicap is indicated for extended core build-ups

Indications

- Core buildups
- Base under amalgams
- Temporary restorations
- Primary teeth fillings
- Small fillings in posterior teeth

Item	Packing	Item no.
Aplicap™	50 Capsules/ Blister Pack	3M-37010
Maxicap™		3M-56210





Restorative

Glass Ionomers

Self-Cure GI

Ketac™ Molar Aplicap™ & Easymix / 3M ESPE

Advantages

- High compressive strength and surface hardness
- Minimal abrasion
- Packable
- Fluoride release

Indications

- Protective lining under composite or amalgam restorations
- Restorations of primary teeth
- Semi-permanent in single- or multi- surface restorations of primary and permanent teeth
- Semi-permanent Class III & V restorations of primary and permanent teeth
- Cervical fillings without cosmetic priorities
- Core build-ups
- Fissure sealings



Ketac™ Molar Aplicap™



Ketac™ Molar Easymix

■ Ketac™ Molar Aplicap™ Refills

Item	Packing	Item no.
A1	50 Capsules/ Pack	3M-56430
A3		3M-56420

■ Ketac™ Molar Easymix MID & ART Kit

Content	Item no.
- 1 x Powder A3 (12.5 g) - 1 x Liquid (8.5 mL) - 1 x Spoon - 1 x Mixing Pad - 2 x Technique Guides	3M-56633

Glaslonomer FX-II Improved / Shofu

✓ An Enhanced Direct Restorative Glass Ionomer Designed In Vita Shades With Remarkable Mechanical Properties And Ease Of Handling



NEW

Advantages

- Superior mechanical properties
- Enhanced aesthetics with high translucency and optimal polishability
- Exceptional handling characteristics
- Continuous release and recharge of fluoride ions, biocompatible and radiopaque with low solubility
- Packable and easy to use

Indications

- Restoration of pits and fissures, cervical erosion and abrasive
- Intermediate restoration for heavy stress bearing areas
- Final restoration for non-stress bearing Class I and II in adult dentition
- Created for Atraumatic Restorative Treatment (ART) and Minimal Intervention (MI) dentistry
- Geriatric and Pediatric restorations

■ Complete Set

Shade	Packing	Item no.
A2	- 1 x Powder (15 gm) - 1 x Liquid (8.4 mL)	SF-3245
A3		SF-3246
A3.5		SF-3247
B2		SF-3248

Glaslonomer FX-II Capsule / Shofu

✓ The Proven Direct Restorative Glass Ionomer Now Available In Capsule Delivery System To Achieve Simple And Predictable Application Every Time



NEW

Advantages

- Chemically adheres to tooth structure with continuous release of fluoride ions
- Vita shaded with fluorescence close to natural teeth for better aesthetics
- Packable consistency with high compressive strength and excellent biocompatibility
- Assures a consistent mix every time with adequate working time
- Fast set at intra-oral temperature

■ Capsule Refills

Shade	Packing	Item no.
A2	50 x Capsules (0.55 gm)	SF-3270
A3		SF-3271
A3.5		SF-3272
B2		SF-3273



- Anaesthetics
- Articulating
- Bur and Abrasives
- Cosmetics
- Diagnostic
- Disposables
- Endodontics
- Equipment
- Handpieces
- Implant
- Infection Control
- Instruments
- Isolation
- Laboratory
- Oral Surgery
- Orthodontics
- Patient Education and Gift
- Pedodontics
- Periodontal
- Post and Core
- Preventive
- Prosthodontics
- Restorative
- Uniform

Glass Ionomers

Self-Cure GI

Zirconomer Improved / Shofu

✓ *Zirconia-Reinforced Glass Ionomer*

✓ *Ideal for restoration of selective class I & II permanent restorations, as structural base in sandwich restoration, core build-up, pediatric and geriatric restoration, long-term temporary replacement of fractured cusps, fractured amalgam restoration, ART as well as restoration of root surfaces where overdentures rest*

Indications

- Reinforced with special zirconia fillers to match the strength and durability of amalgam
- Sustained high fluoride release for anti-cariogenic benefits especially in cases with his caries risk
- Packable and condensable like amalgam without the hazard of mercury or the risk of corrosion, expansion and thermal conductivity
- High flexural modulus and compressive strength ensures longevity in stress bearing areas
- Chemically bonds to enamel/dentin with co-efficient of thermal expansion similar to tooth structure for greater marginal integrity
- Easy mixing and handling characteristics with adequate working time and snap-set reaction
- Excellent resistance to abrasion and erosion



Zirconomer 1-1 Set

Content	Item no.
- 1 x Powder (12 gm) - 1 x Liquid (5 mL) - 1 x Mixing Pad - Spoon & Spatula	SF-3130

GI Conditioner & Primer

Ketac™ Conditioner Refill / 3M ESPE

/ 3M ESPE

✓ *A Conditioner For The Pretreatment Of Dentin, Ketac™ Conditioner Pretreatment Solution Removes The Dentin Smear Layer Associated With Tooth Preparation Procedures*



Packing	Item no.
10 mL/bottle	3M-37470

Ketac™ Nano Glass Ionomer Primer / 3M ESPE

/ 3M ESPE

✓ *Ketac™ Nano Primer Is An One Part, Visible Light-Cure Liquid Specifically Designed For Use With Ketac™ Nano Restorative*



- Comprised of the Vitrebond™ copolymer, HEMA, water, and photoinitiators. The primer is acidic in nature. Its function is to modify the smear layer and adequately wet the tooth surface to facilitate adhesion of Ketac Nano restorative to the hard tissue.

Packing	Item no.
6.5 mL/bottle	3M-3304P

Base & Liner

GI Base & Liner

Vitrebond™ Light Cure Glass Ionomer Liner/Base / 3M ESPE

✓ *A Light-Cure Glass Ionomer Liner/Base That Combines Light-Cure Convenience With True Glass Ionomer Fluoride Release And Bond Strength*

Advantages

- Fast technique and excellent handling characteristics allow for quick and easy placement
- Forms a strong bond to dentin, sealing the tooth to reduce the incidence of microleakage which helps reduce post-operative sensitivity
- Produces zones of inhibition against bacteria commonly found in the oral cavity in vitro
- Sustained fluoride release to help protect the tooth long after the restoration is complete
- Provides excellent adhesion to 3M dental adhesive systems, whether using composite or amalgam

Indications

- Under direct and indirect restorations, including composite, amalgam, ceramic and metal restorations



Intro Kit

Content	Item no.
- 1 x Powder (9 g) - 1 x Liquid (5.5 mL) - 1 x Spoon - Mixing Pads	3M-7510

Refills

Item	Packing	Item no.
Powder	9 g/bottle	3M-7512P
Liquid	5.5 mL/bottle	3M-7512L

Anaesthetics
Articulating
Bur and Abrasives
Cosmetics
Diagnostic
Disposables
Endodontics
Equipment
Handpieces
Implant
Infection Control
Instruments
Isolation
Laboratory
Oral Surgery
Orthodontics
Patient Education and Gift
Pedodontics
Periodontal
Post and Core
Preventive
Prosthodontics
Restorative
Uniform



Restorative

Base & Liner

GI Base & Liner

Vitrebond™ Plus Light Cure Glass Ionomer Liner/Base / 3M ESPE

✓ *This Convenient Paste-Liquid Formula In Our Unique Clicker™ Dispenser Offers The Protection Of A Conventional Glass Ionomer Liner Without The Mess Of Measuring*

NEW



Advantages

- Faster, easier dispensing and the right mix every time
- Conventional paste-liquid formula offers the protection of a conventional glass ionomer liner without the mess of measuring
- Award-winning Clicker dispenser assures a consistent mix
- Helps protect against post-op sensitivity: offers a strong bond to dentin, sealing the tooth structure to protect against micro-leakage
- Low microleakage comparable to original Vitrebond™ Liner/Base
- Fluoride release of a true glass ionomer
- Excellent handling characteristics for easy placement
- Cures in 20 seconds - a faster cure time than powder-liquid formulas
- Delivers excellent adhesion with any dental adhesive system
- Radiopaque
- No additional hardware required

Indications

- Liner/Base under direct and indirect restorations, including composite, amalgam, ceramic and metal restorations

Content	Item no.
- 1 x Clicker™ Dispenser (10 g)	3M-7571
- 1 x Mixing Pad	

Calcium Hydroxide Base & Liner

Ultra-Blend® Plus / Ultradent



Features & Advantages

- No mixing necessary
- Will not dissolve over time
- Radiopaque
- Controlled, precise syringe delivery
- Light cured when clinician is ready

Ultra-Blend plus Kit

Content	Item no.
- 2 x Syringes Ultra-Blend plus Dentin (1.2 mL)	UL-415
- 2 x Syringes Ultra-Blend plus Opaque White (1.2 mL)	
- 20 x Black Micro Tips	
- 20 x Black Mini Tips	

Ultra-Blend plus Refills

Shade	Packing	Item no.
Dentin	4 x Syringes (1.2 mL)	UL-416
Opaque White		UL-417

Resin Base & Liner

Ionosit-Baseliner / DMG

✓ *Resin Ionomer Base Liner*



Advantages

- Unique compensation of composite polymerisation shrinkage through controlled expansion
- Continuous fluoride release
- Contains bactericidal zinc ions
- No complicated handmixing
- Protects the pulp
- Seals the dentinal tubuli
- Prevents post-operative sensitivities

Speciality - Easy and fast application with the Saferinge

Indications

- Underfilling material for composite, amalgam and ceramic restorations and for blocking out undercuts

Package

Content	Item no.
2 x Syringes (1.5 g)	DMG-213002
3 x Saferinge (0.33 g)	DMG-220091

Economy Package

Content	Item no.
20 x Saferinge (0.33 g)	DMG-209538

Amalgam



HENRY SCHEIN® Cyber Tech / Henry Schein

- Excellent handling characteristics
- Uniform consistency
- High copper content for greater strength, less tarnish and corrosion, and less creep

- Firm but smooth carving and burnishing surface
- Good marginal sealing

Non Gamma 2 Amalgam with High Copper

Items	Packing	Item no.
Yellow Cap	1 spill	HSPL-900-6197
Orange Cap	2 spills	HSPL-900-6198
Green Cap	3 spills	HSPL-900-6199



Opaquers

Beautifil® Opaquers

/ Shofu

- Light-cure
- Fluoride releasing
- Flowable opaquer to mask intense discolorations

Refills

Item	Packing	Item no.
Universal Opaquer (UO)	1 x Syringe (2 g)	SF-1363
Light Opaquer (LO)		SF-1364



Temporary Fillings

Cavit™ Temporary Filling Material / 3M ESPE

✓ A Self-Cured Temporary Sealing Compound For The Temporary Restoration Of Cavities

- Proven filling material for temporary fillings
- Simple to apply with filling instrument
- Quick and void-free curing in a moist environment
- Slight expansion of the filling material ensures a well-sealed margin

- **Cavit (High Surface Hardness)**
Temporary filling of occlusion-headed-loaded restorations
- **Cavit-W (Reduced Final Hardness, Increased Adhesion)**
Temporary filling after endodontic treatments
- **Cavit-G (Removable Completely without Burs)**
Inlay preparation



Item	Packing	Item no.
Cavit		3M-44030
Cavit-W	28 g/bottle	3M-44130
Cavit-G		3M-44313

UltraTemp® Regular / Ultradent

✓ Polycarboxylate, Non-Eugenol, Temporary Luting/Filling Material

- Provides easy clean up
- Seals to prevent bacterial leakage
- Has antimicrobial properties
- Mixes and delivers in one action
- Provides excellent retention
- Is kind to pulp
- Is non-eugenol



UltraTemp Regular (4-5 mins Set)

Content	Item no.
- 1 x Dual Barrel PermaFlo DC Syringe	UL-5916
- 20 x Mixing Tips	

Temporary Luting + Filling 2 In 1

Accessories

Capsule Activators & Appliers

Activator & Appliers

/ 3M ESPE

Aplicap™ & Maxicap™ Activators

- Easy activation of the powder and liquid to prepare the Aplicap/Maxicap capsule for mixing



Item	Packing	Item no.
Aplicap™ Activator	Each	3M-73040
Maxicap™ Activator		3M-73060

Aplicap™ & Maxicap™ Appliers

- Easy delivery of Aplicap/Maxicap capsule material to the tooth preparation or to the restoration



Item	Packing	Item no.
Aplicap™ Applier	Each	3M-73050
Maxicap™ Applier		3M-73070

Dappen Dishes

HENRY SCHEIN® Glass Dappen Dishes / Henry Schein



Packing: Each

Color	Item no.	Color	Item no.
Clear	HS-100-9211	Blue	HS-100-0259
Amber	HS-100-5501	Green	HS-100-6529

Plastic Dappen Dishes / Henry Schein




- Disposable
- Assorted Colors

Packing	Item no.
1,000's/box	HS-900-1878

Mixing Pads

Mixing Pads / 3M ESPE

- Allows for easy mixing of 3M™ ESPE™ Dental materials
- Individual mixing sheets are simply discarded after use, for fast and easy clean up
- Tear-off, disposable coated sheets allow for easy hand-mixing of dental materials



Item	Size (mm)	Packing	Item no.
Extra Large	150 x 240	34 sheets/pad, 10's/pack	3M-70180
Large	125 x 180	24 sheets/pad, 10's/pack	3M-70190
Small	87 x 155	20 sheets/pad, 10's/pack	3M-70128
Universal	70 x 80	60 sheets/pad, 10's/pack	3M-70144







Restorative

Accessories

Tips & Syringes




HENRY SCHEIN®
Prebent Multipurpose Dispensing Tip / Henry Schein



- Prebent to the proper angle with a blunt end for comfort
- Designed for placing etchants, sealants, cements and other flowable materials
- Universal, color-coded, Luer-lock design

Volume	Color	Description	Packing & Item no. (HS-)	
			Standard (20's/pack)	Bulk (100's/pack)
18 ga.	Pink	Cavity Liners Bases & Composites	900-4585	900-4581
19 ga.	Black	Composites & Cements	900-4582	900-4577
20 ga.	Yellow	Composites & Sealants	900-4584	900-4579
22 ga.	Gray	Sealants Etchants & Cements	900-4586	900-4580
25 ga.	Blue	Etchants	900-4583	900-4578

Filtek™ Flow Dispensing Tips / 3M ESPE



- Single-use tips allow for easy dispensing of Filtek™ Z350XT Flowable restorative from syringes

Packing	Item no.
20's/pack	3M-3700T

Scotchbond™ Dispensing Tips / 3M ESPE



- Disposable syringe tips for dispensing Scotchbond™ Phosphoric Etchant

Packing	Item no.
25's/pack	3M-41274

Beautiful Flow Needle Tips / Shofu




Packing	Item no.
50's/pack	SF-1484

Micro® Tips / Ultradent

- Provides pinpoint precision
- Narrow cannulae accurately deliver materials


Volume	Color	Application	Packing	Item no.
22 ga	Black	Ultra-Blend plus, EtchArrest, OpalDam, UltraSeal XT plus, and Ultracare	100's/pack	UL-1085
25 ga	Blue	Ultra-Etch		UL-127



Metal Dento-Infusor® Tips / Ultradent

- Places hemostatic agents precisely and effectively removes superficial coagulum
- Blunt, bent cannula with padded brush enables gentle pressure on the sulcus
- Ultradent's first tip, the "MDI" remains paramount for successful tissue management

Volume	Color	Application	Packing	Item no.
19 ga	Yellow	Astringent, Astringent X, and ViscoStat Regular, Wintermint, and Clear	100's/pack	UL-125



Barrier Sleeve / Shofu



Packing	Item no.
50 pcs	SF-2070
200 pcs	SF-2071

Applicator Brush Tips / 3M ESPE

- Single-use brush tips for use with 3M™ ESPE™ Applicator Handles

Packing	Item no.
60's/pack	3M-1919B



Applicator Handles / 3M ESPE



- A set of three color-coded applicator handles for use with 3M™ ESPE™ adhesive systems

Packing	Item no.
- 3 x Handpieces, colored	3M-1919
- 120 x Disposable Brush Tips	



Crown Forms

CoForm™ / Directa

- Aids pressure to push material into cavities and etched tubules
- Does not adhere to composite so is easy to remove without drag
- Conforms easily to patient's dentition to give a more natural restoration
- Reduces wastage

Packing	Item no.
64's/set (11 Mesial/Distal, 21 Mesial/Distal, 12 Mesial/Distal, 22 Mesial/Distal)	DIR-6040-34

Composite Capsule Dispensers

HENRY SCHEIN® Multi-Use Dispenser / Henry Schein



- Heavy-duty auto-clavable plastic
- Ergonomic design

Packing	Item no.
Each	HS-900-0790

Restorative Dispenser for Capsules / 3M ESPE

- Provides quick and easy delivery of 3M ESPE restorative materials packaged in capsules
- Easy-to-use capsules and dispenser allow for direct, single-patient delivery



Packing	Item no.
Each	3M-5707SD

Composite Wetting Resin / Ultradent




- Composite Wetting Resin is a 45% filled, light-cured, liquid resin
- It is significantly superior to single-component adhesives, which contain solvents and inhibit composite polymerization
- Facilitates composite adaptation
- Does not act as a detriment to polymerization
- Moistens dry composite during contouring

Mini Refill	
Content	Item no.
2 x Syringes (1.2 mL)	UL-3059

Consepsis® Scrub / Ultradent

✓ Chlorhexidine Anti-bacterial Slurry



- Reduces post-op sensitivity
- Does not compromise bond strength
- Comes in a non-splatter formula
- Scrubs, cleans, and disinfects in one operation

Refill	
Content	Item no.
4 x Syringes (1.2 mL)	UL-732



Anaesthetics
Articulating
Bur and Abrasives
Cosmetics
Diagnostic
Disposables
Endodontics
Equipment
Handpieces
Implant
Infection Control
Instruments
Isolation
Laboratory
Oral Surgery
Orthodontics
Patient Education and Gift
Pedodontics
Periodontal
Post and Core
Preventive
Prosthodontics
Restorative
Uniform