



03

Bur & Abrasives

Diamond Burs	p.13 - p.17	Polishing Paste	p.25
Carbide Burs	p.18 - p.19	Polishing Cups, Points, Wheels & Rubbers	p.26 - p.27
Surgical Burs	p.19	Polishing Discs	p.28
Acrylic Burs	p.19 - p.20	Cutting Disks	p.28
Finishing Abrasives	p.21	Bur Kits	p.29 - p.32
Polishing Abrasives	p.22 - p.24	Bur Blocks & Accessories	p.33
Abrasive Strips	p.25		
Polishing Strips	p.25		

Bur & Abrasives



Bur & Abrasives

Diamond Burs




HENRY SCHEIN®
Maxima® Diamonds / Henry Schein

- Special Triton™ coating for strengthened diamond coating
- 60% diamond grain exposure for fast cutting

- SUPER FINE**
- FINE**
- MEDIUM**
- COARSE**
- SUPER COARSE**

*COLOR CODED for Easy Identification



Round

(HS-)				
SHAPE	801	801	801	801
Size (1/10mm)	010	014	016	018
Length (mm)	1.0	1.4	1.6	1.8
Total Length (mm)	19	19	19	19
SUPER FINE				
FINE			900-7084	900-7088
MEDIUM	900-7077	900-7081	900-7085	900-7089
COARSE		900-7080	900-7083	900-7087
SUPER COARSE		900-7082	900-7086	900-7090

Interproximal

(HS-)	
SHAPE	392
Size (1/10mm)	016
Length (mm)	5.0
Total Length (mm)	19
SUPER FINE	900-7075
FINE	900-7073
MEDIUM	900-7074

Inverted Cone

(HS-)				
SHAPE	805	805	805	805
Size (1/10mm)	009	010	012	016
Length (mm)	0.9	1.0	1.5	1.5
Total Length (mm)	19	19	19	19
SUPER FINE				
FINE				
MEDIUM	900-7106	900-7107	900-7109	900-7110
COARSE			900-7108	
SUPER COARSE				

Depth Marker

(HS-)		
SHAPE	834	834
Size (1/10mm)	016	021
Length (mm)	4.5	4.7
Total Length (mm)	20	20
SUPER FINE		
FINE		
MEDIUM	900-7135	900-7136
COARSE		
SUPER COARSE		

Pear

(HS-)	
SHAPE	830
Size (1/10mm)	010
Length (mm)	3.0
Total Length (mm)	19
SUPER FINE	
FINE	
MEDIUM	900-7129
COARSE	900-7128
SUPER COARSE	900-7131

Flame

(HS-)			
SHAPE	862	862	862
Size (1/10mm)	009	010	018
Length (mm)	8.0	8.0	8.5
Total Length (mm)	22	22	22.5
SUPER FINE	900-7329	900-7333	900-7338
FINE	900-7327	900-7331	900-7335
MEDIUM	900-7328	900-7332	900-7336
COARSE			
SUPER COARSE			

(HS-)	
SHAPE	889
Size (1/10mm)	009
Length (mm)	4.0
Total Length (mm)	23
SUPER FINE	900-7416
FINE	900-7414
MEDIUM	900-7415
COARSE	
SUPER COARSE	900-7413

Flat End Cylinder

(HS-)		
SHAPE	835	837
Size (1/10mm)	010	012
Length (mm)	4.0	8.0
Total Length (mm)	19	22
SUPER FINE		
FINE		900-7158
MEDIUM	900-7138	900-7156
COARSE	900-7137	900-7155
SUPER COARSE	900-7139	900-7160

Flat End Taper

(HS-)		
SHAPE	847	847
Size (1/10mm)	014	016
Length (mm)	8.0	8.5
Total Length (mm)	22	22.5
SUPER FINE		
FINE	900-7189	
MEDIUM	900-7190	900-7193
COARSE	900-7188	900-7192
SUPER COARSE	900-7191	900-7195

(HS-)		
SHAPE	848	862
Size (1/10mm)	014	018
Length (mm)	10.0	10.0
Total Length (mm)	24	24
SUPER FINE		
FINE		900-7205
MEDIUM	900-7203	900-7206
COARSE	900-7202	900-7204
SUPER COARSE		900-7211

Round End Taper

(HS-)			
SHAPE	850	856	856L
Size (1/10mm)	016	018	018
Length (mm)	10.0	8.0	9.0
Total Length (mm)	24	22	22
SUPER FINE			
FINE	900-7229	900-7260	900-7276
MEDIUM	900-7230	900-7261	900-7277
COARSE	900-7228	900-7259	900-7275
SUPER COARSE	900-7231	900-7262	900-7278

Needle

(HS-)		
SHAPE	859	859
Size (1/10mm)	010	014
Length (mm)	10.0	10.0
Total Length (mm)	24	24
SUPER FINE	900-7298	900-7306
FINE	900-7296	900-7304
MEDIUM	900-7297	900-7305
COARSE	900-7295	
SUPER COARSE		



Bur & Abrasives

Diamond Burs

Shofu Diamond Point FG / Shofu



- Made by hi-tech robots in a computer automated environment untouched by human hands
- Guarantees an even distribution of natural diamond particles onto a single stainless steel shank
- Multi-coating technique with a special surface treatment for durability
- Available in 4 different grits:
 - ✓ Regular Grit (Unbanded),
 - ✓ Fine Grit (Red banded),
 - ✓ Superfine Grit (Yellow banded) and
 - ✓ Coarse Grit (Green banded)

Regular Grit (Unbanded Shank)

✓ Used for cavity cutting and tooth preparation

Benefits & Features

- Able to maintain cutting efficiency with less clogging
- Gentle on hand piece cartridge
- 84 different shapes, some of which are in short shank (categorized under 1000 class) and in Robot Points

Packing	Item no.
Each	SF-DP (See the table)



Regular ▶	401M	401	301	301M	201S	201	101S	101	101L	202	103	301CR	202S	202CR	201CR	101CR	103CR	206CR	207CR	201SR	201R	104R	101R
Short ▶	1401	1301			1201			1101	1101L	1202	1103										1201R	1104R	1101R
SuperShort ▶	3401	3301								3202											3201R	3104R	
Robot Points ▶																		0825	0826				
ISO# ▶	170/010	170/011	170/016	170/013	171/014	171/016	173/014	173/016	173/018	170/021	172/022	544/018	170/017	544/021	545/018	553/019	553/023	544/027	544/031	222/014	222/016	193/016	199/016
L(mm) ▶	3	4	4.5	4	7	7	10	10	10	5	9	4	5	5	7	10	10	3.5	3.5	7	7	10	10
Regular ▶	101SR	102R	103R	105R	106RD	107RD	201Ra	104Ra	101Ra	304	204	104	104L	411	411M	311	211	214	311RM	212R	213R	114S	114
Short ▶		1102R	1103R								1204	1104				1311	1211						
SuperShort ▶		3102R			0836	0839					3204	3104											
Robot Points ▶		0833																					
ISO# ▶	199/013	223/017	223/022	223/022	198/021	198/023	218/016	220/016	220/016	164/016	165/013	166/016	166/018	109/011	109/010	109/013	110/014	156/012	139/013	288/010	289/012	158/010	158/012
L(mm) ▶	9	9	9	11.5	9	9	7	10	10	6	8	10	10	4	3	4.5	7	5	4.5	6	8	8	8
Regular ▶	316	319	215S	215	216	115	116	421M	421	321M	321	120	421SS	421S	321S	435	335	150	440SS	440S	440	340	340S
Short ▶				1215	1216								1421SS	1421S									
SuperShort ▶				3215									3421SS										
Robot Points ▶																		0899	3440SS		3440	1340	
ISO# ▶	247/013	243/017	249/012	249/014	249/017	250/013	250/015	012/010	012/015	012/016	012/018	226/023	019/010	019/014	019/016	032/013	032/015	031/035	001/010	001/013	001/015	001/018	002/016
L(mm) ▶	4	3	7.5	7.5	7.5	10	10	1	1.5	1.4	1.7	7	2.2	3.5	3.5	1.5	1.7	5.7	0.9	1.2	1.3	1.6	3.5
Regular ▶	461R	462R	364R	165R	265R	144	330	230	130	325	111	246	117	145	148	149	248	121	122				
Short ▶																							
SuperShort ▶	3461R																						
Robot Points ▶														0893				0896	0897				
ISO# ▶	237/009	237/011	237/015	237/025	237/021	237/032	041/025	041/043	041/050	304/028	068/045	162/033	254/032	277/029	039/032	038/039	005/018	552/026	552/036				
L(mm) ▶	2	1.7	4.5	6.5	5	5	0.5	0.6	0.7	0.8	1.7	4	5	4.9	5	4	7.5	3.6	4				





Bur & Abrasives

Anaesthetics
Articulating
Bur and Abrasives
Cosmetics
Diagnostic
Disposables
Endodontics
Equipment
Handpieces
Implant
Infection Control
Instruments
Isolation
Laboratory
Oral Surgery
Orthodontics
Patient Education and Gift
Pedodontics
Periodontal
Post and Core
Preventive
Prosthodontics
Restorative
Uniform

Diamond Burs

Fine Grit (Red banded Shank)

✓ Used for fine finishing of tooth preparation

Benefits & Features

- Designed utilizing fine diamond particles to ensure a smooth uniform surface
- 12 shapes, some in short shank (categorized under 1000 class)

Regular ▶	F201R	F104R	F102R	F105R	F204	F104	F215	F216	F440R	F414SS	F414S	F414
Short ▶	F1201R				F1204		F1215					
SuperShort ▶												
Robot Points ▶												
ISO# ▶	222/015	193/015	223/016	223/021	165/012	166/015	249/013	249/016	138/007	150/010	150/012	150/014
L(mm) ▶	7	10	9	11.5	8	10	7.5	7.5	0.7	9	9	9



Packing	Item no.
Each	SF-DP (See the table)


SuperFine Grit (Yellow banded Shank)

✓ Used for super finishing crown preparation, composites, porcelain and glass ionomers

Benefits & Features

- Utilizing ultra-fine diamond particles to ensure a smooth surface
- 27 shapes, some in short shank (categorized under 1000 class)

Regular ▶	SF201	SF101	SF301CR	SF202CR	SF201CR	SF101CR	SF103CR	SF206CR	SF207CR	SF104R	SF102R	SF102RD	SF106RD	SF107RD	SF201Ra	SF206Ra
Short ▶																
SuperShort ▶																
Robot Points ▶																
ISO# ▶	171/015	173/015	544/018	544/021	545/018	553/019	553/023	544/027	544/031	193/015	223/016	198/018	198/021	198/023	218/016	218/012
L(mm) ▶	7	10	4	5	7	10	10	3.5	3.5	10	9	9	9	9	7	7



Packing	Item no.
Each	SF-DP (See the table)

Regular ▶	SF204	SF215	SF217	SF416	SF114	SF318	SF319	SF219	SF150	SF440	SF265R	SF313R	SF145
Short ▶													
SuperShort ▶													
Robot Points ▶													
ISO# ▶	165/010	249/012	249/013	249/010	158/013	254/023	243/014	243/019	031/035	001/015	237/020	284/016	277/029
L(mm) ▶	8	7.5	6	4	8	2	3.4	5	5.7	1.3	5	3.3	4.9



Packing	Item no.
Each	SF-DP (See the table)

Coarse Grit (Green banded Shank)

✓ Used for fast gross cutting

Benefits & Features

- Utilizing diamond particles which are larger than regular type
- More durable and available in 8 shapes

Regular ▶	C101	C202	C102R	C104LL	C165R	C101X	C103X	C105RX
Short ▶								
SuperShort ▶								
Robot Points ▶								
ISO# ▶	173/017	170/022	222/018	160/022	237/025	173/017	173/023	223/023
L(mm) ▶	10	5	9	10	7	10	10	12





Bur & Abrasives

Diamond Burs

Strauss Diamond Point FG / Strauss

- Strauss uses the vast technical knowledge and expertise to offer 100% controlled top quality diamond burs
- Each and every one of Strauss's burs is microscopically tested as part of our strict quality controls



Icon	Description	Colour of Ring	μ
	Extra Extra Fine Burs (XXF)	Purple	15μ
	Super Fine Burs (SF)	Yellow	30μ
	Fine Burs (F)	Red	50μ
	Medium Burs (Neutral) (M)	White	107-120μ
	Coarse Burs (C)	Green	150-180μ
	Super Coarse Burs (SC)	Black	180-250μ
Packing		Item no.	
6's/box		ST-(Cat. No.)-(XXF/SF/F/M/C/SC) (See the table)	

Round

CAT. No.	A0	A1	A2	A3	A4	A5
Head Size	009	012	014	018	022	032
Total Length	19	19	19	19	19	19
SUPER FINE						
FINE						
MEDIUM						
COARSE						

Inverted Cone

CAT. No.	B0	B1	B3
Head Size	009	012	014
Total Length	19	19	19
MEDIUM			
COARSE			

Egg/Football

CAT. No.	K1	K2	K2R	K11
Head Size	016	023	023	018
Total Length	20	20	19	18.5
SUPER FINE				
FINE				
MEDIUM				
COARSE				
SUPER Coarse				

Double Inverted Cone

CAT. No.	M51
Head Size	009
Total Length	18.5
MEDIUM	
COARSE	

Barrel

CAT. No.	M4
Head Size	033
Total Length	21
MEDIUM	
COARSE	

Round Wheel

CAT. No.	G5
Head Size	054
Total Length	19
COARSE	

Gross Reduction

CAT. No.	BR4	BR1
Head Size	018	018
Total Length	20.5	19.5
COARSE		

Pear

CAT. No.	AB2	AB1
Head Size	014	016
Total Length	17	19
MEDIUM		
COARSE		

Flat End Cylinder

CAT. No.	D0	D1	D31	D2	D3	D4
Head Size	008	010	012	016	012	014
Total Length	19	19	19.5	19	21	21
MEDIUM						
COARSE						

Flat End Taper

CAT. No.	Z31	Z32	Z3	F11	F1	F2	F21
Head Size	014	016	018	018	016	018	021
Total Length	18.5	18.5	19.5	21	22.5	22.5	22
MEDIUM							
COARSE							

Pointed Cylinder

CAT. No.	E7	PR3
Head Size	012	014
Total Length	21	23
MEDIUM		
COARSE		

Round End Cylinder

CAT. No.	D4R
Head Size	014
Total Length	21
MEDIUM	
COARSE	

Round End Taper

CAT. No.	C1	C2	C3	C11	M31	M3	F1R	F1RSP	F11R	F2R	C4
Head Size	014	018	025	012	016	018	021	016	018	018	021
Total Length	21	19	21	21	21	21	22.5	22	21.5	22.5	22
MEDIUM											
COARSE											
SUPER Coarse											

Torpedo Cylinder

CAT. No.	S11	PR1
Head Size	012	012
Total Length	21	23.5
FINE		
MEDIUM		
COARSE		

Needle

CAT. No.	PR15	E1	E2	Z11	E31	E3
Head Size	012	012	014	012	014	016
Total Length	20	19	21	22.5	23	23.5
SUPER FINE						
FINE						
MEDIUM						
COARSE						
SUPER Coarse						

Flame

CAT. No.	E4	E51	E5	E6	E9
Head Size	013	012	014	016	023
Total Length	19	21	21	22	21
XXF					
SUPER FINE					
FINE					
MEDIUM					

Gingival Curettage

CAT. No.	S1	S4	S2	T3	T1
Head Size	016	014	014	018	020
Total Length	22	19	21	20	20.5
MEDIUM					
COARSE					




Bur & Abrasives

Anaesthetics
Articulating
Bur and Abrasives
Cosmetics
Diagnostic
Disposables
Endodontics
Equipment
Handpieces
Implant
Infection Control
Instruments
Isolation
Laboratory
Oral Surgery
Orthodontics
Patient Education and Gift
Pedodontics
Periodontal
Post and Core
Preventive
Prosthodontics
Restorative
Uniform

Carbide Burs


HENRY SCHEIN®
Metal Cutting Carbide Burs / Henry Schein



10's/pack (HSPL-)

Bur No.	1558 Turbo	1557 Turbo	331 Turbo	330 Turbo	245 Turbo	557 Turbo
Item No.	900-4174	900-4175	900-4176	900-4177	900-4178	900-4179

Endo Safe-End Access Bur
/ Henry Schein
- Opens up the pulp chamber after initial penetration



length	Packing	Item no.
21mm	5's/pack	HSPL-900-7038

HENRY SCHEIN® Maxima® Carbide Burs / Henry Schein

High-Speed Operative Burs

	Round			Inverted Cone		Amalgam Prep	Pear			Flat Fissure (Plain Cut)	
10's/pack (HSPL-)											
US No.	1/2	2	4	33-1/2	35	245	330	330 SS	331	56	57
FG (19mm)	006	010	014	006	010	009	008	008	010	008	010
Item No.	900-4508	900-4509	900-4510	900-4511	900-4512	900-4515	900-4516	900-6922	900-4517	900-4513	900-4514

	Flat Fissure (Cross Cut)			Taper Fissure (Cross Cut)	Domed Fissure		
10's/pack (HSPL-)							
US No.	556	557	558	700	1157	1158	1557
FG (19mm)	008	010	012	010	010	012	010
Item No.	900-4518	900-4519	900-4520	900-4521	900-4522	900-4523	900-4524

Low-Speed Operative Burs

US No.	4	8
RA (22mm)	014	023
Item No.	900-6923	900-6924

HENRY SCHEIN® 12 Blade Trimming & Finishing Burs / Henry Schein

- Perfect for occlusal, lingual, labial and interproximal finishing

	Round		Flame			Short Flame	Long Taper			Egg		
5's/pack (HSPL-)												
US No.	7006	7008	7104	7106	7108	274	7204	7205	7206	7404	7406	7408
Fig. No.	H-41	H-41	H-46	H-46	H-46	H-274	H-378	H-378	H-378	H-379	H-379	H-379
Diameter	018	023	014	018	023	016	014	016	018	014	018	023
Head L (mm)	-	-	3.5	3.7	4.5	3.5	8.5	8.5	8.5	3.3	3.7	3.9
Item No.	900-4355	900-4356	900-4357	900-4358	900-4359	900-4386	900-4360	900-4361	900-4362	900-4363	900-4364	900-4365

	Cone	Taper Pointed				Taper Round End		Bullet			Needle		
5's/pack (HSPL-)													
US No.	7606	7610	7611	7612	7613	7642	7664	7801	7802	7803	7901	7902	7903
Fig. No.	H-50	H-132	H-133	H-134	H-135	H-375R	H-375R	H-247	H-247	H-247	H-246	H-246	H-246
Diameter	018	008	010	014	014	010	014	009	010	012	009	010	012
Head L (mm)	5.1	3.0	4.2	6.0	9.0	7.0	8.2	4.0	3.4	3.4	3.8	3.8	3.9
Item No.	900-4366	900-4367	900-4368	900-4369	900-4370	900-4384	900-4385	900-4371	900-4372	900-4373	900-4374	900-4375	900-4376

Fissure Carbide Bur HP/FG / Shofu

- The special three-blade design is highly abrasive with superior cutting efficiency to provide precise replication of intricate surface details
- Made of tungsten carbide, these carbide burs assure remarkable durability
- **The Fissure Carbide in HP Shank** is used for carving of indirect composites (eg. Ceramage) as well as acrylic dentures
- **The Fissure Carbide in FG Shank** for use with Lab Air Z is preferred by ceramists for anatomical carving of grooves and fissures in porcelain restorations



Bur No.	HP	FG
Packing	6's/box	3's/box
Item No.	SF FC-HP	SF FC-FG




Bur & Abrasives

Carbide Burs

Strauss Carbide Dental Burs / Strauss

- CNC (Computer Numerical Control) precision manufacturing produces a top quality rotary instrument product that exceeds industry standards for quality and precision
- Highest quality raw materials: One-piece fine-grain tungsten carbide design ground to shape resulting in zero vibration, superior durability, optimal cutting efficiency, and the ability to withstand repeated sterilization without rusting
- The most concentric bur ever produced - consistent in size resulting in less chatter, superior control, efficient cutting and finish
- Task and material specific blade configurations produce optimal cutting efficient finishing and outstanding durability
- Comprehensive range of ideal anatomical shapes and sizes - the right instrument for every situation




(ST-)

Item No.	FG6	FG330	FG1157	FG1158
Head Size	018	008	010	012
Total Length	-	1.6	4.2	4.2
Short Shank 16mm		■		
Standard Shank 19mm	■	■	■	■
Surgical Length 24mm	■			■

Packing	Item no.
10's/pack	ST-(See the table)

Surgical Burs



HENRY SCHEIN®
Surgery Burs
- FG Carbide 25mm
/ Henry Schein

US No.	2	57	170	558
Packing	5's/pack			
Item No.	HS-100-1344	HS-100-4394	HS-100-1226	HS-100-0282

Endo Surgical Burs / Henry Schein

- Transversing cutting tip & 6 helical cutting edges
- Ideal for sectioning wisdom tooth prior to extraction

US No.	Packing	Item no.
151	5's/pack	HS-900-7039

US No.	Packing	Item no.
150L	5's/pack	HS-900-7040


Oral Surgery Sterile Burs - HP Carbide / Henry Schein

- For efficient bone removal and sectioning
- Single-use, gamma sterilized; individually wrapped

Round		Taper Fissure (Cross Cut)	
US No.	8	US No.	702
Diameter	023	Diameter	016
Packing	10's/pack	Packing	10's/pack
Shank 2 (51mm)	HSPL-112-5597	Shank 1 (44.5mm)	HSPL-112-5592
		Shank 2 (51mm)	HSPL-112-5598

- **Shank 1 (44.5mm)**
- Fits most standard handpieces
- **Shank 2 (51mm)**
- Fits Hall® handpieces - for medium bur guards

Acrylic Burs



HENRY SCHEIN®
Laboratory Carbide Burs, Regular Cut
/ Henry Schein

- For trimming, carving & finishing of acrylic
- Non-clogging
- All with 3/32" HP shank

US No.	84T	61A	54D	41P	51A
Packing	Each				
Item No.	HS-100-9680	HS-100-6866	HS-101-2194	HS-101-2733	HS-100-7866



Bur & Abrasives

Acrylic Burs

Acrypoint - HP / Shofu

✓ Designed for adjusting, finishing and polishing of dentures and other acrylic appliances to a high lustre



- Available in 3 grits:

- ✓ Coarse (Dark-Grey) for adjustment
- ✓ Medium (Brown) for finishing
- ✓ Fine (Light-Grey) for polishing



	Packing	BP1	PC2
ISO No.		100	050
Coarse (Dark-Gray)	6's/box	SF-0426	
	12's/box		SF-0423
Medium (Brown)	6's/box	SF-0427	
	12's/box		SF-0424
Fine (Light-Gray)	6's/box	SF-0428	
	12's/box		SF-0425

Lab Carbide Bur - HP / Shofu

Enhanced Cutting Efficiency with Remarkable Durability



✓ New range of carbide burs in different shapes to suit a variety of applications such as fast chairside or bench side adjustment of acrylic resins, chromecobalt alloys, chrome-nickel alloys and precious metal alloys

- Available in 30 shapes:

- For gross reduction of plaster
 - C - Coarse
- For trimming and initial finishing
 - NC - Normal Coarse
 - N - Normal
 - B - Basic Plain Cut
- For anatomical contouring and finishing
 - NF - Normal Fine Cut
 - DF - Diamond Fine Cut

Bur No.	1N	3N	5N	7C	7NC	7N	7B	7NF	7DF	8C	8NC	9N
Bur No.	13N	13B	14N	14B	14NF	14DF	15NF	15DF	16N	17N	17B	
Bur No.	21N	21DF	23NF	24DF	37N	38N	39NF					
								Packing	Item no.			
								Each	SF-LB-(Bur No.)			



Anaesthetics
Articulating
Bur and Abrasives
Cosmetics
Diagnostic
Disposables
Endodontics
Equipment
Handpieces
Implant
Infection Control
Instruments
Isolation
Laboratory
Oral Surgery
Orthodontics
Patient Education and Gift
Pedodontics
Periodontal
Post and Core
Preventive
Prosthodontics
Restorative
Uniform

Bur & Abrasives

- Anaesthetics
- Articulating
- Bur and Abrasives**
- Cosmetics
- Diagnostic
- Disposables
- Endodontics
- Equipment
- Handpieces
- Implant
- Infection Control
- Instruments
- Isolation
- Laboratory
- Oral Surgery
- Orthodontics
- Patient Education and Gift
- Pedodontics
- Periodontal
- Post and Core
- Preventive
- Prosthodontics
- Restorative
- Uniform

Finishing Abrasives

Finishing Stones

Dura-Green Stones / Shofu



- Made of high quality silicon carbide grit, Dura Green stones are durable concentric with excellent cutting efficiency
- Dura Green Stones are used for fast contouring and finishing of porcelain, composite, precious alloys and amalgam

- Available in HP and CA shank and selected number of shapes in FG shank

Packing	Item no.
12's/box	See the table

Shape	CN1	CN2	CN3	CY2	CY3	CY4	FL1	FL2	FL3
HP	SF-0002			SF-0006	SF-0007	SF-0008	SF-0009	SF-0010	SF-0011
CA	SF-0052	SF-0053	SF-0054		SF-0057		SF-0059	SF-0060	SF-0061
FG	SF-0101			SF-0103				SF-0104	

Shape	FL4	FL5	IC1	IC2	IC3	IC4	IC5	IC6	IC7
HP	SF-0012	SF-0013	SF-0014	SF-0015	SF-0016	SF-0017	SF-0018	SF-0019	SF-0020
CA	SF-0062		SF-0064		SF-0066		SF-0068		
FG				SF-0105					

Shape	IC8	IC9	KN1	KN3	KN4	KN5	KN6	KN7	PC1
HP	SF-0021	SF-0022	SF-0023	SF-0025	SF-0026	SF-0027	SF-0028	SF-0029	SF-0030
CA									SF-0080

Shape	PC2	RD1	RD2	RD3	RD4	RE1	RE2	TC1	TC2
HP	SF-0031	SF-0032	SF-0033	SF-0034	SF-0035	SF-0036	SF-0037	SF-0039	SF-0040
CA	SF-0081	SF-0082	SF-0083	SF-0084	SF-0085	SF-0086		SF-0089	SF-0090
FG		SF-0107	SF-0108					SF-0109	

Shape	TC3	TC4	WH1	WH2	WH3	WH4	WH5	WH6
HP	SF-0041	SF-0042	SF-0043	SF-0044	SF-0045	SF-0046	SF-0047	SF-0048
CA		SF-0092			SF-0095	SF-0096		SF-0098
FG			SF-0110					

Dura-White Stones / Shofu



- Ideal for finishing enamel, glass ionomer cement, composite, porcelain and prepared tooth abutment
- These micro-grained aluminium oxide abrasives are perfectly balanced to provide vibration-free performance

- Available in HP / CA / FG shank

Packing	Item no.
12's/box	See the table

Shape	CN1	CY2	FL1	FL2	FL3	IC2	IC4	KN1
HP	SF-0201		SF-0202	SF-0203	SF-0204		SF-0205	
CA	SF-0221		SF-0222	SF-0223	SF-0224		SF-0225	
FG	SF-0241	SF-0243	SF-0240	SF-0244		SF-0245		SF-0246

Shape	KN3	RD1	RD2	RE1	TC1	TC4	WH1	WH4
HP	SF-0206	SF-0207	SF-0208	SF-0209	SF-0210	SF-0211		SF-0212
CA		SF-0227	SF-0228	SF-0229	SF-0230			SF-0232
FG		SF-0247	SF-0248		SF-0249		SF-0250	





Bur & Abrasives



Polishing Abrasives

Metal & Amalgam Polishing

Brownie, Greenie & SuperGreenie / Shofu



- These popular silicone polishers are impregnated with high quality polishing ingredients for fast, brilliant polishing of amalgam, gold and cast alloys
- Autoclavable or can be disinfected in cold sterilization and doesn't require any polishing paste or pumice

- **Brownies** - for achieving a smooth surface
- **Greenies** - for a lustrous polish
- **Supergreenies (with yellow band)** - for a final superpolish

- Mini-points are recommended for polishing occlusal surfaces and Cups for proximal and occlusal areas

Packing	Item no.
12's/box	See the table



	Shape	Cup	Mini-Point	KN7	PC2	WH6	Floppie	Disk	Unmounted Wheel
Brownie	HP			SF-0271	SF-0272	SF-0273	SF-0458		SF-0453
	CA	SF-0401	SF-0403	SF-0281	SF-0282	SF-0283		SF-0405	
	FG	SF-0411	SF-0413						
Greenie	HP			SF-0274	SF-0275	SF-0276	SF-0459		SF-0454
	CA	SF-0402	SF-0404	SF-0284	SF-0285	SF-0286		SF-0406	
	FG	SF-0412	SF-0414						
Supergreenie (Super-polish)	HP			SF-0277	SF-0278	SF-0279			
	CA	SF-402A	SF-404B	SF-0287	SF-0288	SF-0289			SF-454A
	FG	SF-412A	SF-414B						

Composite Polishing

CompoMaster Coarse & CompoMaster / Shofu



- Diamond impregnated polishers designed with a unique blend and high density of natural diamond particles to provide an exceptionally smooth surface on any composite material

- **CompoMaster Coarse (Yellow Polisher)** - for finishing & basic polishing
- **CompoMaster (Yellow & White Bands)** - for superpolishing

- Provides a simple, fast and efficient super polishing for intra oral, chair side or laboratory use

Packing	Item no.
3's/box	See the table



	Shape	Bullet	Knife	Cup	Mini-Point	KN7	WH6	Unmounted Wheel
CompoMaster (Yellow & White Bands)	CA	SF-0111	SF-0112	SF-0113	SF-0114			
	HP	SF-0116			SF-0119	SF-0117	SF-0118	
CompoMaster Coarse (Yellow Polisher)	CA	SF-111C	SF-112C	SF-113C	SF-114C			
	HP	SF-116C			SF-119C	SF-117C	SF-118C	SF-131C



Bur & Abrasives

Polishing Abrasives

Composite Polishing

One-Gloss / Shofu



- ✓ **2-in-1 finishing and polishing of all tooth coloured direct restorations**
- ✓ **Excellent for:**
 - a) Polishing of enamel after tooth scaling
 - b) Removal of orthodontic cement and polishing after debanding and bracket removal
 - c) Surface stain removal
- Easy to use, replaceable and disposable snap-on system
- For polishing Composite, apply heavy pressure for finishing and feather light pressure for polishing
- Saves time and cost as it both finishes and polishes
- Removes excess resin or cement without marring the enamel

Silicone Midi / Shofu
For finishing and polishing composite restorations

Packing	Item no.
25's/box	SF-SILICONE MIDI

One-Gloss Complete Set

Packing	Item no.
- 20 pcs each of Unmounted Cups, Midi Points and Invented Cones	SF-0180
- 3 pcs of Mandrels CA shank	

One-Gloss Individual Pack

Shape	Cup	Midi-Point	Inverted Cone	Mini-Point	Mandrel
Packing	50 pcs	50 pcs	50 pcs	12 pcs	6 pcs
Item No.	SF-0181	SF-0182	SF-0183	SF-0186	SF-0435

Composite & Composite Fine / Shofu



- These polishers are perfect for finishing glass ionomer and hybrid composites without scratching or chipping the surrounding enamel
- **CompoSite Polishers** - for finishing & pre-polishing, aluminum oxide impregnated
- **CompoSite Fine Polishers (White Band)** - for a lustrous polish, zirconium dioxide impregnated
- Cups and Disks are best for polishing subgingival and interproximal surfaces, while Pear-Point and Midi-Point are ideal for polishing occlusal and lingual surfaces

Packing	Item no.
12's/box	See the table

Shape	Bullet	Knife	Wheel	Disk	Cup	Pear-Point	Inverted Cone	Midi-Point	
CompoSite Polishers	CA	SF-0191	SF-0192	SF-0193	SF-0194	SF-0195	SF-0196	SF-0197	SF-0198
	FG					SF-196B			SF-0199
CompoSite Fine Polishers (White Banded)	CA	SF-0291	SF-0292	SF-0293	SF-0294	SF-0295	SF-0296	SF-0297	SF-0298
	FG					SF-296B			SF-0299

Super-Snap System / Shofu

✓ *Super-Snap is a comprehensive polishing system for all direct aesthetic restorative material*



- Super-snap disks have soft, safe, plastic snap-on mounts for mandrels
- Without metal centres or exposed mandrels to avoid discoloration or gouging of composite restorations
- Improved mandrel design for maximum control
- The disks are available in 4 different colour coded grits in 12mm and 8mm diameter
- Mini disks provide perfect control while polishing cervical areas

Super-Snap Refills

Super-Snap® Disks	Contouring		Finishing				Polishing		Super Polishing	
	BLACK (Coarse)		VIOLET (Medium)		GREEN (Fine)		RED (Superfine)			
(50's/box)	L506	L507	L508	L509	L528	L522	L501	L519	L502	L521
Super-Snap® Mandrel	CA-0439		Super-Snap® Polystrips™		BLACK / VIOLET (Coarse / Medium)		GREEN / RED (Fine / Superfine)			
					L525		L526			


Super-Snap X-Treme Refills

Super-Snap® X-Treme Disks	Polishing		Super Polishing	
	GREEN (Fine)		RED (Superfine)	
(50's/box)	Standard	Mini	Standard	Mini
	L511	L531	L512	L532

Super-Snap X-Treme



Bur & Abrasives



Polishing Abrasives

Porcelain Polishing

Ceramaster Coarse & Ceramaster / Shofu



- Diamond impregnated polishers designed with a unique blend of silicone and even distribution of diamond particles make these instruments the best choice for polishing porcelain and ceramic polymer to a glazed-like surface
- A simple 2 step finishing and polishing system, ideal for intra-oral, chairside or laboratory use

- **CeraMaster Coarse (Blue Polisher)** - for finishing/polishing and adjustment margins of alumina copings
- **CeraMaster (Blue and White Bands)** - for superpolishing

- Durable and cost-effective

Packing	Item no.
3's/box	See the table

Shape		Bullet	Knife	Cup	Mini-Point	KN7	WH6	Unmounted Wheel
CeraMaster (Blue & White Band)	CA	SF-0121	SF-0122	SF-0123	SF-0124			
	HP	SF-0126			SF-0129	SF-0127	SF-0128	
CeraMaster Coarse (Blue Band)	CA	SF-121C	SF-122C	SF-123C	SF-124C			
	HP	SF-126C			SF-129C	SF-127C	SF-128C	SF-132C

Ceramisté / Shofu



- Designed for polishing porcelain and enamel to achieve a glazed-like polish

- **Standard Polishers** - for pre-polishing
- **Ultra (Yellow Band)** - for polishing
- **Ultra II (White Band)** - for superpolishing to a glazed-like surface

Packing	Item no.
12's/box	See the table

Shape		Bullet	Cup	Knife	Wheel	Midi-Point	KN7	PC2	WH6	Unmounted Wheel
Standard	HP						SF-0261	SF-0262	SF-0263	SF-0451
	CA	SF-0251	SF-251A	SF-0252	SF-0253	SF-253A				
	FG		SF-251B			SF-253B				
Ultra (Yellow Band)	HP						SF-0264	SF-0265	SF-0266	SF-0452
	CA	SF-0254	SF-254A	SF-0255	SF-0256	SF-256A				
	FG		SF-254B			SF-256B				
Ultra II (White Band)	HP						SF-0267	SF-0268	SF-0269	SF-452A
	CA	SF-0257	SF-257A	SF-0258	SF-0259	SF-259A				
	FG		SF-257B			SF-259B				

SoftCut - HP / Shofu

- ✓ *Created for contouring and polishing porcelain an indirect composite restorations*
- ✓ *Recommended for polishing porcelain surfaces prior to glazing to lower the firing temperature and enhance aesthetics*



Packing	Item no.
12's/box	See the table

Shape	PC2	WH6	Unmounted Wheel
SoftCut PA (Green Banded)		SF-0564	SF-0595
SoftCut PB (Blue Banded)	SF-0563	SF-0565	SF-0596




Bur & Abrasives

Abrasive Strips

Sof-Lex™ Finishing Strips / 3M ESPE

✓ *Soft-Lex™ Finishing Strips are designed to produce a smooth restoration surface to help prevent interproximal plaque buildup for superior finishing and polishing results*



Indications

- Finishing and polishing composite, ceramic and glass ionomer dental restorations

Product Benefits

- Easy-to-use system allows you to create smooth proximal surfaces
- Strips are gapped at their centers for easy interproximal insertion
- Precision-coated with an aluminum-oxide grit to produce a smooth restoration surface that helps prevent interproximal plaque buildup
- Unique narrow strip is available for easier use in short interproximal areas

Grit Sizes	Coarse / Medium	Coarse / Medium-Narrow	Fine / Superfine
Size	7" x 5/32"	7" x 5/64"	7" x 5/32"
Packing	150's/pack	100's/pack	120's/pack
Item No.	3M-1954	3M-1954N	3M-1956

Polishing Strips



HENRY SCHEIN® Maestro™ Finishing and Polishing Strips

/ Henry Schein

Made of
Aluminium
Oxide



Size	Style	Packing	Fine / Superfine
Coarse / Medium	Regular	100's/pack	HS-112-5696
	Narrow		HS-112-5697
Fine / Superfine	Regular	100's/pack	HS-112-5698
	Narrow		HS-112-5699

Polishing Paste



HENRY SCHEIN® Diamond Polishing Paste

/ Henry Schein

- Actual diamond particles produce a high-gloss lustrous finish on porcelain, composites, gold and other metals



- High viscosity for easy handling

Packing	Item no.
2gm/syringe	HS-101-4116

Dura-Polish Dia / Shofu

- Dura-Polish DIA is a diamond impregnated polishing wax for super polishing of indirect composite and porcelain



Packing	Item no.
5gm paste/jar	SF-0554

Ultradent® Diamond Polish Mint

/ Ultradent



- High-grade white microcrystalline diamond particles
- Unsurpassed esthetic polish
- Ideal for porcelain or composite restorations
- For use with Jiffy brushes and felt wheels

Microns	Packing	Item No.
0.5µm	2 x Syringes (1.2mL)	UL-5540
1.0µm		UL-5541

Brinell L High-Lustre Polishing Paste

/ Renfert

- Impregnated with diamond powder for precious metal bonding alloys and porcelain
- The finest diamond particles achieve the optimum high-lustre shine
- Sparing application via the measuring syringe
- Efficient take-up using felt wheels
- Excellent high-lustre using cotton buffs



Packing	Item no.
- 1 x Syringe (5mL) with 1 Cotton Buff - 1 x Felt Wheel - 1 x Mandrel	RF-519-0001



Bur & Abrasives

Polishing Cups, Points, Wheels & Rubbers

Anaesthetics
 Articulating
 Bur and Abrasives
 Cosmetics
 Diagnostic
 Disposables
 Endodontics
 Equipment
 Handpieces
 Implant
 Infection Control
 Instruments
 Isolation
 Laboratory
 Oral Surgery
 Orthodontics
 Patient Education and Gift
 Pedodontics
 Periodontal
 Post and Core
 Preventive
 Prosthodontics
 Restorative
 Uniform



HENRY SCHEIN®



Maestro™ Amalgam/Gold Polishers

/ Henry Schein

Amalgam/Gold Polishing Kit

Packing	Item no.
- 12 pcs (1 each Mini-Point, Point, Cup, Bullet in 3 different grits)	HS-112-5745

Refill

Packing	Item no.
12's/pack	See the table

Amalgam/Gold Polishers Shaping - Coarse



Amalgam/Gold Polishers Smoothing - Medium



Amalgam/Gold Polishers High Shine - Fine



RA	Shape	Mini-Point	Point	Cup	Bullet	Mini-Point	Point	Cup	Bullet	Mini-Point	Point	Cup	Bullet
	Item no. (HS-)		112-5732	112-5733	112-5734	112-5735	112-5736	112-5737	112-5738	112-5739	112-5740	112-5741	112-5743



Maestro™ Porcelain Polishers

/ Henry Schein

Porcelain Polishing Kit

Packing	Item no.
- 12 pcs (1 each Small Point, Point, Cup, Disc in 3 different grits. Polishing of ceramic, metal and zirconia)	HS-112-5731

Refill

Packing	Item no.
3's/pack	See the table

Porcelain Polishers Shaping - Coarse



Porcelain Polishers Smoothing - Medium



Porcelain Polishers High Shine - Fine



Shape	Wheel	Small Point	Point	Cup	Disc	Wheel	Small Point	Point	Cup	Disc	Wheel	Small Point	Point	Cup	Disc	
	Item no. (HS-)		112-5715	112-5716	112-5717	112-5719	112-5720	112-5721	112-5722	112-5723	112-5724	112-5725	112-5726	112-5727	112-5728	112-5729



Maestro™ Two Step Composite Polishing System

/ Henry Schein

Two Step Polishing Kit

Packing	Item no.
- 8 pcs (1 each Small Point, Point, Cup, and Knife Edge in both coarse and fine grits)	HS-112-5768

Refill

Packing	Item no.
3's/pack	See the table

Coarse Pre-Polishers



Fine Pre-Polishers



Shape	Small Point	Point	Cup	Knife Edge	Small Point	Point	Cup	Knife Edge
Item no. (HS-)	112-5760	112-5761	112-5762	112-5763	112-5764	112-5765	112-5766	112-5767



Bur & Abrasives

Polishing Cups, Points, Wheels & Rubbers

HENRY SCHEIN®
Denture Adjustment Polishing System
/ Henry Schein

HP	Shape	0666	0667	0668
	Item No. (HS-)	900-2122	900-2123	900-2124

■ Denture Adjustment Polishing Kit

Packing	Item no.
- 5 pcs HP Polishers in green (coarse), gray (medium), and beige (fine) (3 x #0666, 1 x #0667, 1 x #0668 & 1 x TC cutter)	HS-900-2104

■ Refill

Packing	Item no.
12's/pack	See the table

Prophy Cups Webbed & Ribbed
/ Henry Schein

Description	Latch Type - Firm White (#S20-02)
Packing	36's/pack
Item No.	HS-100-5906

Maestro™ Composite Discs and Strips / Henry Schein

■ Complete Disc/Strip Kit

Packing	Item no.
- 200 Discs (25 of each: 3/8" & 1/2" in Superfine, Fine, Medium & Coarse) - 50 Strips (25 of each: Regular Coarse/Medium & Regular Fine/Superfine) - 2 RA Mandrels	HS-112-5686

■ Disc Refill

Packing: 100's/pack

Size	Color	Item no.
3/8"	Coarse	White HS-112-5690
	Medium	Green HS-112-5689
	Fine	Blue HS-112-5688
	Superfine	Brown HS-112-5687
1/2"	Coarse	White HS-112-5695
	Medium	Green HS-112-5694
	Fine	Blue HS-112-5693
	Superfine	Brown HS-112-5692

■ Mandrels RA

Packing	Item no.
3's/pack	HS-112-5704

Sof-Lex™ Spiral Finishing and Polishing Wheels / 3M ESPE

✓ *Polish anterior and posterior restorations with one versatile tool, rather than switch between points, cups, discs and brushes*

3M ESPE Sof-Lex
Introductory Kit - Spiral Wheels

CONTENTS: 20 RA

Ideal for use with:

- Composites
- Resin-modified glass ionomers
- Bisacrylic temporary materials
- Lava™ Ultimate CAD/CAM Restorative
- Precious and semi-precious metal

Spiral Polishing Wheels

Spiral Finish Wheels

■ Sof-Lex™ Spiral Finishing Wheels, Fine

Packing	Item no.
- 12 x Spiral Finishing Wheels (Beige)	3M-5080

■ Sof-Lex™ Spiral Polishing Wheels

Packing	Item no.
- 12 x Spiral Polishing Wheels (White)	3M-5081

■ Sof-Lex™ Introductory Kit - Spiral Wheels

Packing	Item no.
- 6 x Spiral Finishing Wheels (Beige) - 6 x Spiral Polishing Wheels (White) - 1 x RA Mandrel - 1 x Step-by-Step Card	3M-5082I


■ Sof-Lex™ System Kit - Spiral Wheels, Discs & Strips

Packing	Item no.
- 10 x Spiral Finishing Wheels (Beige) - 10 x Spiral Polishing Wheels (White) - 30 x Sof-Lex™ XT Polishing Discs 3/4" (C, M, F, SF) - 30 x Sof-Lex™ XT Polishing Discs 1/2" (C, M, F, SF) - 15 x Sof-Lex™ Finish Strips (Coarse/Medium) - 1 x RA Mandrel	3M-5082K





Bur & Abrasives



Polishing Discs

Sof-Lex™ Finishing and Polishing System / 3M ESPE

- Easy-to-use system allows you to create high-luster restorations that offer lifelike esthetics
- Discs are color coded to allow for easy identification of grit sequence
- Pop-on disc technology makes changing discs easier to save you time
- Pop-on discs conform to the contours of the teeth for easier finishing and polishing
- **Extra-thin (XT) discs** allow for easier access to interproximal surfaces while offering a slightly stiffer disc

■ Sof-Lex™ Contouring and Polishing Discs



■ Sof-Lex™ Extra-Thin Contouring and Polishing Discs



◆ Intro Kit

Packing	Item no.
- 30 each of Discs in Coarse, Medium, Fine and Superfine Grits in Both 9.5mm and 12.7mm Diameters - 1 x RA Mandrel	3M-4851

◆ Refills (packing: 30's/box)

Discs	Coarse	Medium	Fine	Superfine
Ø9.5mm	3M-4850C	3M-4850M	3M-4850F	3M-4850SF
Ø12.7mm	3M-4851C	3M-4851M	3M-4851F	3M-4851SF

◆ Intro Kit

Packing	Item no.
- 30 each of XT Discs in Coarse, Medium, Fine and Superfine Grits in Both 9.5mm and 12.7mm Diameters - 1 x RA Mandrel	3M-4931

◆ Refills (packing: 30's/box)



XT Discs	Coarse	Medium	Fine	Superfine
Ø9.5mm	3M-4930C	3M-4930M	3M-4930F	3M-4930SF
Ø12.7mm	3M-4931C	3M-4931M	3M-4931F	3M-4931SF

Cutting Discs

Separating Discs for Porcelain, BIFLEX, 22mm

/ Renfert

- Highly flexible separating disc for porcelain with coating on both sides
- Fine diamond coating for fine separation and contouring



Packing	Item no.
Each	RF-27-1000





Bur & Abrasives

Bur Kits

Porcelain Adjustment Kit

/ Shofu

✓ For chairside adjustment and polishing of porcelain restorations to a glaze like finish



- Saves time and avoids the tedious re-finishing process

- Available in HP Shank

	Packing	Item no.
	- 3 pcs x Dura White Stones for Adjusting (1 each of CN1, FL3, RE1)	SF-0301
	- 9 pcs x Ceramiste Points (1 each of KN7, PC2, WH6) in Standard, Ultra and Ultra II for polishing	

Composite Finishing Kit

/ Shofu

✓ Designed for contouring and finishing aesthetic restorations



- Available in CA / FG Shank

	Packing	Item no.
FG	- 6 pcs x Dura Green Stones for Adjusting (3 each of CN1, FL2)	SF-0302
CA	- 6 pcs x White Stones for Finishing (3 each of CN1, FL2)	SF-0306

Enamel Adjustment Kit

/ Shofu

✓ Designed for cosmetic contouring of natural dentition or porcelain and super polishing the adjusted surfaces



- Available in CA Shank

	Packing	Item no.
	- 3 pcs x Dura White Stones for Adjusting & Finishing (1 each of TC1, FL2, RD1)	SF-0307
	- 9 pcs x Ceramiste Points (1 each of KN7, PC2, WH6) in Standard, Ultra and Ultra II for polishing	

Amalgam Polishing Kit

/ Shofu

✓ An ideal assortment of abrasives impregnated with polishing compounds to provide fast and brilliant results on Amalgam restorations without pumice or tin oxide



- Available in CA / FG Shank

	Packing	Item no.
FG	- 3 each Brownie Cups & Mini-Points for Pre-polishing	SF-0309
	- 2 each of Greenie Cups & Mini-Points for Polishing	
CA	- 1 each Supergreenie Cup & Mini-Point for Super Polishing	SF-0308

Composite Polishing Kit

/ Shofu

✓ A comprehensive assortment for finishing and polishing of composites and glass ionomers



- Available in CA Shank

	Packing	Item no.
	- 2 pcs x Dura Green Stones CA for reduction of excess material (1 each of CN1 & FL2)	SF-0310
	- 2 pcs x Dura White Stones CA for finishing (1 each of CN1 & FL2)	
	- 6 pcs x Composite Polishers CA for Pre-Polishing (2 each of Bullet, Knife, Wheel)	
	- 2 pcs x Composite Fine Polishers CA for Super-Polishing (2 each of Cup, Disk) without damaging enamel	

Acrylic Contouring & Finishing Kit

/ Shofu

✓ Designed for adjustment and polishing of acrylic dentures and appliances chairside or in the laboratory



- Available in HP Shank

	Packing	Item no.
	- 4 pcs x Cross-Cut Carbide Burs for smooth and fast adjustments of acrylic materials	SF-0329
	- 8 pcs x AcryPoints for finishing and polishing (Coarse grit x 4, Medium grit x 2, Fine grit x 2)	



Bur & Abrasives

Bur Kits

Super-Snap System / Shofu

✓ *Super-Snap is a comprehensive polishing system for all direct aesthetic restorative material*

- Super-snap disks have soft, safe, plastic snap-on mounts for mandrels
- Without metal centres or exposed mandrels to avoid discoloration or gouging of composite restorations
- Improved mandrel design for maximum control
- The disks are available in 4 different colour coded grits in 12mm and 8mm diameter
- Mini disks provide perfect control while polishing cervical areas



Super-Snap

Super-Snap® Disks	Contouring		Finishing				Polishing		Super Polishing	
	BLACK (Coarse)		VIOLET (Medium)				GREEN (Fine)		RED (Superfine)	
	L506	L507	L508	L509	L528	L522	L501	L519	L502	L521

Super-Snap® Mandrel		Super-Snap® Polystrips™	BLACK / VIOLET (Coarse / Medium)		GREEN / RED (Fine / Superfine)	
			L525		L526	
	CA-0439		L525		L526	

◆ Super-Snap Rainbow Technique Kit

Packing	Item no.
- 2 pcs x Dura White Stones for finishing	SF-0500
- 1 pc x Composite Fine Midi Point for gross polishing	
- 180 pcs x Super-Snap Disks	
- 4 x Super-Snap Mandrels CA	
- 40 pcs x Super-Snap Polystrips	

◆ Super-Snap Mini Kit

Packing	Item no.
- 48 x Super-Snap Disks	SF-0505
- 1 x Composite Fine-Point	
- 1 x Super-Snap Mandrel CA	

Super-Snap X-Treme

Super-Snap® X-Treme Disks	Polishing		Super Polishing	
	GREEN (Fine)		RED (Superfine)	
	L511	L531	L512	L532



◆ Super-Snap X-Treme Standard Disk Kit

Packing	Item no.
- 100 x Standard X-Treme Disks (50 x Green Disks; 50 x Red Disks)	SF-0506

◆ Super-Snap X-Treme Mini Disk Kit

Packing	Item no.
- 100 x Mini X-Treme Disks (50 x Green Disks; 50 x Red Disks)	SF-0507

SuperBuff Kit / Shofu

✓ *Superpolishing of direct and indirect composite restorations*

- Aluminium oxide impregnated felt disk with an easy to use snap-on system
- Just wet with water to activate it before use
- Disposable polishers prevent cross-contamination
- Thick disk design to super-polish pits and fissures
- Fast and efficient in polishing multiple surfaces in a single bond



Packing	Item no.
- 25 x Polishers	SF-0535
- 2 x Mandrels CA	

TPE Diamond Kit / Shofu

✓ *Specially developed for preparing porcelain shoulder margins*

- Bevelled at cutting edge to prevent laceration of adjacent tissue or undercut the axial walls when preparing the shoulder

- TPE burs come in 2 grits:

- **Coarse grit:** recommended for rapid reduction of the shoulder margin into the gingival sulcus
- **Fine grit (red-banded shank):** recommended for smooth finishing of shoulder margin, particularly effective when finishing subgingivally
- **Short Shank TPE Point (in Coarse grit only):** recommended when accessibility is limited



Packing	Item no.
9 x Selected Points & Instructions per Kit	SF-0688

	Regular Shank				Short Shank		
Length (mm)	1.0	1.2	1.4	1.6	1.2	1.4	
COARSE	680C	682C	684C	686C	682S	684S	
FINE	680F	682F	684F	686F			



Anaesthetics
Articulating
Bur and Abrasives
Cosmetics
Diagnostic
Disposables
Endodontics
Equipment
Handpieces
Implant
Infection Control
Instruments
Isolation
Laboratory
Oral Surgery
Orthodontics
Patient Education and Gift
Pedodontics
Periodontal
Post and Core
Preventive
Prosthodontics
Restorative
Uniform



Bur & Abrasives

Bur Kits

Crown & Bridge Preparation Kit

/ Shofu

✓ *Ideal for speedy, simple and accurate tooth preparation for a wide range of Crown and Bridge restorations, Veneers, Inlays Onlays*



- Precise depth orientation grooves can be achieved easily with the specially designed depth cutting diamond points
- 2 special "Tissue-Protective End-Cutting Burs" are developed by Dr R.E. Goldstein for shoulder preparation and finishing under gingival margins without injuring the soft tissue
- These points are hygienically organized in an aluminum bur block for convenient sterilization autoclaving and easy replacement

Packing	Item no.
- 14 pcs x FG Diamond Points in 3 different grits: • Standard (no marking on shank) • Fine (with red band) • Super-Fine (with Yellow band)	SF-C&B

Ceramage Polishing Kit

/ Shofu


✓ *Specially created to achieve a lustrous super polished surface on Ceramic Polymer materials such as Ceramage, Solidex and indirect composites*



- **Recommended Polishing Sequence**
- Adjust small parts on the occlusal surface with Robot Carbide Fissure Bur
- Eliminate scratches on the surface and adjust the surface with Shofu Silicone Point HR2 and Shofu Pivot Brush SC
- Prepolish with Shofu Pivot Brush using DURA-POLISH
- Superpolish with Shofu Felt Wheel using DURA-POLISH DIA

Packing	Item no.
- 1 x Robot Carbide Fissure Bur - 4 x Shofu Silicone Points, HR2 • 1 each of #11 (SF-SP-11) • 2 pieces #28 (SF-SP-28) - 1 x Shofu Pivot Brush (SF-5950) - 3 x Shofu Felt Wheels • 2 x #1602 (Ø16mm, 2mm thickness) • 1 x #2405 (Ø24mm, 5mm thickness)	- 1 x Shofu Pivot Brush SC - 2 x Shofu Mandrels #13S (SF-SP-13S) (20g/jar) - 1 x DURA-POLISH DIA (5g/jar)
	SF-POLISH & FINISH


Strauss Composite Kit / Strauss



Profile	Item No.	Application
	E4XXF	Very fine polishing of composite excess at margins
	E5XXF	Very fine polishing of composite excess at margins
	K2XXF	Very fine polishing of palatal / lingual surface
	A3XXF	Very fine polishing of the marginal ridge
	A4XXF	Very fine polishing of the marginal ridge
	E4SF	Finishing and polishing of composite excess restoration margins
	E5SF	Finishing and polishing of composite restoration margins
	T5SF	Finishing of composite excess at margins
	K2SF	Fine polishing of palatal / lingual surface
	A4SF	Finishing and polishing of the marginal ridge
	PR2XXF	Very fine polishing of the occlusal surface
	PR2SF	Fine finishing of the occlusal surface

Packing	Item no.
Kit	ST-COMSET1

Strauss Porcelain Veneer Kit / Strauss



Profile	Item No.	Application
	PR12M	Producing 0.5mm depth cut
	PR13M	Producing 0.3mm depth cut
	A0F	Defining the finish line
	PR15SF	Finishing and polishing porcelain margins
	PR15XXF	Fine polishing of porcelain margins
	K2SF	Occlusal adjustment of porcelain at the incisal or palatal side
	K2XXF	Finishing and polishing of porcelain at the incisal or palatal side
	PR11C	Occlusal and incisal edge reduction
	F1RM	Labial and interproximal reduction and creation of the finish line around the preparation

Packing	Item no.
Kit	ST-PLSET1



Bur Kits

Bur Kits / Coltene & Whaledent

■ Diatech Inlay & Crown Preparation Kit

- Everything depends on the geometry of the preparation being correct
- Make use of this ingenious selection of 8 Multilayer FG preparation diamonds and 6 Multilayer FG finishing diamonds to perfect your inlay-, onlay- and crown preparations systematically and with the correct geometry
- The outcome: Minimal loss of tooth structure and optimum retention
- The best preconditions for producing long lasting restorations, for your patients' benefit

Diatch Inlay & Crown Preparation Kit

Diatch Universal Shapening & Finishing Kit



Diatch Composite Polishing Kit

Fig. No.	836R	836R	878	878	879	850	850	859	836R	836R	878	879	368	859
ISO Ø 1/10 mm	010	014	010	012	014	012	016	010	010	014	010	014	023	010
Head L mm/Grit	6ML	6ML	8ML	8ML	10ML	10ML	10ML	10ML	6F	6F	8F	10F	5F	10F
Grit Size	125µm	125µm	125µm	125µm	125µm	125µm	125µm	125µm	45µm	45µm	45µm	45µm	45µm	45µm
	Packing												Item no.	
	14 x FG Multilayer Diamonds												CW-250006AA	

■ Diatch Universal Shapening & Finishing Kit

- Comprises 8 Multilayer FG diamonds for contouring and 6 Multilayer FG diamonds for finishing amalgam, composite, compomer, porcelain, lassionomer and gold restorations
- Diatch fine grit Multilayer diamonds are highly efficient and will not harm the surface of the restoration nor the adjacent tooth structure

Packing	Item no.
14 x FG Multilayer Diamonds	CW-250007AA

Fig. No.	905	905	368	368	863	862	801	801	368	368	863	862	801	801
ISO Ø 1/10 mm	023	023	016	023	012	012	018	023	016	023	012	012	018	023
Head L mm/Grit	2.7F	2.9F	3.5XF	5XF	10XF	8XF	XF	XF	3.5UF	5UF	10UF	8UF	UF	UF
Grit Size	45µm	45µm	25µm	25µm	25µm	25µm	25µm	25µm	15µm	15µm	15µm	15µm	15µm	15µm

■ Diatch Crown & Bridge Separation Kit

- Time saving
- These 3 cross-cut burs and 1 tungsten carbide bur plus the coarse Multilayer diamond cylinder separate even the hardest of alloys in the shortest possible time
- Even if these instruments are used repeatedly, their efficiency hardly drops

Packing	Item no.
- 1 x FG Multilayer Diamond - 4 x FG Tungsten Carbide Burs	CW-250012AA

Fig. No.	CX21R	CX21R	CX23R	C31R	880
ISO Ø 1/10 mm	010	012	012	010	012
Head L mm/Grit	4	4	4	4	6ML
Grit Size / Type of Blades	cross	cross	cross	Reg.	125µm

■ Diatch Composite Polishing Kit

- Smooth surfaces offer the best protection against plaque accumulation and discolouration, especially important for composite materials
- The DIATECH silicone polishers with medically-tested silicone create a smooth, shiny surface on composites – even on tricky fine hybrid composites
- The kit comprises 8 light green RA instruments for prepolishing and 8 dark green RA instruments for high lustre polishing

Fig. No.	2101RA	2103RA	2111RA	2127RA	2201RA	2203RA	2211RA	2227RA	
ISO Ø 1/10 mm	030	050	060	100	030	050	060	100	
	Prepolishing				High Shine Polishing				
No. of Piece	x 2	x 2	x 2	x 2	x 2	x 2	x 2	x 2	
	Packing								Item no.
	16 x RA Silicone Polishers								CW-250014AA



Bur & Abrasives

Bur Block & Accessories

Bur Blocks



HENRY SCHEIN® Anodized Aluminum Bur Blocks

/ Henry Schein

- Anodized aluminum bur blocks allow an organized, hygienic and ergonomic approach to every task
- Pivoting cover keeps burs safely in place when not in use




- Autoclavable
- Dimensions: 60 x 31 x 26mm

Description	Packing	Item no.
15-Hole: FG	Each	HS-112-5480
15-Hole: 10 FG & 5 RA		HS-112-5481

Autoclavable Bur Caddy

/ Henry Schein



- Holds 8 FG and 4 CA burs
- Autoclavable to 325°F
- Cover/handle swivels 360° with stops at each 45°

Color	Packing	Item no.
Beige	Each	HS-100-5329
Blue		HS-100-8158

Directa RotoStand™

/ Directa

- For 18 FG burs and 44 HP/RA burs

- Simplifies handling and organization of burs
- The pivoted bur platform rotates to facilitate quick and easy selection of burs
- Can be cold sterilized
- Close fitting lid to protect burs from dust and contamination



Description	Packing	Item no.
60 FG Burs	Each	DIR-4122-85
18 FG Burs & 44 HP/RA Burs		DIR-4122-87

Bur Station II

/ Shofu

- ✓ *Specially designed for easy sterilization and storage of FG and CA burs*

- Open for easy placement for chair-side use and closed for disinfection
- Holds 11 FG and 5 CA burs
- Can be sterilized by autoclaving to 135°C (270°F), ultrasonic cleaning and cold sterilization
- Dimension: 42 x 23 x 31mm

Handy & Convenient



Description	Packing	Item no.
Blue	Each	SF-0797
Yellow		SF-0796

Mandrels

Mandrels (HP/CA)

/ Shofu

- Stainless steel Mandrel (HP/CA) secures all unmounted wheels and disks for optimal control



Shank	Packing	Item no.
CA	12's/box	SF-CA-MANDREL
HP		SF-HP-MANDREL

RA Mandrels Refill

/ 3M ESPE

- ✓ *A latch-attachment, pop-on mandrel for use with Sof-Lex™ Contouring and Polishing Discs and Sof-Lex™ Extra-Thin Contouring and Polishing Discs*

- Low-profile mandrel design reduces risk of marring the restoration for easier use and improved results



Packing	Item no.
3's/box	3M-1983RA

Others

Mesh Bur and Diamond Holder / Miltex

- Securely holds FG, RA, HP burs, diamonds, and other small items, preventing their loss
- Stainless steel design for use in all types of autoclaves
- Holds crowns and bridges that are placed in an ultrasonic cleaner for cement removal



Size	Packing	Item no.
Small	Each	HS-953-2143
Large		HS-100-3082

One-Gloss Mandrels CA

/ Shofu



Packing	Item no.
6's/box	SF-0435

Super-Snap® Mandrels CA

/ Shofu



Packing	Item no.
6's/box	SF-0439

Diamond Dresser / Shofu

- ✓ *Diamond plated instrument for cleaning, shaping and reshaping of abrasive stones and silicone points*



Packing	Item no.
Each	SF-0781

