

3M™ Bair Hugger™ Therapy
Temperature Management Unit
Model 750



Clinical Warming With Control

**Performance Precision. Ease of Use.
It's everything you want from your
forced-air warming unit...and more.**



The standard of care in temperature management technology

Forced-air warming was introduced to the world in 1987 with the launch of Bair Hugger therapy. Since that time, Bair Hugger warming units and blankets have safely and effectively warmed more than 160 million patients worldwide.

Performance

- High airflow for greater heat transfer to the patient
- Hose end temperature sensing ensures consistent temperatures

Precision

- Microprocessor technology regulates temperature settings within $\pm 1.5^{\circ}\text{C}$ for accurate temperature delivery
- Three independent sensors continuously monitor system performance

Ease of Use

- Advanced user interface provides finger-tip control of temperature settings for one-touch temperature selection
- No need to open the unit for over-temperature tests or calibration



Ordering Information

| Specifications | Bair Hugger Model 750 |
|-------------------------------|--|
| Dimensions | 12.5"h x 10"w x 13.5"d (31.5h x 25w x 34d cm) |
| Weight | 16.3 lb (7.4 kg) |
| Operating Temperatures | High: $43^{\circ} \pm 1.5^{\circ}\text{C}$ ($109.4^{\circ} \pm 2.7^{\circ}\text{F}$) Med: $38^{\circ} \pm 1.5^{\circ}\text{C}$ ($100.4^{\circ} \pm 2.7^{\circ}\text{F}$) Low: $32^{\circ} \pm 1.5^{\circ}\text{C}$ ($89.6^{\circ} \pm 2.7^{\circ}\text{F}$) Ambient: Room Temperature |
| Fan Speeds and Relative Noise | High Fan Setting: 41 CFM (17.5 L/s) 53 dBA |
| Leakage Current | Meets regulatory standards for leakage current |
| Alarm Systems | Audible and visual alarms; Over-temperature $\leq 56^{\circ}\text{C}$, $53^{\circ}\text{C} \pm 3^{\circ}\text{C}$ typical |
| Filter | High-efficiency 0.2 μm filter |
| Device Ratings | 110-120 VAC, 50/60 Hz, 11.7 Amperes 220-240 VAC, 50/60 Hz, 7.2 Amperes 100 VAC, 50/60 Hz, 15 Amperes |
| Compatibility | Use only with Bair Hugger blankets, 3M™ 241™ fluid warming and 3M™ Bair Paws™ gowns |
| Accessories | Rolling cart, accessory handle and teathered clip also available |

To order, please contact your 3M Patient Warming representative or call Customer Service at 1-800-733-7775.



**Arizant Healthcare Inc.,
a 3M company**
10393 West 70th St.
Eden Prairie, MN 55344 USA
Phone 800-733-7775
Fax 800-775-0002
www.bairhugger.com