



# ACCESSORIES

for use with the  
**LIFEPAK® 12**  
 DEFIBRILLATOR/MONITOR



**MANUAL DEFIBRILLATION**  
 1 Push EPR. Ready indicator will go to heart position or apply combination electrodes.  
 2 Select ENERGY.  
 3 Push CHARGE. Stand clear.  
 Push SHOCK to deliver energy.

**AED OPERATION**  
 • Push EPR.  
 • Push ANALYZE.  
 • Push SHOCK when directed to deliver energy.

**PACER OPERATION**  
 • Push PACER to turn pacer on.  
 • Push RATE button and adjust up or down as needed.  
 • Push CURRENT button and adjust to capture.

**Medtronic**  
**LIFEPAK 12**  
 BIPHASIC DEFIBRILLATOR/MONITOR

12-LEAD  
 TRANSMIT  
 CODE SUMMARY  
 PRINT

DANGER (EXPLOSION HAZARD) - DO NOT USE IN THE PRESENCE OF FLAMMABLE GASES.  
 WARNING (HARMFUL ELECTRICAL OUTPUT) - FOR USE ONLY BY QUALIFIED PERSONNEL.

Genuine accessories from Physio-Control



# POWER OPTIONS

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR

Thank you for choosing Physio-Control as your partner in helping you save lives and improve patient care. You already know the LIFEPAK reputation for quality and reliability in our defibrillators and AEDs. When it comes to LIFEPAK accessories, we know you won't accept any less. Inside this catalog you'll find accessories designed to meet your needs and reliability standards. All the accessories and disposables in this catalog are designed to work with the LIFEPAK 12 Defibrillator/Monitor. If you have any questions about accessories or disposables, please contact your Physio-Control representative.



### FASTPAK® Battery

1.2 amp hour capacity, rechargeable nickel-cadmium.

**11141-000044**



### FASTPAK 2 Battery

1.2 amp hour capacity, rechargeable nickel-cadmium, with fuel gauge.

**11141-000025**



### LIFEPAK SLA Battery

Rechargeable, sealed, lead-acid, 2.5 amp hour capacity.

**11141-000028**



### LIFEPAK NiCd Battery

1.7 amp hour capacity, rechargeable nickel-cadmium, with fuel gauge ...**11141-000027**

### LIFEPAK NiCd Battery

2.4 amp hour capacity, rechargeable nickel-cadmium, with fuel gauge. Must be used with an AC or DC Power Adapter that is compatible with 2.4 amp hour batteries. Same size and weight as 1.7 amp hour battery (see above) .....**11141-000026**

# POWER OPTIONS

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR



## Mobile Battery Service Station (MBSS)

DC and AC operation for mobile transport applications, for use with any of the above batteries. Includes AC power cord, mounting bracket, and manual. (Batteries and DC cable not included.)

**99996-00009**



**AC Power Adapter.....11400-00038**

**Bracket for Vertical or Horizontal Mounting of the MBSS Charger .....99996-00010**

**12-Volt DC Cable for MBSS Charger.....99996-00011**



## Battery Support System 2

For use with any of the above batteries. Includes power cord and operating instructions (batteries not included) .....**99407-00002**

**Wall Mount Bracket for Battery Support System 2.....11210-00002**



## Power Adapter Extension Cable for use with LIFEPAK 12 Defibrillator/Monitor

Optional use with AC power adapter.

**11110-000051**

# ECG MONITORING ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR



**12-Lead ECG Patient Cable**

Main 5-foot trunk (right angle connector) ..... **11110-000111**

Main 8-Foot Trunk (Right-Angle Connector) ..... **11110-000110**



**4-Wire Limb Lead Attachment**

**11110-000102**



**6-Wire Precordial Lead Attachment**

**11110-000103**



**3-Wire ECG Cable**

Right-angle connector

**11110-000029**



**5-Wire ECG Cable**

Right-angle connector, 4-wire limb plus 1 chest lead, labeled "V1" on the LIFEPAK 12 defibrillator reports.

**11110-000066**



**LIFE•PATCH® ECG Electrodes**

Adult, pregelled (package of 3 electrodes)

**11100-000001**



**LIFE•PATCH ECG Electrodes**

Adult, pregelled (package of 4 electrodes)

**11100-000002**



**Box of ECG Printer Paper**

50mm x 30m (3 rolls per box)

**11240-000013**

**Box of Strip Chart Recorder Paper**

100mm x 22m (2 rolls per box)

**11240-000016**

# THERAPY DELIVERY ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR



## Standard Hard Paddles

1 pair, for use with the LIFEPAK 12 defibrillator/monitor ..... **11130-00001**

## External Sterilizable Paddles

1 pair, for use with the LIFEPAK 12 defibrillator/monitor ..... **11134-00003**



## Pediatric Paddle, External

1 each (2 required); Slips onto standard adult or external sterilizable paddle.

**11133-00001**



## SIGNAGEL® Electrode Gel

One 8.5 oz. tube for use with hard paddles. Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel. Recommended for ECG, defibrillation, biofeedback, and EMG.

**21300-005847**



## Internal Paddle Handles with Discharge Control

1 pair, for use with the LIFEPAK 12 defibrillator/monitor.

**11131-00001**

## QUIK-COMBO® Therapy Cable

For use with the LIFEPAK 12 defibrillator/monitor.

**11110-000040**

## FAST-PATCH® Adapter Cable

For use with the QUIK-COMBO therapy cable for the LIFEPAK 12 defibrillator/monitor.

**11110-000052**



## Internal Paddles

1." size (6.25" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000010**

## Internal Paddles

2" size (8.75" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000022**

## Internal Paddles

1.5" size (6.00" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000011**

## Internal Paddles

2.5" size (5.50" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000013**

## Internal Paddles

1.5" size (9.00" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000021**

## Internal Paddles

2.5" size (8.50" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000019**

## Internal Paddles

1.5" size (14.25" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000024**

## Internal Paddles

3.5" size (5.00" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000014**

## Internal Paddles

2" size (5.75" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000012**

## Internal Paddles

3.5" size (8.00" shaft length), 1 pair, requires Internal Paddle Handles.

**11131-000023**

# THERAPY DELIVERY ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR

## EDGE SYSTEM™ ELECTRODES FOR PACING/DEFIBRILLATION/ECG WITH QUIK-COMBO CONNECTORS

18-month minimum shelf life remaining at time of shipment from Physio-Control except where noted.



EDGE System Electrodes with QUIK-COMBO Connector .....11996-000091



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK™ Preconnect System

11996-000017



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector

11996-000090



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only on manual defibrillator/monitors. 12-month minimum shelf life at time of shipment.

11996-000093

## OTHER EDGE SYSTEM ELECTRODES FOR PACING/DEFIBRILLATION/ECG

18-month minimum shelf life at time of shipment.

EDGE System Electrodes with FAST-PATCH Connector

Must be used with QUIK-COMBO to FAST-PATCH

Adapter Cable. (11110-000052) ..... 11996-000092

# SpO<sub>2</sub> MONITORING ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR

## MASIMO® SET® LNCS SENSORS AND CABLES



### Masimo SET LNC-4 LNCS Patient Cable

4-foot reusable connector cable .....11171-00024

### Masimo SET LNC-10 LNCS Patient Cable

10-foot reusable connector cable .....11171-00016

### Masimo SET LNC-14 LNCS Patient Cable

14-foot reusable connector cable .....11171-00025



### Masimo SET LNCS DCIP Pediatric Reusable Sensor

Multiuse sensor for patients 10-50 kg.

11171-00018



### Masimo SET LNCS DCI Adult Reusable Sensor

Multiuse sensor for patients >30 kg .....11171-00017



### Masimo SET LNC-4-Ext LNCS Extension Cable

4-foot extension patient cable (May be used with devices equipped with Nellcor conventional SpO<sub>2</sub>)

11171-00027



### Masimo SET LNCS Pdtx Pediatric Adhesive Sensor

Single patient use sensor for patients 10-50 kg (20/box).

11171-00020



### Masimo SET LNCS Adtx Adult Adhesive Sensor

Single patient use sensor for patients >30 kg (20/box).

11171-00019

# SpO<sub>2</sub> MONITORING ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR

## MASIMO SET LNCS SENSORS AND CABLES (continued)



### Masimo SET LNCS Neo-L Neonatal Adhesive Sensor

Single patient use sensor for patients <3 kg or >40 kg (20/box).

11171-000021



### Masimo SET LNCS NeoPt-L Neonatal Adhesive Sensor

Single patient use sensor for patients <1 kg (20/box).

11171-000022



### Masimo SET LNCS Inf-L Infant Adhesive Sensor

Single patient use sensor for patients 3–20 kg (20/box).

11171-000026

## MASIMO SET LNOP SENSORS AND CABLES



### Masimo SET PC-04 LNOP Patient Cable

4-foot reusable connector cable ..... 11171-000006

### Masimo SET PC-08 LNOP Patient Cable

8-foot reusable connector cable ..... 11171-000008

### Masimo SET PC-12 LNOP Patient Cable

12-foot reusable connector cable ..... 11171-000009



### Masimo SET LNOP DCI Adult Reusable Sensor

Multiuse sensor for patients >30 kg

11171-000007



### Masimo SET LNOP DCIP Reusable Pediatric Sensor

Multiuse sensor for patients 10-50 kg

11171-000010



## MASIMO SET LNOP SENSORS AND CABLES (continued )



### Masimo SET LNOP Adt Adult Adhesive Sensor

Single patient use sensor for patients >30 kg (20/box).

11171-000011



### Masimo SET LNOP Pdt Pediatric Adhesive Sensor

Single patient use sensor for patients 10–50 kg, (20/box).

11171-000012



### Masimo SET LNOP Neo Neonatal Adhesive Sensor

Single patient use sensor for patients <10 kg (20/box).

11171-000013



### Masimo SET LNOP Neo-Pt Neonatal Adhesive Sensor

Single patient use sensor for patients <1 kg (20/box).

11171-000014

---

## Adapter Cables

### MNC-1 Adapter Cable (10 ft)

Allows LIFEPAK 12 defibrillator/monitor with Masimo SpO<sub>2</sub> to connect to Nellcor sensors.

11996-000183

### MNC-1 Adapter Cable (4 ft)

Allows LIFEPAK 12 defibrillator/monitor with Masimo SpO<sub>2</sub> to connect to Nellcor sensors.

11996-000198

# SpO<sub>2</sub> MONITORING ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR

## NELLCOR® SpO2 SENSORS AND CABLES



### DS100A DURASENSOR®

Adult reusable, includes 1 reusable clip,  
1-year warranty.

**11996-00060**



### Oxisensor II Disposable Pediatric Sensor

(box of 24)

**11996-000116**

### DURA-Y® Multisite Sensor

Reusable, includes 1 reusable sensor, 6-month warranty.....**11996-000106**

### Oxiband® Adult/Neonatal Sensor

Reusable, 3-month warranty (sensor), includes 50 adhesive bandage wraps.....**11996-000061**



### Oxisensor II Disposable Infant Sensor

(box of 24)

**11996-000115**

### Oxiband Pediatric/Infant Sensor

Reusable, 3-month warranty (sensor), includes 50 adhesive bandage wraps.....**11996-000062**

### Disposable Adhesive Foam Wrap

For OXI-A/N (2 bags of 50).....**11996-000110**

### Disposable Adhesive Bandage Wrap

For OXI-A/N (2 bags of 50).....**11996-000048**

### Disposable Adhesive Foam Wrap

For OXI-P/I (2 bags of 50).....**11996-000108**

### Disposable Adhesive Bandage Wrap

For OXI-P/I (2 bags of 50).....**11996-000049**

**Disposable Posey Wrap for Oxiband (bag of 12) .....11996-000107**

### Oxisensor II Disposable Neonatal Sensor

(box of 24) .....**11996-000117**

### Oxisensor II Disposable Adult Sensor Assorted Pack

Includes 2xD25 and 1xD25L .....**11996-000052**

### Oxisensor II Disposable Pediatric Sensor Assorted Pack

Includes 2xD20 and 1xL20 .....**11996-000051**

### Oxisensor II Disposable Neonatal Sensor Assorted Pack

Includes 2xN25 and 1xL20 .....**11996-000050**



### Oxisensor® II Disposable Adult Sensor

(box of 24)

**11996-000113**

Long cable (box of 24)

**11996-000114**



### DEC-4 Cable Extension

Includes 2xN25 and 1xL20

**11110-000042**

### DEC-8 Cable Extension

8-foot extension cable

**11110-000176**

# END TIDAL CO<sub>2</sub> (ETCO<sub>2</sub>) MONITORING ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR

## ORIDION® INTUBATED FILTERLINES



### FilterLine® SET Adult/Pediatric

Includes airway adapter. FilterLine (79") for short-term intubated patients, 24 hours typical (pack of 10).

11996-000299



### FilterLine H SET Adult/Pediatric

Includes airway adapter, High-humidity FilterLine (79") for long-term intubated patients, 72 hours typical (box of 25).

11996-000080

### FilterLine SET Long Adult/Pediatric

Includes airway adapter, Long-length FilterLine (185") for short-term intubated patients, 24 hours typical (box of 25).

11996-000164

### FilterLine SET Adult/Pediatric

Includes airway adapter. FilterLine (79") for short-term intubated patients, 24 hours typical (box of 25).

11996-000081



### FilterLine H SET Infant/Neonatal

Includes airway adapter, High-humidity FilterLine (79") for long-term intubated patients, 72 hours typical (box of 25).

11996-000001

## ORIDION NON-INTUBATED FILTERLINES



### Smart CapnoLine® Plus with O<sub>2</sub> Delivery Adult/Intermediate

For patients above 44 lbs, O<sub>2</sub> delivery and CO<sub>2</sub> sampling via mouth or nose, 12 hours equals typical length of use (pack of 10).

11996-000298



### Smart CapnoLine with O<sub>2</sub> Delivery Pediatric

For patients less than 44 lbs, O<sub>2</sub> delivery and CO<sub>2</sub> sampling via mouth or nose, 12 hours equals typical length of use (box of 25).

11996-000128

### Smart CapnoLine® Plus with O<sub>2</sub> Delivery Adult/Intermediate

For patients above 44 lbs, O<sub>2</sub> delivery and CO<sub>2</sub> sampling via mouth or nose, 12 hours equals typical length of use (box of 25) ..... 11996-000163

### Smart CapnoLine Plus with O<sub>2</sub> Delivery Adult/Intermediate

For patients above 44 lbs, O<sub>2</sub> delivery and CO<sub>2</sub> sampling via mouth or nose, 12 hours equals typical length of use (case of 100) ..... 11996-000167

### Smart CapnoLine Plus Long with O<sub>2</sub> Delivery Adult/Intermediate

For patients above 44 lbs, long-length nasal cannula (165") for O<sub>2</sub> delivery and CO<sub>2</sub> sampling via mouth or nose, 12 hours equals typical length of use (box of 25) ..... 11996-000165



### Smart CapnoLine Pediatric

For patients less than 44 lbs, CO<sub>2</sub> sampling via the mouth or the nose, 12 hours equals typical length of use (box of 25).

11996-000120

### Smart CapnoLine Plus Adult/Intermediate

For patients above 44 lbs, CO<sub>2</sub> sampling via the mouth or the nose, 12 hours equals typical length of use (case of 100)..... 11996-000166

### Smart CapnoLine Plus Adult/Intermediate

For patients above 44 lbs, CO<sub>2</sub> sampling via the mouth or the nose, 12 hours equals typical length of use (pack of 10) ..... 11996-000297

### Smart CapnoLine Plus Adult/Intermediate

For patients above 44 lbs, CO<sub>2</sub> sampling via the mouth or the nose, 12 hours equals typical length of use (box of 25)..... 11996-000162

# NIBP MONITORING ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR



- Adult Cuff**, 14x37.5cm, reusable..... **11996-00024**
- Large Adult Cuff**, 16x42cm, reusable ..... **11996-00025**
- Small Adult Cuff**, 12x30cm, reusable ..... **11996-00023**
- Child Cuff (Pediatric)**, 9x27cm, reusable..... **11996-00022**
- Small Child Cuff**, 7x21cm, reusable ..... **11996-00021**
- Thigh Cuff**, 20x54cm, reusable..... **11996-00026**
- Infant Cuff**, 6x18cm, reusable ..... **11996-00020**



- Small Adult Cuff**, 12x30cm, disposable ..... **11996-00029**
- Adult Cuff**, 14x35cm, disposable ..... **11996-00030**
- Large Adult Cuff**, 16x40cm, disposable ..... **11996-00031**
- X-Large Adult Cuff**, 20x50cm, disposable..... **11996-00032**
- Child Cuff (Pediatric)**, 9x22.5cm, disposable ..... **11996-00028**
- Infant Cuff**, 6x15cm, disposable ..... **11996-00027**



**NIBP Hose**  
Coiled, 2–9 feet  
**11996-000125**



**NIBP Hose**  
9 feet  
**11996-000033**  
  
12 feet  
**11996-000089**

# CARRYING CASES AND ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR



## Basic Carrying Case

Includes shoulder strap, right pouch, left pouch, and front cover.

**11260-000030**

## Carrying Case for use with LIFEPAK 12 Defibrillator/Monitor with AC Power Adapter

Includes shoulder strap, right pouch, left pouch and front cover.

**11260-000032**

## Carrying Case for use with LIFEPAK 12 Defibrillator/Monitor with Voice Recorder

Includes shoulder strap, right pouch, left pouch and front cover.

**11260-000033**



## Back Pouch

Ideal for additional accessory storage. (Not recommended for use with SLA batteries.)

**11260-000029**



## Top Pouch

Storage for sensors and electrodes; insert in place of standard paddles.

**11220-000028**

## Replacement Shoulder Strap

**11260-000010**

## Right Pouch Replacement

(Note: Included with basic case.)

**21300-007203**

## Left Pouch Replacement

(Note: Included with basic case.)

**21300-007201**

## Carrying Case Base and Side Supports Replacement

**21300-006361**

## Front Cover

Soft plastic cover for display and keypads. Cannot be used with Removable Acrylic Screen Shield. (Note: Included with basic case.)

**11220-000033**

## Removable Acrylic Screen Shield

Hard acrylic shield for display protection (cannot be used with front cover.)

**11998-000063**

## Therapy Connector Guard

Recommended for protection of the therapy connector on the defibrillator/monitor from direct impacts during use. Permanently attaches to the lower front edge of the defibrillator/monitor.

**21330-001112**

## Therapy Connector Shield

Recommended for protection of the therapy connector, as an alternative to the Therapy Connector Guard. Attaches to the defibrillator/monitor with VELCRO®.

**21330-006463**

# COMMUNICATION ACCESSORIES

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR

## Internal Modem (PC Card Modem)

For cellular and land transmission, requires modem door assembly.

**11150-000015**



## Modem Door Assembly

**11150-000009**

## LIFEPAK 12 defibrillator/monitor to External Modem Cable

6 feet

**11150-000007**

## LIFEPAK 12 defibrillator/monitor to External Modem Cable

10 feet

**11150-000008**

## External Modem

For land transmission only, requires external modem cable.

**11150-000010**



## LIFEPAK 12 defibrillator/monitor Configuration Transfer Cable

**11110-000055**



## LIFEPAK 12 defibrillator/monitor to PC Serial Port Cable

Use with LIFENET® System

**11230-000005**



## Analog ECG Output Cable

**11110-000044**

## MOUNTING OPTIONS AND ROLL STAND

### LIFESTAND™ Mobile Roll Stand for the LIFEPAK 12 Defibrillator/Monitor

**11996-000157**

### Optional IV Pole for the LIFESTAND Mobile Roll Stand

**11996-000158**

### Surface Mount Bracket

Designed to securely mount the LIFEPAK 12 defibrillator/monitor in an ambulance. Bracket is lightweight and suitable for use in aircraft.

**11996-000309**



# TRAINING TOOLS, TESTERS AND MATERIALS

for the LIFEPAK 12 DEFIBRILLATOR/MONITOR



## QUIK-COMBO 12-Lead Patient Simulator

Simulates patient heart rhythms. Allows users to practice and demonstrate monitoring with the 12-lead or 5-wire ECG cable and defibrillating with the QUIK-COMBO therapy cable. Simulates: asystole, bradycardia, normal sinus rhythm, ventricular fibrillation (VF), ventricular tachycardia, and acute MI. Also simulates patient motion. Includes QUIK-COMBO cable, batteries, and operating instructions .....**11200-00005**

## QUIK-COMBO 3-Lead Patient Simulator

Simulates patient heart rhythms. Allows users to practice and demonstrate monitoring with the 3-wire ECG cable and defibrillating with the QUIK-COMBO therapy cable. Simulates: asystole, bradycardia, normal sinus rhythm, ventricular fibrillation (VF), and ventricular tachycardia. Also simulates patient motion. Includes QUIK-COMBO cable, batteries, and operating instructions.....**11201-00002**

**Operating Instructions: LIFEPAK 12 Defibrillator/Monitor** .....**26500-002481**

**Operating Instructions: Battery Support System 2** .....**26500-000168**

## Inservice Video: LIFEPAK 12 Defibrillator/Monitor

VHS format.....**26500-000213**

**Service Manual on CD-ROM: LIFEPAK 12 Defibrillator/Monitor and Battery Support System 2** .....**21500-006152**

**Service Manual (Paper Version): LIFEPAK 12 Defibrillator/Monitor and Battery Support System 2** .....**26500-000234**



## Clip-on Training Electrodes

For use with QUIK-COMBO patient simulators. Allows users to practice defibrillation electrode placement. Includes instructions for use and five (5) pair of reusable training electrodes.

**11250-000052**



## Defibrillator Checker

Tests integrity of energy delivery through the standard hard paddles. Neon light indicates energy has been delivered.

**11998-000060**

## Replacement AED Training Electrodes

Includes five (5) pair of reusable training electrodes. Adhesive electrode portion only.

**11101-000003**

## QUIK-COMBO Therapy Cable Tester

Tests integrity of the QUIK-COMBO therapy cable.

**11998-000160**



## Test Load

Use to perform therapy cable performance checks. Connects directly into the QUIK-COMBO therapy cable on the defibrillator.

**11998-000022**



© 2008 Physio-Control, Inc. LIFEPAK, LIFE•PATCH, FASTPAK, FAST-PATCH, LIFENET, and QUIK-COMBO are registered trademarks of Medtronic, Inc. EDGE System and REDI-PAK are trademarks of Medtronic, Inc. FilterLine and Smart CapnoLine are registered trademarks of Oridion Systems, Ltd. LIFESTAND is a trademark of MODO, Inc. Durasensor, DURA-Y, Oxiband and Oxisensor are registered trademarks of Mallinckrodt, Ltd. MASIMO SET is a registered trademark of MASIMO CORPORATION. Dura-Y and Oxiband are registered trademarks of Nellcor Puritan Bennett Inc. SIGNAGEL is a registered trademark of Parker Laboratories. VELCRO is a registered trademark of Velcro Industries B.V. Specifications subject to change without notice. Not all products available outside the U.S. Please consult your LIFEPAK Operating Instructions or Declaration of Conformity for a list of approved accessories in your region.

MIN 3207694-001 / CAT 26500-002680

