



Resuscitation supplies

PHILIPS

Changes to this edition

April 2013

New product added

Multi Patient Disposable Cuffs – Multi Care

| | | |
|--------------|---|----|
| 989803183311 | Infant (orange) | 41 |
| 989803183321 | Pediatric (green) | 41 |
| 989803183331 | Small Adult (light blue) | 41 |
| 989803183341 | Adult (navy blue) | 41 |
| 989803183351 | Adult XL (navy blue with edge striping) | 41 |
| 989803183361 | Large Adult (burgundy) | 41 |
| 989803183371 | Thigh (grey) | 41 |

Single Care Disposable Cuffs

| | | |
|--------------|---|----|
| 989803182281 | Pediatric (green) | 42 |
| 989803182291 | Small Adult (light blue) | 42 |
| 989803182301 | Adult (navy blue) | 42 |
| 989803182311 | Adult XL (navy blue with edge striping) | 42 |
| 989803182321 | Large Adult (burgundy) | 42 |

Product revised

Multifunction Electrode Pads

| | | |
|--------------|---|---|
| M3501A | Adult/Child with barrel connector | 7 |
| M3504A | Infant with barrel connector | 7 |
| M3713A | Adult/Child Plus with plug style connector | 7 |
| M3716A | Adult/Child Radiolucent with plug style connector | 7 |
| M3717A | Infant Plus with plug style connector | 7 |
| M3718A | Adult/Child Radiotransparent (Reduced Skin Irritation) with plug style connector | 7 |
| M3719A | Infant Radiotransparent (Reduced Skin Irritation) with plug style connector | 7 |
| 989803166021 | Adult/Child Preconnect | 7 |

Multi Patient Traditional Cuffs – Single Hose

| | | |
|--------|-------------|----|
| 40401A | Infant | 36 |
| 40401B | Pediatric | 36 |
| 40401C | Adult | 36 |
| 40401D | Large Adult | 36 |
| 40401E | Thigh | 36 |

Multi Patient Comfort Cuffs – Single Hose

| | | |
|--------|---|----|
| M1572A | Pediatric, color coded: green tubing | 37 |
| M1573A | Small Adult, color coded: light blue tubing | 37 |
| M1574A | Adult, color coded: navy blue tubing | 37 |
| M1575A | Large Adult, color coded: burgundy tubing | 37 |

Multi Patient Easy Care Cuffs – Single Hose

| | | |
|---------|----------------------------------|----|
| M4552B | Infant, orange | 38 |
| M4552B5 | Infant, orange | 38 |
| M4553B | Pediatric, green | 38 |
| M4553B5 | Pediatric, green | 38 |
| M4554B | Small Adult, light blue | 39 |
| M4554B5 | Small Adult, light blue | 39 |
| M4555B | Adult, navy blue | 39 |
| M4555B5 | Adult, navy blue | 39 |
| M4556B | Adult Extra-Long, navy blue | 39 |
| M4556B5 | Adult Extra-Long, navy blue | 39 |
| M4557B | Large Adult, burgundy | 39 |
| M4557B5 | Large Adult, burgundy | 39 |
| M4558B | Large Adult Extra-Long, burgundy | 39 |
| M4558B5 | Large Adult Extra-Long, burgundy | 40 |
| M4559B | Thigh, grey | 40 |
| M4559B5 | Thigh, grey | 40 |

Gentle Care Disposable Cuffs - Single Hose

| | | |
|--------|----------------------------------|----|
| M4572B | Infant, orange | 43 |
| M4573B | Pediatric, green | 43 |
| M4574B | Small Adult, light blue | 43 |
| M4575B | Adult, navy blue | 43 |
| M4576B | Adult Extra-Long, navy blue | 43 |
| M4577B | Large Adult, burgundy | 43 |
| M4578B | Large Adult Extra-Long, burgundy | 44 |
| M4579B | Thigh, grey | 44 |

Blood Pressure Air Hoses

| | | |
|--------|---------------|----|
| M1598B | NIBP Air Hose | 45 |
| M1599B | NIBP Air Hose | 45 |

Product removed**Multifunction Electrode Pads**

| | | |
|--------|--|---|
| M3502A | Adult Multifunction Pads (IEC – barrel connector) | 8 |
| M3503A | Infant/Child Multifunction Pads (IEC - barrel connector) | 8 |

Disposable Blood Pressure Kits

| | | |
|--------|---|----|
| M1567A | Single line blood pressure kit | 26 |
| M1568A | Dual line blood pressure kits | 26 |
| M1634A | Reusable adapter cable for M1567A and M1568A kits | 26 |
| M2271A | Holder for disposable pressure transducer | 26 |
| M2272A | IV Pole Mount | 26 |

Dear Customer:

If you're looking for high-quality medical supplies to optimize the performance of your medical devices, you'll find them right here. Whether your focus is on quality patient care, reliability, or cost containment, you can turn to a single source – Philips Medical Supplies. Explore our complete family of supplies, from electrodes to SpO₂ sensors, for every area of care.

With Philips, you can buy your medical equipment and supplies from a single source, simply and easily. Use this guide to find the Philips supplies you need – or go to our website to purchase supplies: <http://www.philips.com/healthcarestore>.

Using this guide

Latex Cautions

Philips Medical Supplies currently distributes over 600 products, of which less than 9% contain natural rubber latex or dry natural rubber that could contact humans. Because these substances may cause allergic reactions, products in this guide that contain them are identified with the following symbol:

LATEX

Caution: Products that are identified with this symbol contain dry natural rubber or natural rubber latex that can come into contact with humans and may cause allergic reactions.

To the best of our knowledge, products not identified with this symbol are not known to contain natural rubber latex or dry natural rubber. The latex symbol indicates **only** which **medical supplies** contain natural rubber latex or dry natural rubber that can come into contact with humans and may cause allergic reactions. If you want to know whether a Philips **medical device** contains natural rubber latex or dry natural rubber, refer to the product labeling.

The FDA requires that products be labeled with human contacting latex content information when appropriate. Always refer to product labeling for accurate information on human contacting dry natural rubber or natural rubber latex.

CE Marking

Products labeled with this symbol conform to the European Union (EU) directives as mandated by EU member state legislation and are available for sale in Europe to countries that accept CE marking with restrictions as may be noted. Medical devices labeled with this marking means the product is in conformity with the European

Medical Directive (MDD): 93/42/EEC Council Directive of 14 June 1993 on the approximation of the laws of the Member States concerning medical devices.

In order for a product to be able to be labeled with the CE marking, product quality and manufacturing quality systems must be documented to comply with certain standards. Compliance with the European directives and regulations ensures a standard level of performance for device function, safety, and health protection for users, patients, and others.

The CE symbol provided with some products in the Philips Medical Supplies Reference Guide and online store is intended only as an indication of availability in EU countries and is not intended as a legal product marking.

Countries that accept CE marked products include Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Liechtenstein, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Russia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom. A CE mark does not necessarily indicate that a product is marketed or distributed in all of the countries listed above.

Non-medical products that comply with other harmonized directives may also bear the CE mark.

Note: Some products identified in this catalogue may not be licensed in accordance with Canadian law.

Ordering Information

To order the medical supplies in this guide or to obtain more information or additional literature, contact your Philips Medical Supplies representative near you or visit our website: <http://www.philips.com/healthcarestore>

Online Ordering

Online ordering is a simple process:

1. Log in using a secure password.
2. Enter the product number and quantity in the shopping cart.
3. Enter your PO or credit card number at check-out.
4. Your order is immediately transferred into our fulfillment system, where you can track it anytime.

If you are new to the web, we can guide you through the process on the phone in minutes. Many new users tell us that the ordering process is so intuitive that they do not need help, even the first time.

You can search for supplies in the online catalog three easy ways: by application category, Philips instrument, or keyword.

If you currently transmit via MMS/EDI or auto fax, the interactive system provides real-time information and services not available through one-way transmissions.

Online ordering is available to U.S. customers at the following url: <http://www.philips.com/newhealthcarestore>

Resuscitation Supplies

ALS Monitors-Defibrillators

| | |
|--|----|
| ALS External Multifunction Electrode Pads and Cables | 7 |
| ALS External and Internal Paddle Sets | 9 |
| ALS Data Management Solutions | 11 |
| Q-CPR | 12 |
| ALS Accessories | 13 |
| ALS Batteries | 14 |
| Cases/Pouches | 15 |
| Capnography | 16 |
| ECG Cables and Leads | 18 |
| Patient Monitoring ECG Cables and Leads | 20 |
| ECG Electrodes | 25 |
| Paper | 26 |
| Invasive Blood Pressure | 27 |
| NIBP Cuffs and Air Hoses | 36 |
| SpO ₂ Sensors and Cables | 46 |
| Temperature Probes | 49 |

AED Supplies

| | |
|-------------------------------|----|
| External AED Pads | 52 |
| AED Accessories and Kits | 54 |
| AED Data Management Solutions | 55 |
| AED Training Materials | 56 |
| AED Batteries | 59 |
| AED Cases/Mounts/Cabinets | 60 |

ALS Monitors – Defibrillators


ALS External Multifunction Electrode Pads and Cables

Multifunction Electrode Pads


| Product # | Description | Packaging | Use With: |
|------------------------------|--|--------------|---|
| REVISÉ M3501A CE | Adult/Child with barrel connector | 10 sets/case | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |
| REVISÉ M3504A CE | Infant with barrel connector | 5 sets/case | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |
| REVISÉ M3713A CE | Adult/Child Plus with plug style connector | 10 sets/case | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |
| REVISÉ M3716A CE | Adult/Child Radiolucent with plug style connector | 10 sets/case | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |
| REVISÉ M3717A CE | Infant Plus with plug style connector | 5 sets/case | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |
| REVISÉ M3718A CE | Adult/Child Radiotransparent (Reduced Skin Irritation) with plug style connector | 10 sets/case | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |
| REVISÉ M3719A CE | Infant Radiotransparent (Reduced Skin Irritation) with plug style connector | 5 sets/case | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |
| REVISÉ 989803166021 CE | Adult/Child Preconnect | 10 sets/case | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |




M3713A



M3717A



M3718A



989803166021

External Multifunction Pads Cables, Adapters, Test Loads


| Product # | Description | Packaging | Use With: |
|--------------------|--|-------------------------|---|
| 989803171271 CE | Test Plug for M3508A Cable | 1 pack = 1 test plug | HeartStart XL+ Defibrillator/Monitor 861290 |
| M3507A CE | Hands-free Pads Cable (Barrel connector) | 1 cable | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |
| M3508A CE | Hands-free Pads Cable (Plug connector) | 1 cable | HeartStart MRx: M3535A, M3536A; HeartStart XL+: 861290; HeartStart XL: M4735A |
| 05-10000 CE | Medtronic™ Pads Adapter | 1 adapter | Medtronic Physio-Control LifePak 500, 9, 10c, 11, 12, 20 |
| 05-10100 CE | Zoll™ Pads Adapter | 1 adapter | Zoll 1200, 1400, 1600 and M-Series |
| 05-10200 CE | Defibrillator Pads Adapter | 1 adapter | Fits Philips hands-free cable M3507A to connect to: M3535A, M3536A, M4735A |
| M1781A CE | 50-ohm test load for M3507A cable | 1 | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M3725A CE | HeartStart 50-ohm test load for M3507A. Can also be used in combination with 05-10200 adapter. | 1 | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M4748A CE | HeartStart extension cable Length: 4" (10cm) | 1 cable | HeartStart XL: M4735A; 989803171271 Test Plug |



M3507A



M3508A



05-10000



05-10100



05-10200

ALS External and Internal Paddle Sets

External Anterior/Anterior Paddle Sets

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|---|
| M3543A CE | HeartStart MRx Water Resistant External Paddles with PCI | 1 set | HeartStart MRx: M3535A, M3536A |
| M4745A CE | Sterilizable External Paddles without PCI | 1 set | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M4746A CE | HeartStart XL External Paddles with PCI | 1 set | HeartStart XL: M4735A |

Replacement Electrodes for External Paddle Sets

| Product # | Description | Packaging | Use With: |
|--------------|--|-------------|--|
| M1789A CE | Rectangular Paddle Electrode Replacement | 1 electrode | For use with M4746A External Paddle Set with PCI |
| M4759A | Rectangular Paddle Electrode Replacement | 1 electrode | M3543A Paddle Set |

Switchless Reusable Internal Paddles (Discharge Switch on Defibrillator)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|---|
| M1741A CE | Adult Internal Paddles: 3", 7.5cm, switchless | 1 pair | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M1742A CE | Adult Internal Paddles: 2.4", 6.0cm, switchless | 1 pair | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M1743A CE | Pediatric Internal Paddles: 1.8", 4.5cm, switchless | 1 pair | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M1744A CE | Infant Internal Paddles: 1.1", 2.8cm, switchless | 1 pair | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |



M1741A-M1744A

Switched Reusable Internal Paddles (Discharge Switch on Handle)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|---|
| M4741A CE | Adult Internal Paddles: 3", 7.5cm, switched | 1 pair | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M4742A CE | Adult Internal Paddles: 2.4", 6.0cm, switched | 1 pair | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M4743A CE | Pediatric Internal Paddles: 1.8", 4.5cm, switched | 1 pair | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M4744A CE | Infant Internal Paddles: 1.1", 2.8cm, switched | 1 pair | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |

Internal Paddle Adapters

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|---|
| M4740A CE | HeartStart Switchless Internal Paddle Adapter | 1 adapter | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M4748A CE | HeartStart Extension Cable | 1 adapter | HeartStart XL: M4735A |

Sync Cables

| Product # | Description | Packaging | Use With: |
|--------------|---------------------------------------|-----------|---|
| M1783A CE | Sync Cable, 12-pin. Length: 8', 2.5m | 1 cable | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M5526A CE | Sync Cable, 12-pin. Length: 25', 7.8m | 1 cable | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |



M4741A-M4744A



M4740A

ALS Data Management Solutions

HeartStart Data Cards and Software

| Product # | Description | Packaging | Use With: |
|--------------------|-------------------------------|-------------------------|---|
| 989803147711 CE | HeartStart XL Data Card | 1 | HeartStart XL: M4735A |
| M3544A CE | HeartStart MRx Data Card Tray | 1 | HeartStart MRx: M3535A, M3536A |
| 989803171261 CE | USB Data Drive | 1 pack = 1 USB drive | HeartStart XL+ Defibrillator/Monitor 861290 |

Q-CPR

Q-CPR Compression Sensor and Accessories

| Product # | Description | Packaging | Use With: |
|-----------------------|---|-----------|---|
| M4761A CE | Q-CPR Compression Sensor | 1 sensor | HeartStart MRx: M3535A, M3536A |
| M4762A CE | Q-CPR Compression Sensor Adhesive Pads | 10 pads | HeartStart MRx: M3535A, M3536A |
| M4763A LATEX CE | Pads/CPR Cable | 1 cable | HeartStart MRx: M3535A, M3536A |
| 989803162401 CE | CPR Meter – next generation Q-CPR | 1 sensor | HeartStart MRx: M3535A, M3536A Software Rev. F or higher |
| 989803163291 CE | CPR Meter Adhesive Pads – next generation Q-CPR | 10 pads | HeartStart MRx: M3535A, M3536A Software Rev. F or higher |
| 989803158661 CE | Pads/CPR Cable – next generation Q-CPR | 1 cable | HeartStart MRx: M3535A, M3536A Software Rev. F or higher |












M4761A and M4762A



989803162401

ALS Accessories

HeartStart MRx Accessories

| Product # | Description | Packaging | Use With: |
|---|--|--|---|
| 453564063841 | MRx NIBP Calibration Kit | 1 kit | HeartStart MRx: M3535A, M3536A |
| 453564063851 | MRx etCO ₂ Calibration Kit | 1 kit | HeartStart MRx: M3535A, M3536A |
| M3537A | Bed Rail/Roll Stand Hook for Defibrillator | 1 hook | HeartStart MRx: M3535A, M3536A |
|  | | | |
| M3549A | MRx Wide Bed Rail/Roll Stand Hook | 1 hook | HeartStart MRx: M3535A, M3536A |
| M4737A | HeartStart MRx Monitor/Defibrillator Display Cover | 1 cover | HeartStart MRx: M3535A, M3536A |
|  | | | |
| M5521A | Green Handle and label set | 1 handle | HeartStart MRx: M3535A, M3536A |
|  | | | |
| M5522A | Blue Handle and label set | 1 handle | HeartStart MRx: M3535A, M3536A |
|  | | | |
| M5523A | Yellow Handle and label set | 1 handle | HeartStart MRx: M3535A, M3536A |
|  | | | |
| M5524A | Rose Handle and label set | 1 handle | HeartStart MRx: M3535A, M3536A |
|  | | | |
| M5525A | Gray Handle and label set | 1 handle | HeartStart MRx: M3535A, M3536A |
|  | | | |
| M5528A | Mounting Bracket | 1 bracket | HeartStart MRx: M3535A, M3536A |
|  | | | |
| M5529A | DC Power Module | 1 | HeartStart MRx: M3535A, M3536A |
|  | | | |
| 989803171701 | Bed Rail/Roll Stand Hook | 1 pack = 1 bed rail/ roll stand mount | HeartStart XL+ Defibrillator/Monitor 861290 |
|  | | | |
| 989803176541 | HeartStart MRx Quick Disconnect DC Power Cable | 1 pack = 1 complete set of cables | 861593, 861494, 861495, M3536A |



M5521A



M5528A



M5529A

ALS Batteries

HeartStart ALS Batteries

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--------------------------------|
| M3516A CE | Sealed Lead Acid Battery | 1 battery | HeartStart XL: M4735A |
| M3538A CE | Lithium Ion Battery with Capacity Gauge | 1 battery | HeartStart MRx: M3535A, M3536A |
| M3506A CE | Battery Charger Adapter Tray | 1 tray | HeartStart MRx: M3516A ALS |
| M4747A CE | HeartStart XL Battery Charger Kit Consists of: M3516A battery, M3506A charger tray and M3517A AC Power Module | 1 kit | HeartStart XL: M4735A |



M3516A



M3538A

Cases/Pouches

Cases and Pouches

| Product # | Description | Packaging | Use With: |
|--------------------|---|---|---|
| M3541A CE | MRx Complete Carrying Case (includes 3 accessory pouches, shoulder strap and sleeve) | 1 case | HeartStart MRx: M3535A, M3536A |
| M4751A | HeartStart XL Accessory Pouch | 1 pouch | HeartStart XL: M4735A |
| 989803172501 | Black MRx Carry Case | 1 box = 1 case | HeartStart MRx M3535A, M3536A |
| 989803171281 CE | Accessory Storage Bags | 1 pack = 1 Therapy Accessory Pouch & 1 Monitoring Accessory Pouch | HeartStart XL+ Defibrillator/Monitor 861290 |
| 989803171291 CE | Cable Retention Straps | 1 pack = 3 cable wraps | HeartStart XL+ Defibrillator/Monitor 861290 |
| 989803180871 | MRx Black Complete Carrying Case | 1 pack = 1 complete carrying case | M3536A |



M3541A

Capnography

IntelliVue Microstream® CO₂ Solution

| Product # | Description | Packaging | Use With: |
|--------------|--|----------------|---|
| 989803177951 | Smart CapnoLine H O ₂ Adult | 25 Filterlines | |
| CE | | | |
| 989803177961 | Smart CapnoLine H O ₂ Adult Long | 25 Filterlines | |
| CE | | | |
| 989803177971 | Smart CapnoLine H O ₂ Pediatric | 25 Filterlines | |
| CE | | | |
| 989803177981 | Smart CapnoLine H O ₂ Pediatric Long | 25 Filterlines | |
| CE | | | |
| 989803178001 | CapnoLine H O ₂ Infant/Neonatal | 25 Filterlines | |
| CE | | | |
| 989803178011 | CapnoLine H Infant/Neonatal Long | 25 Filterlines | |
| CE | | | |
| 989803178021 | Nasal FilterLine Infant/Neonatal | 25 Filterlines | |
| CE | | | |
| 989803178031 | Smart CapnoLine Guard | 25 Filterlines | |
| CE | | | |
| 989803178041 | Smart CapnoLine Guard O ₂ | 25 Filterlines | |
| CE | | | |
| 989803178051 | Smart CapnoLine Guard O ₂ Long | 25 Filterlines | |
| CE | | | |
| 989803178071 | Hook and Loop Strap | 1 box = 25 | |
| CE | | | |
| 989803179101 | Nasal Filterline O ₂ Adult | 25 Filterlines | |
| CE | | | |
| 989803179111 | Nasal Filterline O ₂ Adult Long | 25 Filterlines | |
| CE | | | |
| 989803179121 | Nasal Filterline O ₂ Pediatric | 25 Filterlines | |
| CE | | | |
| M1920A | FilterLine Set, Adult/Pedi, Intubated | 25 sets | HeartStart MRx: M3535A Opt A06, M3536A Opt A06 |
| CE | | | |
| M1921A | FilterLine H Set, Adult/Pediatric, Intubated (yellow) | 25 sets | HeartStart MRx: M3535A Opt A06, M3536A Opt A06 |
| CE | | | |



M1920A



M1921A

IntelliVue Microstream® CO₂ Solution (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| M1923A CE | FilterLine H Set, Infant/Neonatal, Intubated (yellow) | 25 sets | HeartStart MRx: M3535A Opt A06, M3536A Opt A06 |
| M2520A CE | Smart CapnoLine O ₂ , Pediatric, Non-intubated dual purpose | 25 sets | HeartStart MRx: M3535A Opt A06, M3536A Opt A06 |
| M2522A CE | Smart CapnoLine O ₂ plus, Adult/intermed, Non-intubated dual purpose | 25 sets | HeartStart MRx: M3535A Opt A06, M3536A Opt A06 |
| M2524A CE | Smart CapnoLine, Pediatric, Non-intubated single purpose | 25 sets | HeartStart MRx: M3535A Opt A06, M3536A Opt A06 |
| M2526A CE | Smart CapnoLine plus, Adult/intermediate, Non-intubated single purpose | 25 sets | HeartStart MRx: M3535A Opt A06, M3536A Opt A06 |



M1923A



M2520A



M2522A



M2524A



ECG Cables and Leads

10-Lead ECG Cables and Leads (3+7)

| Product # | Description | Packaging | Use With: |
|--------------|---|------------|--------------------------------|
| M3525A CE | 10-lead ECG patient trunk cable for 3 and 7 electrode combination sets AAMI/IEC 12-pin connector Length: 2.7m | 1 cable | HeartStart MRx: M3535A, M3536A |
| M3526A | 3-electrode limb set, snap, AAMI | 1 lead set | HeartStart MRx: M3535A, M3536A |
| M3527A | 7-lead chest lead set, snap, AAMI | 1 lead set | HeartStart MRx: M3535A, M3536A |
| M3528A CE | 3-electrode limb lead set, snap, IEC | 1 lead set | HeartStart MRx: M3535A, M3536A |
| M3529A CE | 7-electrode chest lead set, snap, IEC | 1 lead set | HeartStart MRx: M3535A, M3536A |
| M5530A CE | Combiner Plug for 3-Wire Lead Set | 1 lead set | HeartStart MRx: M3535A, M3536A |



M3525A

10-Lead ECG Cables and Leads (5+5)

| Product # | Description | Packaging | Use With: |
|--------------------|--|------------|--------------------------------|
| 989803147691 CE | 10-lead ECG patient trunk cable (5+5) AAMI/IEC 12-pin connector Length: 1.3m | 1 cable | HeartStart MRx: M3535A, M3536A |
| 989803158061 | 5-lead set, limb, AAMI | 1 lead set | HeartStart MRx: M3535A, M3536A |
| 989803158071 | 5-lead set, chest, AAMI | 1 lead set | HeartStart MRx: M3535A, M3536A |
| 989803158081 CE | 5-lead set, limb, IEC | 1 lead set | HeartStart MRx: M3535A, M3536A |
| 989803158091 CE | 5-lead set, chest, IEC | 1 lead set | HeartStart MRx: M3535A, M3536A |



989803158061 and 989803158071



989803158081 and 989803158091

Patient Monitoring ECG Cables and Leads

IntelliVue 10-Lead (5+5) ECG Monitoring – AAMI/IEC

(for conventional diagnostic 12-lead ECG, use one limb 5-lead set and one chest 5-lead set)

| Product # | Description | Packaging | Use With: |
|--------------|---|------------|--|
| M1663A CE | 10-Lead ECG Patient Trunk (5+5) – AAMI/IEC 12-pin input connector, AAMI & IEC labels included Length: 6.5' (2m) | 1 cable | HeartStart MRx: M3535A, M3536A |
| M1949A CE | 10-Lead ECG Patient Trunk (5+5) – AAMI/IEC 12-pin input connector, AAMI & IEC labels included Length: 8.9' (2.7m) | 1 cable | HeartStart MRx: M3535A, M3536A, 860338, 860335 |
| 989803176161 | Cbl 5-Lead, snap, shld, AAMI, limb rgd (5.2') | 1 lead set | 861493, 861494, 861495, M3536A, M3535A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M6, M3536MC |
| 989803176171 | Cbl 5-Lead, snap, shld, AAMI, chest, rgd (3.9') | 1 lead set | 861493, 861494, 861495, M3536A, M3535A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M6, M3536MC |
| 989803176181 | Cbl 5-Lead, snap, shld, IEC, limb, rgd (5.2') | 1 lead set | 861493, 861494, 861495, M3536A, M3535A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M6, M3536MC |
| 989803176191 | Cbl 5-Lead, snap, shld, IEC, chest, rgd | 1 lead set | 861493, 861494, 861495, M3536A, M3535A, M3536M, M3536M1, M3536M2, M3536M3, M3536M4, M3536M6, M3536MC |




M1663A



M1949A


5-Lead Sets for 10-Lead Monitoring – AAMI Color Coding

| Product # | Description | Packaging | Use With: |
|--------------------|--|--------------------------|---|
| M1968A | ICU grabber, limb, AAMI | 1 lead set | Use with M1663A or M1949A 10-lead patient trunk cable for HeartStart MRx: M3535A, M3536A, 860335, 860338 |
| M1976A | ICU grabber, chest, AAMI | 1 lead set | |
| M1644A | ICU snap, limb, AAMI | 1 lead set | |
| M1602A | ICU snap, chest, AAMI | 1 lead set | |
| M1973A | OR grabber, limb, AAMI | 1 lead set | |
| M1979A | OR grabber, chest, AAMI | 1 lead set | |
| 989803173131 CE | 5-Lead Set, Disposable, Bedside, AAMI, round | 1 pack = 20 lead sets | Compatible Defibrillators M3535A (MRx monitor/ defib), M3536A (MRx monitor/defib), M4735A (HeartStart XL) |



M1968A

M1976A




M1973A




M1979A

5-Lead Sets for 10-Lead Monitoring – IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------------|---|-----------------------|--|
| M1604A CE | ICU snap, chest, IEC | 1 lead set | Use with M1663A or M1949A 10-lead patient trunk cable for HeartStart MRx: M3535A, M3536A, 860335, 860338 |
| M1645A CE | ICU snap, limb, IEC | 1 lead set | |
| M1971A CE | ICU grabber, limb, IEC | 1 lead set | |
| M1974A CE | OR grabber, limb, IEC | 1 lead set | |
| M1978A CE | ICU grabber, chest, IEC | 1 lead set | |
| M1984A CE | OR grabber, chest, IEC | 1 lead set | |
| 989803174211 CE | 5-Lead Set, Disposable, Bedside, IEC, round | 1 pack = 20 lead sets | |



M1971A



M1974A



M1978A



M1984A

3-Lead Monitoring – AAMI/IEC

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|---|
| M1669A CE | 3-Lead ECG Patient Trunk Cable, AAMI/IEC 12-pin connector. AAMI and IEC labels included. Length: 8.9' (2.7m) | 1 cable | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |



M1669A

3-Lead Sets for 3-Lead Monitoring – AAMI/IEC Color Coding

| Product # | Description | Packaging | Use With: |
|--------------------|--|--------------------------|--|
| M1671A | ICU grabber, AAMI | 1 lead set | Use with M1669A 3-Lead ECG Patient Trunk Cable, 860335, 860338 |
| M1673A | ICU snap, AAMI | 1 lead set | |
| M1675A | OR grabber, AAMI | 1 lead set | |
| 989803173121 CE | 3-Lead Set, Disposable, Bedside, AAMI, round | 1 pack = 20 lead sets | Compatible Defibrillators M3535A (MRx Monitor/defib), M3536A (MRx Monitor/defib), M4735A (HeartStart XL) |
| M1672A CE | ICU grabber, IEC | 1 lead set | Use with M1669A 3-Lead ECG Patient Trunk Cable, 860335, 860338 |
| M1674A CE | ICU snap, IEC | 1 lead set | |
| M1678A CE | OR grabber, IEC | 1 lead set | |
| 989803174201 CE | 3-Lead Set, Disposable, Bedside, IEC, round | 1 pack = 20 lead sets | Compatible Defibrillators M3535A (MRx Monitor/defib), M3536A (MRx Monitor/defib), M4735A (HeartStart XL) |



M1671A



M1672A



M1675A



M1678A



989803174201

5-Lead Monitoring – AAMI/IEC

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|---|
| M1668A CE | 5-Lead ECG Patient Trunk Cable, AAMI/IEC 12-pin connector. AAMI and IEC labels included. Length: 8.9' (2.7m) | 1 cable | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A, 860335, 860338 |




M1668A

5-Lead Sets for 5-Lead Monitoring – AAMI/IEC Color Coding


| Product # | Description | Packaging | Use With: |
|--------------------|--|--------------------------|---|
| M1644A | ICU snap, limb, AAMI | 1 lead set | Use with M1668A 5-Lead ECG Patient Trunk Cable for HeartStart MRx: M3535A, M3536A; 860335, 860338; HeartStart XL: M4735A |
| M1645A CE | ICU snap, limb, IEC | 1 lead set | 860338, 860335 |
| M1968A | ICU grabber, limb, AAMI | 1 lead set | Use with M1668A 5-Lead ECG Patient Trunk Cable for HeartStart MRx: M3535A, M3536A; 860335, 860338; HeartStart XL: M4735A |
| M1971A CE | ICU grabber, limb, IEC | 1 lead set | 860338, 860335 |
| M1973A CE | OR grabber, limb, AAMI | 1 lead set | Use with M1668A 5-Lead ECG Patient Trunk Cable for HeartStart MRx: M3535A, M3536A; 860335, 860338; HeartStart XL: M4735A |
| M1974A CE | OR grabber, limb, IEC | 1 lead set | 860338, 860335 |
| 989803173131 CE | 5-Lead Set, Disposable, Bedside, AAMI, round | 1 pack = 20 lead sets | Compatible Bedside Monitors M1001B, M1002A, M3000A, M3001A, M8105A (MP5), 863064 (VM6), 863065 (VM6), 863066 (VM8), 863068 (VM6), Compatible Defibrillators M3535A (MRx monitor/ defib), M3536A (MRx monitor/defib), M4735A (HeartStart XL) |
| 989803174211 CE | 5-Lead Set, Disposable, Bedside, IEC, round | 1 pack = 20 lead sets | Compatible Defibrillators M3535A (MRx monitor/ defib), M3536A (MRx monitor/defib), M4735A (HeartStart XL) |



M1644A




M1645A



M1968A



M1971A



M1973A



M1974A



989803174211

ECG Electrodes

| Product # | Description | Packaging | Use With: |
|--------------------|---|----------------------|---|
| 989803148821 CE | Small Adult, radiolucent foam, solid gel ECG electrode Size: 1.4" x 1.6" (available in U.S. and Japan) | 30/pouch 600/case | HeartStart MRx: M3535A, M3536A |
| M2202A CE | Adult radiotranslucent foam solid gel electrode Size: 2" (50mm) diameter | 5/pouch 300/case | HeartStart MRx: M3535A, M3536A; HeartStart XL: M4735A |
| M4606A CE | ECG skin preparation paper | 100 preps | All ECG Electrodes |
| M4612A CE | Adult solid gel ECG electrode Shape: teardrop | 5/pouch 300/case | HeartStart MRx: M3535A, M3536A |
| M4613A CE | Adult solid gel ECG electrode Shape: teardrop | 30/pouch 300/case | HeartStart MRx: M3535A, M3536A |



989803148821



M2202A



M4606A



M4612A/M4613A

Paper

Recording Paper

| Product # | Description | Packaging | Use With: |
|--------------------|------------------------------|-----------|---|
| 989803138171 CE | MRx Wide Printer Paper | 10 rolls | M3535A, M3536A, M3536J Opt B04 (75mm printer) |
| 989803138181 CE | MRx Wide Printer Paper | 80 rolls | M3535A, M3536A, M3536J Opt B04 (75mm printer) |
| 40457C CE | 1-Channel Chem/Thermal Paper | 10 rolls | M4735A, M3535A, M3536A, M3636J |
| 40457D CE | 1-Channel Chem/Thermal Paper | 80 rolls | M4735A, M3535A, M3536A, M3536J |



989803138181



40457C

Invasive Blood Pressure

TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits – Single Line

| Product # | Description | Packaging | Use With: |
|--------------------------------|--|-----------|---|
| 989803177901 CE 0123 | Single DPT TP4 60" (152cm) with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Dash Tubing | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803181211 CE 0123 | Single DPT TP4 60" (152cm) with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Premium Stripe Tubing | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179771 CE 0123 | Single DPT TP4 72" (183cm) with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Dash Tubing | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803181221 CE 0123 | Single DPT TP4 72" (183cm) with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Premium Stripe Tubing | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179721 CE 0123 | Single DPT TP4 84" (213cm), with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Dash Tubing | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803181231 CE 0123 | Single DPT TP4 84" (213cm), with Drip Chamber, Two 3-Way Stopcocks, 3ml/hr Flush Device, Red Premium Stripe Tubing | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |



98903177901



989803181211

TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits – Single Line *(continued)*

| Product # | Description | Packaging | Use With: |
|--------------------------------|---|-----------|---|
| 989803179871 CE 0123 | Single DPT TP4 9" (23cm) Patient Mount with Drip Chamber, One 3-Way Stopcock, 3ml/hr Flush Device | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803181141 CE 0123 | Single DPT TP4 24" (61cm) Compartmental Pressure, One 4-Way Stopcock | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |



9809803179871



989803181141

TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits with SAFESET™ Blood Sampling/Conservation System – Single Line

| Product # | Description | Packaging | Use With: |
|-------------------------------|--|-----------|--|
| 989803179741 CE0123 | Single DPT TP4 84" (213cm) w/SAFESET, 2 Cannula Sampling Ports, In-line Reservoir | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179751 CE0123 | Single DPT TP4 84" (213cm) w/SAFESET, 2 Luer Sampling Ports, In-line Reservoir | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179861 CE0123 | Single DPT TP4 24" (61cm) w/SAFESET, Patient Mount, 1 Cannula Sampling Port, In-Line Reservoir | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179761 CE0123 | Single DPT TP4 60" (152cm) w/SAFESET, 1 Cannula Sampling Port, In-Line Reservoir | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179781 CE0123 | Single DPT TP4 60" (152cm) w/SAFESET, 2 Cannula Sampling Ports, In-Line Reservoir | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179791 CE0123 | Single DPT TP4 60" (152cm) w/SAFESET, 2 Luer Sampling Ports, In-Line Reservoir | 20 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |



989803179751



98980179761



989803179781

TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits – Multi-Lines

| Product # | Description | Packaging | Use With: |
|-------------------------|--|-----------|---|
| 989803177911 CE 0123 | Dual DPT TP4 72" (183cm) with Drip Chamber, Four 3-Way Stopcocks, Two 3ml/hr Flush Devices, Red/Blue Dash Tubing | 10 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179731 CE 0123 | Triple DPT TP4 60" (152cm) with Drip Chamber, Six 3-Way Stopcocks, Three 3ml/hr Flush Devices, Red/Blue/Yellow Dash Tubing | 10 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M8102A, M8105A |
| 989803181241 CE 0123 | Dual DPT TP4 72" (183cm) with Drip Chamber, Four 3-Way Stopcocks, Two 3ml/hr Flush Devices, Red/Blue Premium Stripe Tubing | 10 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803181251 CE 0123 | Triple DPT TP4 60" (152cm) with Drip Chamber, Six 3-Way Stopcocks, Three 3ml/hr Flush Devices, Red/Blue/Yellow Premium Stripe Tubing | 10 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M8102A, M8105A |



989803177911



989803179731



989803181241




989803181251

TRANSPAC® 4 Disposable Pressure Transducer Monitoring Kits with SAFESET™ Blood Sampling/Conservation System – Multi- Lines

| Product # | Description | Packaging | Use With: |
|------------------------|---|-----------|---|
| 989803179801 CE0123 | Dual DPT TP4 84" (213cm) w/SAFESET, (Red/Blue Dash Tubing), 2 Cannula Sampling Ports, In-Line Reservoir | 10 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179811 CE0123 | Dual DPT TP4 84" (213cm) w/SAFESET, (Red/Blue Dash Tubing), 2 Luer Sampling Ports, In-Line Reservoir | 10 kits | M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |



989803179801



989803179811

SAFESET™ Disposable Blood Sampling/Conservation System Only

| Product # | Description | Packaging | Use With: |
|------------------------|---|-----------|---|
| 989803179891 CE0123 | SAFESET kit w/ 84" (213cm) Tubing, Two In-Line Cannula Sampling Ports, 10ml In-Line Reservoir | 20 kits | TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803180851 CE0123 | SAFESET kit w/ 60" (152cm) Tubing, One In-line Cannula Sampling Ports, 10ml In-Line Reservoir | 20 kits | TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |



989803180851

TRANSPAC® 4 Reusable Transducer Cables

| Product # | Description | Packaging | Use With: |
|------------------------|--|-----------|---|
| 989803177921 CE0123 | Transpac 4 Cable for Philips Monitors 10' (3m) | 1 cable | TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179941 CE0123 | Transpac 4 Cable for Philips Monitors 15' (4.5m) | 1 cable | TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |
| 989803179951 CE0123 | Transpac 4 Cable for Philips Monitors 15' (4.5m), bulk package | 20 cables | TRANSPAC 4 Disposable Transducer Kits with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |



989803179941

TRANSPAC® 4 Disposable Transducer Kit Accessories

| Product # | Description | Packaging | Use With: |
|------------------------|---|------------------------------|--|
| 989803177931 CE0123 | 3-Position Transducer Mount, Reusable | 28 mounts | TRANSPAC 4 Disposable Transducer Kits |
| 989803177941 CE0123 | Transducer Holder for IV Pole, Reusable | 5 holders | TRANSPAC 4 Disposable Transducer Kits |
| 989803179901 CE0123 | SAFESET Mounting Bracket, Reusable | 30 brackets | TRANSPAC 4 Disposable Transducer Kits and SAFESET System |
| 989803179911 CE0123 | Single Transducer Mount, Reusable | 30 mounts | TRANSPAC 4 Disposable Transducer Kits |
| 989803179921 CE0123 | SAFESET Shielded Blunt Cannula | 50 blunt cannulas | TRANSPAC 4 Disposable Transducer Kits and SAFESET System |
| 989803179931 CE0123 | SAFESET Shielded Blunt Cannula w/ Blood Tube Holder | 25 cannula with tube holders | TRANSPAC 4 Disposable Transducer Kits and SAFESET System |
| 989803180861 CE0123 | Transducer Simulator, Reusable | 1 simulator | TRANSPAC 4 Disposable Transducer Kits and Cables with instruments: M1006B, M3001A, M3002A, M3012A, M3014A, M3015A, M3015B, M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M8102A, M8105A |



989803177931



989803177941



989803179901



989803179921

Memscap® Reusable Invasive Blood Pressure Solution

| Product # | Description | Packaging | Use With: |
|----------------|--------------------------------|-----------------|---|
| CPJ84046 CE | Pressure Transducer Holder | 1 holder | CPJ840J6 transducer with HeartStart MRx: M3535A, M3536A, M3015A |
| CPJ84447 CE | Pressure Transducer Pole Mount | 1 mount | CPJ840J6 transducer with HeartStart MRx: M3535A, M3536A, M3015A |
| CPJ840J6 CE | Reusable Pressure Transducer | 1 transducer | HeartStart MRx: M3535A, M3536A, M3015A |
| CPJ84022 CE | Disposable Domes | 50 domes | CPJ840J6 transducer with HeartStart MRx: M3535A, M3536A, M3015A |
| CPJ84026 CE | Monitoring Line Set | 10 sets per bag | 863065, 863066, 863068, 863087, 863088, M3536M, M3002A, M8102A, M3535A, M8105AT, M8105A, M3536A, M3015A, M3015B |



CPJ840J6

NIBP Cuffs and Air Hoses


Multi-Patient Traditional Cuffs – Single Hose

| Product # | Description | Packaging | Use With: |
|-------------------------|---|-----------|--|
| REVISED 40401A CE | Infant Bladder width = 2.4" (6cm) Bladder length = 4.7" (12cm) Limb circumference = 10 – 19cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED 40401B CE | Pediatric Bladder width = 3.5" (8.9cm) Bladder length = 7.1" (18cm) Limb circumference = 18 – 26cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED 40401C CE | Adult Bladder width = 4.8" (12.1cm) Bladder length = 9.1" (23cm) Limb circumference = 25 – 35cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED 40401D CE | Large Adult Bladder width = 6" (15.2cm) Bladder length = 13" (33cm) Limb circumference = 33 – 47cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED 40401E CE | Thigh Bladder width = 7" (17.8cm) Bladder length = 14.2" (36cm) Limb circumference = 46 – 66cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |



Multi-Patient Comfort Cuffs – Single Hose

| Product # | Description | Packaging | Use With: |
|-------------------------|---|-----------|--|
| REVISED M1572A CE | Pediatric, color coded: green tubing Bladder width: 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M1573A CE | Small Adult, color coded: light blue tubing Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference: 20.5 – 28cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M1574A CE | Adult, color coded: navy blue tubing Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference: 27 – 35cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M1575A CE | Large Adult, color coded: burgundy tubing Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference: 34 – 43cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |



Multi-Patient Easy Care Cuffs – Single Hose

| Product # | Description | Packaging | Use With: |
|--------------------------|--|-----------|--|
| REVISÉD M4552B CE | Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10 – 15cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4552B5 CE | Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10 – 15cm | 5 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4553B CE | Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4553B5 CE | Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm | 5 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |



Multi-Patient Easy Care Cuffs – Single Hose (continued)

| Product # | Description | Packaging | Use With: |
|--------------------------|---|-----------|--|
| REVISÉD M4554B CE | Small Adult, light blue Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5 – 28.5cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4554B5 CE | Small Adult, light blue Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5 – 28.5cm | 5 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4555B CE | Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5 – 36cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4555B5 CE | Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5 – 36cm | 5 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4556B CE | Adult Extra-Long, navy blue Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5 – 36cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4556B5 CE | Adult Extra-Long, navy blue Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5 – 36cm | 5 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4557B CE | Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35 – 45cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4557B5 CE | Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35 – 45cm | 5 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISÉD M4558B CE | Large Adult Extra-Long, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35 – 45cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |

Multi-Patient Easy Care Cuffs – Single Hose (continued)

| Product # | Description | Packaging | Use With: |
|--------------------------|--|-----------|--|
| REVISED M4558B5 CE | Large Adult Extra-Long, burgundy Bladder width: 6.7" (17cm) Bladder length: 15" (38cm) Limb circumference = 35 – 45cm | 5 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M4559B CE | Thigh, grey Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm) Limb circumference = 44 – 56cm | 1 cuff | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M4559B5 CE | Thigh, grey Bladder width: 8.3" (21cm) Bladder length: 18.9" (48cm) Limb circumference = 44 – 56cm | 5 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |









Multi Patient Disposable Cuffs – Multi Care

| Product # | Description | Packaging | Use With: |
|----------------------------------|--|-----------|--|
| NEW 989803183311 CE | Infant, orange Bladder width = 2.2" (5.5cm) Bladder length = 6.1" (15.4cm) Limb circumference = 10 – 15cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803183321 CE | Pediatric, green Bladder width = 3.1" (8cm) Bladder length = 8.5" (21.5cm) Limb circumference = 14 – 21.5cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803183331 CE | Small Adult, light blue Bladder width: 4.1" (10.5cm) Bladder length: 9.4" (24cm) Limb circumference = 20.5 – 28cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803183341 CE | Adult, navy blue Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference = 27 – 35cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803183351 CE | Adult XL, navy blue with edge striping Bladder width: 5.1" (13cm) Bladder length: 11" (28cm) Limb circumference = 27 – 35cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803183361 CE | Large Adult, burgundy Bladder width: 6.3" (16cm) Bladder length: 13.8" (35cm) Limb circumference = 34 – 43cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803183371 CE | Thigh, grey Bladder width: 7.9" (20cm) Bladder length: 16.9" (43cm) Limb circumference: 42 – 54cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |

Single Care Disposable Cuffs

| Product # | Description | Packaging | Use With: |
|----------------------------------|--|-----------|--|
| NEW 989803182281 CE | Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm | 20 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803182291 CE | Small Adult, light blue Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5 – 28.5cm | 20 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803182301 CE | Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5 – 36cm | 20 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803182311 CE | Adult XL, navy blue with edge striping Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5 – 36cm | 20 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| NEW 989803182321 CE | Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length 15" (38cm) Limb circumference = 35 – 45cm | 20 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |

Gentle Care Disposable Cuffs – Single Hose

| Product # | Description | Packaging | Use With: |
|--|---|-----------|--|
| REVISED M4572B  | Infant, orange Bladder width 2.2" (5.5cm) Bladder length: 5.9" (15cm) Limb circumference: 10 – 15cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M4573B  | Pediatric, green Bladder width 3.1" (8cm) Bladder length: 8.5" (21.5cm) Limb circumference: 14 – 21.5cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M4574B  | Small Adult, light blue Bladder width: 4.2" (10.6cm) Bladder length 9.4" (24cm) Limb circumference = 20.5 – 28.5cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M4575B  | Adult, navy blue Bladder width: 5.3" (13.5cm) Bladder length 11.8" (30cm) Limb circumference = 27.5 – 36cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M4576B  | Adult Extra-Long, navy blue Bladder width: 5.3" (13.5cm) Bladder length: 11.8" (30cm) Limb circumference = 27.5 – 36cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M4577B  | Large Adult, burgundy Bladder width: 6.7" (17cm) Bladder length 15" (38cm) Limb circumference = 35 – 45cm | 10 cuffs | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |

Gentle Care Disposable Cuffs – Single Hose *(continued)*

| | | | |
|--|---|-----------------|---|
| <p>REVISED</p> <p>M4578B</p> <p>CE</p> | <p>Large Adult Extra-Long, burgundy</p> <p>Bladder width: 6.7" (17cm)</p> <p>Bladder length: 15" (38cm)</p> <p>Limb circumference = 35 – 45cm</p> | <p>10 cuffs</p> | <p>HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495</p> |
| <p>REVISED</p> <p>M4579B</p> <p>CE</p> | <p>Thigh, grey</p> <p>Bladder width: 8.3" (21cm)</p> <p>Bladder length 18.9" (48cm)</p> <p>Limb circumference = 44 – 56cm</p> | <p>10 cuffs</p> | <p>HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495</p> |



Blood Pressure Air Hoses

| Product # | Description | Packaging | Use With: |
|-------------------------|---------------------------------------|-----------|---|
| REVISED M1598B CE | NIBP Air Hose Length: 4.92' (1.5m) | 1 cable | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |
| REVISED M1599B CE | NIBP Air Hose Length: 9.84' (3.0m) | 1 cable | HeartStart MRx: M3535A, M3536A, M3536M, M3536M2, M3536M3, M3536M4, M3536M5, M3536M6, M3536MC, 861290, 861493, 861494, 861495 |



M1599B

SpO₂ Sensors and Cables

Philips Disposable SpO₂ Sensors

| Product # | Description | Packaging | Use With: |
|--------------|---|------------|---|
| M1131A CE | Adult/Pediatric Finger Clip Disposable Sensor (single patient use) Length = 45cm Patient size: >44lb (>20kg) | 20 sensors | HeartStart Defibrillators: M3535A, M3536A, M4735A, 860335, 860338, 863276, 863278 |



M1131A

Philips Reusable SpO₂ Sensors


| Product # | Description | Packaging | Use With: |
|---------------|--|-----------|---|
| M1191B CE | Adult Finger Glove Sensor Length = 2m Patient size: >110lb (>50kg) 8-pin connector | 1 sensor | HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278 |
| M1191BL CE | Adult Finger Glove Sensor Length = 3m Patient size: >110lb (>50kg) 8-pin connector | 1 sensor | HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278 |
| M1191T CE | Adult Finger Glove Sensor Length = 0.45m Patient size: >110lb (>50kg) 9-pin (D-sub) connector | 1 sensor | HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278 |
| M1192A CE | Pediatric/Small-Adult Glove Sensor Length = 1.5m Patient size: 33lb to 110lb (15kg to 50kg) 8-pin connector | 1 sensor | HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278 |



M1191B



M1191T



M1192A

Philips Reusable SpO₂ Sensors (continued)

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|---|
| M1192T CE | Pediatric/Small-Adult Glove Sensor Length = 0.45m Patient size: 33 – 110lb (15 – 50kg) 9 pin (D-sub) connector | 1 sensor | HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278 |
| M1194A CE | Pediatric/Adult Ear Clip Sensor Length = 1.5m Patient size: >88lb (>40kg) 8 pin connector | 1 sensor | HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278 |
| M1195A CE | Reusable Infant Finger Glove Sensor Length = 1.5m Patient size: 8.8 – 33lb (4 – 15kg) 8-pin connector | 1 sensor | HeartStart MRx only: M3535A, M3536A, 863276, 863278 |
| M1196A CE | Pediatric/Adult Finger Clip Sensor Length = 3m Patient size: >88lb (>40kg) 8 pin connector | 1 sensor | HeartStart Defibrillators: M3535A, M3536A, M4735A, 860335, 860338, 863276, 863278 |
| M1196S CE | Reusable Clip Adult SpO ₂ Sensor Length = 2m Patient size: > 88lb (>40kg) 8-pin connector | 1 sensor | HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278 |
| M1196T CE | Pediatric/Adult Finger Clip Sensor Length 90cm Patient size: >88lb (>40kg) 9-pin (D-sub) connector | 1 sensor | HeartStart Defibrillators: M3535A, M3536A, M4735A, 863276, 863278 |
| M1627A CE | Wristband Velcro straps for SpO ₂ glove sensors | 10 straps | M1191B/BL/T sensors, all ALS Defibrillators |



M1194A



M1195A



M1196A

SpO₂ Adapter Cables for Philips Sensors

| Product # | Description | Packaging | Use With: |
|--------------|---|-----------|--|
| M1941A CE | SpO ₂ Extension Cable Adds additional length to 8-pin sensors with 8-pin sockets Length: 6.6' (2m) | 1 cable | To add additional length to M1191B, M1192A, M1194A, M1195A sensors only, connecting to HeartStart Defibrillators: M3535A, M3536A, M4735A, 860335, 860338, 863276, 863278 |
| M1943A CE | SpO ₂ Adapter Cable Adapts 9-pin sensors to 8-pin sockets Length: 3.3' (1m) | 1 cable | Adapter cable for M1191T, M1192T and M1196T sensors only, connecting to HeartStart Defibrillators: M3535A, M3536A, M4735A, 860335, 860338, 863276, 863278 |



M1941A



M1943A


Temperature Probes

Disposable Temperature Probes

| Product # | Description | Packaging | Use With: |
|--------------|---|-------------------|--|
| M1837A CE | Esophageal/Rectal Temperature Probe, 9 FR Length: 36" (91cm), sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| 21090A CE | Esophageal/Rectal Temperature Probe, 12 FR Length: 36" (91cm) sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| 21091A CE | Disposable sterile skin surface temperature probe Length: 30" (76cm) Shelf life: at least 12 months | 1 box = 20 probes | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| 21093A CE | Esophageal/Stethoscope Temperature Probe, 12 FR Length: 20" (51cm) sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| 21094A CE | Esophageal/Stethoscope Temperature Probe, 18 FR Length: 20" (51cm) sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| 21095A CE | Esophageal/Stethoscope Temperature Probe, 24 FR Length: 20" (51cm) sterilized, continuous monitoring Shelf life: at least 12 months | 20 probes | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |



21090A



21091A



21094A

Disposable Temperature Probes *(continued)*

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| 21096A CE | Foley Catheter Temperature Probe, 16FR Length: 17" (43cm) sterilized, continuous monitoring Shelf life: at least 12 months | 10 probes | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| 21097A CE | Foley Catheter Temperature Probe, 18FR Length: 17" (43cm) sterilized, continuous monitoring Shelf life: at least 12 months | 10 probes | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| M2255A CE | Foley Catheter Temperature Probe, 14FR Length: 17" (43cm) sterilized, continuous monitoring Shelf life: at least 12 months | 10 probes | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |



21096A



21097A



M2255A

Disposable Temperature Probe Extension Cables

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| 21082A CE | Disposable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 10' (3m) | 1 cable | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| 21082B CE | Disposable temperature probe extension cable Accepts 1/8-inch mini phone plug. (2-pin input connector) Length: 5' (1.5m) | 1 cable | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |



21082A

Reusable Temperature Probes for Continuous Monitoring (meets EN12470 Standard)

| Product # | Description | Packaging | Use With: |
|--------------|--|-----------|--|
| 21075A CE | Esophageal/Rectal Temperature Probe, 12FR (autoclavable) 2-prong plug, reusable, Length: 10' (3m) | 1 probe | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| 21076A CE | Esophageal/Rectal Temperature Probe, 10FR (autoclavable) 2-prong plug, reusable, continuous monitoring Length: 10' (3m) | 1 probe | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, M3015A |
| 21078A CE | Skin Surface Temperature Probe (autoclavable) 0.37" (9.5mm) diameter 2-prong plug, reusable tape-on, continuous monitoring Length: 10' (3m) | 1 probe | HeartStart MRx: M3535A Opt A06, M3536A Opt A06, 860335, 860338, M3015A |



21075A



21076A



21078A

AED Supplies

External AED Pads

External AED Pads and Adapters

| Product # | Description | Packaging | Use With: |
|--------------------|-------------------------------------|-----------|---|
| M5071A CE | Adult SMART Pads Cartridge | 1 set | HeartStart HS1: M5066A, M5067A, M5068A |
| M5072A CE | Infant/Child SMART PADS Cartridge | 1 set | HeartStart HS1: M5066A, M5067A, M5068A |
| 989803139261 CE | SMART Pads II | 1 set | HeartStart FRx: 861304, 861305 |
| 989803158211 CE | HeartStart Adult Defibrillator Pads | 1 set | ForeRunner: E, S, EM FR series: 940010XX, 940020XX, 940030XX FR2 series: M3840A, M3841A, M3860A, M3861A |
| 989803158221 CE | HeartStart Adult Defibrillator Pads | 5 sets | ForeRunner: E, S, EM FR series: 940010XX, 940020XX, 940030XX FR2 series: M3840A, M3841A, M3860A, M3861A |



M5071A



M5072A



989803139261

External AED Pads and Adapters *(continued)*


| Product # | Description | Packaging | Use With: |
|----------------|--|-----------|---------------------------------|
| M3870A CE | Infant/Child Reduced-Energy Defibrillator Pads | 1 set | HeartStart FR2+: M3860A, M3861A |
| 05-10000 CE | Medtronic™ Pads Adapter | 1 adapter | 989803158211, 989803158221 |
| 05-10100 CE | Zoll™ Pads Adapter | 1 adapter | 989803158211, 989803158221 |



M3870A



05-10000



05-10100


AED Accessories and Kits

AED Fast Response Kit

| Product # | Description | Packaging | Use With: |
|-----------|-------------------|--|---|
| 68-PCHAT | Fast Response Kit | 1 kit: 2 pairs/gloves, 1 CPR mask, 1 scissors, 1 prep razor, 1 paper towel, 1 zippered pouch | HeartStart HS1: M5066A, M5067A, M5068A ForeRunner AED: E, S, EM HeartStart FR: 940010XX, 940020XX, 940030XX HeartStart FR2: M3860A, M3861A, M3840A, M3841A HeartStart FRx: 861304, 861305 |

AED Accessories

| Product # | Description | Packaging | Use With: |
|--------------------|--|--------------------------------|---|
| 989803139311 CE | Infant/Child Key. Under 8 years old (<55lb or <25kg) | 1 key | HeartStart FRx: 861304, 861305 |
| M3873A CE | FR2+ ECG Assessment Module, AAMI | 1 module | HeartStart FR2+: M3860A, M3840A |
| M3874A CE | FR2+ ECG Assessment Module, IEC | 1 set electrodes User guide | HeartStart FR2+: M3860A, M3840A |
| 989803137771 | Replacement ECG Monitoring Electrodes for the ECG Assessment Module | 17 packs of 3 electrodes | M3873A FR2+ ECG Assessment Module – AAMI M3874A FR2+ ECG Assessment Module – IEC |
| 861476 | AED Awareness Posters 11" (28cm) w, 17" (43cm) h | 1 package = 4 posters | 861304 (FRx), M5066A (861282-OnSite), M3860A (861264-FR2), M3861A (861265-FR2) |
| 861478 | AED Signage Bundle. AED Wall Sign (989803170921), AED Awareness Placard (989803170901) and AED Awareness Poster 4-pack (861476 – English). Wall Sign: 9" (23cm) h, 6.1" (15cm) d, Awareness Placard: 10.25" (26cm) w, 4. Available in US only | 1 bundle | FRx (861304), FR2 (M3860A/861264, M3861A/861265), OnSite (M5066A/861282) |
| 861477 | AED Wall Mount and Signage Bundle. Wall Mount: 10.5" (27cm) w, 8" (20cm) h, 6.9" (17cm) d; Wall Sign: 9" (23cm) h, 6.1" (15cm) d; Awareness Placard: 10.25" (26cm) w, 4.5" (11cm) h; Awareness Posters: 11" (28cm) w, 17" (43cm) h. Available in US only | 1 bundle | FRx (861304), FR2 (M3860A/861264, M3861A/861265), OnSite (M5066A/861282) |



989803139311



M3873A

AED Data Management Solutions

AED Data Cards and Data Management Accessories

| Product # | Description | Packaging | Use With: |
|--------------------|--|----------------------|--|
| M3524A | Data Card Reader | 1 card | CodeRunner SW: M3831A, M3833A; HeartStart FR2 series data card/tray: M3854A; HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A; |
| M3834A CE | HeartStart Event Review Software – Single License | 1 CD with license | HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A HeartStart HS1: M5066A, M5067A HeartStart FRx: 861304, 861305 |
| M3853A CE | AED Tray for Data Card | 1 tray | HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A |
| M3854A CE | Data Card and Tray | 1 card and tray | |
| SDCF05 CE | Data Card Adapter | 1 adapter | |
| 989803141811 CE | HeartStart Event Review Software – Organization-wide License PC | 1 CD with license | HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A HeartStart HS1: M5066A, M5067A HeartStart FRx: 861304, 861305 |
| ACT-IR CE | Infrared Data Cable | 1 cable | HeartStart HS1: M5066A, M5067A, M5068A, HeartStart FRx: 861304, 861305 HeartStart Event Review Pro SW: 861276 |

AED Training Materials

AED Training Materials

| Product # | Description | Packaging | Use With: |
|--------------------|---|--|---|
| M5088A | Internal Manikin Adapter | 1 set | HeartStart Trainer: M5085A and M5086A (for CPR/AED manikins) HeartStart HS1: M5073A Adult Training Pads Cartridge and M5074A Infant/Child Training Pads Cartridge |
| M5089A CE | External Manikin Adapter | 1 set | HeartStart Trainer: M5085A and M5086A (for CPR/AED manikins) HeartStart HS1: M5073A Adult Training Pads Cartridge and M5074A Infant/Child Training Pads Cartridge |
| M5093A CE | Replacement Adult Training Pads with Replaceable Seal | 1 set | HeartStart HS1: M5073A (Adult Training Pads Cartridge) |
| M5094A CE | Replacement Infant/Child Training Pads | 1 set | HeartStart HS1: M5074A (Infant/Child Training Pads Cartridge) |
| M3855A CE | Charger for FR2 Training & Admin Pack | 1 charger | M3864A (Training and Admin Pack for FR2 series) |
| M3864A CE | Training & Administration Pack for FR2 series Defibrillator | 1 set | HeartStart FR2 and FR2+: M3860A, M3861A, M3840A, M3841A |
| 861306 CE | HeartStart FRx Trainer | Trainer, case, 1 set pads, 1 external manikin adapter 1 IFU | 989803139311 (Infant/child key) M5088A (Internal manikin adapter) M5089A (External Manikin Adapter) 989803139271 (HeartStart Training Pads II) 989803139291 (Replacement training pads II) M5090A (Adult pad placement guide) 989803139281 (Infant/child pad placement guide) |
| 989803139271 CE | Training Pads II | 1 cartridge | HeartStart FRx: 861304, 861305 HeartStart FRx Trainer: 861306 |



M3864A




861306




989803139271

AED Training Materials (continued)

| Product # | Description | Packaging | Use With: |
|-------------------------|--|---|---|
| 989803139281 CE | Infant/Child Pad Placement Guide | 1 guide | HeartStart FRx Trainer: 861306 HeartStart Trainer: M5085A AED Trainer 2: M3752A |
| 989803139321 | HeartStart FRx Defibrillator Training Toolkit | 1 kit (CD & DVD) | HeartStart FRx: 861304, 861305 |
| 989803150291 | HeartStart FR2 Defibrillator Training Toolkit | 1 kit (CD & DVD) | HeartStart FR2+: M3860A, M3861A, M3840A, M3841A |
| 07-10900 CE | Defibrillator Training Pads | 1 set | M3752A M3840A, M3841A, M3860A and M3861A when equipped with Training and Administration Pack (M3864A) |
| 07-11000 CE | Defibrillator Trainer 2 Carrying Case | 1 case | M3752A (AED Trainer 2) |
| M3752A CE | AED Trainer 2 | FR2 AED trainer, 1-set of 07-10900, quick ref. card | 07-10900, 07-11000, M3753A, M3754A, M3755A, M3756A, M3871A FR2 Infant/Child Training Pads |
| M3753A CE | Remote Control for AED Trainer 2 | 1 remote | M3752A (AED Trainer 2) M3756A (AED Little Anne Training System) |
| M3754A CE | Programming Kit for AED Trainer 2 | CD-Rom training SW, serial port connector cable | M3752A (AED Trainer 2) M3756A (AED Little Anne Training System) |
| M3755A CE | AED Little Anne Training Pads | 1 set | M3752A (AED Trainer 2) M3756A (AED Little Anne Training System) |
| M3756A CE (LATEX) | AED Little Anne Training System | Contains: AED Trainer 2, 1 remote, 1 AED Little Anne training manikin | M3755A Training pads for AED Little Anne |
| M3860-97800 | Quick Reference Card, FR2 Defibrillator, English | 1 card | FR2: M3860A, M3861A, M3868A, M3840A |
| M3871A | FR2 Infant/Child Training Pads | 1 set | FR2: M3860A, M3861A, M3840A, M3841A when used with M3864A (Training and Admin Pack) M3752A (AED Trainer 2) |



07-10900



M3871A

AED Training Materials (continued)

| Product # | Description | Packaging | Use With: |
|--------------|--------------------------------------|---|---|
| M5066-89100 | Training Toolkit | 1 CD & DVD | HeartStart HS1: M5066A, M5067A |
| M5066-97800 | Quick Reference Guide, OnSite | 1 guide | HeartStart HS1: M5066A, M5067A, M5068A |
| M5073A CE | Adult Training Pads Cartridge | 1 cartridge | HeartStart HS1: M5066A, M5067A, M5068A |
| M5074A CE | Infant/Child Training Pads Cartridge | 1 cartridge | HeartStart HS1: M5066A, M5067A, M5068A |
| M5085A CE | HeartStart Trainer | 1 case, one set of reusable adult training pads, one external manikin adapter | M5073A, M5074A, M5085-91900, M5087A, M5090A, 989803139281 |
| 989803139291 | Replacement Training Pads II | 1 set | HeartStart FRx Trainer: 861306 |
| M5090A CE | Adult Pad Placement Guide | 1 guide/per package | M5066A, M5068A |
| 989803173221 | HS1 Trainer, Cantonese | 1 guide | M5085A |
| 989803173231 | HS1 Trainer, Traditional Chinese | 1 guide | M5085A |
| 989803173241 | HS1 Trainer, Simplified Chinese | 1 guide | M5085A |
| 989803173341 | HS1 Trainer, Korean | 1 guide | M5085A |
| 989803173461 | FRx Trainer, Cantonese | 1 guide | 861306 |
| 989803173471 | FRx Trainer, Traditional Chinese | 1 guide | 861306 |
| 989803173481 | FRx Trainer, Simplified Chinese | 1 guide | 861306 |
| 989803173561 | FRx Trainer, Korean | 1 guide | 861306 |



M5073A



M5074A



M5085A

AED Batteries

External AED Batteries

| Product # | Description | Packaging | Use With: |
|--------------------|--|--------------------|--|
| 989803136291 CE | FR2 Series Aviation Battery | 1 battery | HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A |
| 989803136301 CE | ForeRunner Aviation Battery | 1 battery | ForeRunner AED: E, S, EM HeartStart FR AED: 940010XX, 940020XX, 94030XX |
| 989803139301 CE | FRx Aviation Battery | 1 battery | HeartStart FRx: 861304, 861305 |
| BT1 CE | ForeRunner Disposable LitMnO ₂ Battery Pack | 1 battery | ForeRunner AED: E, S, EM HeartStart FR AED: 940010XX, 940020XX, 94030XX |
| M3848A CE | FR2+ Rechargeable Li Ion Battery | 1 battery | HeartStart FR2+ Series: M3840A, M3841A, M3860A, M3861A, M3849A |
| M3849A CE | Charger for the Rechargeable Li Ion Battery | 1 charger | HeartStart FR2+ Series: M3840A, M3841A, M3860A, M3861A, M3848A |
| M3863A CE | Long-Life LiMnO ₂ Battery | 1 battery | HeartStart FR2+ Series: M3840A, M3841A, M3860A, M3861A |
| M5070A CE | HeartStart Four-Year Battery | 1 battery | HeartStart HS1: M5066A, M5067A, M5068A FRx: 861304, 861305 |
| 989803167281 CE | Lithium Ion Battery Module | 1 Pack = 1 Battery | HeartStart XL+ Defibrillator/Monitor 861290 |



M3848A



M3863A



M5070A

AED Cases/Mounts/Cabinets

AED Wall Mounts/Accessories

| Product # | Description | Packaging | Use With: |
|--------------|------------------------------|-------------|--|
| M3857A CE | Wall Mount Bracket | 1 mount | HeartStart FRx: 861304, 861305 HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A HeartStart HS1: M5066A, M5067A ForeRunner: E, S, EM |
| M3859A CE | Secure-Pull Seal | 1 box of 10 | HeartStart FRx: 861304, 861305 HeartStart HS1: M5066A ForeRunner: E, S, EM HeartStart FR series: 940010XX, 940020XX, 94030XX HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A Wall Bracket: M3857A (Wall Mount Bracket) Cases: M3868A, M3869A, M5075A, M5076A, YC Cabinets: PFE7023D, PFE7024D |
| 989803170891 | AED Wall Mount Bracket | 1 unit | FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282 |
| 989803170921 | AED Wall Sign, red | 1 sign | FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282 |
| 989803170931 | AED Wall Sign, green | 1 sign | FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282 |
| 989803170901 | AED Awareness Placard, red | 1 placard | FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282 |
| 989803170911 | AED Awareness Placard, green | 1 placard | FRx 861304, FR2 M3860A/861264, M3861A/861265, OnSite M5066A/861282 |



M3857A



M3859A

AED Cases/Wall Cabinets

| Product # | Description | Packaging | Use With: |
|--------------------|--|-----------|---|
| 989803133171 | HeartStart Temperature Control Carrying Case | 1 case | HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A |
| 989803136531 CE | Defibrillator Cabinet (basic) | 1 cabinet | HeartStart FRx: 861304, 861305 HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A HeartStart HS1: M5066A, M5067A ForeRunner: E, S, EM |
| 989803139251 CE | FRx Carrying Case | 1 case | HeartStart FRx: 861304, 861305 |
| HC | ForeRunner Semi-Rigid Carrying Case | 1 case | ForeRunner: E, S, EM HeartStart FR series: 940010XX, 940020XX, 94030XX |
| M3868A | Carrying Case for FR2 series | 1 case | HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A |
| M3869A CE | Vinyl Carrying Case for FR2 Series | 1 case | HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A |
| M5075A CE | Standard Carry Case for HeartStart HS1 | 1 case | HeartStart HS1: M5066A, M5067A, M5068A |
| M5076A CE | Slim carrying case for HeartStart HS1 | 1 case | HeartStart HS1: M5066A, M5067A, M5068A |
| M5087A | Replacement Carry Case, HeartStart Trainer | 1 case | HeartStart Trainer: M5085A |



M3868A



M5076A

AED Cases/Wall Cabinets *(continued)*


| Product # | Description | Packaging | Use With: |
|----------------|---|-----------|--|
| PFE7023D CE | Defibrillator Cabinet (semi-recessed) | 1 cabinet | HeartStart FRx: 861304, 861305 HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A HeartStart HS1: M5066A, M5067A ForeRunner: E, S, EM |
| PFE7024D CE | Defibrillator Cabinet (surface mount) | 1 cabinet | HeartStart FRx: 861304, 861305 HeartStart FR2 series: M3840A, M3841A, M3860A, M3861A HeartStart HS1: M5066A, M5067A ForeRunner: E, S, EM |
| YC | Carrying Case, Plastic Waterproof Shell | 1 case | HeartStart FRx: 861304, 861305 ForeRunner: E, S, EM HeartStart FR series: 940010XX, 940020XX, 94030XX HeartStart FR2 series: M3860A, M3861A, M3840A, M3841A, M3854A HeartStart HS1: M5066A, M5067A |

Index

| | | | | | |
|--------------------|----|--------------------|----|--------------------|----|
| 05-10000..... | 8 | 989803137771 | 54 | 989803173341 | 58 |
| 05-10000..... | 53 | 989803138171 | 26 | 989803173461 | 58 |
| 05-10100..... | 8 | 989803138181 | 26 | 989803173471 | 58 |
| 05-10100..... | 53 | 989803139251 | 61 | 989803173481 | 58 |
| 05-10200..... | 8 | 989803139261 | 52 | 989803173561 | 58 |
| | | 989803139271 | 56 | 989803174201 | 23 |
| 07-10900..... | 57 | 989803139281 | 57 | 989803174211..... | 22 |
| 07-11000..... | 57 | 989803139291 | 58 | 989803174211..... | 24 |
| | | 989803139301 | 59 | 989803176161..... | 20 |
| 68-PCHAT | 54 | 989803139311 | 54 | 989803176171..... | 20 |
| | | 989803139321 | 57 | 989803176181..... | 20 |
| 21075A | 51 | 989803141811 | 55 | 989803176191..... | 20 |
| 21076A | 51 | 989803147691 | 19 | 989803176541 | 13 |
| 21078A | 51 | 989803147711 | 11 | 989803177901 | 27 |
| 21082A | 50 | 989803148821 | 25 | 989803177911 | 30 |
| 21082B | 50 | 989803150291 | 57 | 989803177921 | 33 |
| 21090A | 49 | 989803158061 | 19 | 989803177931 | 34 |
| 21091A | 49 | 989803158071 | 19 | 989803177941 | 34 |
| 21093A | 49 | 989803158081 | 19 | 989803177951 | 16 |
| 21094A | 49 | 989803158091 | 19 | 989803177961 | 16 |
| 21095A | 49 | 989803158211 | 52 | 989803177971 | 16 |
| 21096A | 50 | 989803158221 | 52 | 989803177981 | 16 |
| 21097A | 50 | 989803158661 | 12 | 989803178001 | 16 |
| | | 989803162401 | 12 | 989803178011 | 16 |
| 40401A | 36 | 989803163291 | 12 | 989803178021 | 16 |
| 40401B | 36 | 989803166021 | 7 | 989803178031 | 16 |
| 40401C | 36 | 989803167281 | 59 | 989803178041 | 16 |
| 40401D | 36 | 989803170891 | 60 | 989803178051 | 16 |
| 40401E | 36 | 989803170901 | 60 | 989803178071 | 16 |
| 40457C | 26 | 989803170911 | 60 | 989803179101 | 16 |
| 40457D | 26 | 989803170921 | 60 | 989803179111..... | 16 |
| | | 989803170931 | 60 | 989803179121 | 16 |
| 861306 | 56 | 989803171261..... | 11 | 989803179721 | 27 |
| 861476..... | 54 | 989803171271 | 8 | 989803179731 | 30 |
| 861477..... | 54 | 989803171281 | 15 | 989803179741 | 29 |
| 861478..... | 54 | 989803171291 | 15 | 989803179751 | 29 |
| | | 989803171701 | 13 | 989803179761 | 29 |
| 453564063841 | 13 | 989803172501 | 15 | 989803179771 | 27 |
| 453564063851 | 13 | 989803173121 | 23 | 989803179781 | 29 |
| | | 989803173131 | 21 | 989803179791 | 29 |
| 989803133171 | 61 | 989803173131 | 24 | 989803179801 | 31 |
| 989803136291 | 59 | 989803173221 | 58 | 989803179811 | 31 |
| 989803136301 | 59 | 989803173231 | 58 | 989803179861 | 29 |
| 989803136531 | 61 | 989803173241 | 58 | 989803179871 | 28 |

| | | | | | |
|--------------|-----|--------|-----|--------|-----|
| 989803179891 | .32 | M1191T | .46 | M1968A | .21 |
| 989803179901 | .34 | M1192A | .46 | M1968A | .24 |
| 989803179911 | .34 | M1192T | .47 | M1971A | .22 |
| 989803179921 | .34 | M1194A | .47 | M1971A | .24 |
| 989803179931 | .34 | M1195A | .47 | M1973A | .21 |
| 989803179941 | .33 | M1196A | .47 | M1973A | .24 |
| 989803179951 | .33 | M1196S | .47 | M1974A | .22 |
| 989803180851 | .32 | M1196T | .47 | M1974A | .24 |
| 989803180861 | .34 | M1572A | .37 | M1976A | .21 |
| 989803180871 | .15 | M1573A | .37 | M1978A | .22 |
| 989803181141 | .28 | M1574A | .37 | M1979A | .21 |
| 989803181211 | .27 | M1575A | .37 | M1984A | .22 |
| 989803181221 | .27 | M1598B | .45 | M2202A | .25 |
| 989803181231 | .27 | M1599B | .45 | M2255A | .50 |
| 989803181241 | .30 | M1602A | .21 | M2520A | .17 |
| 989803181251 | .30 | M1604A | .22 | M2522A | .17 |
| 989803182281 | .42 | M1627A | .47 | M2524A | .17 |
| 989803182291 | .42 | M1644A | .21 | M2526A | .17 |
| 989803182301 | .42 | M1644A | .24 | M3501A | .7 |
| 989803182311 | .42 | M1645A | .22 | M3504A | .7 |
| 989803182321 | .42 | M1645A | .24 | M3506A | .14 |
| 989803183311 | .41 | M1663A | .20 | M3507A | .8 |
| 989803183321 | .41 | M1668A | .24 | M3508A | .8 |
| 989803183331 | .41 | M1669A | .23 | M3516A | .14 |
| 989803183341 | .41 | M1671A | .23 | M3524A | .55 |
| 989803183351 | .41 | M1672A | .23 | M3525A | .18 |
| 989803183361 | .41 | M1673A | .23 | M3526A | .18 |
| 989803183371 | .41 | M1674A | .23 | M3527A | .18 |
| | | M1675A | .23 | M3528A | .18 |
| ACT-IR | .55 | M1678A | .23 | M3529A | .18 |
| | | M1741A | .9 | M3537A | .13 |
| BT1 | .59 | M1742A | .9 | M3538A | .14 |
| | | M1743A | .9 | M3541A | .15 |
| CPJ840J6 | .35 | M1744A | .9 | M3543A | .9 |
| CPJ84022 | .35 | M1781A | .8 | M3544A | .11 |
| CPJ84026 | .35 | M1783A | .10 | M3549A | .13 |
| CPJ84046 | .35 | M1789A | .9 | M3713A | .7 |
| CPJ84447 | .35 | M1837A | .49 | M3716A | .7 |
| | | M1920A | .16 | M3717A | .7 |
| HC | .61 | M1921A | .16 | M3718A | .7 |
| | | M1923A | .17 | M3719A | .7 |
| M1131A | .46 | M1941A | .48 | M3725A | .8 |
| M1191B | .46 | M1943A | .48 | M3752A | .57 |
| M1191BL | .46 | M1949A | .20 | M3753A | .57 |

| | | | | | |
|-------------|-----|-------------|-----|----------|-----|
| M3754A | .57 | M4606A | .25 | PFE7023D | .62 |
| M3755A | .57 | M4612A | .25 | PFE7024D | .62 |
| M3756A | .57 | M4613A | .25 | | |
| M3834A | .55 | M4737A | .13 | SDCF05 | .55 |
| M3848A | .59 | M4740A | .10 | | |
| M3849A | .59 | M4741A | .10 | YC | .62 |
| M3853A | .55 | M4742A | .10 | | |
| M3854A | .55 | M4743A | .10 | | |
| M3855A | .56 | M4744A | .10 | | |
| M3857A | .60 | M4745A | .9 | | |
| M3859A | .60 | M4746A | .9 | | |
| M3860-97800 | .57 | M4747A | .14 | | |
| M3863A | .59 | M4748A | .8 | | |
| M3864A | .56 | M4748A | .10 | | |
| M3868A | .61 | M4751A | .15 | | |
| M3869A | .61 | M4759A | .9 | | |
| M3870A | .53 | M4761A | .12 | | |
| M3871A | .57 | M4762A | .12 | | |
| M3873A | .54 | M4763A | .12 | | |
| M3874A | .54 | M5066-89100 | .58 | | |
| M4552B | .38 | M5066-97800 | .58 | | |
| M4552B5 | .38 | M5070A | .59 | | |
| M4553B | .38 | M5071A | .52 | | |
| M4553B5 | .38 | M5072A | .52 | | |
| M4554B | .39 | M5073A | .58 | | |
| M4554B5 | .39 | M5074A | .58 | | |
| M4555B | .39 | M5075A | .61 | | |
| M4555B5 | .39 | M5076A | .61 | | |
| M4556B | .39 | M5085A | .58 | | |
| M4556B5 | .39 | M5087A | .61 | | |
| M4557B | .39 | M5088A | .56 | | |
| M4557B5 | .39 | M5089A | .56 | | |
| M4558B | .39 | M5090A | .58 | | |
| M4558B5 | .40 | M5093A | .56 | | |
| M4559B | .40 | M5094A | .56 | | |
| M4559B5 | .40 | M5521A | .13 | | |
| M4572B | .43 | M5522A | .13 | | |
| M4573B | .43 | M5523A | .13 | | |
| M4574B | .43 | M5524A | .13 | | |
| M4575B | .43 | M5525A | .13 | | |
| M4576B | .43 | M5526A | .10 | | |
| M4577B | .43 | M5528A | .13 | | |
| M4578B | .44 | M5529A | .13 | | |
| M4579B | .44 | M5530A | .18 | | |



© 2013 Koninklijke Philips N.V. All rights are reserved.
Philips Healthcare is part of Royal Philips

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands
4522 962 95491 * JUN 2013

How to reach us
www.philips.com/healthcare
healthcare@philips.com

Product information
www.philips.com/healthcarestore