

3Dpro-zir[®]

— More than nature —

Color: A1 Height: 18mm

CE 0197

Aidite[®]

ZIRCONIA

Aidite[®]

Recommend Indication: Veneer, Inlay, Anterior Crown, Posterior Crown, Full Crown Bridge, Screw Retained Full Crown Bridge, Implant



Colors: A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, Hollywood Bleach

750 Mpa



1050 Mpa



57% Translucency



43% Translucency



High Strength 1050MPa

Layerless natural transition

High Transparent 57%

Optimum Esthetics

- All the benefits of porcelain layered to zirconia without the cost with no graying
- Multi-layer translucency for optimum esthetics
- Perfect color intensity for darker preps and implant abutments

Optimum Strength

- High flexural strength of > 1,050 MPa's
- Maximum fracture toughness • Indicated for all tooth positions and bridgework



Bio Link Medical Surgical Equipment & Instruments Trading

Address: Office 301& 302, 3rd Floor Jumeirah Terrace Building - Jumeirah 1 - Dubai

Phone: 04 359 6242, Website: www.biolinkdubai.com, Email: info@biolinkdubai.com