

BISCO's Intraoral Repair Kit



BISCO's Intraoral Repair Kit provides all the materials necessary to repair the broken or damaged restorations intraorally. Once the operative site is properly prepared, an esthetic and durable repair can be successfully achieved utilizing a resin composite material. The Intraoral Repair Kit is a cost effective option as a practical restorative solution. The Intraoral Repair Kit is complete with BISCO's Z-PRIME™ Plus, Porcelain Primer, Porcelain Bonding Resin, Dual-Cured Opaquer, Porcelain Etchant (9.5% HF) and Barrier Gel.



INDICATIONS FOR USE

- Repair of Porcelain & Lithium Disilicate Restorations
- Repair of Porcelain-Fused-to-Metal or Zirconia/Alumina Restorations
- Repair of Direct or Indirect Composite Restorations



BISCO products are proudly made in the U.S.A.

INTRAORAL REPAIR KIT

B-22300K

1 Bottle Porcelain Primer (3ml), 1 Syringe 9.5% Porcelain Etchant (5g), 1 Syringe Barrier Gel (5g), 1 Bottle Z-PRIME™ Plus (2ml), 1 Bottle Porcelain Bonding Resin (6ml), 1 Bottle Opaquer Catalyst (3ml), 1 Bottle Opaquer Base Universal (3ml), Accessories, Instructions/MSDS



Bio Link Medical Surgical Equipment & Instruments Trading
 Address: Office 301 & 302, 3rd Floor Jumeirah Terrace Building -
 Jumeirah 1 - Dubai Phone: 04 359 6242, Website: www.biolinkdubai.com,
 Email: info@biolinkdubai.com