


RETREATMENT

HyFlex™ Remover

NEW



The essential piece for your retreatment

remover.coltene.com

 **COLTENE**

- One file to remove the canal obturation
- Respects the root canal anatomy
- Efficient without any solvent
- Safe thanks to non-cutting tip



ORDER INFORMATION



30/.07

Maximum torque: **2.5 N.cm**

Speed of rotation: **400-800 rpm**

HyFlex Remover, 5 pcs

19 mm


23 mm

30/.07 Remover File



60023649

60023648



◀ For best results complete the procedure with HyFlex EDM shaping files. You are also free to use a different file system such as HyFlex CM or system of your choice.



For more information about HyFlex, please take a look at the brochure or visit our website hyflex.coltene.com.

© COLTENE – www.coltene.com



Coltene/Whaledent Ltd.
The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF / UK
T +44 1444 235486
F +44 1444 870640
info.uk@coltene.com

Coltene/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten / Switzerland
T +41 71 757 5300
F +41 71 757 5301
info.ch@coltene.com

Coltene/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223 / USA
T +1 330 916 8800
F +1 330 916 7077
info.us@coltene.com

COLTENE