



LP — 20

**H Y G I E N E**

**A U S T R I A**

## **Respirator Mask, FFP2**

1. Self and third party protection
2. Great comfort
3. High breathability due to high quality non-woven materials



1. Wash your hands using soap and water or hand sanitizer.



2. Pull the loops of the mask back to fit them around your ears.



3. Pull the upper and lower ends of the mask over nose and chin so that no open gaps remain between face and mask.



4. Firmly press the nose clip in your mask across nose and cheeks to fit.



5. Avoid touching the mask while wearing.

# Mask Specifications

Manufacturer	- Hygiene Austria LP GmbH
Model	- Respirator Mask, FFP2
Type	- HA PP02
Device Category	- PPE-Category III
Classification	- FFP2 NR
Quantity/Units	- 1 pc., 2 pc., 20 pc., 40 pc., 350 pc., 400 pc. (in various types of packaging)
Colour	- White
Instructions and Durability	- Read user instructions - Shelf life: See Mask - Non-reusable, disposable - It is recommended to remove and dispose of the mask every 4 hours or when humid
Storage Conditions	- Store above -20°C and below 40°C in a dry place without direct sunlight - Store below 80% relative humidity
Size	- 16 x 11 cm folded (slight deviations possible)
Weight	- ca. 5,3 g (slight deviations possible)
Certificates	- Declaration of Conformity (EU 2016/425) - EN 149:2001+A1:2009 - Oekotex 100
Product specifications	- Made of multiple layers of non-woven material (PP) - extra strong and flexible double core nose clip