

CARRYALL & CUSHION LEATHER

DON'T THINK TWICE,
JUST BRING IT ALONG!



movinglife[®]

ATTO Accessories are created specifically to enhance, support and work with the ATTO mobility scooters in all user modes; while driving, while wheeling along as trolley and when folded for storage or transportation.

SCAN TO LEARN
HOW TO COMBINE
OUR ACCESSORIES

SEATING & SUPPORT



ATTO /ATTO
SPORT Armrests



Seat cushion



Seat cushion
leather



Luxury Swivel
Seat Leather



Belt & Support
Strap



Luxury Orthopedic
Backrest

STORAGE



Essentials
Pouch



Front & rear
basket



Carryall
& Cushion



Carryall
& Cushion leather



Pet Carrier &
Cushion Black



Flight Kit



Oxygen Tank
Holder and Case

TECHNICAL



Cane / Crutches
Holders




Mobile Phone
Holder



Rear Light Lifting
Handle



ATTO Left Hand
Throttle Kit



Standard Flight
Battery




Launch - H2 2023

Long distance
battery - not flight
approved



Cable kit

How to Use CARRYALL & CUSHION LEATHER



01

Slip the cushion over the seat (use handles to pull if needed)




02

Simply pull hook on elastic strap to fold under and connect to ATTO chassis

NOTE!

Take off before folding ATTO



03

One-hand "pull to open" secure funnel access, so nothing falls out. Quick "no zip" flap, to ensure everything stays inside when closing.



04

Attach to ATTO Seat Cushion dedicated upper zip