

# LUXURY SWIVEL SEAT LEATHER

EASY ACCESS EVERY TIME

**movinglife**<sup>®</sup>



**ATTO Accessories** are created specifically to enhance, support and work with the ATTO mobility scooters in all user modes; while driving, while wheeling along as trolley and when folded for storage or transportation.

SCAN TO LEARN HOW TO COMBINE OUR ACCESSORIES

## SEATING & SUPPORT



ATTO /ATTO SPORT Armrests



Seat cushion



Seat cushion leather



Luxury Swivel Seat Leather



Belt & Support Strap



Luxury Orthopedic Backrest

## STORAGE



Essentials Pouch



Front & rear basket



Carryall & Cushion



Carryall & Cushion leather



Pet Carrier & Cushion Black



Flight Kit



Oxygen Tank Holder and Case

## TECHNICAL



Cane / Crutches Holders




Mobile Phone Holder



Rear Light Lifting Handle



ATTO Left Hand Throttle Kit



Standard Flight Battery




Launch - H2 2023  
Long distance battery - not flight approved



Cable kit


# How to Use

## LUXURY SWIVEL SEAT LEATHER




01

Fit the cushion over the bottom of the seat and pass the connecting strap through the gap in the seat




02

Fasten the connection strap to the bottom of the cushion



03

The product is ready to use. Must be removed before folding



04

Attach to ATTO Carryall & Cushion dedicated upper zip