




DESS

> DENTAL SMART SOLUTIONS



Quick Reference
Guide 2023

PRODUCT LIST



	HEALING ABUTMENTS	04		DESSLoc®	17		C-Base®	31		CASTABLES	43
	SCAN ABUTMENTS	06		MULTI-UNIT ABUTMENTS	19		ABUTMENTS	33		TI-BASE CASTABLES	45
	TRANSFERS	09		TI-BASES	23		SCREWS	34		AURUM CASTABLES	49
	ANALOGUES	12		AURUMBase®	27		PRE-MILLED BLANKS & BLOCKS	38		HOLDERS	51
	TEMPORARY ABUTMENTS	16		ELLIPTIBase®	30		CoCr BASES	41		TOOLS	52

Healing Abutments

Our healing abutments have been manufactured in titanium Grade 5 ELI. They are identified with the platform and gingival height and have the same connection as the implant screw.
- Manufactured in titanium Grade 5 ELI.



A) Single Crown
B) For Bridges & structures


		UniGrip®	(3.0) ^A	NP(3.5) ^A	NP(3.5) ^B	RP(4.3/5.0) ^A	RP(4.3/5.0) ^B	WP(5.5) ^A	WP(5.5) ^B	
Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	3,0mm	20.068/3*	20.039/3	20.047/3	20.040/3	20.048/3	20.069/3*	20.070/3*	
		4,0mm	20.068/4*	20.039/4	20.047/4	20.040/4	20.048/4	20.069/4*	20.070/4*	
		5,0mm	20.068/5*	20.039/5	20.047/5	20.040/5	20.048/5	20.069/5*	20.070/5*	
	Nobel Replace Select™ Trilobe	3,0mm	UniGrip®	NP (3.5)	RP (4.3)	WP (5.0)				
		4,0mm	20.004/3	20.005/3	20.006/3					
		5,0mm	20.004/4	20.005/4	20.006/4					
	Nobel Brånemark® External Hex Universal	3,0mm	UniGrip®	NP (3.5)	RP (4.1)	WP (5.0)				
		4,0mm	20.001/3	20.002/3	20.003/3					
		5,0mm	20.001/4	20.002/4	20.003/4					
	Multi-Unit® Multiunit Connection	3,0mm	UniGrip®	NP (3.5)	RP (4.1)	WP (5.0)				
		4,0mm	20.001/5	20.002/5	20.003/5					
		5,0mm	20.001/5	20.002/5	20.003/5					
Torx® Ball		20.100	20.054							
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)					
		3,0mm	20.017/3	20.018/3	20.019/3					
		4,0mm	20.017/4	20.018/4	20.019/4					
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)					
		3,0mm	20.011/3	20.012/3	20.013/3					
		4,0mm	20.011/4	20.012/4	20.013/4					
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)					
		3,0mm	20.014/3	20.015/3	20.016/3					
		4,0mm	20.014/4	20.015/4	20.016/4					
	Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0			
			3,0mm		20.065/3	20.066/1				
			4,0mm		20.065/4	20.066/4				
Biohorizons® Implants Internal Internal Hex BH		Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)				
		3,0mm		20.087/3	20.088/3	20.089/3				
		4,0mm		20.087/4	20.088/4	20.089/4				
Biohorizons® Implants External External Hex BH		Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)					
		3,0mm	20.020/3*	20.021/3*	20.022/3*					
		4,0mm	20.020/4*	20.021/4*	20.022/4*					

		Hex ø 1,27mm	EV (3.0)	EV (3.6)	EV (4.2)	EV (4.8)	EV (5.4)	
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	3,5mm	ø4 20.406/3	ø4 20.059/3	ø5 20.060/3	ø5 20.061/3	ø5 20.407/3	
		4,5mm	ø4 20.406/4	ø4 20.059/4	ø5 20.060/4	ø6,5 20.061/4	ø6,5 20.407/4	
		6,5mm	ø4 20.406/6	ø4 20.059/6	ø5 20.060/6	ø6,5 20.061/6	ø6,5 20.407/6	
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)			
		3,0mm		20.024/3	20.025/3			
		4,0mm		20.024/4	20.025/4			
		5,0mm		20.024/5	20.025/5			
		Anatomic vers.						
		4,0mm		ø4,5 20.424/4				
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5		
		3,0mm	20.036/3	20.037/3	20.038/3			
		4,0mm	20.036/4	20.037/4	20.038/4			
Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm							
	1,5mm		20.046/1					
	2,5mm		20.046/3					
MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8			
	3,0mm		20.095/3	20.096/3				
	4,0mm		20.095/4	20.096/4				
	5,0mm		20.095/5	20.096/5				
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)			
		2,0mm		20.009/2	20.045/2			
		3,0mm		20.009/3	20.045/3			
	Straumann® Bone Level Conical BL	On SynOcta®			20.010	20.049		
		4,0mm		20.009/4	20.045/4			
		5,0mm						
	Straumann® BLX Conical BLX	Torx®	NP (3.3)	RP (4.1)	NP/RP			
		3,0mm		20.041/3	20.042/3			
		4,0mm		20.041/4	20.042/4			
		5,0mm		20.041/5	20.042/3			
		Torx®	RB (3.5)	WB (4.5)	RB/WB			
		0,8mm ^{A/B}		20.463/1				
1,5mm ^A			20.461/2					
1,5mm ^B			20.462/2					
2,5mm ^A		20.463/2						
2,5mm ^B		20.461/3						
3,5mm ^A		20.462/3						
3,5mm ^B		20.461/4						
		20.462/4						

^A For Single Crown Restoration
^B For Bridges and Structures Restoration

Scan Abutments

Take a digital impression with our desktop, intraoral or Ti-based scanbodies. Our comprehensive range, combined with our free libraries for ExoCad, 3Shape and Dental Wings, allows you to design with unique precision.



Brand	Model	Scan Body	Part Numbers					Engagement
			3.0	NP (3.5)	RP (4.3)	WP (5.5)	6.0	
Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	50.068	50.041	50.042	50.069		
		Internal Hex	52.068	52.041	52.042	52.069		
		53.100-P10	53.100-P10	53.100-P10	53.100-P10	53.100-P10	non-engaging	
		53.200-P10	53.200-P10	53.200-P10	53.200-P10	53.200-P10	engaging	
	Nobel Replace Select™ Trilobe	UniGrip	50.004	50.005	50.006	50.097		
		Internal Hex	52.004	52.005	52.006	52.097		
		53.100-P10	53.100-P10	53.100-P10	53.100-P10	53.100-P10	non-engaging	
		53.200-P10	53.200-P10	53.200-P10	53.200-P10	53.200-P10	engaging	
	Nobel Brånemark® External Hex Universal	UniGrip	50.001	50.002	50.003			
		Internal Hex	52.001	52.002	52.003			
		53.100-P10	53.100-P10	53.100-P10	53.100-P10		non-engaging	
		53.200-P10	53.200-P10	53.200-P10	53.200-P10		engaging	
Multi-Unit® Multiunit Connection	UniGrip®	50.007	50.047					
	Internal Hex	52.007	52.047					
	53.107	53.154						
	53.107	53.155						
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	50.017	50.018	50.019			
		Internal Hex	52.017	52.018	52.019			
		53.100-P10	53.100-P10	53.100-P10		non-engaging		
		53.200-P10	53.200-P10	53.200-P10		engaging		
	3i Osseotite® External Hex USA	Hex ø 1,27mm	50.011	50.012	50.013			
		Internal Hex	52.011	52.012	52.013			
		53.100-P10	53.100-P10	53.100-P10		non-engaging		
		53.200-P10	53.200-P10	53.200-P10		engaging		
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	50.014	50.015	50.016			
		Internal Hex	52.014	52.015	52.016			
		53.100-P10	53.100-P10	53.100-P10		non-engaging		
		53.200-P10	53.200-P10	53.200-P10		engaging		
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	50.059	50.060	50.061			
		On Uniabutment®	52.406	52.059	52.060	52.061	52.407	50.058
	Internal Hex	On Uniabutment®	52.058					52.058
		53.100-P10	53.100-P10	53.100-P10	53.148-P10	53.148-P10	non-engaging	
		53.200-P10	53.200-P10	53.200-P10	53.248-P10	53.248-P10	engaging	
		53.200-P10	53.200-P10	53.200-P10	53.248-P10	53.248-P10	engaging	

Brand	Model	Scan Body	Part Numbers				Engagement
			NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)	5.5	
Astra Tech Osseospeed™ Implants	Internal Hex Conic	Hex ø 1,27mm	50.023	50.024	50.025		
		Internal Hex	52.023	52.024	52.025		
		20°	50.036				
		45°	50.037				
	Internal Hex FD	Hex ø 1,22mm	50.038	50.039	50.040	50.071	
		Internal Hex	52.038	52.039	52.040	52.071	
		53.100-P10	53.100-P10	53.100-P10	53.100-P10	53.100-P10	non-engaging
		53.200-P10	53.200-P10	53.200-P10	53.200-P10	53.200-P10	engaging
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm				50.049	
		Internal Hex				52.049	
		53.100-P10	53.100-P10	53.100-P10	53.100-P10	53.100-P10	non-engaging
		53.200-P10	53.200-P10	53.200-P10	53.200-P10	53.200-P10	engaging
MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	50.056	50.017	50.018			
	Internal Hex	52.056	52.017	52.018	52.085		
	53.100-P10	53.100-P10	53.100-P10	53.100-P10		non-engaging	
	53.200-P10	53.200-P10	53.200-P10	53.200-P10		engaging	
MIS® C1 Implants Conic MC	Hex ø 1,27mm	52.056	52.017	52.018			
	Internal Hex	53.100-P10	53.100-P10	53.100-P10		non-engaging	
	53.200-P10	53.200-P10	53.200-P10	53.200-P10		engaging	
	53.200-P10	53.200-P10	53.200-P10	53.200-P10		engaging	
MIS® V3 Implants Conic MV	Hex ø 1,27mm	52.084	52.082	52.085			
	Internal Hex	53.100-P10	53.100-P10			non-engaging	
	53.200-P10	53.200-P10	53.200-P10			engaging	
	53.200-P10	53.200-P10	53.200-P10			engaging	
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	50.070	50.009	50.045		
		On Abutment		50.010	50.048		
		Internal Hex	52.070	52.009	52.045		
		On Abutment		52.010	52.048		
	Straumann® Bone Level Conical BL	Hex ø 1,27mm	53.100-P10	53.110-P10	53.148-P10		non-engaging
		Internal Hex	53.200-P10	53.210-P10	53.248-P10		engaging
		Torx®	50.043	50.044			
		On Abutment SRA			50.098		
	Straumann® BLX Conical BLX	Hex ø 1,27mm	52.043	52.044			
		Internal Hex			52.098		
		53.100-P10	53.100-P10				non-engaging
		53.200-P10	53.200-P10				engaging
Neodent Grand Morse® NEO GM	Torx®	52.062	52.062				
	Internal Hex	53.100-P10	53.100-P10			non-engaging	
	53.200-P10	53.200-P10				engaging	
	53.200-P10	53.200-P10				engaging	

Scan Abutments

Group	Implant	Hex	3.3	3.8	4.3	5.0		
Straumann® Group	Anthogyr Axiom® Implants Conic Axio	Hex ø 1,18mm	52.062					
			53.100-P10				non-engaging	
			53.200-P10				engaging	
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	50.064	50.065	50.066	50.067		
			52.064	52.065	52.066	52.067		
			53.100-P10	53.100-P10	53.100-P10	53.100-P10	non-engaging	
			53.200-P10	53.200-P10	53.200-P10	53.200-P10	engaging	
			NP (3.3) RP (3.8) WP (5.0)					
			52.074	52.075	52.076			
	Conelog® Implants Internal Con	Hex ø 1,27mm	53.100-P10	53.100-P10	53.100-P10			
			53.200-P10	53.200-P10	53.200-P10			
			NP (3.7) RP (4.2) WP (5.2)					
			50.020	50.021	50.022			
			52.020	52.021	52.022			
			53.100-P10	53.100-P10	53.100-P10		non-engaging	
BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	52.099	52.087	52.088	52.089			
			53.100-P10	53.100-P10	53.100-P10	non-engaging		
			53.200-P10	53.200-P10	53.200-P10	engaging		
		NP (3.7) RP (4.2) WP (5.2)						
		50.020	50.021	50.022				
		52.020	52.021	52.022				
BioHorizons® Implants External External Hex BH	Hex ø 1,27mm	53.100-P10	53.100-P10	53.100-P10		non-engaging		
		53.200-P10	53.200-P10	53.200-P10		engaging		
		MINI STD						
		50.072	50.073					
		52.072	52.073					
		53.100-P10	53.100-P10			non-engaging		
Hiossen/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	53.200-P10				engaging	
			53.100-P10				non-engaging	
			53.200-P10				engaging	
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm	50.057					
			52.057					
			53.100-P10				non-engaging	
			53.200-P10				engaging	
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm	50.080					
			52.080					
			53.100-P10				non-engaging	
			53.200-P10				engaging	
GlobalD®	GlobalD® Implants Conic Ink	Hex ø 1,22mm	52.062					
			53.100-P10				non-engaging	
			53.200-P10				engaging	
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28	50.090					
			52.090					
			53.100-P10				non-engaging	
			53.200-P10				engaging	

Group	Implant	Hex	3.25/3.75	4.1	4.5	5.5				
Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27	52.091				52.092	52.093	52.094	
			53.100-P10				53.100-P10	53.100-P10	53.100-P10	non-engaging
			53.200-P10				53.200-P10	53.200-P10	53.200-P10	engaging
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22	NP		RP					
			52.078							
			53.100-P10						non-engaging	
53.200-P10						engaging				
BTI®	BTI® Implants Internal BT	Hex ø 1,22	4.1							
			52.079							
			53.100-P10						non-engaging	
53.200-P10						engaging				
Phibo®	Phibo® TSH® External Hex PH	Hex ø 1,27	NP		SP		WP			
			52.030				52.031	52.031		
			53.100-P10				53.100-P10	53.100-P10		non-engaging
53.200-P10				53.200-P10	53.200-P10		engaging			

Transfers

With a body manufactured in stainless steel AISI 303, and the open and closed tray screws included manufactured in titanium, the DESS Transfer has been designed with a flat face to improve retention in the model and grooves to easily splint them together. The Transfer screw can be tightened using different method-specific system screwdriver, spatula, sharpened object or manually.



Group	Implant	Hex	3.0	NP (3.5)	RP (4.3)	WP (5.5)			
Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0 NP (3.5) RP (4.3) WP (5.5)						
			17.068				17.044	17.045	17.069
			NP (3.5) RP (4.3) WP (5.0) 6.0						
			17.004				17.005	17.006	17.097
Nobel Biocare®	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5) RP (4.1) WP (5.0)						
			17.041				17.042	17.043	
			RP (4.8) WP (6.5)						
			17.007				17.050		
Nobel Biocare®	Nobel Brånemark® External Hex Universal	UniGrip®	5mm						
			5mm 4 un.						
			17.507-P4				For closed tray only		
			8mm						
8mm 4 un.				17.807-P4					
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5) RP (4.5) WP (5.7)						
			17.017				17.018	17.019	
			NP (3.4) RP (4.1) WP (5.0)						
			17.011				17.012	17.013	

Transfers

Zimmer® 3i® Biomet®	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)	
			17.014	17.015	17.016	
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6) 17.059	EV (4.2) 17.060	EV (4.8) 17.061
			17.058			
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm		NP (3.0) 17.023	RP (3.5/4.0) 17.024	WP (4.5/5.0) 17.025
		20°	17.036			
		45°	17.037			
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5
			17.038	17.039	17.040	17.071
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm GH 1,5mm GH 2,5mm				17.049 17.250
	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	
			17.056	17.017	17.018	
MIS® C1 Implants Conic MC	Hex ø 1,27mm	NP (yellow)	SP (purple)	WP (green)	Multi-unit 4.8 17.085	
		17.081	17.082	17.083		
MIS® V3 Implants Conic MV	Hex ø 1,27mm	NP (blue)	SP (purple)			
		17.084	17.082			
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
			17.070	17.009	17.053	
		On Abutment Torx®		17.010	17.048	
			NP (3.3)	RP (4.1)		
Straumann® Bone Level Conical BL			17.046	17.047		
Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)			
		17.046	17.046			

Straumann® Group	Neodent Grand Morse® NEO GM	NEO			GM	
					17.035	
	Anthogyr Axiom® Implants Conic Axio	Hex ø 1,18mm			17.077	
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0
			17.064	17.065	17.066	17.067
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)	
			17.074	17.075	17.076	
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
			17.087	17.088	17.089	
BioHorizons® Implants External External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)		
		17.020	17.021	17.022		
Hiossen® Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD		
			17.072	17.073		
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm			17.057	
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm			17.057	
GlobalID®	GlobalID® Implants Conic Ink	Hex ø 1,22mm			17.086	
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28mm			17.090	


Transfers

Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27mm	3.25/3.75	4.1	4.5	5.5
			17.091	17.092	17.093	17.094
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22mm		NP	RP	
			17.078			
BTI®	BTI® Implants Internal BT	Hex ø 1,22mm	4.1			
			17.079			
Phibo®	Phibo® TSH® External Hex PH		NP	SP	WP	Std. 1un.
		17.030	17.031	17.032		
DESSLoc®	DESSLoc®					17.054

Analogues

Manufactured in surgical stainless steel (metal finish) or titanium Grade 5 ELI colour coded (anodised) to match the implant platform.

- Self-locking
- Excellent stability in stone models
- Digital analogue includes a screw to secure the analogue to the bottom of the model



Traditional: xx.xxx
Digital: xx.xxx/D

Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	Digital 1un.	3.0	NP (3.5)	RP (4.3)	WP (5.5)
		Std. 1un.	14.068/D	14.041/D	14.042/D	14.069/D
		Std. 10un.	14.068	14.041	14.042	14.069
	Nobel Replace Select™ Trilobe	Digital 1un.	NP (3.5)	RP (4.3)	WP (5.0)	6.0
		Std. 1un.	14.004/D	14.005/D	14.006/D	14.097/D
		Std. 10un.	14.004-P10	14.005-P10	14.006	14.097
	Nobel Brånemark® External Hex Universal	Digital 1un.	NP (3.5)	RP (4.1)	WP (5.0)	
		Std. 1un.	14.001/D	14.002/D	14.003/D	
		Std. 10un.	14.001	14.002	14.003	
	Multi-Unit® Multiunit Connection	Digital 1un.	RP (4.8)	WP (6.5)		
		Std. 1un.	14.007/D	14.047/D		
		Std. 10un.	14.007	14.047		
Std. 10un.		14.007-P10	14.047-P10			

Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	Digital 1un.	3.0	NP (3.5)	RP (4.3)	WP (5.5)			
		Std. 1un.	14.068/D	14.041/D	14.042/D	14.069/D			
		Std. 10un.	14.068	14.041	14.042	14.069			
	Nobel Replace Select™ Trilobe	Digital 1un.	NP (3.5)	RP (4.3)	WP (5.0)	6.0			
		Std. 1un.	14.004/D	14.005/D	14.006/D	14.097/D			
		Std. 10un.	14.004-P10	14.005-P10	14.006	14.097			
	Nobel Brånemark® External Hex Universal	Digital 1un.	NP (3.5)	RP (4.1)	WP (5.0)				
		Std. 1un.	14.001/D	14.002/D	14.003/D				
		Std. 10un.	14.001	14.002	14.003				
	Multi-Unit® Multiunit Connection	Digital 1un.	RP (4.8)	WP (6.5)					
		Std. 1un.	14.007/D	14.047/D					
		Std. 10un.	14.007	14.047					
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Digital 1un.	NP (3.5)	RP (4.5)	WP (5.7)				
		Std. 1un.	14.017/D	14.018/D	14.019/D				
		Std. 10un.	14.017	14.018	14.019				
	3i Osseotite® External Hex USA	Digital 1un.	NP (3.4)	RP (4.1)	WP (5.0)				
		Std. 1un.	14.011/D	14.012/D	14.013/D				
		Std. 10un.	14.011	14.012	14.013				
	3i Certain® Implants Internal Hex "Click"	Digital 1un.	NP (3.4)	RP (4.1)	WP (5.0)				
		Std. 1un.	14.014/D	14.015/D	14.016/D				
		Std. 10un.	14.014	14.015	14.016				
	Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Digital 1un.	Of Uniabutment®	EV (3.0)	EV (3.6)	EV (4.2)	EV (4.8)	EV (5.4)
			Std. 1un.	14.058/D	14.406/D	14.059/D	14.060/D	14.061/D	14.407/D
			Std. 10un.	14.058	14.406	14.059	14.060	14.061	14.407
Astra Tech Osseospeed™ Implants Internal Hex Conic		Digital 1un.	Of Uniabutment®	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)			
		Digital 1un. 20°	14.036/D	14.023/D	14.024/D	14.025/D			
		Digital 1un. 45°	14.037/D						
		Std. 1un.		14.023	14.024	14.025			
		Std. 10un.		14.023-P10	14.024-P10	14.025-P10			
		Std. 1un. 20°	14.036						
		Std. 10un. 20°	14.036-P10						
		Std. 1un. 45°	14.037						
		Std. 10un. 45°	14.037-P10						
Dentsply Friadent® Xive® Implants Internal Hex FD	Digital 1un.	NP (3.4)	RP (3.8)	WP (4.5)	5.5				
	Std. 1un.	14.038/D	14.039/D	14.040/D	14.071/D				
	Std. 10un.	14.038	14.039	14.040	14.071				
Dentsply Ankylos® C/X Implants Internal ANK	Digital 1un. ø 3,5				14.048/D				
	Digital 1un. ø 4,5				14.049/D				
	Std. 1un. ø 3,5				14.048				
	Std. 10un. ø 3,5				14.048-P10				
	Std. 10un. ø 4,5				14.049				

Analogues

Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	Digital 1un.	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	
		Digital 1un.	14.056/D	14.095/D	14.096/D	14.085/D	<i>mis</i> Multi-Unit
		Std. 1un.	14.056	14.095	14.096		
		Std. 10un.	14.056-P10	14.095-P10	14.096-P10		
						14.085	<i>mis</i>
	MIS® C1 Implants Conic MC	Digital 1un.	NP (yellow)	SP (purple)	WP (green)	Multi-unit 4.8	
		Digital 1un.	14.081/D	14.082/D	14.083/D	14.085/D	
		Std. 1un.	14.081	14.082	14.083		
	MIS® V3 Implants Conic MV	Std. 1un.				14.085	
		Digital 1un.	NP (blue)	SP (purple)	Multi-unit 4.8		
Digital 1un.		14.084/D	14.082/D		14.085/D		
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Digital 1un.	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
		Digital 1un.		14.070/D	14.009/D	14.045/D	
		Std. 1un.		14.070	14.009	14.045	On Abutment
		Std. 10un.			14.009-P10		
		Std. 1un.			14.010	14.046	On Abutment
	Straumann® Bone Level Conical BL	Std. 10un.			14.010-P10		
		Digital 1un.	Torx®	NP (3.3)	RP (4.1)	NP/RP	
		Digital 1un.		14.043/D	14.044/D		
		Std. 1un.				14.098/D	On Abutment SRA
		Std. 10un.		14.043-P10	14.044-P10		
	Straumann® BLX Conical BLX	Std. 1un.				14.098	On Abutment SRA
		Digital 1un.	Torx®	RB (3.5)	WB (4.5)		
		Std. 1un.		14.062/D	14.063/D		
	Neodent Grand Morse® NEO GM	Std. 1un.		14.062	14.063		
		Digital 1un.	NEO			GM	
Std. 1un.					14.035/D		
Anthogyr Axiom® Implants Conic Axio	Std. 1un.				14.035		
	Digital 1un.	Hex ø 1,18mm					
	Std. 1un.				14.077/D		
Henry Schein® Group	Camlog® Implants Internal Cam	Digital 1un.	Hex ø 1,27mm	3.3	3.8	4.3	5.0
		Std. 1un.		14.064/D	14.065/D	14.066/D	14.067/D
	Conelog® Implants Internal Con	Std. 1un.		14.064	14.065	14.066	14.067
		Digital 1un.	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)	
		Std. 1un.		14.074/D	14.075/D	14.076/D	
			14.074	14.075	14.076		

Henry Schein® Group	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		Digital 1un.	14.099/D	14.087/D	14.088/D	14.089/D
	Std. 1un.	14.099	14.087	14.088	14.089	
Henry Schein® Group	BioHorizons® Implants External External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)	
		Digital 1un.	14.020/D	14.021/D	14.022/D	
		Std. 1un.	14.020	14.021	14.022	
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD		
		Digital 1un.	14.072/D	14.073/D		
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm				
		Digital 1un.	14.057/D			
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm				
		Digital 1un.	14.080/D			
GlobalD®	GlobalD® Implants Conic Ink	Hex ø 1,22mm				
		Digital 1un.	14.086/D			
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28mm				
		Digital 1un.	14.090/D			
Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27mm	3.25/3.75	4.1	4.5	5.5
		Digital 1un.	14.091/D	14.092/D	14.093/D	14.094/D
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22mm		NP	RP	
		Digital 1un.	14.078/D			
BTT®	BTT® Implants Internal BT	Hex ø 1,22mm	4.1			
		Digital 1un.	14.079/D			
Phibo®	Phibo® TSH® External Hex PH	Hex ø 1,27mm	NP	SP	WP	
		Digital 1un.	14.030/D	14.031/D	14.031/D	
		Std. 1un.	14.030	14.031	14.031	

Temporary Abutments

Manufactured in titanium Grade 5 ELI. Our Temporary abutment comes with a 10,5 mm chimney height that can be cut to the required length. With parallel stretch marks to improve adhesion of the material.

- Screw included




Brand	Model	Hex	Height				Engagement
			3.0	NP (3.5)	RP (4.3)	WP (5.5)	
Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	24.041	24.042		non-engaging	
			25.041	25.042		engaging	
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	non-engaging
			24.004	24.005	24.006		engaging
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		non-engaging
			24.001	24.002	24.003		engaging
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			non-engaging
			24.007	24.030			non-engaging
			hex screw 1,20mm	24.007H120			non-engaging
			hex screw 1,22mm	24.007H122			non-engaging
			hex screw 1,27mm	24.007H127			non-engaging
	Multi-Unit® PEEK Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			non-engaging
24.707			25.754			non-engaging	
hex screw 1,20mm			24.757			non-engaging	
hex screw 1,22mm			24.767			non-engaging	
hex screw 1,27mm			24.777			non-engaging	
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)	non-engaging	
			24.017	24.018	24.019	engaging	
3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		non-engaging	
			24.011	24.012	24.013	engaging	
3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		non-engaging	
			24.014	24.015	24.016	engaging	

Brand	Model	Hex	Height				Engagement
			3.0	NP (3.5)	RP (4.3)	WP (5.5)	
Dentsply® Sirona®	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm	24.048			non-engaging	
			24.250			non-engaging	
MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	non-engaging	
			24.095*	24.096*		engaging	
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	non-engaging	
			24.026	24.045		engaging	
			25.009	25.045		non-engaging	
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP	non-engaging	
			24.043	24.044*		engaging	
			25.043	25.044		non-engaging	
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0	non-engaging
			24.065	24.066		engaging	
BioHorizons® Implants External External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)		non-engaging	
			24.020*	24.021*	24.022*	engaging	
Hiossen® Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD		non-engaging	
			24.072	24.073		engaging	

DESSLoc®

Removable overdenture attachment. With a zirconium nitride coating that produces a stronger abutment and reduces inflammation and plaque adhesion. Available in up to 5 different heights.



Brand	Model	Hex	Height				Engagement
			3.0	NP (3.5)	RP (4.3)	WP (5.5)	
Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	GH 1,0mm	88.039/1	88.040/1			
			GH 2,0mm	88.039/2	88.040/2		
			GH 3,0mm	88.039/3	88.040/3		
			GH 4,0mm	88.039/4	88.040/4		
Nobel Replace Select™ Trilobe	Hex ø 1,27mm	NP (3.5)	RP (4.3)	WP (5.0)	6.0		
			88.004/1	88.005/1			
			88.004/2	88.005/2			
			88.004/3	88.005/3			
GH 4,0mm	88.004/4	88.005/4					

DESSLoc®

			NP (3.5)	RP (4.1)	WP (5.0)		
Nobel Biocare®	Nobel Brånemark® External Hex Universal	GH 1,0mm	88.001/1	88.002/1			
		GH 2,0mm	88.001/2	88.002/2			
		GH 3,0mm	88.001/3	88.002/3			
		GH 4,0mm		88.002/4			
		GH 5,0mm		88.002/5			
Multi-Unit® Multiunit Connection			RP (4.8) 88.007	WP (6.5)			
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	GH 1,0mm	88.017/1	88.018/1	88.019/1		
		GH 2,5mm	88.017/2	88.018/2	88.019/2		
		GH 3,5mm	88.017/3	88.018/3	88.019/3		
		GH 4,5mm	88.017/4	88.018/4	88.019/4		
		GH 5,5mm	88.017/5	88.018/5			
	3i Osseotite® External Hex USA	GH 1,0mm	88.011/1	88.002/1	88.013/1		
		GH 2,0mm	88.011/2	88.002/2	88.013/2		
		GH 3,0mm	88.011/3	88.002/3	88.013/3		
		GH 4,0mm	88.011/4	88.002/4	88.013/4		
		GH 5,0mm	88.011/5	88.002/5	88.013/5		
		GH 6,0mm	88.011/6	88.002/6	88.013/6		
	3i Certain® Implants Internal Hex "Click"	GH 1,0mm	88.014/1	88.015/1	88.016/1		
		GH 2,0mm	88.014/2	88.015/2	88.016/2		
		GH 3,0mm	88.014/3	88.015/3	88.016/3		
		GH 4,0mm	88.014/4	88.015/4	88.016/4		
GH 5,0mm		88.014/5	88.015/5	88.016/5			
GH 6,0mm		88.014/6	88.015/6	88.016/5			
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	GH 1,0mm	88.059/1	88.060/1	88.061/1		
		GH 2,0mm	88.059/2	88.060/2	88.061/2		
		GH 3,0mm	88.059/3	88.060/3	88.061/3		
		GH 4,0mm	88.059/4	88.060/4	88.061/4		
	Astra Tech Osseospeed™ Implants Internal Hex Conic	GH 1,0mm	88.023/1	88.024/1	88.025/1		
		GH 2,0mm	88.023/2	88.024/2	88.025/2		
		GH 3,0mm	88.023/3	88.024/3	88.025/3		
		GH 4,0mm	88.023/4	88.024/4	88.025/4		
		GH 5,0mm			88.025/5		
	Dentsply Friadent® Xive® Implants Internal Hex FD	GH 1,0mm	88.036/1	88.037/1			
		GH 2,0mm	88.036/2	88.037/2			
		GH 3,0mm	88.036/3	88.037/3			
		GH 4,0mm	88.036/4	88.037/4			
	Dentsply Ankylos® C/X Implants Internal ANK	GH 2,0mm		88.046/2			
		GH 3,0mm		88.046/3			
		GH 4,0mm		88.046/4			

			NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	GH 1,0mm	88.017/1	88.018/1		
		GH 2,5mm	88.017/2	88.018/2		
		GH 3,5mm	88.017/3	88.018/3		
		GH 4,5mm	88.017/4	88.018/4		
		GH 5,5mm	88.017/5			
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	GH 1,0mm	NNC (3.5)	RN (4.8)	WN (6.5)	
		GH 2,0mm		88.009/1		
		GH 3,0mm		88.009/2		
		GH 4,0mm		88.009/3		
	Straumann® Bone Level Conical BL	GH 1,0mm	NP (3.3)	RP (4.1)	NP/RP	
		GH 2,0mm		88.042/1		
		GH 3,0mm	88.041/2	88.042/2		
		GH 4,0mm	88.041/3	88.042/3		
Henry Schein® Group	Camlog® Implants Internal Cam	GH 1,0mm	3.3	3.8	4.3	5.0
		GH 2,0mm		88.065/1	88.066/1	
		GH 3,0mm		88.065/2	88.066/2	
		GH 4,0mm		88.065/3	88.066/3	
	BioHorizons® Implants Internal Hex BH	GH 1,0mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		GH 2,5mm		88.087/1	88.088/1	
		GH 3,5mm		88.087/2	88.088/2	
		GH 4,5mm		88.087/3	88.088/3	
		GH 5,5mm		88.087/4	88.088/4	
				88.087/5		

Multi-Units

Available in 14 compatibilities, different gingival heights, and a choice of 0°, 17° and 30° angulation, the DESS® range offers you the most complete range in the market.



		UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	GH 1,5mm 0°		40.007/1	40.008/1		non-engaging
		GH 2,5mm 0°		40.007/2	40.008/2		non-engaging
		GH 3,5mm 0°		40.007/3	40.008/3		non-engaging
		GH 4,5mm 0°			40.008/4		non-engaging
		GH 2,5mm 17°		41.707/2	41.708/2		engaging
		GH 3,5mm 17°		41.707/3	41.708/3		engaging
		GH 4,5mm 30°		43.007/3	43.008/3		engaging
		GH 4,5mm 30°		43.007/4	43.008/4		engaging

Multi-Units

Nobel Biocare®	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		GH 1,2mm 0°	40.004/1	40.005/1	40.006/1		non-engaging
GH 2,5mm 0°	40.004/2	40.005/2	40.006/2		non-engaging		
GH 3,5mm 0°	40.004/3	40.005/3	40.006/3		non-engaging		
GH 4,5mm 0°		40.005/4			non-engaging		
GH 4,5mm 0°		40.005/5			non-engaging		
GH 2,5mm 17°	41.704/2*	41.705/2*			engaging		
GH 3,5mm 17°	41.704/3*	41.705/3*			engaging		
GH 2,5mm 17°		41.705/4*			engaging		
GH 3,5mm 30°		43.005/4*			engaging		
GH 3,5mm 30°		43.005/5*			engaging		
GH 4,5mm 30°		43.005/5R*			engaging		
Nobel Biocare®	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
		GH 1,0mm 0°	40.001/1	40.002/1	40.003/1		non-engaging
GH 2,0mm 0°	40.001/2	40.002/2	40.003/2		non-engaging		
GH 3,0mm 0°	40.001/3	40.002/4	40.003/3		non-engaging		
GH 4,0mm 0°		40.002/4			non-engaging		
GH 5,0mm 0°		40.002/5			engaging		
GH 2,0mm 17°	41.701/2*	41.702/2*			engaging		
GH 3,0mm 17°	41.701/3*	41.702/3*			engaging		
GH 4,0mm 17°		43.702/4*			engaging		
GH 4,0mm 30°		43.002/4*			engaging		
GH 5,0mm 30°		43.002/5*			engaging		
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
		GH 1,0mm 0°		40.019/1	40.020/1		non-engaging
GH 2,0mm 0°		40.019/2	40.020/2		non-engaging		
GH 3,0mm 0°		40.019/3	40.020/3		non-engaging		
GH 4,0mm 0°		40.019/4	40.020/4		non-engaging		
GH 5,0mm 0°		40.019/5	40.020/5		non-engaging		
GH 2,5mm 17°		41.719/2	41.720/2		engaging		
GH 3,5mm 17°		41.719/3	41.720/3		engaging		
GH 3,5mm 30°		43.019/3	43.020/3		engaging		
GH 4,5mm 30°		43.019/4	43.020/4		engaging		
Zimmer® 3i® Biomet®	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
		GH 1,0mm 0°	40.009/1	40.017/1			non-engaging
GH 2,0mm 0°	40.009/2	40.017/2			non-engaging		
GH 3,0mm 0°	40.009/3	40.017/3			non-engaging		
GH 4,0mm 0°		40.017/4			non-engaging		
GH 5,0mm 0°		40.017/5			non-engaging		
GH 2,0mm 17°		41.712/2*			engaging		
GH 3,0mm 17°		41.712/3*			engaging		
GH 4,0mm 17°		41.712/4*			engaging		
GH 2,0mm 30°		43.012/2*			engaging		
GH 3,0mm 30°		43.012/3*			engaging		
GH 4,0mm 30°		43.012/4*			engaging		

Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)			
		GH 1,0mm	40.059/1	40.060/1	40.061/1			non-engaging
GH 2,0mm	40.059/2	40.060/2	40.061/2			non-engaging		
GH 3,0mm	40.059/3	40.060/3	40.061/3			non-engaging		
GH 5,0mm	40.059/5	40.060/5	40.061/5			non-engaging		
Dentsply® Sirona®	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	EV (3.0)	EV (3.6)	EV (4.2)	EV (4.8)	EV (5.4)	
		GH 2,5mm 0°	40.406/1	40.027/1	40.036/1	40.037/1	40.407/1	
GH 3,5mm 0°	40.406/2	40.027/2	40.036/2	40.037/2	40.407/2		non-engaging	
GH 4,5mm 0°	40.406/3	40.027/3	40.036/3	40.037/3	40.407/3		non-engaging	
GH 2,5mm 17°		41.727/3	41.736/3	41.737/3			engaging	
GH 3,5mm 17°		41.727/4	41.736/4	41.737/4			engaging	
GH 3,5mm 30°		43.027/4	43.036/4	43.037/4			engaging	
GH 4,5mm 30°		43.027/5	43.036/5	43.037/5			engaging	
Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)			
		GH 0,5mm 20°		40.013/0	40.015/0			non-engaging
GH 1,0mm 20°		40.013/1	40.015/1				non-engaging	
GH 2,0mm 20°		40.013/2	40.015/2				non-engaging	
GH 4,0mm 20°		40.013/4	40.015/4				non-engaging	
GH 6,0mm 20°		40.013/6	40.015/6				non-engaging	
GH 0,5mm 45°		40.014/0	40.016/0				non-engaging	
GH 1,0mm 45°		40.014/1	40.016/1				non-engaging	
GH 2,0mm 45°		40.014/2	40.016/2				non-engaging	
GH 4,0mm 45°		40.014/4	40.016/4				non-engaging	
GH 6,0mm 45°		40.014/6	40.016/6				non-engaging	
Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)			
		GH 2,5mm 0°		40.024/1	40.025/1			non-engaging
GH 3,5mm 0°		40.024/2	40.025/2				non-engaging	
GH 4,5mm 0°		40.024/3	40.025/3				non-engaging	
GH 2,5mm 17°		41.724/3	41.725/3				engaging	
GH 3,5mm 17°		41.724/4	41.725/4				engaging	
GH 3,5mm 30°		43.024/4	43.025/4				engaging	
GH 4,5mm 30°		43.024/5	43.025/5				engaging	
Straumann® Group	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP			
		GH 1,5mm 0°		40.043/1	40.044/1			non-engaging
GH 2,5mm 0°		40.043/2	40.044/2				non-engaging	
GH 3,5mm 0°		40.043/3	40.044/3				non-engaging	
GH 4,5mm 0°			40.044/4				non-engaging	
GH 2,5mm 17°			41.743/2	41.744/2			engaging	
GH 3,5mm 17°			41.7043/3	41.744/3			engaging	
GH 3,5mm 30°			43.043/3	43.044/3			engaging	
GH 4,5mm 30°			43.043/4	43.044/4			engaging	

Multi-Units

Straumann® Group	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)	RB/WB		
		GH 1,5mm 0°			40.062/1		non-engaging
		GH 2,5mm 0°		40.062/2		non-engaging	
		GH 3,5mm 0°		40.062/3		non-engaging	
		GH 2,5mm 17°		41.762/2		engaging	
		GH 3,5mm 17°		41.762/3		engaging	
		GH 3,5mm 30°		43.062/3		engaging	
		GH 4,5mm 30°		43.062/4		engaging	
Neodent Grand Morse® NEO GM	NEO	GM					
	GH 1,5mm 0°	40.035/1				non-engaging	
	GH 2,5mm 0°	40.035/2				non-engaging	
	GH 3,5mm 0°	40.035/3				non-engaging	
	GH 2,5mm 17°	41.735/2				engaging	
	GH 3,5mm 17°	41.735/3				engaging	
	GH 3,5mm 30°	43.035/3				engaging	
	GH 4,5mm 30°	43.035/4				engaging	
Henry Schein® Group	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
		GH 1,0mm 0°		40.087/1	40.088/1		non-engaging
		GH 2,0mm 0°		40.087/2	40.088/2		non-engaging
		GH 3,0mm 0°		40.087/3	40.088/3		non-engaging
		GH 4,0mm 0°		40.087/4			non-engaging
		GH 5,0mm 0°		40.087/5			non-engaging
		GH 2,5mm 17°		41.787/2*	41.720/2*		engaging
		GH 3,5mm 17°		41.787/3*	41.720/3*		engaging
		GH 3,5mm 30°		43.087/3*	43.020/3*		engaging
GH 4,5mm 30°		43.087/4*	43.020/4*		engaging		
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
		GH 1,5mm 0°	40.072/1	40.073/1			non-engaging
		GH 2,5mm 0°	40.072/2	40.073/2			non-engaging
		GH 3,5mm 0°	40.072/3	40.073/3			non-engaging
		GH 2,5mm 17°	41.772/2*	41.773/2*			engaging
		GH 3,5mm 17°	41.772/3*	41.773/3*			engaging
		GH 3,5mm 30°	43.072/3*	43.073/3*			engaging
GH 4,5mm 30°	43.072/4*	43.073/4*			engaging		
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28					
		GH 1,5mm 0°	40.090/1				non-engaging
		GH 2,5mm 0°	40.090/2				non-engaging
		GH 3,5mm 0°	40.090/3				non-engaging
		GH 2,5mm 17°	41.790/2*				engaging
		GH 3,5mm 17°	41.790/3*				engaging
		GH 2,5mm 30°	43.090/2*				engaging
GH 3,5mm 30°	43.090/3*				engaging		

Ti-Bases



Manufactured in titanium Grade 5 ELI in engaging and non-engaging versions with up to three different gingival heights to best suit your needs and clinical cases. The SelectGrip® patented surface treatment ensures greater adhesion of the structure with the cement.

- Final screw included
- Up to 3 different gingival heights

Straumann® Group	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)	RB/WB		
		GH 1,5mm 0°			40.062/1		non-engaging
		GH 2,5mm 0°		40.062/2		non-engaging	
		GH 3,5mm 0°		40.062/3		non-engaging	
		GH 2,5mm 17°		41.762/2		engaging	
		GH 3,5mm 17°		41.762/3		engaging	
		GH 3,5mm 30°		43.062/3		engaging	
		GH 4,5mm 30°		43.062/4		engaging	
Neodent Grand Morse® NEO GM	NEO	GM					
	GH 1,5mm 0°	40.035/1				non-engaging	
	GH 2,5mm 0°	40.035/2				non-engaging	
	GH 3,5mm 0°	40.035/3				non-engaging	
	GH 2,5mm 17°	41.735/2				engaging	
	GH 3,5mm 17°	41.735/3				engaging	
	GH 3,5mm 30°	43.035/3				engaging	
	GH 4,5mm 30°	43.035/4				engaging	
Henry Schein® Group	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
		GH 1,0mm 0°		40.087/1	40.088/1		non-engaging
		GH 2,0mm 0°		40.087/2	40.088/2		non-engaging
		GH 3,0mm 0°		40.087/3	40.088/3		non-engaging
		GH 4,0mm 0°		40.087/4			non-engaging
		GH 5,0mm 0°		40.087/5			non-engaging
		GH 2,5mm 17°		41.787/2*	41.720/2*		engaging
		GH 3,5mm 17°		41.787/3*	41.720/3*		engaging
		GH 3,5mm 30°		43.087/3*	43.020/3*		engaging
GH 4,5mm 30°		43.087/4*	43.020/4*		engaging		
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
		GH 1,5mm 0°	40.072/1	40.073/1			non-engaging
		GH 2,5mm 0°	40.072/2	40.073/2			non-engaging
		GH 3,5mm 0°	40.072/3	40.073/3			non-engaging
		GH 2,5mm 17°	41.772/2*	41.773/2*			engaging
		GH 3,5mm 17°	41.772/3*	41.773/3*			engaging
		GH 3,5mm 30°	43.072/3*	43.073/3*			engaging
GH 4,5mm 30°	43.072/4*	43.073/4*			engaging		
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28					
		GH 1,5mm 0°	40.090/1				non-engaging
		GH 2,5mm 0°	40.090/2				non-engaging
		GH 3,5mm 0°	40.090/3				non-engaging
		GH 2,5mm 17°	41.790/2*				engaging
		GH 3,5mm 17°	41.790/3*				engaging
		GH 2,5mm 30°	43.090/2*				engaging
GH 3,5mm 30°	43.090/3*				engaging		

Straumann® Group	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)	RB/WB		
		GH 1,5mm 0°			40.062/1		non-engaging
		GH 2,5mm 0°		40.062/2		non-engaging	
		GH 3,5mm 0°		40.062/3		non-engaging	
		GH 2,5mm 17°		41.762/2		engaging	
		GH 3,5mm 17°		41.762/3		engaging	
		GH 3,5mm 30°		43.062/3		engaging	
		GH 4,5mm 30°		43.062/4		engaging	
Neodent Grand Morse® NEO GM	NEO	GM					
	GH 1,5mm 0°	40.035/1				non-engaging	
	GH 2,5mm 0°	40.035/2				non-engaging	
	GH 3,5mm 0°	40.035/3				non-engaging	
	GH 2,5mm 17°	41.735/2				engaging	
	GH 3,5mm 17°	41.735/3				engaging	
	GH 3,5mm 30°	43.035/3				engaging	
	GH 4,5mm 30°	43.035/4				engaging	
Henry Schein® Group	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
		GH 1,0mm 0°		40.087/1	40.088/1		non-engaging
		GH 2,0mm 0°		40.087/2	40.088/2		non-engaging
		GH 3,0mm 0°		40.087/3	40.088/3		non-engaging
		GH 4,0mm 0°		40.087/4			non-engaging
		GH 5,0mm 0°		40.087/5			non-engaging
		GH 2,5mm 17°		41.787/2*	41.720/2*		engaging
		GH 3,5mm 17°		41.787/3*	41.720/3*		engaging
		GH 3,5mm 30°		43.087/3*	43.020/3*		engaging
GH 4,5mm 30°		43.087/4*	43.020/4*		engaging		
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
		GH 1,5mm 0°	40.072/1	40.073/1			non-engaging
		GH 2,5mm 0°	40.072/2	40.073/2			non-engaging
		GH 3,5mm 0°	40.072/3	40.073/3			non-engaging
		GH 2,5mm 17°	41.772/2*	41.773/2*			engaging
		GH 3,5mm 17°	41.772/3*	41.773/3*			engaging
		GH 3,5mm 30°	43.072/3*	43.073/3*			engaging
GH 4,5mm 30°	43.072/4*	43.073/4*			engaging		
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28					
		GH 1,5mm 0°	40.090/1				non-engaging
		GH 2,5mm 0°	40.090/2				non-engaging
		GH 3,5mm 0°	40.090/3				non-engaging
		GH 2,5mm 17°	41.790/2*				engaging
		GH 3,5mm 17°	41.790/3*				engaging
		GH 2,5mm 30°	43.090/2*				engaging
GH 3,5mm 30°	43.090/3*				engaging		

Straumann® Group	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)	RB/WB		
		GH 1,5mm 0°			40.062/1		non-engaging
		GH 2,5mm 0°		40.062/2		non-engaging	
		GH 3,5mm 0°		40.062/3		non-engaging	
		GH 2,5mm 17°		41.762/2		engaging	
		GH 3,5mm 17°		41.762/3		engaging	
		GH 3,5mm 30°		43.062/3		engaging	
		GH 4,5mm 30°		43.062/4		engaging	
Neodent Grand Morse® NEO GM	NEO	GM					
	GH 1,5mm 0°	40.035/1				non-engaging	
	GH 2,5mm 0°	40.035/2				non-engaging	
	GH 3,5mm 0°	40.035/3				non-engaging	
	GH 2,5mm 17°	41.735/2				engaging	
	GH 3,5mm 17°	41.735/3				engaging	
	GH 3,5mm 30°	43.035/3				engaging	
	GH 4,5mm 30°	43.035/4				engaging	
Henry Schein® Group	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
		GH 1,0mm 0°		40.087/1	40.088/1		non-engaging
		GH 2,0mm 0°		40.087/2	40.088/2		non-engaging
		GH 3,0mm 0°		40.087/3	40.088/3		non-engaging
		GH 4,0mm 0°		40.087/4			non-engaging
		GH 5,0mm 0°		40.087/5			non-engaging
		GH 2,5mm 17°		41.787/2*	41.720/2*		engaging
		GH 3,5mm 17°		41.787/3*	41.720/3*		engaging
		GH 3,5mm 30°		43.087/3*	43.020/3*		engaging
GH 4,5mm 30°		43.087/4*	43.020/4*		engaging		
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
		GH 1,5mm 0°	40.072/1	40.073/1			non-engaging
		GH 2,5mm 0°	40.072/2	40.073/2			non-engaging
		GH 3,5mm 0°	40.072/3	40.073/3			non-engaging
		GH 2,5mm 17°	41.772/2*	41.773/2*			engaging
		GH 3,5mm 17°	41.772/3*	41.773/3*			engaging
		GH 3,5mm 30°	43.072/3*	43.073/3*			engaging
GH 4,5mm 30°	43.072/4*	43.073/4*			engaging		
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28					
		GH 1,5mm 0°	40.090/1				non-engaging
		GH 2,5mm 0°	40.090/2				non-engaging
		GH 3,5mm 0°	40.090/3				non-engaging
		GH 2,5mm 17°	41.790/2*				engaging
		GH 3,5mm 17°	41.790/3*				engaging
		GH 2,5mm 30°	43.090/2*				engaging
GH 3,5mm 30°	43.090/3*				engaging		

Ti-Bases

Zimmer © 3i © Biomet®		Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)	
		3i Osseotite® External Hex USA	GH 0,3mm	15.011	15.012	15.013
	GH 1,5mm	15.211	15.212	15.213		non-engaging
	GH 3,0mm	15.311	15.312	15.313		non-engaging
	GH 0,3mm	16.011	16.012	16.013		engaging
	GH 1,5mm	16.211	16.212	16.213		engaging
	GH 3,0mm	16.311	16.312	16.313		engaging
	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
3i Certain® Implants Internal Hex "Click"	GH 0,3mm	15.014	15.015	15.016		non-engaging
	GH 1,5mm	15.214	15.215	15.216		non-engaging
	GH 3,0mm	15.314	15.315	15.316		non-engaging
	GH 0,3mm	16.014	16.015	16.016		engaging
	GH 1,5mm	16.211	16.215	16.216		engaging
	GH 3,0mm	16.314	16.315	16.316		engaging

Dentsply® Sirona®		Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)	EV (5.4)	
		ASTRATECH Implants System™ EV Conic EVO	GH 1,0mm	15.059	15.060	15.061	15.407
GH 3,0mm	15.359		15.360	15.361			non-engaging
GH 1,0mm	16.059		16.060	16.061	16.407		engaging
GH 3,0mm	16.359		16.360	16.361			engaging
	On Uniabutment®	15.058					engaging
	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)			
Astra Tech Osseospeed™ Implants Internal Hex Conic	GH 0,9mm			15.025			non-engaging
	GH 1,0mm	15.023					non-engaging
	GH 1,1mm		15.024				non-engaging
	GH 2,0mm	15.023/2	15.024/2	15.025/2			non-engaging
	GH 3,0mm	15.023/3	15.024/3	15.025/3			non-engaging
	GH 0,9mm			16.025			engaging
	GH 1,0mm	16.023					engaging
	GH 1,1mm		16.024				engaging
	GH 2,0mm	16.023/2	16.024/2	16.025/2			engaging
	GH 3,0mm	16.023/3	16.024/3	16.025/3			engaging
	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5		
Dentsply Friudent® Xive® Implants Internal Hex FD	GH 0,3mm	15.038	15.039	15.040	15.071		non-engaging
	GH 1,5mm	15.238	15.239	15.240	15.271		non-engaging
	GH 3,0mm	15.338	15.339	15.340	15.371		non-engaging
	GH 0,3mm	16.038	16.039	16.040	16.071		engaging
	GH 1,5mm	16.238	16.239	16.240	16.271		engaging
	GH 3,0mm	16.338	16.339	16.340	16.371		engaging
	Hex ø 1,00mm						
Dentsply Ankylos® C/X Implants Internal ANK	GH 1,5mm			15.050			non-engaging
	GH 2,5mm			15.250			non-engaging
	GH 1,5mm			16.050			engaging
	GH 2,5mm			16.250			engaging

Dentsply® Sirona®		Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	
		MIS® Seven Implants Internal Hex MI	GH 0,3mm				15.096*
GH 0,8mm	15.056*		15.095*				non-engaging
GH 1,5mm	15.256*		15.295*	15.296*			non-engaging
GH 3,0mm	15.356*		15.395*	15.096*			non-engaging
GH 0,3mm				16.096*			engaging
GH 0,8mm	16.056*		16.095*				engaging
	GH 1,5mm	16.256*	16.295*	16.296*			engaging
	GH 3,0mm	16.356*	16.395*	16.396*			engaging
	Hex ø 1,27mm	NP (yellow)	SP (purple)	WP (green)	Multi-unit 4.8		
MIS® C1 Implants Conic MC	GH 0,3mm				15.085*		non-engaging
	GH 1,0mm	15.081*	15.082*	15.083*			non-engaging
	GH 0,3mm				16.085*		engaging
	GH 1,0mm	16.081*	16.082*	16.083*			engaging
	Hex ø 1,27mm		NP (blue)	SP (purple)	Multi-unit 4.8		
MIS® V3 Implants Conic MV	GH 0,3mm				15.085*		non-engaging
	GH 1,0mm		15.084*	15.082*			non-engaging
	GH 0,3mm						non-engaging
	GH 1,0mm		16.084*	16.082*			engaging

Straumann® Group		Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
		Straumann® Soft Tissue Level & synOcta® Octagon	GH 0,3mm			
GH 0,4mm				15.046		non-engaging
GH 0,5mm	15.070					non-engaging
GH 0,3mm				16.045		engaging
GH 0,4mm				16.046		engaging
GH 0,5mm	16.070					engaging
	On Abutment		15.010	15.048		non-engaging
	On Abutment		16.010	16.048		engaging
	Torx®	NP (3.3)	RP (4.1)	NP/RP		
Straumann® Bone Level Conical BL	GH 1,2mm	15.043	15.044			non-engaging
	GH 2,0mm	15.043/2	15.044/2			non-engaging
	GH 3,0mm	15.343	15.344			non-engaging
	GH 1,2mm	16.043	16.044			engaging
	GH 2,0mm	16.043/2	16.044/2			engaging
	GH 3,0mm	16.343	16.344			engaging
	On SRA			15.098		non-engaging
	On SRA			16.098		engaging
	Torx®	RB (3.5)	WB (4.5)	RB/WB		
Straumann® BLX Conical BLX	GH 0,8mm		15.063			non-engaging
	GH 1,5mm	15.062/2	15.063/3	15.062/2		non-engaging
	GH 2,5mm	15.062/3		15.062/3		non-engaging
	GH 3,5mm	15.062/4		15.062/4		non-engaging
	GH 0,8mm		16.063			engaging
	GH 1,5mm	16.062/2	16.063/3	16.062/2		engaging
	GH 2,5mm	16.062/3		16.062/3		engaging
	GH 3,5mm	16.062/4		16.062/4		engaging
	NEO	GM				
Neodent Grand Morse® NEO GM	GH 0,8mm	15.035				non-engaging
	GH 1,5mm	15.235				non-engaging
	GH 2,5mm	15.335				non-engaging
	GH 0,8mm	16.035				engaging
	GH 1,5mm	16.235				engaging
GH 2,5mm	16.335				engaging	
	Hex ø 1,18mm					
Anthogyr Axiom® Implants Conic Axio	GH 2,9mm			15.077		non-engaging
	GH 2,9mm			16.677		engaging


Ti-Bases

Group	Implant	Hex ø	GH				Engaging
			3.3	3.8	4.3	5.0	
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm					
		GH 0,3mm		15.065	15.066	15.067	non-engaging
		GH 0,5mm	15.064				non-engaging
		GH 1,5mm		15.265	15.266		non-engaging
		GH 3,0mm		15.065	15.066		non-engaging
		GH 0,3mm		16.065	16.066	16.067	engaging
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)		
		GH 1,4mm	15.074*	15.075*	15.076*		non-engaging
		GH 3,0mm	15.074/3*	15.075/3*	15.076/3*		non-engaging
		GH 1,4mm	16.074*	16.075*	16.076*		engaging
		GH 3,0mm	16.074/3*	16.075/3*	16.076/3*		engaging
		GH 0,3mm					
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
		GH 0,3mm			15.088	15.089	non-engaging
		GH 0,8mm		15.087			non-engaging
		GH 0,3mm			16.088	16.089	engaging
BioHorizons® Implants External External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)			
	GH 0,3mm	15.051*	15.052*	15.053*		non-engaging	
Hiossen®/Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
		GH 1,0mm	15.072	15.073		non-engaging	
		GH 3,0mm	15.372	15.373		non-engaging	
		GH 1,0mm	16.072	16.073		engaging	
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm					
		GH 2,0mm			15.057	non-engaging	
		GH 3,0mm			15.357	non-engaging	
		GH 2,0mm			16.057	engaging	
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm					
					15.080*	non-engaging	
GlobalID®	GlobalID® Implants Conic Ink	Hex ø 1,22mm					
					16.086*	engaging	
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28mm					
		GH 0,9mm			15.090*	non-engaging	
		GH 1,5mm			15.190*	non-engaging	
		GH 3,0mm			15.390*	non-engaging	
		GH 0,9mm			16.090*	engaging	
		GH 1,5mm			16.190*	engaging	
GH 3,0mm			16.390*	engaging			

Group	Implant	Hex ø	GH				Engaging
			3.25/3.75	4.1	4.5	5.5	
Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27mm					
		GH 0,7mm			15.091*	15.094*	non-engaging
		GH 0,9mm			15.091*		non-engaging
		GH 1,1mm	15.091*				non-engaging
		GH 0,7mm			16.091*	16.094*	engaging
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22mm		NP	RP		
		GH 0,7mm					non-engaging
		GH 0,9mm	15.078				engaging
		GH 1,1mm	16.078				engaging
		GH 0,7mm					
BTI®	BTI® Implants Internal BT	Hex ø 1,22mm		4.1	MULTI IM®		
		GH 0,7mm		15.079*	15.065*		non-engaging
		GH 0,9mm		16.079*	16.665*		engaging
		GH 1,1mm					

AURUMBase®

Angulated screw channel solution. A 25° entry angle, a 360° rotation and a short shaft makes the AURUMBase® the perfect solution for those complicated cases. The Torx® screw combined with the Torx® ball screwdriver allows for maximum torque to be applied at a 25° angle.



Group	Implant	Hex ø	GH				Engaging
			3.0	NP (3.5)	RP (4.3)	WP (5.5)	
Nobel Active® & Nobel Replace®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip					
		GH 0,3mm			36.041	36.042	non-engaging
		GH 1,0mm			37.041	37.042	engaging
Nobel Replace Select™	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		GH 0,3mm	36.004	36.005	36.006		non-engaging
		GH 1,0mm	37.004	37.005	37.006		engaging
Nobel Biocare®	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
		GH 0,3mm	36.001	36.002			non-engaging
		GH 1,5mm	36.001/2	36.002/2			non-engaging
		GH 2,5mm	36.001/3	36.002/3			non-engaging
		GH 0,3mm	37.001	37.002			engaging
		GH 1,5mm	37.001/2	37.002/2			engaging
		GH 0,3mm	37.001/3	37.002/3			engaging
Multi-Unit®	Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			
		GH 0,3mm	36.007	36.054			non-engaging
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
		GH 0,3mm	36.017	36.018	36.019	non-engaging	
		GH 1,5mm	37.017	37.018	37.019	engaging	
		GH 0,3mm					
3i Osseotite®	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
		GH 0,3mm	36.011	36.012	36.013	non-engaging	
		GH 1,5mm	36.011/2	36.012/2		non-engaging	
		GH 2,5mm	36.011/3	36.012/3		non-engaging	
		GH 0,3mm	37.011	37.012	37.013	engaging	
		GH 1,5mm	37.011/2	37.012/2		engaging	
GH 0,3mm	37.011/3	37.012/3		engaging			

AURUMBase®


Zimmer® 3i® Biomet®	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
			36.014	36.015	36.016	non-engaging	
			37.014	37.015	37.016	engaging	
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)	
				36.059	36.060	36.061	non-engaging
		On Abutment	36.058				engaging
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)		
				36.024	36.025		non-engaging
				37.024	37.025		engaging
Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5		
		36.038	36.039	36.040		non-engaging	
			37.038	37.039	37.040	engaging	
MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8		
			36.017	36.018		non-engaging	
			37.017	37.018		engaging	
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)		
				36.046	36.047		non-engaging
			37.046	37.047		engaging	
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP		
				36.043	36.044		non-engaging
				37.043	37.044		engaging
	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)	RB/WB		
			GH 1,5mm	36.062	36.063	36.062	non-engaging
		GH 1,5mm	37.062	37.063	37.062	engaging	
Neodent Grand Morse® NEO GM	NEO	GM					
			GH 2,5mm	36.035/3			non-engaging
			GH 2,5mm	37.035/3		engaging	

Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0	
				36.065	36.066		non-engaging
				37.065	37.066		engaging
Henry Schein® Group	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)		
				36.074*	36.075*	36.076*	non-engaging
			37.074*	37.075*	37.076*		engaging
Henry Schein® Group	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
					36.087	36.088	36.89
				37.087	37.088	37.089	engaging
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
				36.072	36.073		non-engaging
			37.072	37.073			engaging
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm					
				36.057			non-engaging
				37.057			engaging
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm	GH 2,2mm				non-engaging
			GH 0,9mm	36.080*			non-engaging
			GH 2,2mm	36.180*			engaging
			GH 0,9mm	37.080*			engaging
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28mm					
				36.057			non-engaging
				37.057			engaging
Phibo®	Phibo® TSH® External Hex PH	GH 0,3mm	NP	SP	WP		
				36.030*	36.031*	36.032*	non-engaging
				36.030/2*	36.031/2*	36.032/2*	non-engaging
				36.030/3*	36.031/3*	36.032/3*	non-engaging
				37.030*	37.031*	37.032*	engaging
				37.030/2*	37.031/2*	37.032/2*	engaging
		37.030/3*	37.031/3*	37.032/3*	engaging		

ELLIPTIBase®

Narrow interdental spaces don't represent a challenge anymore. The ELLIPTIBase®, Has been specially developed for narrow implants (3.0-3.5). The elliptical shape achieves almost the same cementation area as a standard Ti-base with the possibility of achieving a more aesthetic result with the 20° entry angle.

- ZrN coating
- SelectGrip® enhanced surface
- Corner/flat versions




Brand	Implant System	Hex	Width				
			3.0	NP (3.5)	RP (4.3)	WP (5.5)	
Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip Corner	56.068				
		UniGrip Flat	57.068				
Zimmer® 3i Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm Corner	56.017	NP (3.5)	RP (4.5)	WP (5.7)	
		Hex ø 1,27mm Flat	57.017				
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm Corner	56.406	3.0	EV (3.6)	EV (4.2)	EV (4.8)
		Hex ø 1,27mm Flat	57.406				
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx® Corner	56.070	NNC (3.5)	RN (4.8)	WN (6.5)	
		Torx® Flat	57.070				
Straumann® Group	Straumann® Bone Level Conical BL	Torx® Corner	56.043	NP (3.3)	RP (4.1)	NP/RP	
		Torx® Flat	57.043				
Henry Schein® Group	BioHorizons® Implants Internal Hex BH	Hex ø 1,27mm Corner	56.099	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		Hex ø 1,27mm Flat	57.099				

C-Base®

The iconic design of the Sirona Ti-base with a twist. Engaging and non-engaging versions, different gingival heights, gold anodised and SelectGrip® surface treatment create the perfect upgrade.

-ExoCad and 3Shape libraries available for download at www.dessdental.com



Brand	Implant System	Hex	Width				Engaging	
			3.0	NP (3.5)	RP (4.3)	WP (5.5)		
Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip GH 0,8mm		58.041	58.042		non-engaging	
		UniGrip GH 3,0mm		58.041/3	58.042/3		non-engaging	
		UniGrip GH 0,8mm		59.041	59.042		engaging	
		UniGrip GH 3,0mm		59.041/3	59.042/3		engaging	
	Nobel Replace Select™ Trilobe	UniGrip GH 0,4mm	NP (3.5)	58.004		58.006	58.097	6.0 non-engaging
		UniGrip GH 0,5mm			58.005			non-engaging
		UniGrip GH 1,4mm		58.004/2	58.005/2	58.006/2		non-engaging
		UniGrip GH 3,0mm		58.004/3	58.005/3	58.006/3		non-engaging
		UniGrip GH 0,4mm		59.004		59.006	59.097	engaging
		UniGrip GH 0,5mm			59.005			engaging
		UniGrip GH 1,4mm		59.004/2	59.005/2	59.006/2		engaging
		UniGrip GH 3,0mm		59.004/3	59.005/3	59.006/3		engaging
Nobel Brånemark® External Hex Universal	UniGrip GH 1,5mm	NP (3.5)	58.001/2*	58.002/2*			non-engaging	
	UniGrip GH 3,0mm		58.001/3*	58.002/3*			non-engaging	
	UniGrip GH 0,6mm		59.001*	59.002*			engaging	
	UniGrip GH 1,5mm		59.001/2*	59.002/2*			engaging	
	UniGrip GH 3,0mm		59.001/3*	59.002/3*			engaging	
	Hex ø 1,27mm GH 0,5mm	NP (3.5)	58.017	58.018	58.019		non-engaging	
Zimmer® 3i Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm GH 1,5mm		58.017/2	58.018/2	58.019/2		non-engaging
		Hex ø 1,27mm GH 3,0mm		58.017/3	58.018/3	58.019/3		non-engaging
		Hex ø 1,27mm GH 0,5mm		59.017	59.018	59.019		engaging
		Hex ø 1,27mm GH 1,5mm		59.017/2	59.018/2	59.019/2		engaging
		Hex ø 1,27mm GH 3,0mm		59.017/3	59.018/3	59.019/3		engaging
		Hex ø 1,27mm GH 0,5mm		59.017/2	59.018/2	59.019/2		engaging
	3i Osseotite® External Hex USA	Hex ø 1,27mm GH 1,5mm	NP (3.4)	58.011/2	58.012/2	58.013/2		non-engaging
		Hex ø 1,27mm GH 3,0mm		58.011/3	58.012/3	58.013/3		non-engaging
		Hex ø 1,27mm GH 0,6mm		59.011	59.012	59.013		engaging
		Hex ø 1,27mm GH 1,5mm		59.011/2	59.012/2	59.013/2		engaging
		Hex ø 1,27mm GH 3,0mm		59.011/3	59.012/3	59.013/3		engaging
		Hex ø 1,27mm GH 0,4mm		58.014	58.015			non-engaging
3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm GH 0,5mm				58.016		non-engaging	
	Hex ø 1,20mm GH 1,5mm		58.014/2	58.015/2	58.016/2		non-engaging	
	Hex ø 1,20mm GH 3,0mm		58.014/3	58.015/3	58.016/3		non-engaging	
	Hex ø 1,20mm GH 0,4mm		59.014	59.015			engaging	
	Hex ø 1,20mm GH 0,5mm				59.016		engaging	
	Hex ø 1,20mm GH 1,5mm		59.014/2	59.015/2	59.016/2		engaging	
	Hex ø 1,20mm GH 3,0mm		59.014/3	59.015/3	59.016/3		engaging	
	Hex ø 1,20mm GH 0,4mm						engaging	

Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)		
		GH 0,8mm	58.059	58.060	58.061	non-engaging	
		GH 3,0mm	58.059/3	58.060/3	58.061/3	non-engaging	
		GH 0,8mm	59.059	59.060	59.061	engaging	
		GH 3,0mm	59.059/3	59.060/3	59.061/3	engaging	
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)		
		GH 0,8mm		58.024		non-engaging	
		GH 0,9mm			58.025	non-engaging	
		GH 3,0mm		58.024/3	58.025/3	non-engaging	
		GH 0,8mm		59.024		engaging	
		GH 0,9mm			59.025	engaging	
		GH 3,0mm		59.024/3	59.025/3	engaging	
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5	
		GH 0,3mm	58.038	58.039		non-engaging	
		GH 0,6mm			58.040	58.071	non-engaging
		GH 1,5mm	58.038/2	58.039/2	58.040/2	58.071/2	non-engaging
GH 3,0mm		58.038/3	58.039/3	58.040/3	58.071/3	non-engaging	
GH 0,3mm		59.038	59.039			engaging	
GH 0,6mm				59.040	59.071	engaging	
GH 1,5mm		59.038/2	59.039/2	59.040/2	59.071/2	engaging	
Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm						
	GH 1,4mm			58.050	non-engaging		
	GH 2,5mm			58.050/2	non-engaging		
	GH 1,4mm			59.050	engaging		
	GH 2,5mm			59.050/2	engaging		
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)		
		GH 1,0mm		58.009		non-engaging	
		GH 1,2mm			58.045	non-engaging	
		GH 1,0mm		59.009		engaging	
		GH 1,2mm			59.045	engaging	
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP		
		GH 0,7mm	58.043*	58.044		non-engaging	
		GH 3,0mm	58.043/3	58.044/3		non-engaging	
GH 0,7mm		59.043	59.044		engaging		
	GH 3,0mm	59.043/3	59.044/3		engaging		
Henry Schein® Group	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)		
		GH 1,0mm	58.074*	58.075*	58.076*	non-engaging	
		GH 2,1mm	58.074/2*	58.075/2*	58.076/2*	non-engaging	
		GH 1,0mm	59.074*	59.075*	59.076*	engaging	
		GH 2,1mm	59.074/2*	59.075/2*	59.076/2*	engaging	
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
		GH 0,5mm		58.087	58.088	58.089	non-engaging
		GH 3,0mm		58.087/3	58.088/3	58.089/3	non-engaging
GH 0,5mm			59.087	59.088	59.089	engaging	
	GH 3,0mm		59.087/3	59.088/3	59.089/3	engaging	
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD			
		GH 1,0mm	58.072	58.073	non-engaging		
		GH 3,0mm	58.072/3*	58.073/3*	non-engaging		
		GH 1,0mm	59.072	59.073	engaging		
		GH 3,0mm	59.072/3	59.073/3	engaging		

Abutments

Available in 0°, 15° and 25° with SelectGrip® enhanced surface treatment, the DESS abutments offer an excellent alternative.

- Includes final screw



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		0°		13.041	13.042		
		15°		22.041*	22.042*		
			25°			23.042*	
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		0°	13.004	13.005	13.006		
		15°	22.004*	22.005*	22.006*		
			25°		23.005*		
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
0°		13.001	13.002	13.003			
15°		22.001*	22.002*	22.003*			
		25°		23.002*			
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
		0°	13.017	13.018	13.019		
		15°	22.017	22.018	22.019		
			25°		23.017	23.018	
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
		0°	13.011	13.012	13.013		
15°		22.011*	22.012*	22.013*			
		25°		23.012*			
3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)			
	0°	13.014	13.015	13.016			
	15°	22.014*	22.015*	22.016*			
		25°		23.015*			
Dentsply® Sirona®	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,22mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)		
		0°		13.024	13.025		
		15°		22.024*	22.025*		
			25°		23.024*		
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5	
		0°	13.038	13.039	13.040		
		15°	22.038*	22.039*	22.040*		
				25°		23.039*	
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm					
		GH 0,7mm	0°		13.047		
		GH 2,5mm	0°		13.250		
		GH 1,5mm	15°		22.048*		
GH 2,5mm		15°		22.250*			
GH 1,5mm		25°		23.048*			
	GH 2,5mm	25°		23.250*			
MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8		
	0°		13.095*	13.096*			
	15°		22.095*	22.096*			
			25°		23.095*	23.096*	

Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)		
		0°		13.009			*inc.captive screw
		15°		22.009*			*inc.captive screw
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP		
		0°	13.043	13.044			
		15°	22.043	22.044			
Henry Schein® Group	BioHorizons® Implants External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)		
		0°	13.020*	13.021*	13.022*		
		15°	22.020*	22.021*			

Screws

Our screws have the same geometry as the ones of the implant manufacturer. They have the same footprint for the screwdriver tip, seat angle, diameter, length, metric, and more, making them identical to the one supplied by the implant manufacturer.

- Manufactured in titanium Grade 5 ELI
- Available in Standard, DLC*, TIN
- Single unit or 10 packs





Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		Std. Single	19.068	19.027	19.030	19.030	
		Std. 10 pack		19.027-P10	19.030-P10	19.030-P10	
		DLC	19.268*	19.227	19.230	19.230	
		TIN Single	19.668*	19.627*	19.630*	19.630*	
		Torx®	19.468	19.441	19.442		AURUMBase® & ELLIPTIBase®
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		Std. Single	19.004	19.005	19.005	19.005	
		10 pack	19.004-P10	19.005-P10	19.005-P10	19.005-P10	
		DLC	19.204	19.205	19.205	19.205	
		TIN Single	19.604	19.605	19.605	19.605	
		Torx®	19.404	19.405	19.405		AURUMBase® & ELLIPTIBase®
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
		Std. Single	19.027	19.028	19.029		
		10 pack	19.027-P10	19.028-P10	19.029-P10		
		DLC	19.227	19.228	19.229		
		TIN Single	19.627*	19.628*	19.629*		
		Torx®	19.401	19.402			AURUMBase® & ELLIPTIBase®
	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			
		Std. Single	19.006	19.031			
Std. 10 pack		19.006-P10	19.031-P10				
DLC		19.206	19.231				
TIN Single		19.606*	19.631*				
Torx®		19.306	19.331*			AURUMBase® & ELLIPTIBase®	
Hex 1,20mm		19.055					
Hex 1,22mm		19.065					
Hex 1,27mm	19.110						

Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)	
		Std. Single	19.042	19.042	19.042	
		10 pack	19.042-P10	19.042-P10	19.042-P10	
		TIN Single	19.642*	19.642*	19.642*	
		Torx®	19.417	19.417	19.417	AURUMBase® & ELLIPTIBase®
		3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)
	Std. Single	19.034	19.034	19.034		
	10 pack	19.034-P10	19.034-P10	19.034-P10		
	TIN Single	19.634*	19.634*	19.634*		
	Torx®	19.402	19.402	19.402	AURUMBase® & ELLIPTIBase®	
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)	
		Std. Single	19.011	19.011	19.011	
10 pack		19.011-P10	19.011-P10	19.011-P10		
TIN Single		19.611*	19.611*	19.611*		
Torx®		19.438	19.438	19.438	AURUMBase® & ELLIPTIBase®	

ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (3.6)	EV (4.2)	EV (4.8)	EV (5.4)	
	Std. Single	19.071	19.059	19.060	19.061	19.061	
	TIN Single		19.659*	19.660*	19.661*	19.661*	
	Torx®	19.468	19.441	19.460	19.461		AURUMBase® & ELLIPTIBase®
	On Uniabutment®						
	Std. Single						19.058
	TIN Single						19.658*
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)		
		Std. Single	19.015	19.016	19.017		
		Std. 10 pack		19.016-P10	19.017-P10		
		TIN Single	19.615*	19.616*	19.617*		
		Torx®		19.441	19.425		AURUMBase® & ELLIPTIBase®
On Uniabutment®							
Std. Single						19.018*	
Std. 10 pack						19.018-P10	
TIN Single						19.618*	
Dentsply Sirona® Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5		
	Std. Single	19.024	19.024	19.024	19.024		
	Std. 10 pack	19.024-P10	19.024-P10	19.024-P10	19.024-P10		
	TIN Single	19.624*	19.624*	19.624*	19.624*		
	Torx®	19.438	19.438	19.438		AURUMBase® & ELLIPTIBase®	
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm					
Std. Single 1,5mm			19.048				
Std. 10 pack 1,5mm			19.048-P10				
TIN Single 1,5mm			19.648*				
Std. Single 2,5mm			19.050				
C-Base Std Screw			19.702				
MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8		
	Std. Single	19.056*	19.081*	19.081*	19.110		
	TIN Single	19.656*	19.681*	19.681*	19.085*		
	Torx®		19.417	19.417	19.685*		
	MIS® C1 Implants Conic MC	Hex ø 1,27mm	NP (yellow)	SP (purple)	WP (green)		
		Std. Single	19.056*	19.081*	19.081*		
TIN Single		19.656*	19.681*	19.681*			
MIS® V3 Implants Conic MV	Hex ø 1,27mm	NP (blue)	SP (purple)				
	Std. Single	19.056	19.081				
	TIN Single	19.656	19.681				

Screws

Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)	
		Std. Single	19.070	19.026	19.026	
		10 pack		19.026-P10	19.026-P10	
		TIN Single	19.670*	19.626*	19.626*	
		Torx®	19.441	19.446	19.446	
		On SynOcta®				
		Single		19.009	19.009	
10 Pack		19.009-P10	19.009-P10			
			19.609*	19.609*		
Straumann® Group	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP	
		Std. Single	19.032	19.033	19.037	
		10 pack	19.032-P10	19.033-P10	19.037-P10	
		TIN Single	19.632*	19.633*	19.637*	
		Torx®	19.443	19.443		
		On SRA®				
		Std. Single			19.098*	
	10 pack			19.098-P10*		
	TIN Single			19.698*		
	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)	RB/WB	
		Std. Single	19.054	19.054	19.054	
		10 pack	19.700*	19.700*	19.700*	
		Torx®	19.441	19.441	19.441	
	Neodent Grand Morse® NEO GM	NEO	GM			
		Std. Single	19.051			
	Anthogyr Axiom® Implants Conic Axio	Hex ø 1,18mm				
		Std. Single			19.077*	
TIN Single				19.677*		
On Abutment						
	Single			19.062*		
	TIN Single			19.662*		
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0
		Std. Single	19.044	19.044	19.044	19.045
		TIN Single	19.644*	19.644*	19.644*	19.645*
		Torx®		19.444	19.444	
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)	
		Std. Single	19.074	19.074	19.076	
		TIN Single	19.674	19.674	19.676	
	BioHorizons® Implants Internal Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		Std. Single	19.043	19.043	19.043	19.043
		TIN Single	19.643*	19.643*	19.643*	19.643*
	BioHorizons® Implants External External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)	
		Std. Single	19.013*	19.014*	19.014*	
Std. 10 pack			19.014-P10*	19.014-P10*		
TIN Single		19.613*	19.614*	19.614*		

Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD		
		Std. Single	19.072	19.073		
		DLC	19.272	19.273		
		TIN	19.672	19.673		
		Torx®	19.444	19.461		
AURUMBase® & ELLIPTIBase®						
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm				
		Std. Single	19.057			
		TIN	19.657*			
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm				
		Std. Single	19.080*			
		TIN Single	19.680*			
		Torx® 0,2mm	19.444			
		Torx® 2,2mm	19.443			
AURUMBase® AURUMBase®						
GlobalD®	GlobalD® Implants Conic Ink	Hex ø 1,22mm				
		Std. Single	19.086*			
		TIN Single	19.686*			
		On Abutment				
		Std. Single	19.035*			
		TIN Single	19.635*			
		Single	19.036*			
TIN Single	19.636*					
MUA  MUA 						
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28mm				
		Std. Single	19.090*			
		TIN Single	19.690*			
		Torx®	19.405			
		On Abutment				
	Single	19.066*				
	TIN Single	19.666*				
AURUMBase®						
Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27mm	3.25/3.75	4.1	4.5	5.5
		Std. Single	19.091*	19.091*	19.091*	19.091*
		TIN Single	19.691*	19.691*	19.691*	19.691*
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22mm				
		Std. Single	19.078			
		TIN Single	19.678			
		On Abutment		NP	RP	
		Std. Single		19.063	19.064	
TIN Single		19.663	19.664			
BTI®	BTI® TSH® External Hex PH	Hex ø 1,22mm	4.1	MULTI IM®		
		Std. Single	19.079*	19.065*		
		TIN Single	19.679*	19.665*		
PHIBO® Implants Internal BT	Torx	NP	SP	WP		
	Torx®	19.403*	19.402	19.402		

Pre-Milled Blanks

DESS pre-milled blanks are available in titanium Grade 5 ELI (ø 10mm - ø 14mm) and Micro Melt (CrCo ø 14mm). With 28 different compatibilities, it is the most comprehensive range in the market.

- Total millable length of 20 mm
- Gingival transition only 0,3 mm to maximize the design freedom
- Includes final screw



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		ø 10mm Ti	61.068	61.041	61.042	61.069	
		ø 14mm Ti		62.041	62.042	62.069	
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		ø 10mm Ti	61.004	61.005	61.006	61.097	
		ø 14mm Ti		62.005	62.006	62.097	
	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
		ø 10mm Ti	61.001	61.002	61.003		
		ø 14mm CoCr	71.001	71.002	71.003		
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
		ø 10mm Ti	61.017	61.018	61.019		
		ø 14mm Ti	62.017	62.018	62.019		
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
		ø 10mm Ti	61.011	61.012	61.013		
		ø 14mm CoCr	71.011	71.012	71.013		
	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
		ø 10mm Ti	61.014	61.015	61.016		
		ø 14mm Ti		62.015	62.016		
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.06)	EV (3.6)	EV (4.2)	EV (4.8)	EV (5.4)
		ø 10mm Ti	61.406	61.059	61.060	61.060	61.407
		ø 14mm Ti		62.059	62.060	62.060	62.407
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)		
		ø 10mm Ti	61.023	61.024	61.025		
		ø 14mm Ti		62.024	62.035		
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5	
		ø 10mm Ti	61.038	61.039	61.040	61.071	
		ø 14mm Ti		62.039	62.040	62.071	
Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm						
	ø 10mm Ti	61.038					
	ø 14mm CoCr	71.038					

Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8
		ø 10mm Ti	61.056*	61.095*	61.096*	
		ø 14mm Ti				
	MIS® C1 Implants Conic MC	Hex ø 1,27mm	NP (yellow)	SP (purple)	WP (green)	
		ø 10mm Ti	61.081*	61.082	61.083*	
MIS® V3 Implants Conic MV	Hex ø 1,27mm	NP (blue)	SP (purple)			
	ø 10mm Ti	61.084*	61.082			

Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)
		ø 10mm Ti	61.070	61.009	61.045
		ø 14mm Ti		62.009	62.045
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP
		ø 10mm Ti	61.043	61.044	
		ø 14mm Ti	62.043	62.044	
	Straumann® BLX Conical BLX	Torx®	RB (3.5)	WB (4.5)	RB/WB
		ø 10mm Ti	61.062	61.062	61.062
		ø 14mm Ti		62.063	
	Neodent Grand Morse® NEO GM	NEO	GM		
		ø 10mm Ti	61.035		
		ø 14mm Ti	62.035		
Anthogyr Axiom® Implants Conic Axio	Hex ø 1,18mm				
	ø 10mm Ti	61.077*			

Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0
		ø 10mm Ti	61.064	61.065	61.066	61.067
		ø 14mm Ti		62.065	62.066	62.067
	Conelog® Implants Internal Con	Hex ø 1,27mm	NP (3.3)	RP (3.8)	WP (5.0)	
		ø 10mm Ti	61.074*	61.075*	61.076*	
		ø 14mm Ti	62.074*	62.075*	62.076*	
	BioHorizons® Implants Internal Hex BH	Hex ø 1,27mm	3.0	NP (3.5)	RP (4.5)	WP (5.7)
		ø 10mm Ti	61.099	61.087	61.088	61.089
		ø 14mm Ti		62.087	62.088	62.089
	BioHorizons® Implants External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)	
		ø 10mm Ti	61.020*	61.021*	61.022*	
		ø 14mm Ti	62.020*	62.021*	62.022*	
		ø 14mm CoCr	71.020*	71.021*	71.022*	

Pre-Milled Blanks

Hiossen/ Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm		MINI	STD
		ø 10mm	Ti	61.072	61.073
		ø 14mm	Ti	62.072	62.073
		ø 14mm	CoCr	71.072	71.073
Megagen®	Megagen® Anyridge Implants Conic Anyr	Hex ø 1,20mm			
		ø 10mm	Ti	61.057	
		ø 14mm	CoCr	71.057	
Medentis®	Medentis ICX® Implants Conic IC	Hex ø 1,4mm			
		ø 10mm	Ti	61.080*	
		ø 14mm	Ti	62.080*	
		ø 14mm	CoCr	71.080*	
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28			
		ø 10mm	Ti	61.080*	
		ø 14mm	Ti	62.080*	
		ø 14mm	CoCr	71.080*	
Bego®	Bego Semados® Implants Internal Sem	Hex ø 1,27		3.25/3.75	4.1
		ø 10mm	Ti	61.091	61.092
		ø 14mm	Ti		61.093
		ø 14mm	CoCr		5.5
					61.094
					62.094
					71.094
Biotech®	Biotech Kontakt® Implants Conic Bio	Hex ø 1,22			
		ø 10mm	Ti	61.078*	
		ø 14mm	Ti		
		ø 14mm	CoCr	71.078*	

Blocks

Blocks	Blocks					
		Ti		63.020	□15mm	L 20mm
		Ti		63.040	□15mm	L 40mm
		CoCr		73.020	□15mm	L 20mm
		CoCr		73.040	□15mm	L 40mm

CoCr bases



Manufactured in micromelt CoCr alloy with a sanitary POM sleeve, the DESS CoCr bases are perfect for casting techniques. With a higher melting point (1369°C), they will maintain precision of shape in the fusion with other metals.

- Can be used in bone-level restorations or from abutment level (Multiunit)
- Engaging and non-engaging versions available
- Final screw included/met

Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
						30.041	30.042
				31.041	31.042		engaging
Nobel Biocare®	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
			30.004	30.005	30.006		non-engaging
			31.004	31.005	31.006		engaging
Nobel Biocare®	Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)		
			30.001	30.002	30.003		non-engaging
			31.001	31.002	31.003		engaging
				30.102*			non-engaging
Nobel Biocare®	Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)			
			30.007	30.054			non-engaging
			Hex 1,20mm	30.007H120			non-engaging
			Hex 1,22mm	30.007H122			non-engaging
		Hex 1,27mm	30.007H127			non-engaging	
		Torx® Ball	30.007TB			non-engaging	
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
			30.017	30.018	30.019		non-engaging
			31.017	31.018	31.019		engaging
Zimmer® 3i® Biomet®	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
			30.011	30.012	30.013		non-engaging
			31.011	31.012	31.013		engaging
Zimmer® 3i® Biomet®	3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)		
			30.014	30.015	30.016		non-engaging
			31.014	31.015	31.016		engaging

CoCr bases

Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm		EV (3.6)	EV (4.2)	EV (4.8)			
				30.059	30.060	30.061	non-engaging		
		On Uniabutment®	30.058						
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)				
				30.024	30.025		non-engaging		
				31.024	31.025		engaging		
		20°		30.036				On Uniabutment®	
		45°		30.037				On Uniabutment®	
	Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5			
				30.038	30.039	30.040		non-engaging	
			31.038	31.039	31.040		engaging		
Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm								
	GH 1,5mm				30.050		non-engaging		
	GH 2,5mm				30.250		non-engaging		
	GH 1,5mm				31.050		engaging		
	GH 2,5mm				31.250		engaging		
MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8				
			30.056*	30.095*	30.096*		non-engaging		
			31.056*	31.095*	31.096*		engaging		
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)				
				30.009	30.045				
				31.009	31.045				
				30.109*	30.145*		non-engaging	For welding techniques	
				31.109*	31.145*		non-engaging	For welding techniques	
			On Abutment	30.010	30.048		engaging		
		On Abutment	31.010	31.048		engaging			
			30.110*			non-engaging	For welding techniques		
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP				
				30.043	30.044		non-engaging		
			31.043	31.044		engaging			
	On Abutment SRA			30.098*		non-engaging			
	On Abutment SRA			31.098*		engaging			
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0			
				30.065	30.066		non-engaging		
				31.065	31.066		engaging		
Hiossen® Osstem®	Osstem® TS Implants Conic Oss	Hex ø 1,22mm	MINI	STD					
				30.072	30.073		non-engaging		
				31.072	31.073		engaging		
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28mm							
					30.090*		non-engaging		
					31.090*		engaging		

Castables

Unlike other products on the market, DESS castables are machined, which produces a more accurate connection replica, avoiding imperfections in the casting.

- Engaging and non-engaging
- Single and 10 unit packs
- Single unit includes final screw



Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	UniGrip	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
		Single		28.041	28.042		non-engaging
		10 pack		11.041-P10	11.042-P10		
		Single		29.041	29.042		engaging
		10 pack		12.041-P10	12.042-P10		
	Nobel Replace Select™ Trilobe	UniGrip	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
		Single		28.004	28.005	28.006	non-engaging
		10 pack		11.004-P10	11.005-P10		
		Single		29.004	29.005	29.006	engaging
		10 pack		12.004-P10	12.005-P10		
Nobel Brånemark® External Hex Universal	UniGrip	NP (3.5)	RP (4.1)	WP (5.0)			
	Single		28.001	28.002	28.003	non-engaging	
	10 pack			11.002-P10			
	Single		29.001	29.002	29.003	engaging	
	10 pack			12.002-P10			
Multi-Unit® Multiunit Connection	UniGrip	RP (4.8)	WP (6.5)				
	Single		28.007	28.054		non-engaging	
	10 pack		11.007-P10	11.054-P10			
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	Hex ø 1,27mm	NP (3.5)	RP (4.5)	WP (5.7)		
		Single		28.017	28.018	28.019	non-engaging
		10 pack		11.017-P10	11.018-P10		
		Single		29.017	29.018	29.019	engaging
		10 pack		12.017-P10	12.018-P10		
	3i Osseotite® External Hex USA	Hex ø 1,27mm	NP (3.4)	RP (4.1)	WP (5.0)		
		Single		28.011	28.012	28.013	non-engaging
		10 pack			11.012-P10		
		Single		29.011	29.012	29.013	engaging
		10 pack			12.012-P10		
3i Certain® Implants Internal Hex "Click"	Hex ø 1,20mm	NP (3.4)	RP (4.1)	WP (5.0)			
	Single		28.014	28.015	28.016	non-engaging	
	10 pack			11.015-P10			
	Single		29.014	29.015	29.016	engaging	
	10 pack			12.015-P10			

Castables

Brand	Product	Hex ø	Width				Engaging	Angle
			NP	RP	WP	Multi-unit		
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	Hex ø 1,27mm	EV (3.6)	EV (4.2)	EV (4.8)			
		Single	28.059	28.060	28.061		non-engaging	
		Single	28.058	On Uniabutment®				
	Astra Tech Osseospeed™ Implants Internal Hex Conic	Hex ø 1,27mm	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)			
		Single		28.024	28.025		non-engaging	
		10 pack		11.024-P10	11.025-P10			
		Single		29.024	29.025		engaging	
		10 pack		12.024-P10	12.025-P10			
		Single	28.036				non-engaging	20°
		10 pack	11.036-P10	On Uniabutment®			non-engaging	20°
		Single	28.037				non-engaging	45°
		10 pack	11.037-P10				non-engaging	45°
		Dentsply Friadent® Xive® Implants Internal Hex FD	Hex ø 1,22mm	NP (3.4)	RP (3.8)	WP (4.5)	5.5	
	Single		28.038	28.039	28.040		non-engaging	
	10 pack			11.039-P10				
	Single		29.038	29.039	29.040		engaging	
	Dentsply Ankylos® C/X Implants Internal ANK	Hex ø 1,00mm						
		Single		28.050				GH 1,5mm
MIS® Seven Implants Internal Hex MI	Hex ø 1,27mm	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8			
	Single		28.056	28.095	28.096		non-engaging	
	10 pack			11.095-P10	11.096-P10			
	Single		29.056	29.095	29.096			
	10 pack			12.095-P10	12.096-P10			
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	Torx®	NNC (3.5)	RN (4.8)	WN (6.5)			
		Single		28.009	28.045		non-engaging	
		10 pack		11.009-P10				
		Single		29.009	29.045		engaging	
		10 pack		12.009-P10				
		Single		28.010	28.046		non-engaging	
	Straumann® Bone Level Conical BL	Torx®	NP (3.3)	RP (4.1)	NP/RP			
		Single		28.043	28.044		non-engaging	
		10 pack		11.043-P10	11.044-P10			
		Single		29.043	29.044		engaging	
		10 pack		12.043-P10	12.044-P10			

Brand	Product	Hex ø	Width				Engaging
			NP	RP	WP	Multi-unit	
Henry Schein® Group	Camlog® Implants Internal Cam	Hex ø 1,27mm	3.3	3.8	4.3	5.0	
		Single		28.065	28.066		non-engaging
	BioHorizons® Implants External Hex BH	Hex ø 1,27mm	NP (3.7)	RP (4.2)	WP (5.2)		
		Single	28.020	28.021	28.022		non-engaging
	10 pack			11.021-P10			
		Single	29.020	29.021	29.022		engaging
	10 pack			12.021-P10			
Dentium®	Superline™ & Implantium® Conic Dent	Hex ø 1,28					
		Single			28.090		non-engaging
		Single			29.090		engaging

Ti-Base Castables



Manufactured in sanitary POM, the Ti-base sleeves are perfect for casting processes.
 - Available in packs of 5
 - Engaging and non-engaging versions

Brand	Product	Width	Width				Engaging
			NP	RP	WP	Multi-unit	
Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	3.0	NP (3.5)	RP (4.3)	WP (5.5)		
		5 pack	32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
			33.001-P5	33.001-P5	33.001-P5	33.001-P5	engaging
	Nobel Replace Select™ Trilobe	NP (3.5)	RP (4.3)	WP (5.0)	6.0		
		5 pack	32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
			33.001-P5	33.001-P5	33.001-P5	33.001-P5	engaging
	Nobel Brånemark® External Hex Universal	NP (3.5)	RP (4.1)	WP (5.0)			
		5 pack	32.001-P5	32.001-P5	32.003-P5		non-engaging
			33.001-P5	33.001-P5	33.003-P5		engaging
	Multi-Unit® Multiunit Connection	RP (4.8)	WP (6.5)				
5 pack		32.007-P5	32.054-P5			non-engaging	
		32.007/8-P5				engaging	

Ti-Base Castables

Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	5 pack	NP (3.5)	RP (4.5)	WP (5.7)		
			32.001-P5	32.001-P5	32.001-P5	non-engaging	
				33.001-P5	33.001-P5	33.001-P5	engaging
	3i Osseotite® External Hex USA	5 pack	NP (3.4)	RP (4.1)	WP (5.0)		
			32.001-P5	32.001-P5	32.001-P5	non-engaging	
				33.001-P5	33.001-P5	33.001-P5	engaging
3i Certain® Implants Internal Hex "Click"	5 pack	NP (3.4)	RP (4.1)	WP (5.0)			
		32.001-P5	32.001-P5	32.001-P5	non-engaging		
			33.001-P5	33.001-P5	33.001-P5	engaging	
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	5 pack	EV (3.6)	EV (4.2)	EV (4.8)		
			32.001-P5	32.001-P5	32.001-P5	non-engaging	
			33.001-P5	33.001-P5	33.001-P5	engaging	
	Astra Tech Osseospeed™ Implants Internal Hex Conic	5 pack	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)		
			32.001-P5	32.001-P5	32.001-P5	non-engaging	
			33.001-P5	33.001-P5	33.001-P5	engaging	
	Dentsply Friadent® Xive® Implants Internal Hex FD	5 pack	NP (3.4)	RP (3.8)	WP (4.5)	5.5	
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
				33.001-P5	33.001-P5	33.001-P5	engaging
	Dentsply Ankylos® C/X Implants Internal ANK	5 pack				32.001-P5	non-engaging
						33.001-P5	33.001-P5

Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	5 pack	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	
			32.001-P5	32.001-P5	32.001-P5	32.007-P5	non-engaging
				33.001-P5	33.001-P5	33.001-P5	engaging
	MIS® C1 Implants Conic MC	5 pack	NP (yellow)	SP (purple)	WP (green)	Multi-unit 4.8	
			32.001-P5	32.001-P5	32.001-P5	32.007-P5	non-engaging
				33.001-P5	33.001-P5	33.001-P5	engaging
MIS® V3 Implants Conic MV	5 pack		NP (blue)	SP (purple)	Multi-unit 4.8		
			32.001-P5	32.001-P5	32.007-P5	non-engaging	
			33.001-P5	33.001-P5	33.001-P5	engaging	
Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	5 pack	NNC (3.5)	RN (4.8)	WN (6.5)		
			32.001-P5	32.046-P5	32.045-P5	non-engaging	
			33.001-P5	33.009-P5	33.045-P5	non-engaging	
	Straumann® Bone Level Conical BL	5 pack	NP (3.3)	RP (4.1)	NP/RP		
			32.001-P5	32.001-P5		non-engaging	
			33.001-P5	33.001-P5		engaging	
	Straumann® BLX Conical BLX	5 pack			32.098-P5	non-engaging	
					33.098-P5	engaging	
	Straumann® BLX Conical BLX	5 pack	RB (3.5)	WB (4.5)	RB/WB		
			32.001-P5	32.001-P5	32.001-P5	non-engaging	
			33.001-P5	33.001-P5	33.001-P5	engaging	
	Neodent Grand Morse® NEO GM	5 pack			GM		
				32.001-P5	non-engaging		
				33.001-P5	engaging		
Anthogyr Axiom® Implants Conic Axio	5 pack			32.001-P5	non-engaging		
				33.001-P5	engaging		
Henry Schein® Group	Camlog® Implants Internal Cam	5 pack	3.3	3.8	4.3	5.0	
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
				33.001-P5	33.001-P5	33.001-P5	engaging
	Conelog® Implants Internal Con	5 pack	NP (3.3)	RP (3.8)	WP (5.0)		
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	non-engaging
				33.001-P5	33.001-P5	33.001-P5	engaging
	BioHorizons® Implants Internal Internal Hex BH	5 pack	3.0	NP (3.5)	RP (4.5)	WP (5.7)	
				32.001-P5	32.001-P5	32.001-P5	non-engaging
				33.001-P5	33.001-P5	33.001-P5	engaging
	BioHorizons® Implants External External Hex BH	5 pack	NP (3.7)	RP (4.2)	WP (5.2)		
			32.001-P5	32.001-P5	32.001-P5	non-engaging	
				33.001-P5	33.001-P5	33.001-P5	engaging


Ti-Base Castables

			MINI	STD			
Hiossen®/Osstem®	Osstem® TS Implants Conic Oss	5 pack	32.001-P5	32.001-P5	non-engaging		
			33.001-P5	33.001-P5	engaging		
Megagen®	Megagen® Anyridge Implants Conic Anyr	5 pack		32.001-P5	non-engaging		
				33.001-P5	engaging		
Medentis®	Medentis ICX® Implants Conic IC	5 pack		32.001-P5	non-engaging		
				33.001-P5	engaging		
GlobalD®	GlobalD® Implants Conic Ink	5 pack		32.001-P5	non-engaging		
				33.001-P5	engaging		
Dentium®	Superline™ & Implantium® Conic Dent	5 pack		32.001-P5	non-engaging		
				33.001-P5	engaging		
Bego®	Bego Semados® Implants Internal Sem	5 pack	3.25/3.75	4.1	4.5	5.5	non-engaging
			32.001-P5	32.001-P5	32.001-P5	32.001-P5	33.001-P5
Biotech®	Biotech Kontakt® Implants Conic Bio	5 pack		NP	RP		non-engaging
			32.001-P5				
BTI®	BTI® Implants Internal BT	5 pack	4.1				non-engaging
			32.001-P5				
			33.001-P5				

AURUMBase® Castables

Machine manufactured in sanitary POM and available in 0°, 10°, 20°, the AURUMBase® castables, will facilitate the creation of a precise angulated channel when used on traditional workflows.

- Available in packs of 5 units











Nobel Biocare®	Nobel Active® & Nobel Replace® CC Active Hex	5 pack	3.0	NP (3.5)	RP (4.3)	WP (5.5)	
					33.100-P5	33.100-P5	0°
					33.101-P5	33.101-P5	10°
			33.102-P5	33.102-P5	20°		
Nobel Biocare®	Nobel Replace Select™ Trilobe	5 pack	NP (3.5)	RP (4.3)	WP (5.0)	6.0	
					33.100-P5	33.100-P5	0°
					33.101-P5	33.101-P5	10°
					33.102-P5	33.102-P5	20°
Nobel Biocare®	Nobel Brånemark® External Hex Universal	5 pack	NP (3.5)	RP (4.1)	WP (5.0)		
					33.100-P5	33.100-P5	0°
					33.101-P5	33.101-P5	10°
					33.102-P5	33.102-P5	20°
Nobel Biocare®	Multi-Unit® Multiunit Connection	5 pack	RP (4.8)	WP (6.5)			
					33.700-P5	33.700-P5	0°
					33.701-P5	33.701-P5	10°
					33.702-P5	33.702-P5	20°
Zimmer® 3i® Biomet®	Zimmer Screw-Vent® Internal Hex USA	5 pack	NP (3.5)	RP (4.5)	WP (5.7)		
					33.700-P5	33.700-P5	0°
					33.701-P5	33.701-P5	10°
			33.702-P5	33.702-P5	20°		
Zimmer® 3i® Biomet®	3i Osseotite® External Hex USA	5 pack	NP (3.4)	RP (4.1)	WP (5.0)		
					33.100-P5	33.100-P5	0°
					33.101-P5	33.101-P5	10°
			33.102-P5	33.102-P5	20°		
Zimmer® 3i® Biomet®	3i Certain® Implants Internal Hex "Click"	5 pack	NP (3.4)	RP (4.1)	WP (5.0)		
					33.100-P5	33.100-P5	0°
					33.101-P5	33.101-P5	10°
			33.102-P5	33.102-P5	20°		
Dentsply® Sirona®	ASTRATECH Implants System™ EV Conic EVO	5 pack	EV (3.6)	EV (4.2)	EV (4.8)		
					33.100-P5	33.100-P5	0°
					33.101-P5	33.101-P5	10°
			33.102-P5	33.102-P5	20°		
Dentsply® Sirona®	Astra Tech Osseospeed™ Implants Internal Hex Conic	5 pack	NP (3.0)	RP (3.5/4.0)	WP (4.5/5.0)		
					33.100-P5	33.100-P5	0°
					33.101-P5	33.101-P5	10°
			33.102-P5	33.102-P5	20°		
Dentsply® Sirona®	Dentsply Friadent® Xive® Implants Internal Hex FD	5 pack	NP (3.4)	RP (3.8)	WP (4.5)	5.5	
					33.100-P5	33.100-P5	0°
					33.101-P5	33.101-P5	10°
			33.102-P5	33.102-P5	20°		
Dentsply® Sirona®	MIS® Seven Implants Internal Hex MI	5 pack	NP (narrow)	SP (Std)	WP (wide)	Multi-unit 4.8	
					33.100-P5	33.100-P5	0°
					33.101-P5	33.101-P5	10°
			33.102-P5	33.102-P5	20°		

AURUMBase® Castables

Straumann® Group	Straumann® Soft Tissue Level & synOcta® Octagon	5 pack	NNC (3.5)	RN (4.8)	WN (6.5)	
				33.200-P5	33.300-P5	0°
				33.201-P5	33.301-P5	10°
	Straumann® Bone Level Conical BL	5 pack		33.202-P5	33.302-P5	20°
			NP (3.3)	RP (4.1)	NP/RP	
			33.100-P5	33.100-P5		0°
	Straumann® BLX Conical BLX	5 pack		33.101-P5	33.101-P5	10°
				33.102-P5	33.102-P5	20°
			RB (3.5)	WB (4.5)	RB/WB	
	Neodent Grand Morse® NEO GM	5 pack		33.100-P5	33.100-P5	0°
				33.101-P5	33.101-P5	10°
				33.102-P5	33.102-P5	20°
Henry Schein® Group	Camlog® Implants Internal Cam	5 pack	3.3	3.8	4.3	5.0
				33.100-P5	33.100-P5	0°
				33.101-P5	33.101-P5	10°
	Conelog® Implants Internal Con	5 pack		33.102-P5	33.102-P5	20°
			NP (3.3)	RP (3.8)	WP (5.0)	
			33.100-P5	33.100-P5	33.100-P5	0°
	BioHorizons® Implants Internal Hex BH	5 pack		33.101-P5	33.101-P5	10°
				33.102-P5	33.102-P5	20°
			3.0	NP (3.5)	RP (4.5)	WP (5.7)
				33.100-P5	33.100-P5	0°
				33.101-P5	33.101-P5	10°
				33.102-P5	33.102-P5	20°
Hiossen®/ Osstem®	Osstem® TS Implants Conic Oss	5 pack	MINI	STD		
				33.100-P5	33.100-P5	0°
				33.101-P5	33.101-P5	10°
			33.102-P5	33.102-P5	20°	
Megagen®	Megagen® Anyridge Implants Conic Anyr	5 pack				0°
						10°
						20°
Medentis®	Medentis ICX® Implants Conic IC	5 pack				0°
						10°
						20°
Dentium®	Superline™ & Implantium® Conic Dent	5 pack				0°
						10°
						20°
Phibo®	Phibo® TSH® External Hex PH	5 pack	NP	SP	WP	
				33.100-P5		0°
				33.101-P5		10°
			33.102-P5		20°	








HOLDERS

Holders		For standard 98mm disc machines	Any milling machine with a 98 mm disc holder
	HLDSTD		
		For Imes-Icore® machine	Coritec 245i touch / Coritec 250i touch Coritec 350i touch / Coritec 350i touch Loader Coritec 450i / Coritec 550i Coritec 650i / Coritec 650i Loader
	HLDIMES		
		For VHF®	VHF® S1 VHF® S2 VHF® J5
	HLDVHF		
		For VHF®	VHF® R5
	HLDVHFR5		
		For YENADENT®	D6 D15 D43N
	HLDYENA		
	For VHF®	VHF® N4	
HLDVHFN4			
	For VHF®	VHF® Z4	
HLDVHFZ4			
	For machines with support / changer 3R	D14P DC40 D43N	
HLD3R			

Tools

Our screws have the same geometry as the implant manufacturer's ones. Having the same footprint for the screwdriver tip, seat angle, diameter, length, metric, and more makes them identical to the ones supplied by the implant manufacturer.











- Manufactured in titanium Grade 5 ELI
- Available in Standard, DLC*, TIN
- Single unit or packs of 10








Screwdriver	Tip	18mm	20mm	24mm	30mm	35mm	Image
		HEXAGON					
	0,90 mm		DTC20H090	DTC24H090			
	1,00 mm		DTC20H100	DTC24H100		DTC35H100	
	1,20 mm	DTC18H120	DTC20H120	DTC24H120	DTC30H120	DTC35H120	
	1,22 mm		DTC20H122	DTC24H122		DTC35H122	
	1,25 mm		DTC20H125	DTC24H125			
	1,27 mm		DTC20H127	DTC24H127	DTC30H127	DTC35H127	
	1,30 mm		DTC20H130	DTC24H130			
	1,50 mm		DTC20H150	DTC24H150			
	1,70 mm		DTC20H170	DTC24H170			
	1,75 mm		DTC20H175	DTC24H175			
SQ							
			DT20S130	DT24S130			
Torx® 06							
			DT20TC06	DT24TC06	DT30TC06	DT35TC06	
Torx® ball							
			DT20TB06	DT24TB06	DT30TB06	DT35TB06	
UG							
			DT20U	DT24U	DT30U6	DT35U	
NEO							
				DT24NG	DT30NG	DT35NG	
Screwdriver handles	Premium laboratory handle		PHLABGOLD				
			PHLABBLACK				
	Dentist handle		7 mm			MDTMAN700	DT35T06
		9 mm			MDTMAN900	DT35TC06	



DESSLoc® Removal & Insertion tool	DESSLoc® removal & insertion tool	PHLOC	
Ti-Base extracting tool	For Conelog® NP/RP	DCCON	
	For Conelog® WP	DCCON50	
	For Zimmer screw-vent®	DCZIM	
	For Medentis ICX®	DCICX	
Torque wrench tool	Torque wrench tool	Type Straumann 7mm DTSTTW	
		Type Nobel 8mm DTNBTW	
Try Multi-use	Try Multi-use	BOX1	

Tools

Special tools	Multi-unit® RP	DTMURP	
	Multi-unit® WP	DTMUWP	
	UniAbutment®	DTUNI	
	Octagon	DTSYN	
	Torque wrench adaptor	TWADAP	
	DESSLoc®	DTLOC	
	Adaptor from 6 mm to straumann® torque wrench to square connection	DTSTSQ	
	Screwdriver tips adaptor for straumann® torque wrench	DTSTWH	
		DTECTLOC	
		DTSTRLOC	
	Starumann® torque wrench	DTSTCA	

Special tools	Nobel Biocare™ adaptor from 8mm to square connection	DTNBWH	
	Nobel Biocare™ adapter from 8mm to square connection	DTNBSQ	
	Screwdriver tips adaptor from Nobel Biocare™ torque wrench	DTNBWH	
	Nobel Biocare™ contra-angles adapter	DTNBCA	
	Handpiece adaptor for square connection	DTANGSQ	
	Manual driver for Intraoral Scan abutment	DTPEI	
	Manual tool for intraoral Scan abutment	DTPEIPEEK-P2	

DESS HQ



6.200m² Barberà del Vallès, Barcelona Spain



All DESS products are manufactured according to ISO 9001 and ISO 13485 with CE marking for all products
Class IIb

All DESS® AURUMBase® screws and screwdrivers are based on the standard ISO 10664
3i® and 3i Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations,
USA

ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XiVE®, FRIADENT®, ASTRATECH implant system™ EV and
UniAbutment®

Cerec® are registered trademarks of DENTSPLY SIRONA Implants

BioHorizons® is a registered trademark of BioHorizons, INC.

Brånemark System®, NobelReplace® and Replace Select™, NobelActive™, NobelReplace® CC, NobelSpeedy™ and
Multi-unit® are registered trademarks of Nobel Biocare AB,

Neodent®, Grand Morse®, Straumann® and synOcta® are registered trademarks of Straumann Holding AG,
Switzerland

MIS® is a registered trademark of MIS Implant Technologies Ltd.

Zimmer Screw-Vent® and SwissPlus® are registered trademarks of Zimmer Dental INC, USA

Dyna® is a registered trademark of Dyna Dental Engineering BV

Locator® is a registered trademark of ZEST IP Holdings

Camlog® and Conelog® are registered trademarks of CAMLOG Biotechnologies AG

Megagen® is a registered trademark of MEGAGEN IMPLANT CO. LTD

TORX® is a trademark of ACUMENT INTELLECTUAL PROPERTIES, LLC

Bego Semados® is a registered trademark of BEGO GmbH & Co. KG

Anthogyr® is a registered trademark of Anthogyr SAS

Medentis ICX® is a registered trademark of medentis medical GmbH

Osstem®, Hiossen® is a registered trademark of OSSTEM IMPLANT CO. LTD

Dentium® is a trademark of Dentium Co. Ltd.

GOBALD® is a trademark of Global D SASU.

Biotech Kontakt® is a trademark Biotech Dental SAS.

BTI® is a trademark BTI Biotechnology Institute.

PHIBO® is a trademark PHIBO Dental Solutions S.L.

DESS products can be used only by licensed dentists, physicians and labs.

Please refer ever to the instructions of use available on the website

www.dessdental.com

or supplied with the product.

Some products may not be available in your country. Some of the products are cleared for US and Canada.

Please contact your local distributor for more details.



DOC-EN-01-FFR

