



# XTRICATE LYTE

TyTek Medical's XTricate Lyte is a high performance, low cost disposable emergency stretcher. Featuring an independent laboratory tested load capacity of 500 lbs, it has four large handles on each side and two large handles at both the foot and head, it is ready to serve and improve your ability to move patients.

The XTricate Lyte is constructed with heavy-duty polypropylene straps, reinforced polyester thread stitching, and backed by a lightweight non-woven polypropylene fabric. Given XTricate Lyte's durability and value, it should be available in every first responder vehicle, event space, school, sport's field, or any place people gather.

XTricate Lyte is designed for patient carry and is not a drag device.

**Dimensions:**

75" x 30" (Open)

12" x 8" x 1.5" (Packaged)

**Weight:** 12oz

**Capacity:** 500lbs (226kg)



Contact us today to find out more:  
1-513-247-2002  
sales@tytekmedical.com  
www.tytekmedical.com

TyTek Medical Inc. 8904 Beckett Rd, West Chester, Ohio, 45069 USA  
www.tytekmedical.com @BleedingControl  
CAGE: 6TL45 DUNS:078589992 Tax ID: 27-5015102

Rev. 11042024



Scan to visit our website