

**Twenty20
Healthcare
Limited**

**From trauma
wounds to
minor
bleeds,
CollaClot™
starts to
work in
milliseconds
to clot the
blood, stop
bleeding and
improve
healing**











**CollaClot™
is for animal
use only**



“Don't plug the bleed, clot the blood”

COLLACLOT™ Haemostatic Collagen Sponge

CollaClot

-  Formed with a natural collagen derivative
-  For use on major trauma wounds to minor bleeds
-  Versatile for multiple wounds and everyday use
-  Can be removed without tissue damage and minimises re-bleeding
-  Works instantly
-  Quick and simple to apply to any bleeding site
-  Conformable to multiple trauma sites
-  Significantly improves healing
-  Available in two sizes
-  Saves time and money

PLEASE CHECK OUR WEBSITE FOR FURTHER INFORMATION



adapting to change, alternative thinking
for medical equipment and supplies

 info@twenty20hc.com

 www.twenty20hc.com

 023 8098 2020

Rownhams House, Betteridge Drive,
Rownhams, Southampton SO16 8LS

CollaClot™ - Haemostatic Collagen Sponge

The latest innovation in haemostatic agents, which enhances the efficiency of procedures or emergency situations that require a fast and effective control of blood flow.

CollaClot™ is a haemostatic sponge, versatile enough for a variety of animal wounds and appropriate for major and minor bleeds. This product is made from 100% Avian Collagen (which is one of the closest matches to human collagen) and not only stops the bleed but acts as a catalyst to speed up the healing process of a wound by activating the natural blood-clotting cascade.

Collagen is the most prevalent protein within the body. Its fibre-like structure is used to make connective tissue and is a major component of bone, skin, muscle, tendons and cartilage. It plays a critical role in the regulation phases of wound healing and contributes to the mechanical strength and elasticity of tissues.

Unlike other products, CollaClot™ works with the animal's blood to help activate the clotting process, rather than absorbing the blood. Some competitor products absorb the blood and form a plug, which doesn't conform well to the shape of the injury, while others have no clot attached and the blood is only absorbed into the product, forming a light, sticky gel. Due to the fact that this product has to be removed, it counteracts the natural processes of clot formation. CollaClot™, on the other hand, creates a resilient, flexible block that seals blood vessels and resists further bleeding on movements, proving its superiority for general first aid, trauma wounds and emergency bleeding.



**10cm x 7.5cm size
- for use on larger
animals including
equine.**

**Using rough side
face-down to the
wound increases
blood absorption**



**3cm x 3cm size - for smaller
animals (including
domestic pets)**

Product Description	Product Code
CollaClot - 3cm x 3cm	CCV002
CollaClot - 10cm x 7.5cm	CCV003