

Blood Pressure Monitor  
**UA-767F**



Blood Pressure Monitor  
**UA-767S**



# Blood Pressure Monitor UA-767F / UA-767S

## UA-767F

For Family



## UA-767S

Simple



The UA-767F is the 3rd generation of the UA-767 series with new design and innovative features focused on user-friendliness for family use. It has independent memories for 4 users (60 x 4) and newly added %IHB function, which will help users to check the frequency of Irregular Heartbeat.

The UA-767S is also the successor of the UA-767 series for single user focused on "Simple operation". Ease of use as it is, new feature %IHB is installed, also the number of memory has been upgraded to 60.

## About %IHB



### < Features >

The frequency of IHB (Irregular Heart Beat) detection is divided into different grades and shown by icons.

0-24% : No Indication  
25-49% : Grade 1  
50-74% : Grade 2  
75-100% : Grade 3

### < Benefits >

You can see how often IHB has been detected (indicating the level of risk). If the grade is high, you should consult a doctor immediately. This %IHB function is not to be used for diagnosis, but for monitoring purposes only.

## Specifications

	UA-767F	UA-767S
Measurement method	Oscillometric measurement	
Display type	Digital LCD, Pressure/Pulse displayed simultaneously	
Measurement range	Pressure: 0 - 299 mmHg Systolic pressure: 60 - 279 mmHg Diastolic pressure: 40 - 200 mmHg Pulse: 40 - 180 beats / minute	
Measurement accuracy	Pressure: ±3 mmHg Pulse: ±5%	
Power supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233) (Not included)	
Number of measurements	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries) With pressure value 180 mmHg, room temperature 23°C.	
Memory storage	60 x 4 users	60
Memory data average	Maximum 60 data average of systolic, diastolic, and pulse values	
Cuff	SlimFit cuff with plastic D-ring / Arm circumference: 22-32cm	
Operating conditions	+10 to +40 °C / 15 to 85 %RH / 800 to 1060hPa	
Transport/Storage conditions	-20°C to +60°C / 10%RH to 95%RH	
Dimensions	Approx. 140 [W] x 60 [H] x 105 [D] mm	
Weight	Approx. 255g, excluding the batteries	Approx. 245g, excluding the batteries

**UA-767F**

%IHB Indicator to show frequency of Irregular Heartbeat

Cuff Fit Error Indicator

Movement Error Indicator

WHO Classification Indicator

**UA-767S**

**4xM60** 4 user memory (60 x 4)

**31** Clock & Calendar

**AVG.** Average reading for each 4 users (Last 60 measurements)

**M60** 60 memories

**AVG.** Average reading of last 60 measurements

### Standard accessories

- Blood pressure monitor
- SlimFit cuff (22-32cm)
- Carrying case
- Batteries for temporary use x 4
- Instruction manual

### Options

- SlimFit large cuff (31-45cm): CUFF-F-LA
- SlimFit wide range cuff (22-42cm) : CUFF-I
- AC adapter TB-233C: INPUT 100-240V 50/60Hz C plug
- AC adapter TB-233BF: INPUT 240V 50/60Hz BF plug



**A&D Company, Limited**  
3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013 Japan  
Telephone: [81](3)5391-6132 Fax: [81](3)5391-6148  
1-243, Asahi, Kitamoto-shi, Saitama, 364-8585 Japan  
Telephone: [81](0485)93-1111 Fax: [81](0485)93-1119  
<http://www.aandd.jp>

**A&D INSTRUMENTS LTD.**  
Unit 24/26 Blacklands Way Abingdon Business Park,  
Abingdon, Oxon OX14 1DY UNITED KINGDOM  
Telephone: [44](1235) 550420 Fax: [44](1235) 550485  
(Authorized Representative Established in the European Community)