

## 3-193 Whitehouse Gravity Axis Marker For Pre-Op Reference Marking

Designed in cooperation with Dr Geoff Whitehouse, MBBS,  
FRANZCO, Forster Eye Surgery, Australia



### LRI, TORIC

#### Two-Steps Toric Marking

Marker blades dent in the corneal epithelium so the mark placement is visible even if the ink is washed off.



Gravity weight system ensures excellent visualization and balance



- Accurately mark the reference axis to avoid the effect of cyclotorsion
- Marks the limbus at 3, 6 and 9 o'clock - the 0°, 90°, and 180° on the eye while the patient is sitting upright, looking straight ahead