

## **DISPOSABLE DIAMOND DUSTED ILM ELEVATOR\***

A challenging maneuver associated with the removal of the internal limiting membrane (ILM) is to pinch and grab it prior to peeling.

RUMEX introduced a **NEW DISPOSABLE DIAMOND DUSTED ILM ELEVATOR** in 23 and 25 Ga, which is specifically designed to consistently create a precise edge to facilitate the ILM removal with forceps.

12-7523, 23 Ga 12-7525, 25 Ga









**RETRACTABLE VERSION** helps to adjust the length of the tip and provides an easy insertion through the trocar cannula





SOFT SILICONE TIP safe for retinal surface DIAMOND DUSTED tip finish provides an extreme grip

\*Not for sale in Europe

