


IOL Implantation with RUMEX*



STEP 1. CHOOSE THE IOL AND DELIVERY SYSTEM

PRELOADED IOL + DELIVERY SYSTEM



AquaFree Yellow Aspheric Preloaded

FEATURES:

- Fast injection
- Single-use option

OR


NON-PRELOADED IOLS



AquaFree Yellow Aspheric



Hydro-Sense Yellow



Hydro-Sense



Hydro-4
suitable for MICS

DELIVERY SYSTEMS




Single-use cartridge CAT-22 (2.2 mm) or CAT-24 (2.4 mm) + Reusable IOL injector 16-2853

FEATURES:

- Money-wise solution
- Flexibility in choosing the IOL

OR



VE2200 Disposable IOL injector

FEATURES:

- Single-use option
- Flexibility in choosing the IOL

STEP 2. ADD THE OVD

DISPERSIVE

SUPREME
HPMC 2%

BEST FOR:

- Coating the surgical instruments
- Endothelium protection

AND/OR

COHESIVE

SMARTVISC
Sodium Hyaluronate 1,6%
SMARTVISC PLUS
Sodium Hyaluronate 3%

BEST FOR:

- Keeping the empty capsular bag open for IOL insertion
- Flattening the anterior capsule to facilitate capsulorhexis creation

