



# Cleaner Compatibility

There are symbols placed after the values listed under the PRODUCT category below to denote effectivity claims made by governing bodies or manufacturers for COVID-19.

\* United States Environmental Protection Agency (EPA) website (<https://www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2>)

\*\* CDC guideline on disinfection and sterilization in healthcare facilities (<https://www.cdc.gov/infectioncontrol/pdf/guidelines/disinfection-guidelines-H.pdf>)

\*\*\* Cleaning agent manufacturer claims effectiveness

<b>CLINICAL AREA</b>	:	Diagnostic Cardiology
<b>DEVICE</b>	:	MAC 7
<b>PRODUCT</b>	:	Alcohol (Isopropyl) 70%** , Amuchina Disinfectant Spray, CAVIWIPES*** , Cleanisept Wipes*** , Clinell Sporidical Wipes, Clinell Universal Sanitizing Wipes*** , Clorox Healthcare™ Hydrogen Peroxide Cleaner Disinfectant Wipes* , DAX Surface Disinfection, Ethanol (ethyl alcohol) 96% (v/v)** , Hydrogen Peroxide 1-5%, Hydrogen Peroxide 20% (v/v), Incidin Oxy Wipe S*** , Mikrobac® Tissues*** , Mikrozyd® AF Wipes, Mikrozyd® Sensitive Wipes*** , OXIVIR TB WIPES* , PDI Easy Screen Cleaning® , PDI Sani Cloth Chlor, Phenol [max 2% by volume], Sani-Cloth® Active*** , Soap Water, Sodium Hypochlorite (NaOCl) 5% solution), Sodium Hypochlorite 5.25 % (Bleach) Diluted 10:1, Trionic D Cleaning and Disinfection Wipes, Wet Wipe, Wetwipe Chlorine Disinfection Wipes*** , Wetwipe Ethanol, Wetwipe Universal

Clinical Area	Device	Product	ALL(28)	COMPATIBLE(28)	CAUTION(0)	NON-COMPATIBLE(0)
Diagnostic Cardiology	MAC 7	<ul style="list-style-type: none"> <li>✓ Alcohol (Isopropyl) 70%**</li> <li>✓ Amuchina Disinfectant Spray</li> <li>✓ CAVIWIPES***</li> <li>✓ Cleanisept Wipes***</li> <li>✓ Clinell Sporidical Wipes</li> <li>✓ Clinell Universal Sanitizing Wipes***</li> <li>✓ Clorox Healthcare™ Hydrogen Peroxide Cleaner Disinfectant Wipes*</li> <li>✓ DAX Surface Disinfection</li> <li>✓ Ethanol (ethyl alcohol) 96% (v/v)**</li> <li>✓ Hydrogen Peroxide 1-5%</li> <li>✓ Hydrogen Peroxide 20% (v/v)</li> <li>✓ Incidin Oxy Wipe S***</li> <li>✓ Mikrobac® Tissues***</li> <li>✓ Mikrozyd® AF Wipes</li> <li>✓ Mikrozyd® Sensitive Wipes***</li> <li>✓ OXIVIR TB WIPES*</li> <li>✓ PDI Easy Screen Cleaning®</li> <li>✓ PDI Sani Cloth Chlor</li> <li>✓ Phenol [max 2% by volume]</li> <li>✓ Sani-Cloth® Active***</li> <li>✓ Soap Water</li> <li>✓ Sodium Hypochlorite (NaOCl) 5% solution)</li> <li>✓ Sodium Hypochlorite 5.25 % (Bleach) Diluted 10:1</li> <li>✓ Trionic D Cleaning and Disinfection Wipes</li> <li>✓ Wet Wipe</li> </ul>				