



### CS-4 TABLET FRIABILITY TESTER

CS series of instruments are used for detecting friability/ abrasion of tablet.

#### Application Standards

National standard( Chinese codex ) /2020

Corporate standard(Tablet Friability tester )Q/12XQ0189-2010

#### Features

Single line and double cylinder,synchroized operation,auto at timing point.

The cylinder adopts high quality hyaline organic glass.

Rotation speed and number of turns can be setted at will.

Operation is simple,the control precision is very high.

Automation: auto-test, auto-diagnose, auto-alarm.

#### Specifications

Cylinder number: 1

Cylinder radius:  $\Phi 286\text{mm}$

Cylinder depth: 39mm

Tablet falling height: 156mm

Rotation speed: 25rpm

Rotation speed accuracy:  $\pm 1\text{rpm}$

Rotation count: 100r

Rotation count accuracy:  $\pm 1\text{r}$

Power: 220V/ 50Hz/ 60W or 110V/ 60Hz/ 50W

Dimension: 300×200×330 mm