


CJY-300E TABLET FRIABILITY TESTER



FEATURES

- Fully in line with the "Chinese Pharmacopoeia 2015 Edition" standard.
- Higher precision integrated circuit design to ensure experimental accuracy.
- High-definition touch screen with simple interface and minimal operation mode.
- Chinese and English dual language operation interface, can meet the needs of more users.
- The speed and direction of rotation can be adjusted and the number of turns can be set.
- Smooth operation, low energy consumption and high reliability.
- Weak current working voltage is safer.
- Has a fault prompt function.

PARAMETERS

Product model	CJY-300E
Drum dimensions	Internal diameter about 286mm, depth: 39mm
Default rotating speed	25 RPM
Rotating speed setting range	20~50RPM
Number of laps	100 laps
Lap setting range	1 – 1,000 laps
Rotating direction	Clockwise/Anti-clockwise
Number of drums	2 no
Power supply	DC12V , 3A