

Klean



Nitrile Gloves

POWDER FREE | KLH-AC001

100 GLOVES
by weight

COLORS

PINK BLUE VIOLET WHITE GREEN COBALT BLACK PURPLE ORANGE




EN 455



Hand size	XS	S	M	L	XL
Min. length	220	230	240	250	260
Width	76	84	94	105	113

Sizes in mm ± 10



Caution: Components used in making gloves may cause allergic reactions in some users. If allergic symptoms occur, discontinue use of this product and immediately consult a physician.

Storage: Avoid storing at temperatures above 35°C (95°F). Shield gloves from direct sunlight, fluorescent light, X-rays, moisture and ozone.

Usage: These gloves are only to be used for the recommended purposes as covered by the CE ISO EN markings. Please note that gloves are only tested under laboratory conditions and relating only to the tested specimens taken from the palm only and or the chemical(s) tested. Actual conditions in the workplace might differ from the test environment. Users are advised to check that the gloves are suitable for the intended use and operate with caution.

Warning: Before usage, inspect the gloves for any defects or imperfections. Usage and conditions where gloves are subject to abrasion, wear and tear, excessive temperatures and/or chemicals for which the gloves have not been certified, may also cause defects. Users should refrain from using defect gloves and dispose of them immediately.

Disposal: These gloves are designed for single use only. Do not reuse. Dispose of used gloves properly in a protective container.