

THE BOXER ORTHOSIS

EZYWRAP designed the Boxer Orthosis to be a superior alternative to traditional casting for patients with boxer fractures. Designed with the patient in mind, The Boxer Orthosis can be applied single-handedly, with stay relief to accommodate the ulnar styloid. The Boxer Orthosis includes a customizable stay for adjustment by an individual with expertise for tailored care.

Our Boxer Orthosis eliminates ulnar styloid pressure, supports the wrist, forearm, and fingers at the preferred angle for treating fractures of the 4th & 5th metacarpals, and can be modified to address the 2nd & 3rd metacarpal.

Approved **PDAC** L3807 & 3809



INDICATIONS FOR USE:

- Fractures of 4th and 5th metacarpals




P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA   

FEATURES & BENEFITS



1. Ability for single-hand application of a comfortable alternative to casting or splinting.
2. Bendable stays for wrist extension and finger flexion allowing for customization.
3. Adjustable compression to accommodate swelling and skin care during recovery.
4. Relief for Ulnar Styloid with custom stays designed for support and comfort.

PRODUCT ORDERING

Product Name	Ref No.
The Boxer Orthosis	950-XXY

PRODUCT SIZING

Measure circumference of palmar crease

Size	XX	Min		Max	
PED	05	5"	13cm	6"	16cm
XS	10	6"	16cm	7"	18cm
REG	35	7"	18cm	10"	25cm
XL	50	10"	25cm	>10"	>25cm

Note: For Y specify L for left and R for right

R&ID# *282529*
REV A 6/23

FOR THE USE OF:
CODE:

MADE IN THE USA   

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

