

TALON WRIST LACING ORTHOSIS

EZYWRAP® designed the Talon Wrist Lacing Orthosis to be soft, breathable, and ready directly off-the-shelf. Our Wrist Orthosis is designed to contour to the hand using three removable, malleable stays with easy application and adjustable lacing. The unique radial closure design allows the orthosis edges to overlap for a semi-universal fit. Available in both Short and Long for additional comfort.

Approved **PDAC** L3908



INDICATIONS FOR USE:

- Carpal Tunnel Syndrome
- Wrist Sprains and Strains
- Wrist Immobilization
- Control of Post-Op Effusion





P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Soft and breathable materials for compression and comfort
2. Overlapping radial closure edges with lacing for semi-universal sizing and ulnar control
3. Three malleable and removable stays for a supportive custom fit
4. Non-neoprene hypoallergenic material and soft edges for user comfort

PRODUCT ORDERING

Product Name	Model No.
Talon Wrist Lacing Orthosis - Short	1112-XXY
Talon Wrist Lacing Orthosis - Long	1113-XXY

PRODUCT SIZING

Measure the circumference of wrist at narrowest

Size	XX	Min	Max	Min	Max
PED	05	-	4.0"	-	10cm
XS	10	4.0"	6.0"	10cm	15cm
REG	35	5.5"	9.0"	14cm	23cm
XL	50	8.0"	11.0"	20cm	28cm

Sizing Note: For Y specify L for left or R for right

R&ID# *285457*
REV A 2/22

FOR THE USE OF:
CODE:

MADE IN THE USA



FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

