

ULTIMATE THUMB LACING ORTHOSIS

EZYWRAP® designed the Ultimate Thumb Lacing Orthosis to be soft, breathable and ready directly off-the-shelf. Our Wrist Orthosis is designed to contour to the hand using three removable, malleable stays with easy application and adjustable lacing. The unique design allows the orthosis edges to overlap for a semi-universal fit. Available in both Long and Short for additional control.

Approved **PDAC** L3807 & L3809



INDICATIONS FOR USE:

- Thumb Injuries
- CMC Injuries
- Metacarpal Injuries




P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Soft and breathable materials and binding for compression and comfort
2. Anatomically designed for thumb positioning in abduction and opposition
3. Lacing pull tabs with overlapping edges for semi-universal sizing
4. Malleable and removable 6 stay wrist-thumb system for a supportive custom fit
5. Non-neoprene hypoallergenic material

PRODUCT ORDERING

Product Name	Ref No.
Ultimate Thumb Lacing Orthosis - Short	1127-XXY
Ultimate Thumb Lacing Orthosis - Long	1128-XXY

PRODUCT SIZING

Measure from center wrist at narrowest

Size	XX	Min	Max	Min	Max
XS	15	4.0"	6.0"	10cm	15cm
REG	35	5.5"	9.0"	14cm	23cm
XL	50	8.0"	11.0"	20cm	28cm

Sizing Note: For Y specify L for left or R for right

R&ID# *284084*
REV A 6/21

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

