

RAPID WRIST ORTHOSIS

EZYWRAP designed the Rapid Wrist Orthosis to be soft, breathable, and ready directly off-the-shelf. Our Wrist Orthosis is designed to contour to the hand using three removable, malleable stays with easy application and hook and loop closures for rapid application. The unique design allows the orthosis edges to overlap for a semi-universal fit. Available in both Long and Short for additional comfort.

Approved **PDAC** L3908



INDICATIONS FOR USE:

- Carpal Tunnel Syndrome
- Wrist Sprains and Strains
- Wrist immobilization
- Control of Post-op Effusion





P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Soft and breathable materials for compression and comfort
2. Hook and Loop closure for rapid application
3. Overlapping edges for semi-universal sizing
4. Three malleable and removable stays for a supportive custom fit
5. Non-neoprene hypoallergenic material

PRODUCT ORDERING

Product Name	Model No.
Rapid Wrist Orthosis - Short	1132-XXY
Rapid Wrist Orthosis - Long	1133-XXY

PRODUCT SIZING

Measure the circumference of wrist at narrowest

Size	XX	Min	Max	Min	Max
PED	05	-	4.0"	-	10cm
XS	10	4.0"	6.0"	10cm	15cm
REG	35	5.5"	9.0"	14cm	23cm
XL	50	8.0"	11.0"	20cm	28cm

Sizing Note: For Y specify L for left or R for right

R&ID# *285180*
REV A 8/22

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

