

RAPID THUMB ORTHOSIS

EZYWRAP® designed the Rapid Thumb Orthosis to be soft, breathable, and ready directly off-the-shelf. Our Orthosis is designed to contour to the hand and thumb using six removable, malleable stays with easy application and quick release hook and loop closures. The unique design allows the orthosis edges to overlap for a semi-universal fit. Available in both Long and Short for additional comfort.

Approved **PDAC** L3807, L3809



INDICATIONS FOR USE:

- Thumb Injuries
- CMC Injuries
- Metacarpal Injuries




P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Breathable soft materials and binding for compression and comfort
2. Anatomically designed positioning in abduction and opposition
3. Hook and Loop closure for rapid application
4. Overlapping edges for semi-universal sizing
5. Malleable and removable 6 stay system for a custom fit
6. Non-neoprene hypoallergenic material

PRODUCT ORDERING

Product Name	Model No.
Rapid Thumb Orthosis - Short	1137-XXY
Rapid Thumb Orthosis - Long	1138-XXY

PRODUCT SIZING

Measure the circumference of the wrist at narrowest

Size	XX	Min	Max	Min	Max
XS	10	4.0"	6.0"	10cm	15cm
REG	35	5.5"	9.0"	14cm	23cm
XL	50	8.0"	11.0"	20cm	28cm

Sizing Note: For Y specify L for left or R for right

R&ID# *285181*
REV A 8/22

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

