

UNIVERSAL WRIST ORTHOSIS

EZYWRAP designed the Universal Wrist Orthosis to be soft, breathable, and ready directly off-the-shelf. It is made with soft and breathable materials to encourage patient compliance. Our Universal Wrist uses a unique strapping system allowing the user to customize the orthosis for any anatomy. This brace has soft edges and a comfortable tongue to promote patient compliance and long-term wearability. Available in both Short and Long for additional control.

Approved **PDAC** L3908



INDICATIONS FOR USE:

- Carpal Tunnel Syndrome
- Wrist Sprains and Strains
- Wrist Immobilization
- Control of Post-Op Effusion





© 2022 PROFESSIONAL PRODUCTS INC.

P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Universal sizing for stocking optimization
2. Soft and breathable materials for comfort and compliance
3. Three malleable and removable stays for a supportive custom fit
4. Non-neoprene hypoallergenic material and soft edges for user comfort

PRODUCT ORDERING

Product Name	Model No.
Universal Wrist Orthosis - Short	1142-XXY
Universal Wrist Orthosis - Long	1143-XXY

PRODUCT SIZING

Measure the circumference of wrist at narrowest

Size	XX	Min	Max	Min	Max
UNI	00	6.0"	11.0"	15cm	28cm

Sizing Note: For Y specify L for left or R for right

R&ID# *285460*
REV A 2/22

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

