

UNIVERSAL THUMB ORTHOSIS

EZYWRAP® designed the Universal Thumb Orthosis to be soft, breathable and ready directly off-the-shelf or with expertise fitting. Our Orthosis is designed to contour to the hand and thumb using seven removable, malleable stays with easy application and adjustable straps. Single-piece construction gives the brace a refined low profile.

This brace has soft edges and a comfortable stockinette gusset to promote patient compliance and long-term wearability. Available in both short and long for additional control.

Approved **PDAC** L3807 & L3809



INDICATIONS FOR USE:

- Thumb Injuries
- CMC Injuries
- Metacarpal Injuries




P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Universal sizing for stocking optimization
2. Refined design with single-piece construction
3. Anatomically designed for thumb positioning in abduction and opposition
4. Malleable and removable 7 stay wrist-thumb system for a supportive, custom fit
5. Non-neoprene hypoallergenic material and soft edges for user comfort

PRODUCT ORDERING

Product Name	Ref No.
Universal Thumb Orthosis - Short	1147-XXY
Universal Thumb Orthosis - Long	1148-XXY

PRODUCT SIZING

Measure from center wrist at narrowest

Size	XX	Min	Max	Min	Max
UNI	00	6.0"	11.0"	15cm	28cm

Sizing Note: For Y specify L for left or R for right

R&ID# *285557*
REV A 3/22

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

