

# SIDEWINDER WRIST ORTHOSIS

**ezyWRAP**® designed the Sidewinder Wrist Orthosis to be soft, breathable, and ready directly off-the-shelf. It is made with soft and breathable materials to encourage patient compliance. Our Sidewinder includes a unique side lacing system allowing the user to customize the orthosis for compression and stabilization of their forearm anatomy. This brace has soft edges and a comfortable stockinette liner to promote patient compliance and long-term wearability. Available in both short and long for additional comfort.

Approved **PDAC** L3908



## INDICATIONS FOR USE:

- Carpal Tunnel Syndrome
- Wrist Sprains and Strains
- Wrist Immobilization
- Control of Post-Op Effusion





**P:** (800) 234-9004  
**F:** (800) 274-9006  
**E:** SALES@EZYWRAP.COM  
**W:** EZYWRAP.COM

FOR THE USE OF:  
CODE:

MADE IN THE USA   

# FEATURES & BENEFITS



Patent #10,045,876

1. Soft and breathable materials for comfort and compliance
2. Lacing allows customization for varying forearm circumference
3. Three malleable and removable stays for a supportive custom fit
4. Non-neoprene hypoallergenic material and soft edges for user comfort

## PRODUCT ORDERING

Product Name	Model No.
Sidewinder Wrist Orthosis - Short	1152-XXY
Sidewinder Wrist Orthosis - Long	1153-XXY

## PRODUCT SIZING

Measure the circumference of wrist at narrowest

Size	XX	Min	Max	Min	Max
UNI	00	4.0"	8.0"	10cm	20cm

**Sizing Note:** For Y specify L for left or R for right

R&ID# \*285459\*  
REV A 6/23

FOR THE USE OF:  
CODE:

MADE IN THE USA   

**FOR ORDERING:**

**P:** (800) 234-9004

**F:** (800) 274-9006

**E:** SALES@EZYWRAP.COM

**W:** EZYWRAP.COM

*Scan for application  
and accessories.*

