

TRI-PANEL KNEE IMMOBILIZER

EZYWRAP designed the Tri-Panel Knee Immobilizer to be the best in class knee orthosis. Our Tri-Panel Knee Immobilizer is adjustable and suitable for a wide range of anatomies with moldable removable stays allowing for tailored care. All while being easy to apply and comfortable to wear directly off-the-shelf.

Our knee immobilizer is FirstIce® Cold Therapy ready with patented Direct Stick™ Packs as an optional accessory. Our patented Cold Therapy provides up to 4 hours of continuous cold to help patients mitigate pain and **recover quicker**.

Approved **PDAC** L1830



INDICATIONS FOR USE:

- Post-operative rehabilitation
- Knee ligament injuries

*|FirstIce® FEATURES:

- Patented technology that lasts five times longer than gel.
- Direct Stick™ packs eliminate the need for pockets.
- Non-toxic, soft and malleable to conform to the body.



FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM


W: EZYWRAP.COM

FOR THE USE OF:

CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Universal dexterity and sizing for optimized inventory.
2. Tri-Panel suits a wide range of anatomies.
3. Adjustable, removable stays.
4. Easy Application and Cold Therapy Ready.

Note: Foam model is not cold therapy ready

PRODUCT ORDERING

Product Name	Ref No.
Tri-Panel Knee Immobilizer (Foam)	1278-00YY
Tri-Panel Knee Immobilizer	1279-00YY

OPTIONAL COLD THERAPY ACCESSORIES

Product Name	Size	Ref No.
FirstIce® Direct Stick™ Cold Therapy Pack	06 x 07"	100-0607
	06 x 10"	100-0610
	06 x 12"	100-0612
Patella Strap	Universal	120-00

Note: Foam model is not cold therapy ready

PRODUCT SIZING

Measure leg circumference 6" above mid patella

Length	YY	Min	Max
12" 30cm	12	15" 43cm	20" 51cm
14" 36cm	14	17" 43cm	24" 66cm
16" 41cm	16	18" 46cm	28" 71cm
18" 46cm	18	18" 46cm	28" 71cm
20" 51cm	20	19" 48cm	28" 71cm
22" 56cm	22	19" 48cm	28" 71cm
24" 61cm	24	20" 51cm	28" 71cm
26" 66cm	26	20" 51cm	28" 71cm

Note: Specify YY for length

R&ID# *282619*
REV A 4/21

FOR THE USE OF:
CODE:

MADE IN THE USA



FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

