

SIZED KNEE IMMOBILIZER

EZYWRAP designed the Sized Knee Immobilizer to be the best in class knee orthosis. Our Sized Knee Immobilizer is adjustable and suitable for a wide range of anatomies with moldable stays allowing for tailored care. All while being easy to apply and comfortable to wear directly off-the-shelf.

Our knee immobilizer is FirstIce® Cold Therapy ready with patented Direct Stick™ Packs as an optional accessory. Our patented Cold Therapy provides up to 4 hours of continuous cold to help patients mitigate pain and *recover quicker*.

Approved **PDAC** L1830



INDICATIONS FOR USE:

- Post-operative rehabilitation
- Knee ligament injuries

*FirstIce® FEATURES:

- Patented technology that lasts five times longer than gel.
- Direct Stick™ packs eliminate the need for pockets.
- Non-toxic, soft and malleable to conform to the body.



FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM


W: EZYWRAP.COM

FOR THE USE OF:

CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Specifically sized user support and comfort.
2. Universal dexterity for optimized inventory.
3. Moldable stays.
4. Easy Application and Cold Therapy Ready.

Note: Foam model is not cold therapy ready.

PRODUCT ORDERING

Product Name	Ref No.
Sized Knee Immobilizer (Foam)	1288-XXYY
Sized Knee Immobilizer	1289-XXYY

OPTIONAL COLD THERAPY ACCESSORIES

Product Name	Size	Ref No.
FirstIce® Direct Stick™ Cold Therapy Pack	06 x 07"	100-0607
	06 x 10"	100-0610
	06 x 12"	100-0612
Patella Strap	Universal	120-00

Note: Foam model is not cold therapy ready.

PRODUCT SIZING

Measure leg circumference 6" above mid patella

Size	XX	Min	Max	YY	Length
PED	05	12" 30cm	14" 36cm	10*	10" 25cm
XS	10	16" 26cm	18" 41cm	12*	12" 30cm
S	20	18" 41cm	20" 46cm	14	14" 36cm
M	30	20" 46cm	22" 51cm	16	16" 41cm
L	40	22" 51cm	23" 55cm	18	18" 46cm
XL	50	23" 55cm	26" 58cm	20	20" 51cm
XXL	60	26" 58cm	28" 65cm	22	22" 56cm
XXXL	70	28" 71cm	30" 79cm	24	24" 61cm
				26	26" 66cm

Sizing Note: Specify XX for size and YY for length; Only in Pediatric*

R&ID# *282623*
REV A 9/22

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

