

ACCOMMODATOR

eZYWRAP designed the Accommodator to be the best-in-class shoulder and arm immobilization solution for physicians and patients and to be useable directly off-the-shelf. The Accommodator is a self-adjusting shoulder orthosis and results in an exceedingly comfortable and secure fit.

Our straps are custom manufactured using orthopedic foam with removable butterfly closures for strap length customization. We've taken the time to engineer the Velcro® to be strong in shear and easily removable in peel, allowing for strength without struggle.

Model No. 2307 Approved **PDAC** L3670



* Immobilizer shown

INDICATIONS FOR USE:

- Post operative support
- Positioning / Injuries requiring immobilization




P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Posterior gusset allowing for sling to seat on the arm and elbow.
2. Foam straps with subtle sewn seam design for patient comfort.
3. Butterfly closures straps for adjustment.
4. Soft thumb loop eliminates slippage for patient comfort.
5. Velcro® optimized to be strong in shear and ease in peel for strength without struggle.

PRODUCT ORDERING

Product Name	Ref No.
Accommodator Arm Sling	2306-XX
Accommodator Shoulder Immobilizer	2307-XX

PRODUCT SIZING

Measure length of lateral epicondyle to knuckle

Size	XX	Min	Max	Min	Max
XS	10	8"	10"	20cm	25cm
S	20	10"	12"	25cm	30cm
M	30	12"	14"	30cm	35cm
L	40	14"	16"	35cm	40cm
XL	50	16"	18"	40cm	46cm

R&ID# *282549*
REV A 4/21

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

