

TAPERED SHOULDER IMMOBILIZER AND ARM SLING

EZYWRAP designed the Tapered Shoulder Immobilizer to be a comfortable and dependable off-the-shelf solution for physicians and patients. This shoulder immobilizer has a tapered arm envelope for patient comfort and a snug fit.

Our straps are custom manufactured using orthopedic foam with removable butterfly closures for customizable strap length. We've taken the time to engineer the Velcro® to be strong in shear and easily removable in peel, allowing for strength without struggle.

Model No. 2351 Approved **PDAC** **L3670**



* Immobilizer shown

INDICATIONS FOR USE:

- Post operative support
- Injuries requiring immobilization




P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



* Immobilizer shown

1. Foam straps with subtle sewn seam design for patient comfort.
2. Butterfly closures straps for adjustment.
3. Soft thumb loop eliminates slippage for patient comfort.
4. Velcro® optimized to be strong in shear and ease in peel for strength without struggle.

PRODUCT ORDERING

Product Name	Ref No.
Tapered Arm Sling	2349-XX
Tapered Shoulder Immobilizer	2351-XX

PRODUCT SIZING

Measure length of lateral epicondyle to knuckle

Size	XX	Min	Max	Min	Max
XS	10	8"	10"	20cm	25cm
S	20	10"	12"	25cm	30cm
M	30	12"	14"	30cm	35cm
L	40	14"	16"	35cm	40cm
XL	50	16"	18"	40cm	46cm

R&ID# *282628*
REV A 4/21

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

