

TORSO SHOULDER IMMOBILIZER

EZYWRAP® designed the Torso Shoulder Immobilizer to be a comfortable and dependable off-the-shelf solution for physicians and patients. The Torso Shoulder Immobilizer eliminates shoulder straps and immobilizes the shoulder with a torso band and adjustable bicep cuff that can be fixed anywhere on the device.

The torso band and wrist cuff are manufactured using orthopedic foam with a hook-and-loop design engineered to be strong in shear and removable in peel for strength without struggle.

Approved PDAC L3660



INDICATIONS FOR USE:

- Post-operative rehabilitation of the shoulder
- Injuries requiring immobilization




P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM


FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Torso band with elastic insert eliminates shoulder straps.
2. Hook-and-loop cuff wide design supports wrist and hand.
3. Universal sizing and dexterity.
4. Designed for ease of application.



PRODUCT ORDERING

Product Name	Ref No.
Torso Shoulder Immobilizer	2355-XX

PRODUCT SIZING

Measure circumference of chest including affected arm

Size	XX	Max	
PED	05	32"	81cm
S/M	25	44"	112cm
XL/L	45	53"	137cm
UNI	00	53"	137cm

R&ID# *282648*
REV A 6/21

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

