

CLAVICLE SUPPORT

EZYWRAP designed the Clavicle Support to be perfect for patients with clavicle fractures. It is made of comfortable padded material with wide straps for ease of applications and patient comfort. Hook and loop closure enables effortless adjustment, while a movable D-Ring ensures a comfortable custom fit.

Cold Therapy ready with our patented FirstIce® Direct Stick™ technology as an optional accessory. Our specially formulated premium Cold Therapy line provides up to 4 hours of continuous cold.

Approved **PDAC** L3650



INDICATIONS FOR USE:

- Clavicle fractures

*|FirstIce® FEATURES:

- Patented technology that lasts five times longer than gel.
- Direct Stick™ packs eliminate the need for pockets.
- Non-toxic, soft and malleable to conform to the body.



© 2021 PROFESSIONAL PRODUCTS INC.

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM


W: EZYWRAP.COM

FOR THE USE OF:

CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Comfortable breathable material with wide straps for proper immobilization.
2. Adjustable hook and loop closure for sizing on left and right sides.
3. Customizable vertical D-Ring buckle for comfort and fit.
4. Optional Cold Therapy using patented Direct Stick™ for integrated pain relief.

PRODUCT ORDERING

Product Name	Ref No.
Clavicle Support	2415-XX

OPTIONAL COLD THERAPY ACCESSORIES

Product Name	Size	Ref No.
FirstIce® Direct Stick™ Cold Therapy Pack	6 x 7"	100-0607

PRODUCT SIZING

Measure from center of shoulder blades around shoulder, back to center

Size	XX	Min	Max	Min	Max
XS	10	15"	38"	21cm	53cm
S	20	19"	48"	28cm	71cm
M	30	24"	61"	34cm	86cm
L	40	30"	76"	38cm	97cm
XL	50	34"	86"	40cm	102cm

R&ID# *282464*
REV A 4/21

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

