



DONUT PATELLA STABILIZER



EZYWRAP designed our Donut Patella Stabilizer to be exceedingly comfortable and of remarkably high quality. The silicone donut patellar pad is hook-and-loop compatible and repositionable to customize for a patient's anatomy. We use a breathable, hypoallergenic non-neoprene material that keeps the patient cool during exercise and comfortable with anti-migration elements to minimize slippage during extended use.



INDICATIONS FOR USE:

- Chondromalacia
- Patellofemoral Malalignment
- Acute Injuries (Patella Dislocation and Patellar Tendon Rupture)
- Chronic Conditions (Patellar Tendinopathy, Patellar Instability, and Patellar Tracking Disorder)




© 2023 PROFESSIONAL PRODUCTS INC.

P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA   

FEATURES & BENEFITS



1. Repositionable hook-and-loop compatible silicone donut patellar buttress for customizable support
2. Soft and breathable materials for compression and comfort
3. Anti-migration elements to ensure effective support
4. Non-neoprene hypoallergenic material with soft edges for user comfort

PRODUCT SIZING

Measure circumference of leg 6" above mid-patella with leg fully extended

Size	XX	Min	Max
XS	10	13" 33cm	16" 41cm
S	20	16" 41cm	19" 48cm
M	30	19" 48cm	21" 53cm
L	40	21" 53cm	24" 61cm
XL	50	24" 61cm	27" 69cm
XXL	60	27" 69cm	30" 76cm
XXXL	70	30" 76cm	33" 84cm

PRODUCT ORDERING

Product Name	Ref No.
Donut Patella Stabilizer	470-XX

R&ID# *286999*
REV A 3/23

FOR THE USE OF:
CODE:

MADE IN THE USA   

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZWRAP.COM

W: EZYWRAP.COM

Scan for application
and accessories.

