

# DUAL AXIS HINGED KNEE

(PULL-UP)

**EZYWRAP**® designed the Dual Axis Hinged Knee to be fully adjustable, exceedingly comfortable, and highly durable for active individuals. We use a breathable hypoallergenic non-neoprene material that keeps the patient cool during exercise and comfortable with anti-migration elements to minimize slippage during extended use.

Our Knee Orthosis has hidden, customizable dual-axis hinges that are bendable and shapable for customization by an individual with expertise. Additionally, the hinges are removable for ease of maintenance and patient cleaning. All while being designed to be Pull-Up for active applications.

Approved **PDAC L1810 & L1812**



## INDICATIONS FOR USE:

- Knee Ligament Injuries
- Meniscus Injuries
- Knee Sprains and Strains




**P:** (800) 234-9004  
**F:** (800) 274-9006  
**E:** SALES@EZYWRAP.COM  
**W:** EZYWRAP.COM

FOR THE USE OF:  
CODE:

MADE IN THE USA     

# FEATURES & BENEFITS

503 Depicted



1. Shapable Dual Axis hinges for customization
2. Internal hinges for protection during activity
3. Removable hinges for ease of maintenance and cleaning
4. Open popliteal and removable straps
5. Anti-migration elements to ensure effective support
6. Pull-Up design for active applications
7. Soft and breathable non-neoprene hypoallergenic material

## PRODUCT SIZING

Measure circumference of leg 6" above mid-patella with leg fully extended

Size	XX	Min		Max	
XS	10	12"	31cm	16"	41cm
S	20	15"	39cm	18"	46cm
M	30	17"	44cm	20"	51cm
L	40	20"	51cm	23"	59cm
XL	50	23"	59cm	26"	66cm
XXL	60	26"	66cm	29"	74cm
XXXL	70	29"	74cm	32"	81cm

## PRODUCT ORDERING

Product Name	Ref No.
Dual Axis Hinged Knee (Pull-Up) – Short	503-XX
Dual Axis Hinged Knee (Pull-Up) – Long	504-XX

R&ID# \*284181\*  
REV A 6/21

FOR THE USE OF:  
CODE:

MADE IN THE USA



FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

Scan for application  
and accessories.

