

# RALLY HINGED KNEE

(PULL-UP)

**ezyWRAP**® designed the Rally Hinged Knee Orthosis to be a fully adjustable hinged knee orthosis that is both exceedingly comfortable and remarkably high quality. We use a breathable non-neoprene that keeps the patient cool and comfortable with anti-migration elements to ensure effective support.

The Rally has removable dual-axis hinges for ease of maintenance and patient cleaning. Our hinges are customizable with range-of-motion stops allowing for positioning in increments of 10 degrees by an individual with expertise or directly off-the-shelf. All around a Pull-Up design for active applications.

Approved **PDAC** L1832 & L1833



## INDICATIONS FOR USE:

- Post-operative Rehabilitation
- Meniscus Injuries
- Knee Ligament Injuries
- Knee Sprains and Strains




**P:** (800) 234-9004  
**F:** (800) 274-9006  
**E:** SALES@EZYWRAP.COM  
**W:** EZYWRAP.COM

FOR THE USE OF:  
CODE:

MADE IN THE USA     

# FEATURES & BENEFITS

513 Depicted



1. Soft and breathable materials for compression and comfort.
2. Dual Axis ROM with stops for custom hinge settings.
3. Removable hinges for ease of maintenance and cleaning.
4. Anti-migration elements to ensure effective support.
5. Pull-Up design for active applications.
6. Non-neoprene hypoallergenic material.

## PRODUCT SIZING

Measure circumference of leg 6" above mid-patella with leg fully extended

Size	XX	Min	Max
XS	10	12" 31cm	16" 41cm
S	20	15" 39cm	18" 46cm
M	30	17" 44cm	20" 51cm
L	40	20" 51cm	23" 59cm
XL	50	23" 59cm	26" 66cm
XXL	60	26" 66cm	29" 74cm
XXXL	70	29" 74cm	31" 81cm

## PRODUCT ORDERING

Product Name	Ref No.
Rally Hinged Knee (Pull-Up) – Short	513-XX
Rally Hinged Knee (Pull-Up) – Long	514-XX

R&ID# \*283038\*  
REV A 4/21

FOR THE USE OF:  
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZWRAP.COM

W: EZYWRAP.COM

Scan for application  
and accessories.

