



UNIVERSAL KNEE THERAPY SYSTEM



ezyWRAP® designed our Universal Knee Wraps as the premiere hot or cold therapy knee system. This unique, universal single-piece compressive wrap with three separate straps offer maximal mobility and fit in a safe, hypo-allergenic wrap that includes a comfortable ultra-soft interface.

Our patented First Ice® and First Gel® packs allow for accurate, targeted, therapeutic placement of our patented Direct Stick™ packs on our convenient universal wrap. As a result, our joint-targeted knee therapy system is simple to apply, reapply and rotate between user and freezer without any special considerations.



7650 - Universal Knee Compression Wrap Pictured Above

INDICATIONS FOR USE:

- Post-Operative pain and inflammation
- Post-Injury pain and inflammation
- Rehabilitation pain and inflammation
- Acute knee injuries (muscle sprains, strains, contusions, ligament tear, meniscus tear)
- Chronic knee conditions (arthritis, tendonitis, bursitis, instability, ligament tightness)



© 2023 PROFESSIONAL PRODUCTS INC.


P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA   

FEATURES & BENEFITS

*7650 Depicted Below



Patent No. 9,039,747 - and Other Patents Pending

1. First Ice® and First Gel® Direct Stick™ Hot or Cold Therapy application directly to the wrap
2. Three independent compressive straps for a tailored fit
3. Single-piece construction with quick assist closure tabs
4. Universal sizing and dexterity, fits left and right to optimize inventory
5. Optional extension strap for a wide variety of anatomies
6. Non-neoprene hypoallergenic material

PRODUCT SIZING

Measure circumference of leg 6" above mid-patella with leg fully extended

Size	Ref No.	Max	Max
Universal Knee Wrap	7630-XXX	25"	63cm
Universal Knee Compression Wrap	7650-XXX	25"	63cm
Extender Strap	152-00-01	32"	81cm
Bundle Options	Pack Size	Packs	XXX
Wrap Only	-	0	000
First Ice® Bundle (Wrap, Packs, & Tote-Bag)	6" x 12"	4	102
First Gel® Bundle (Wrap, Packs, & Tote-Bag)	6" x 12"	4	202

Note: All options are available in single and 6-pack case quantities, wraps only in 24 packs

FIRST ICE® & FIRST GEL® (PACK ONLY)

Note: DFP decreases pack freeze time and facilitates freezer storage

Product Name	Pack Size	Storage	Packs	Ref No.
First Ice® Direct Stick™ Cold Therapy Packs	6" x 12"	Zip-Lock	2	130-0612
		Direct Freeze Package (DFP)	4	144-0612
First Gel® Direct Stick™ Hot / Cold Therapy Packs	6" x 12"	Zip-Lock	4	230-0612

R&ID# *286727*
REV A 5/23

FOR THE USE OF:
CODE:

MADE IN THE USA



FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZYWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

