

# SPEED LACER ANKLE ORTHOSIS

**EZYWRAP** designed the Speed Lacer Ankle Orthosis to maximize ankle control with a focus on speed of application. Our figure 8 design with calcaneal heel-lock enables maximum inversion and eversion control. Laces allow the brace to open fully for any anatomy. Additionally, we have included the ability to add optional stays for desired stability.

Our brace is breathable, durable, and hook-and-loop compatible without painful underfoot seams. The hook-and-loop is engineered to be strong in shear and easily removable in peel, allowing for strength without struggle.

Approved **PDAC** **L1902**



## INDICATIONS FOR USE:

- Ankle sprains
- Strains or post-injury
- Post-operatively




**P:** (800) 234-9004  
**F:** (800) 274-9006  
**E:** SALES@EZYWRAP.COM  
**W:** EZYWRAP.COM

FOR THE USE OF:  
CODE:

MADE IN THE USA   

# FEATURES & BENEFITS



1. Figure 8 design and calcaneal heel-locking for maximum inversion and eversion control.
2. Constructed of breathable, durable and hook-and-loop compatible material.
3. Maximum comfort from the seamless design, micro thin webbing, and achilles pad.
4. Fully adjustable speed tabs allow for independently customizable tension on all laces.
5. Easy slip on design for ease of application.
6. Hypoallergenic material and soft edges for user comfort.

## PRODUCT ORDERING

Product Name	Ref No.
Speed Lacer Ankle Orthosis	803-XX

## OPTIONAL ACCESSORIES

Product Name	Size	Ref No.
Removable Plastic Stays	Universal	150-00

## PRODUCT SIZING

Measure using shoe size

Size	XX	US Men	US Women	EU Men	EU Women
XXS	05	5-6	6.5-7.5	37-40	37-38
XS	10	6-7.5	7.5-9.0	40-41	38-40
S	20	7.5-9	9.0-10.5	41-43	40-41
M	30	9-11	10.5-12.5	43-45	41-43
L	40	11-13	12.5-14.5	45-47	43-45
XL	50	13-15	14.5-16.5	47-49	45-47
XXL	60	15+	16.5+	49+	47

R&ID# \*282556\*  
REV A 5/23

FOR THE USE OF:  
CODE:

MADE IN THE USA



FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZWRAP.COM

W: EZYWRAP.COM

Scan for application  
and accessories.

