

WAHOO THUMB ORTHOSIS

eZYWRAP® designed the Wahoo Thumb Orthosis to be a comfortable hand finger orthosis ready off-the-shelf or configurable by an individual with expertise. The 3-stay system can be configured and customized to fit the individual patient and provide tailored support. Our design is semi-universal for inventory optimization and single strapped for easy one-handed application.

Approved **PDAC** L3923 & 3924



INDICATIONS FOR USE:

- CMC Stabilization
- Thumb Injuries
- Arthritis




© 2021 PROFESSIONAL PRODUCTS INC.

P: (800) 234-9004
F: (800) 274-9006
E: SALES@EZYWRAP.COM
W: EZYWRAP.COM

FOR THE USE OF:
CODE:

MADE IN THE USA     

FEATURES & BENEFITS



1. Soft and breathable materials for compression and comfort.
2. Stays are customizable, removable, stackable, and bendable.
3. Dual strap design for easy application.
4. Universal sizing and dexterity to reduce inventory.
5. Non-neoprene hypoallergenic material.

PRODUCT ORDERING

Product Name	Ref No.
Wahoo Thumb Orthosis	939-XX

PRODUCT SIZING

Measure circumference of thumb at knuckle and wrist at narrowest

Size	XX	Thumb at Knuckle	Wrist at Narrowest
XS	10	2.1 – 2.5" 5 – 6cm	4.5 – 5.0" 11 – 13cm
REG	35	2.5 – 3.4" 6 – 8cm	5.0 – 8.5" 13 – 22cm
XL	50	3.4 – 3.6" 8 – 9cm	5.0 – 8.5" 13 – 22cm

R&ID# *283036*
REV A 4/21

FOR THE USE OF:
CODE:

MADE IN THE USA     

FOR ORDERING:

P: (800) 234-9004

F: (800) 274-9006

E: SALES@EZWRAP.COM

W: EZYWRAP.COM

*Scan for application
and accessories.*

