


CORE ITEM



Surgical Mask

3-Ply Face Mask - Level 3



SPECIFICATIONS

CLIN	600, 6000, 6001, 6002		
Manufacturer	Henan Huibo Medical Co.		
Model Number	SM-L3	SKU	VZ-HU-FM3
Protection Class	Level 3 (ASTM F2100-19)	BFE	99%
Performance	<ul style="list-style-type: none"> - Sub-micron particulates filtration efficient at 0.1 micron: $\geq 98\%$ - Differential Pressure: ≤ 4.1 mmH₂O/cm² - Bacterial Filtration Efficiency (BFE): $\geq 99\%$ - Resistance to Penetration by Synthetic Blood: No penetration at 80, 120 & 160 mmHg pressure - Flammability: Class 1 		
Material	3-Ply: 2 Non-Woven Layers w/ Hypoallergenic layer inside; Melt Blown Filter in the Middle		
Size	Blue	Size	Medium/Large (Universal)
Packaging	100 masks/box; 20 boxes/case		



PRICING

CLIN #	Quantity	Unit Price
600	1-1,000,000	\$0.07
6000	1,000,001-5,000,000	\$0.07
6001	5,000,001-10,000,000	\$0.07
6002	10,000,001-15,000,000	\$0.07