

TECHNICAL DATA RETRACTORS

SAGA DENTAL SUPPLY AS

PART NR	PRODUCT NAME	SIZE	MATERIAL	COLOUR	STERILIZING
05079	TERIOR - POSTERIOR RETRACTO	-	PC	CLEAR	Autoclav 121C
05049	COMBEE RETR.	LARGE	PC	CLEAR	Autoclav 121C
05048	COMBEE RETR.	SMALL	PC	CLEAR	Autoclav 121C
05047	K&K OCCL. RETRACTOR	LARGE	PC	CLEAR	Autoclav 121C
05046	K&K OCCL. RETRACTOR	SMALL	PC	CLEAR	Autoclav 121C
05059	LC - RETR	MEDIUM	PC	CLEAR	Autoclav 121C
05060	LC - RETR	MEDIUM	PC	ASS. COLOUR	Cold sterilize
05071	LC - RETR	MEDIUM	PP	ASS. COLOUR	Cold sterilize
05070	LC - RETR	MEDIUM	PP	OPAQUE	Autoclav 121C
05064	LIP RING	LARGE	PE	OPAQUE	Cold sterilize
05065	LIP RING	LARGE	PE	ASS. COLOUR	Cold sterilize
05082	LIP RING	SMALL	PE	OPAQUE	Cold sterilize
05052	MAXEE RETR.	LARGE	PC	CLEAR	Autoclav 121C
05069	MAXEE RETR.	LARGE	PP	OPAQUE	Autoclav 121C
05051	MAXEE RETR.	MEDIUM	PC	CLEAR	Autoclav 121C
05050	MAXEE RETR.	SMALL	PC	CLEAR	Autoclav 121C
05068	MAXEE RETR.	SMALL	PP	OPAQUE	Autoclav 121C
05097	OCCL. 45 ANGLED FORK RETR.	SMALL	PC	CLEAR/OPAQUE	Autoclav 121C
05098	OCCL. 45 ANGLED FORK RETR.	LARGE	PC	CLEAR/OPAQUE	Autoclav 121C
05073	OCCLUSAL RETRACTOR	LONG	PC	CLEAR	Autoclav 121C
05072	OCCLUSAL RETRACTOR	SHORT	PC	CLEAR	Autoclav 121C
05061	POSTERIOR RETR.	-	PC	CLEAR	Autoclav 121C
05058	SCANDEE	LARGE	PC	ASS. COLOUR	Cold sterilize
05055	SCANDEE	LARGE	PC	CLEAR	Autoclav 121C
05054	SCANDEE	MEDIUM	PC	CLEAR	Autoclav 121C
05057	SCANDEE	MEDIUM	PC	ASS. COLOUR	Cold sterilize
05067	SCANDEE	MEDIUM	PP	ASS. COLOUR	Cold sterilize
05066	SCANDEE	MEDIUM	PP	OPAQUE	Autoclav 121C
05053	SCANDEE	SMALL	PC	CLEAR	Autoclav 121C
05056	SCANDEE	SMALL	PC	ASS. COLOUR	Cold sterilize
05081	SINGLE SIDED	LARGE	PC	CLEAR	Autoclav 121C
05080	SINGLE SIDED	SMALL	PC	CLEAR	Autoclav 121C
05063	TCL - RETR.	LARGE	PC	CLEAR	Autoclav 121C
05075	TCL - RETR.	LARGE	PP	ASS. COLOUR	Cold sterilize
05062	TCL - RETR.	SMALL	PC	CLEAR	Autoclav 121C
05074	TCL - RETR.	SMALL	PP	ASS. COLOUR	Cold sterilize