# NEW Android

### **SMART OPERATION SYSTEM**




# DISPLAY





### **DISPLAY**



## **PARAMETER**

Main machine					
Electrical power	2400W				
Energy supply	100-240V/50-60HZ				
Wavelength	808/808nm+755nm+1064nm				
Laser	Coherent bars (100w each, peak power 120W, uniform laser)				
Cooling level	5 cooling level				
Cooling system	Water cooling+TEC 450W cooling(exclusive)+Wind cooling+sapphire				
Water tank cooling temperature	MAX 35 °C, above 35, the device will stop working for safety				
Screen size of the device	12 inch android control system (screen resolution 1024*768)				
Screen size of the handle	2.4inch				
Machine memory	8 G				
Functions	Hair removal +Skin rejuvenation				
Warning system	Tank temperature; Handle crystal temperature; Water flow rate				
Language choice	English+Spanish+German+French(Could change base on markets)				
Machine size (L*W*H)	50.2cm×39cm×110cm				
Package size (L*W*H)	65cm×49cm×120cm				
NG/GW	70KG/97KG				
Shield material	Aviation aluminum box				
Device warranty	Whole machine 2 year+Life-time after sale service.				
Handle shooting warranty	20 millions				

### **PARAMETER**





Different Laser Power Options						
Laser power	500W (5barsX100W)	600W (6barsX100W)	800W (8barsX100W)	1000W (10barsX100W)	1200W (12barsX100W)	
Spot size	15*15mm	15*20mm	15*30mm	15*40mm	15*40mm	
Max.energy	130J/CM <sup>2</sup>	160J/CM <sup>2</sup>	190J/CM <sup>2</sup>	270J/CM <sup>2</sup>	310J/CM <sup>2</sup>	
Frequency	1~10HZ	1~10HZ	1~10HZ	1~6HZ	1~6HZ	
Energy density	10~140J/CM <sup>2</sup>	10~140J/CM <sup>2</sup>	10~140J/CM <sup>2</sup>	10~60J/CM <sup>2</sup>	10~60J/CM <sup>2</sup>	
Duration	10~300MS	10~300MS	10~300MS	10~200MS	10~200MS	

### **ADVANTAGE**


#### Intelligent customer management system

It can update and down load clients' recording, so data can be Shared across multiple salons.







Also can make a appointment on your next treatment.





It can update and down load clients' recording, so data can be Shared across multiple salons.





5G Instant after-sales Service System

Save machine operation +clinical+sample problem video, and in the future, 5G will be used to identify and solve machine problems online.





### **ADVANTAGE**

#### Android control system

Energy adjust, skin type and body part can switch smooth.





#### Happy Depilator

Can listen music when do hair removal treatment, enjoy a happy experience.





#### Intelligent customer analysis system


The skin analyzer can analyze and compare the results before and after hair removal to help analyze the skin type and maintenance after hair removal.





#### Uniform laser

Uniform spot more concentrated light, no scattering, stronger energy, hair removal effective. Narrow wave, no sensory depilation.





### **ADVANTAGE**

#### 2.4 inch touch screen

The energy can be adjusted easily and, the temperature of hand can be displayed instantly, the control is simpler, the customer feels more intimate.







50



One of the first published articles describing *laser hair removal* was authored by the group at Massachusetts General Hospital in 1998.

It has been proven to be the safety & efficacy revolution for hair removal treatment


People have more ways to get rid of the unwanted hair.

#### LASER THEORY

**Blend** (808nm+755nm+1064nm) lasers — The ultimate technology for dark and tanned skin. 1 SHOT=3 WAVELENGTHS, safe& effective and treatments can be performed on all patients throughout the year.

The 3 main anatomical targets include the Bulge, Bulb and Papilla.











### Package Inspection and preparation



























0000







