




# RHINOLOGY

NASAL SPLINTS, NASAL PACKING,  
& MORE




# THERMOPLAST

## N A S A L   S P L I N T S

### T H E R M O P L A S T   9 X 1 2 S H E E T S

Rigid, opaque, lightweight, low temperature thermoplastic splinting material. Simply soften with hot water and mold over taped nose. These nasal splints are available in large 9"x12" sheets for convenience to cut custom nasal splints to fit the specific anatomy of the patient.


- Stabilization of the nose after surgery
- Reduction of edema after nasal surgery or trauma



PART #	COLOR	THICKNESS
DM-1022-9X12	IVORY	1/16"
DM-1023-9X12	IVORY	3/32"

### T H E R M O P L A S T   3 X 3 S Q U A R E S

PART #	COLOR	THICKNESS	QTY
DM-1033-IV-3X3-5	IVORY	1/16"	5
DM-1033-IV-3X3-10	IVORY	1/16"	10
DM-1033-IV-3X3-5	IVORY	3/32"	5
DM-1033-IV-3X3-10	IVORY	3/32"	10






# THERMOPLAST

## N A S A L   S P L I N T S

### T H E R M O P L A S T   P R E C U T S P L I N T S



These nasal splints are precut trapezoidal shaped splints made from 1/16" thermoplast.

\*Available in blush or ivory

- Stabilization of the nose after surgery
- Reduction of edema after nasal surgery or trauma
- Dorsal pad reduces dead space and skin wrinkling
- Dorsal pad reduces possibility of polybeak

#### Kit (5/box) Components:

- Instruction sheet
- One alcohol wipe
- One skin prep
- One dorsal pad
- Pre-cut Protecto Strips
- One Thermoplastic Nasal Splint

#### Size Measurements

Petite (P) 58x30mm

Standard (S) 65x35mm

Large (L) 69x43mm

PART #	COLOR	QTY
DM-116-TRP-BLUSH-05-SIZE	BLUSH	5
DM-116-TRP-BLUSH-10-SIZE	BLUSH	10
DM-116-TRP-IVORY-05-SIZE	IVORY	5
DM-116-TRP-IVORY-10-SIZE	IVORY	10

#### Non-Kit (10/box) Components:

- Instruction sheet
- One Thermoplastic Nasal Splint
- One dorsal pad

PART #	COLOR	QTY
DM-116-TRP-BLUSH-05-SIZE	BLUSH	5
DM-116-TRP-BLUSH-10-SIZE	BLUSH	10
DM-116-TRP-IVORY-05-SIZE	IVORY	5
DM-116-TRP-IVORY-10-SIZE	IVORY	10



# DENVER ALUMINUM

## N A S A L   S P L I N T S

### 1 5 0 0   S E R I E S   S P L I N T S



Denver Aluminum Splint: Series 1500 is a two piece external nasal splint designed to provide stabilization after rhinoplasty. This nasal splint is made with soft aluminum and Velcro with an adhesive backing to ensure a 7-10 day wear post surgery. It is easily moldable and provides cushioned pressure to reduce swelling after surgery or trauma.

- Stabilization of the nose after surgery
- Reduction of edema after nasal surgery or trauma
- Dorsal pad reduces dead space and skin wrinkling
- Dorsal pad reduces possibility of polybeak

#### Size Measurements

Mini (M) 50x28mm

Petite (P) 58x30mm

Small/Medium (S) 65x35mm

Large (L) 69x43mm

#### Kit (5/box) Components:

- Instruction sheet, • One alcohol wipe
- One skin prep • One dorsal pad
- Pre-cut Protecto Strips • One Denver Nasal Splint

#### PART #

DM-10-1500-05KM

DM-10-1500-05KP

DM-10-1500-05KS

DM-10-1500-05KL

#### Non-Kit (10/box) Components:

- Instruction sheet • One Denver Nasal Splint
- One dorsal pad

#### PART #

DM-10-1500-05KM

DM-10-1500-05KP

DM-10-1500-05KS

DM-10-1500-05KL



# DENVER ALUMINUM

## N A S A L   S P L I N T S

### 2 0 0 0   S E R I E S   S P L I N T S



Denver Aluminum Splint: Series 2000 is a one piece external nasal splint designed to provide stabilization after rhinoplasty. This nasal splint is made with soft aluminum and 1/8" ventilating foam with an adhesive backing to ensure a 7-10 day wear post surgery. It is easily moldable aluminum and provides cushioned pressure to reduce swelling after surgery or trauma.

- Stabilization of the nose after surgery
- Reduction of edema after nasal surgery or trauma
- Dorsal pad reduces dead space and skin wrinkling
- Dorsal pad reduces possibility of polybeak

#### Size Measurements

- Mini (M) 50x28mm
- Petite (P) 58x30mm
- Small/Medium (S) 65x35mm
- Large (L) 69x43mm

#### Kit (5/box) Components:

- Instruction sheet, • One alcohol wipe
- One skin prep • One dorsal pad
- Pre-cut Protecto Strips • One Denver Nasal Splint

#### PART #

DM-10-2000-05KM  
DM-10-2000-05KP  
DM-10-2000-05KS  
DM-10-2000-05KL

#### Non-Kit (10/box) Components:

- Instruction sheet • One Denver Nasal Splint
- One dorsal pad

#### PART #

DM-10-2000-05KM  
DM-10-2000-05KP  
DM-10-2000-05KS  
DM-10-2000-05KL



# INTERNAL NASAL SPLINTS

## Internal Airway Nasal Splint



Commonly referred as the “Doyle Splint.”

This splint provides support to the septum while maintaining an open airway for natural air flow after surgery. It reduces the risk of the airway closing from accumulation of blood and mucus by allowing the lumen to remain open. Made of soft medical grade silicone to aid with insertion and removal and includes holes for suture

placement. Sterile and single use

DM-SP-79300	1 PAIR
DM-SP-79300-05	5 PAIRS
DM-SP-79300-50	50 PAIRS

## Bivalve Nasal Splint

The Bivalve Nasal Splint provides septal support following surgery. It features a lengthwise slit to facilitate insertion and removal, and suture holes for stabilization. Made of clear medical grade fluoroplastic that is flexible and non-adherent to reduce adhesions between the septum and lateral nasal wall.

Available in two sizes: standard or large and two thicknesses: 0.25 mm or 0.50 mm



STANDARD SIZE: 26.1 MM X 42.6 MM  
OVERSIZE: 39.11 MM X 52 MM

PART #	SIZE	THICKNESS	QTY
DM-SP-27000	STANDARD	THIN (.25MM)	6 PAIRS
DM-SP-27005	STANDARD	THICK (.5MM)	6 PAIRS
DM-SP-27010	OVERSIZED	THIN (.25MM)	6 PAIRS
DM-SP-27015	OVERSIZED	THICK (.5MM)	6 PAIRS



# INTERNAL NASAL SPLINTS

## Nasal Septal Button


One-piece design with 2 flange leaves for maximum security and cleanliness. Designed to reduce or eliminate many symptoms of nasal septal perforations. Easily trimmed to ensure proper fit. Medical-grade silicone. 6mm posts on all sizes.

Sterile, single-use. Qty 1 each

DM-SP-78100	SMALL	3 cm X 4 mm
DM-SP-78105	MEDIUM	5 cm X 4 mm
DM-SP-78107	LARGE	7 cm X 4 mm



# NASAL PACKING



## ABSORBENT™ NASAL PACKS

AbsorbENT nasal and sinus packs utilize three-dimensional interconnected open cells that are soft and pliable.

- Even and fast expansion
- Soft, lint and fiber free
- Rapid wicking properties
- Qty 10/box

PART NO.	SIZE	DIMENSIONS LXWXH (CM)
DM-RH-7401-10	SMALL	4.5 X 1.5 X 2 10
DM-RH-7408-10	SMALL WITH AIRWAY	4.5 X 1.5 X 2 10
DM-RH-7404-10	LARGE	8 X 1.5 X 2 10
DM-RH-7405-10	LARGE WITH AIRWAY	8 X 1.5 X 2 10
DM-RH-7406-10	LARGE - THIN	8 X 1 X 2 10
DM-RH-2310-10	EXTRA LARGE	10 X 1.5 X 2.5 10

