



Contact your local dealer for further information or call toll free 1.877.342.3368.

Miracle Suction

Multifunctional Mirror with Suction Function

A unique and innovative idea!

These mirrors feature multiple suction holes to free up your hands and save you time.

- Fewer instruments in the oral cavity
- Less obstruction for the dentists, assistants, and hygienists
- Increased level of comfort and quality of care for the patients
- Strong tongue retraction
- Anti-fog mirror providing clear view
- Assist with sealant placement
- 8 holes in the mirror head - rapid suction of liquid while preventing suction on cheek and tongue
- Color-coded and comfortable silicone handles for easy identification
- Fits into all suction unit on unit chairs
- Adaptor to fit the unit into HVE (High Volume Evacuation) handpiece to reduce infectious aerosols
- 100% stainless steel and fully autoclavable
- Cleaning brush for the hose to remove debris
- Effective especially in pedodontic treatment, scaling and tooth reduction

Regular Package (#901-001)

Contains:

- 5 autoclavable mirrors with detachable handle
- 2 hoses, 3 HVE adaptors, 2 o-rings and a cleaning brush

Mini Package (#901-002)

Contains:

- 2 autoclavable mirrors with detachable handle
- 1 hose, 1 HVE adaptor, 1 o-ring, and a cleaning brush



Miracle