

Atrium Express Mini 500 and Pneumostat

Mobile chest drains

Atrium Mobile chest drainage products from Getinge are lightweight and ergonomic for easy patient ambulation. The Express Mini 500 and Pneumostat feature a dry seal one-way valve for convenient set-up and maximum patient protection in any position.





Express Mini 500

The Express Mini 500 is a full-featured mobile chest drain. These small, portable drainage collection units with pre-set suction help patients become ambulatory with the convenient, easy-to-wear "get mobile" size.

The Express Mini 500 incorporates a 500 ml collection volume, and a patented automatic pre-set suction regulator at -20 cmH₂O. In addition, the Express Mini design includes a dry seal one-way valve in a lightweight, ergonomic configuration for patient protection.

Pneumostat


The Pneumostat Chest Drain Valve is a disposable, single-use chest drain valve with a small, graduated 30 ml collection capacity integrated with a novel one-way dry seal valve technology.

The Pneumostat Chest Drain Valve is easy-to-use, providing patients excellent protection with a low-cost solution.

Express Mini 500 Dry Seal Chest Drain



Pneumostat Chest Drain Valve



Express Mini 500 Dry Seal Chest Drains

Description	Case qty.	Part no.
Express Mini 500	6	16400
Express Suction Bulb	10	19930

Pneumostat Chest Drain Valve

Description	Case qty.	Part no.
Pneumostat Chest Drain Valve with 30 ml Collection	10	16100

Not made with natural rubber latex.

■ Atrium Express Mini 500 and Pneumostat are manufactured by Atrium Medical Corporation / 40 Continental Blvd.,

Merrimack, NH 03054 / Tel. (603) 880-1433 · Protected by the following international and U.S. patent(s): http://patents.maquet.

com. · Getinge, GETINGE ☆ and Atrium are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates in the United States or other countries · Getinge and Atrium are registered with the U.S. Patent and Trademark Office. ·

Copyright 2018 Atrium Medical Corp. · All rights not expressly granted are reserved. · Refer to Instructions for Use for current indications, warnings, contraindications, and precautions. · Printed in U.S.A. · 05/18 · PN 011071 Rev AA

www.getinge.com

